

Model	SC7132
Weight Capacity	400 lbs.
Maximum Speed <sup>2</sup>	Up to 9.3 mph
Ground Clearance <sup>4</sup>	3.8" at transaxle
Turning Radius <sup>4</sup>	69.3"
Overall Length <sup>4</sup>	57.5"
Overall Width <sup>4</sup>	29.7"
Seat-to-Ground Height	21.9" - 24.9"
Seat-to-Deck Height	17.9" - 20.9"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Front and rear (CTS Suspension)
Range Per Charge <sup>2,3,7</sup>	Up to 40 miles with upgradable dual lithium battery configuration
Total Weight <sup>6</sup>	With standard battery: 242.6 lbs. Without batteries: 220.6 lbs.
Heaviest Piece When Disassembled <sup>4</sup>	N/A (Unit does not disassemble)
Standard Seating <sup>4</sup>	Type: Limited reclining high back with sliders Weight: 54.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>5</sup>	25 volt, 50Ah lithium ion
Battery Charger	8 amp
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 3 years

### HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## INSERT YOURSELF: ANYWHERE



# PRIDE MOBILITY®

800-800-8586 • [pridemobility.com](http://pridemobility.com)



1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) +10%, -5%.
4. Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Includes standard seat weight.
6. Battery weight may vary based on battery manufacturer.
7. With dual battery configuration.

# PURSUIT<sup>2</sup>

SHOWN WITH OPTIONAL  
CANOPY UPGRADE



FDA Class II Medical Device<sup>1</sup>

## Engineered for Adventure



**40 MILE  
RANGE**



**LONGER  
LIFE**



**CRUISING  
SPEED**

# PURSUIT<sup>2</sup>

**NEW**

**OPTIONAL** hard top canopy with durable, powder-coated steel frame, polycarbonate windshield, panoramic rearview mirror, and rear mounting brackets

Memory foam seat with 6" slider for adjustable support

Durable plastic shrouds

Dual rearview mirrors

Tiller adjustment lever

USB charger

Enhanced LED lighting

Front & rear suspension

Pneumatic tires mounted on split rims for off-road thrills

FDA Class II Medical Device<sup>1</sup>

## ENGINEERED FOR ADVENTURE

Built to handle adventures and outdoor terrain, the Pursuit 2 proves the best adventures are yet to come. With a fast cruising speed of 9 mph, a per-charge range up to 40 miles and pneumatic tires that can tackle any road (or backroad), it puts thrill-seekers on the path to happiness, fast.

## FEATURES

- 400-lb. weight capacity
- Maximum speed of up to 9.3 mph
- Up to a 40-mile range with upgradable dual lithium battery configuration
- Durable plastic shrouds
- Memory foam seat with 6" slider for adjustable support and premium, all-day comfort
- Infinitely adjustable tiller with USB charger
- Enhanced LED lighting for clear visibility
- Pneumatic tires for varied terrains and surfaces
- 3-year warranty on lithium battery

## AVAILABLE COLORS

Black

Gray (Matte)

Orange

## ACCESSORIES\*\*

- Wishbone Crutch Holder
- 3" Single/Double Cane or Crutch Holder
- Cane Crutch Holder
- Cell Phone Holder
- Cup Holder\*
- Forearm Crutch Holder
- Lap Belt
- Military Patches
- Oxygen Tank Holder
- Rear Basket\*
- Saddlebag\*
- Safety Flag
- Storage Pod
- Walker Holder
- Weather Cover

\* Cup holder, rear basket and saddlebag can be ordered as a bundle.

\*\* Only one rear accessory can be used at a time.

User-friendly console

Enhanced LED lighting

Front and rear CTS suspension

• **NEW** Canopy and optional canopy accessories:

- Tool clips (pair)
- Jumbo rear basket
- Canopy weather cover