

# Engineered for Powerful Performance and Exceptional Style



The Pursuit® XL is an exciting blend of power and precision that features a 24-volt DC motor, full suspension, 13" low-profile tires and a hydraulic-sealed brake system. Standard luxury touches, such as an infinite adjustable wraparound delta tiller, high-back pillow top seat and LED lighting make it the perfect outdoor scooter.



Full Suspension



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

| [pridemobility.com](http://pridemobility.com) |





## Specifications 4-wheel Model # SC714

Weight Capacity	400 lbs.
Maximum Speed <sup>1</sup>	Up to 7.5 mph
Ground Clearance <sup>2</sup>	5.25" at mid-frame, 2.75" at motor
Turning Radius	82.5"
Overall Length <sup>2</sup>	56"
Overall Width <sup>2</sup>	27.25"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Front and rear
Range Per Charge <sup>1</sup>	Group 24 (70 AH): Up to 22 miles Group 27 (100 AH): Up to 30 miles
Total Weight Without Batteries	254 lbs.
Heaviest Piece When Disassembled	N/A
Standard Seating	Type: High-back pillow top with sliders and headrest Weight: 49 lbs. Material: Black Vinyl Dimensions: 19.5" x 16"
Drivetrain	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Battery Requirements <sup>3</sup>	Size: Group 24 (75 AH) standard Group 27 (100 AH) optional Weight <sup>4</sup> : Group 24: 53 lbs. each Group 27: 60 lbs. each
Battery Charger	8-amp, off-board
Colors	  Red Silver

## Features

- Delta tiller with wraparound handles
- Deluxe reclining high-back seat with headrest, sliders, lumbar support, width, height and adjustable arms
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front and rear suspension
- Front basket
- Full lighting package includes directional signals and lower headlight
- Hydraulic handbrake
- Infinitely-adjustable tiller angle with conveniently located handle
- Large 13" low-profile, pneumatic tires for excellent outdoor performance
- Powerful drivetrain; features 200 amp microprocessor-based controller
- Rearview mirror

## Accessories

- Cell phone holder
- Cup holder
- 50" Lap belt
- Rear basket
- Rear view mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

## Military Patches



<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

**HCPCS Code: K0899**

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC  
OR DOES NOT MEET CRITERIA

(US) 800-800-8586

(CAN) 888-570-1113

pridemobility.com

