

Model	SC134
Weight Capacity	500 lbs.
Maximum Speed ²	Up to 9.3 mph
Ground Clearance ³	3.5" to motor, 4.4" at mid-deck
Turning Radius ³	82.28"
Overall Length ³	55.9"
Overall Width ³	27.05"
Seat-to-Ground Height	24.2" - 25.6"
Seat-to-Deck Height	16.5" - 18.4"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Full suspension
Range Per Charge ^{2,4}	Up to 23.4 miles
Total Weight	201.2 lbs. heaviest piece 253.8 lbs. with seat 298.8 lbs. with batteries 351.4 lbs. with seat and batteries
Heaviest Piece When Disassembled	N/A (Unit does not disassemble)
Standard Seating ³	Type: High back Weight: 52.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ⁵	Type: (2) 12V, sealed lead-acid, non-spillable deep-cycle Size: 75Ah Weight: 48.8 lbs each
Battery Charger	8 amp, off-board
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6 months

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

NOW AVAILABLE IN 4 CAMO OPTIONS



**PRIDE
MOBILITY®**

800-800-8586 • pridemobility.com



1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of +10%, -5%.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4. Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on battery manufacturer.



SHOWN WITH OPTIONAL
CANOPY UPGRADE



MORE ROOM TO ROAM.




**WEIGHT
CAPACITY**



**CRUISING
SPEED**



SUSPENSION



**23.4 MILE
RANGE**

The PX4 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.



NFC-enabled key fob for enhanced security

Memory foam seat with 6" slider for adjustable support and premium, all-day comfort



NEW

OPTIONAL hard top canopy with durable, powder-coated steel frame, polycarbonate windshield, panoramic rearview mirror, and rear mounting brackets

Delta steering wheel with water-resistant LCD console

USB charger port

Standard front basket and dual rearview mirrors (not shown)

Enhanced LED headlight

Full front and rear suspension


Pneumatic tires mounted on split rims for off-road thrills

FDA Class II Medical Device¹

FEATURES

- NFC key fob for enhanced security. Simply place or tap near the receiver to start
- 500 lb. weight capacity
- Impressive top speed of up to 9.3 mph
- Cruises up to 23.4 miles on a single battery charge
- Depth adjustable 18"x18-20" memory foam seat with 6" slider with additional seat size options that include 20"x18-20" and 22"x20-22"
- Delta steering wheel with LCD console
- Dual rear view mirrors and front basket (standard)
- Water-resistant user interface
- Full LED lighting package for maximum visibility
- Pneumatic tires mounted on split rims for excellent handling and traction on every terrain
- Supplemental hand brake for additional stopping power
- Choice of nine unique vinyl accent color options, including 4 camo options, with a 7-year outdoor durability rating

ACCENT COLORS



OPTIONAL ACCESSORIES**


- 3" Single/Double Cane or Crutch Holder
- Cane Crutch Holder
- Cell Phone Holder
- Cup Holder*
- Lap Belt
- Military Patches
- Oxygen Tank Holder
- Rear Basket*
- Saddlebag*
- Safety Flag
- Storage Pod
- Walker Holder
- Weather Cover

* Cup holder, rear basket and saddlebag can be ordered as a bundle.

** Only one rear accessory can be used at a time.



LCD digital display



Enhanced LED headlight



Full suspension

NEW Canopy and optional canopy accessories:

- Tool clips (pair)
- Jumbo rear basket
- Canopy weather cover