


	Revo® 3-wheel	Revo® 4-wheel
Model number:	SC63	SC64
Weight capacity:	300 lbs.	300 lbs.
Maximum speed ¹ :	Up to 4.5 mph	Up to 4.5 mph
Ground clearance ³ :	2.25"	2.25"
Turning radius ³ :	38"	50"
Overall length ³ :	43.25"	44.75"
Overall width ³ :	21.5"	21.5"
Wheels		
Front:	2" x 8" solid	2" x 8" solid
Rear:	3" x 9" solid	3" x 9" solid
Anti-tips:	Rear	Rear
Suspension:	N/A	N/A
Transport weight		
Total w/o batteries ^{3,6} :	106.5 lbs.	125 lbs.
Heaviest piece when disassembled:	48 lbs. (rear section)	51.5 lbs. (front section)
Drivetrain:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	N/A	
Std. seating:	Foldable molded plastic (17"W x 17"D)	
Battery charger:	Off-board, 3 amp	
Per-charge range ^{1,2} :	Up to 20 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)	
Battery & weight size ^{3,4,5} :	U-1 (24.5 lbs.each)	
Warranty:	3-year limited warranty	

Colors:



Red



Blue

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ASNI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ AGM or Gel-Cell type required.

⁵ Battery weight may vary based on battery manufacturer.

⁶ Includes standard seat and weight.

Benefits

- Exclusive black, non-scuffing tires
- Innovative interlocking drop-in battery packs for easy removal/installation
- Easily disassembles into five lightweight, easy-to-manage pieces
- 300 lbs. weight capacity
- State-of-the-art, microprocessor-based motor controller
- Seat offers height adjustment, 360° rotation and flip-up armrests for optimal user comfort
- Adjustable tiller for maximum comfort
- Headlight and taillight for added safety
- Easily-accessible charger port and circuit breaker

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag
- Quad cane holder
- Cane, clip attachment

Silver Star Lifts & Ramps 

- Backpacker® Plus
- Full Platform
- Silver Boom 300
- Backpacker® MV



Backpacker® Plus



Full Platform

These Silver Star® Lifts are compatible with the Revo®.

SLEEK STYLING
AND SUPER PORTABILITY

REVO®

The **Revo**® delivers full-sized scooter style, performance and features in a compact, lightweight design that easily disassembles for transport while remaining highly affordable.

