

# SOPHISTICATED AND STYLISH

## VICTORY®

Versatility knows no bounds with Victory® Scooters. Enjoy simple disassembly with one hand and advanced features designed for convenience and peace of mind. You can feel confident that the Victory's long-lasting batteries will keep you going. Anything is possible with a Victory® Scooter!

### All models feature...

- Feather-touch disassembly
- Standard front basket
- Delta tiller with port to conveniently charge your scooter



FDA Class II Medical Device<sup>1</sup>



SCAN THE QR CODE  
TO LEARN MORE ABOUT  
**VICTORY® SCOOTERS!**



FOLLOW US ON SOCIAL MEDIA

**PRIDE®**  
MOBILITY



FDA  
Class II  
Medical  
Device<sup>1</sup>

**SPEED<sup>2</sup>**  
Up to  
5.2 mph  
—  
**WEIGHT  
CAPACITY**  
400 lbs.  
—  
**COLORS**

### FEATURES

- Available in a 3-wheel model only
- Frosted LED headlights and low-mounted LED pathway light brightens your route
- 18"x17" seat with swivel and sliders
- Integrated storage/cup holder in tiller
- Maximum range up to 13 miles (40AH batteries)<sup>3,4</sup>

### WARRANTY

Lifetime limited on frame; 2-year limited on electronics  
2-year limited on drivetrain; 6-month limited on batteries



FDA  
Class II  
Medical  
Device<sup>1</sup>

**SPEED<sup>2</sup>**  
Up to  
5.3 mph  
—  
**WEIGHT  
CAPACITY**  
400 lbs.  
—  
**COLORS**

### FEATURES

- Available in 3 and 4-wheel models
- LED headlight
- 18"x17" black vinyl seat with 8-point swivel and sliders
- Ambient voltmeter
- Maximum range up to 15.5 miles with 3 and 4-wheel models<sup>3</sup>

### WARRANTY

Lifetime limited on frame; 2-year limited on electronics  
2-year limited on drivetrain; 6-month limited on batteries



FDA  
Class II  
Medical  
Device<sup>1</sup>

**SPEED<sup>2</sup>**  
Up to  
5.3 mph  
—  
**WEIGHT  
CAPACITY**  
300 lbs.  
—  
**COLORS**

### FEATURES

- Available in 3 and 4-wheel models
- LED headlight
- 18"x17" black vinyl seat with 8-point swivel and sliders
- Ambient voltmeter
- Maximum range up to 13 miles<sup>3</sup>

### WARRANTY

Lifetime limited on frame; 2-year limited on electronics  
2-year limited on drivetrain; 6-month limited on batteries

K0899 - POWER MOBILITY DEVICE,  
NOT CODED BY DME PDAC OR DOES  
NOT MEET CRITERIA



FDA  
Class II  
Medical  
Device<sup>1</sup>

**SPEED<sup>2</sup>**  
Up to  
8.6 mph<sup>5</sup>  
—  
**WEIGHT  
CAPACITY**  
400 lbs.  
—  
**COLORS**

### FEATURES

- Available in a 4-wheel model only
- CTS Suspension
- Full LED lighting package includes curb light, rear lights and directional signals
- Puddle light under tiller
- 18"x18" high back seat with sliders
- Maximum range up to 19.7 miles with 200 lb. user weight and up to 13 miles with 400 lb. user weight<sup>3</sup>

### WARRANTY

Lifetime limited on frame; 2-year limited on electronics  
2-year limited on drivetrain; 6-month limited on batteries

1) Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.  
2) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This can be subject to a variance of (+ or -) 10%.  
3) Range varies based on user weight, terrain and battery type.  
4) Maximum range up to 11 miles with U-1 batteries.  
5) Maximum speed is reduced to 8 mph with a user weight over 200 lbs.