

A SMOOTH RIDE for your Active Lifestyle



VICTORY® 10
LX WITH CTS SUSPENSION

CTS SUSPENSION



Front suspension



Rear suspension



The **Victory® 10 LX with CTS Suspension** offers a smooth ride with performance to match. Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. This industry first technology incorporates advanced suspension components for improved absorption over rough surfaces and a more comfortable ride experience over various terrains.

PRIDE® **MOBILITY SCOOTERS**
LIVE YOUR BEST®

| pridemobility.com |



Specifications 4W: S710LX

Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 5.6 mph
Ground Clearance ³	2.4" at transaxle, .9 at anti-tip, 2.9" at front frame
Turning Radius ³	60"
Overall Length ³	47"
Overall Width ³	22.25"
Seat-to-Deck Height	16.5"-18.5"
Seat-to-Ground Height	23.5"-22.5"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	Front and rear
Range Per Charge ^{1,2}	Up to 16 miles
Total Weight Without Batteries ^{3,7}	179 lbs.
Heaviest Piece When Disassembled ³	78 lbs. (front section)
Standard Seating	Type: Reclining high-back Weight: 45.5 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Electronic, regenerative and electromechanical
Battery Requirements ^{3,4,5}	Size: (2) 40 AH Weight: 32 lbs. each
Battery Charger	5-amp, off-board (40 AH)
Colors	 

Candy Apple Red Viper Blue

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).
- AGM or gel cell type required.
- Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Battery weight may vary based on manufacturer.

Features

- Comfort-Trac front and rear suspension
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Delta tiller with wraparound handles
- Easy-to-read LED battery meter
- Exclusive stylish, lightweight, non scuffing, black, low-profile tires
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination with a courtesy light under tiller
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrests
- Storage compartment below tiller
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

Accessories

- Cane/crutch holder (*single, double*)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

Optional Seating

- 22" x 18" black vinyl high-back seat

Military Patches



LIFT FINDER

Pride offers a variety of lifts to **transport your mobility product with ease**. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.