

Zero Turn 10

Weight Capacity	400 lbs.
Maximum Speed ²	Up to 7.3 miles at 400 lbs. Up to 7.5 miles at 200 lbs.
Ground Clearance ⁴	2" at motor; 3" at mid-frame
Turning Radius ⁴	43"
Overall Length ⁴	48"
Overall Width ⁴	24.25"
Seat-to-Ground Height	23" -25"
Seat-to-Deck Height	17.75"-19.25"
Front Tires	9" solid
Rear Tires	10.75" solid
Suspension	Front and rear (CTS Suspension)
Range Per Charge ^{2,3}	Up to 15.9 miles at 400 lbs. Up to 21.2 miles at 200 lbs.
Total Weight ⁶	Without Batteries: 185.8 lbs. With Batteries: 242.3 lbs.
Heaviest Piece When Disassembled ⁴	75 lbs. (front section)
Standard Seating ⁴	Type: Foldable, deluxe contour Weight: 34 lbs. or Foldable, compact Weight: 30 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 17"
Drive System	24-volt DC motor, rear-wheel drive, dual in-line
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{5,7}	Size: (2) 40 Ah Weight: 30.25 lbs. each
Battery Charger	Off-board, 3.5 amp
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

ZERO TURN 10

STABILITY OF A 4-WHEEL WITH THE MANEUVERABILITY OF A 3-WHEEL



PRIDE[®]
MOBILITY

(US) 800.800.8586 • pridemobility.com



** Product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Other conditions apply.

- 1 Pride[®] FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 3 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 5 Battery weight may vary based on battery manufacturer.
- 6 Includes standard seat weight.
- 7 AGM or gel cell type required.

FEATURING **EZ TURN**
TECHNOLOGY

FDA Class II Medical Device¹

NAVIGATE TIGHT CORNERS AND SMALL SPACES EFFORTLESSLY



18"x17"
Black vinyl,
deluxe stadium
seat with swivel
and sliders,
height, width,
depth and angle
adjustable armrests

Ambient voltmeter

Tiller
adjustment
lever

Standard large
front basket

Full LED light
package includes
LED headlights
and rear lights
with sequential
turn signals and
puddle light
under the tiller

REFINEMENT MEETS INNOVATION

With contemporary styling and cutting-edge technology, the Zero Turn 10 is the perfect combination of aggressive performance and maneuverability. Powerful dual motors and two-wheel drive offers aggressive traction, easily tackling rugged terrain on trails. Maneuver smoothly through the halls and doorways of your home with EZ Turn Technology.

Features

- Three-speed control lets you adjust your speed when moving from indoor to outdoor driving
- Delta tiller with wraparound handles
- Easy access tiller-mounted XLR charger port
- USB charger integrated in the tiller for convenient charging of smart phones and portable electronic devices
- Integrated storage/cup holder in the tiller

Available Colors



Black
(matte)**

Pearl
White

Candy
Apple Red

Ocean Blue
(matte)

Accessories*

- 50", 60", 70" lap belt
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Military patches
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Walker holder
- Weather cover
- Wishbone crutch holder

* Only one rear accessory can be used at a time.

** Product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only.

FDA Class II Medical Device¹



Front LED lights



Tiller adjustment lever



Pride's patent-pending EZ Turn Technology provides 4-wheel stability with 3-wheel maneuverability