

THE ULTIMATE 9" SCOOTER
WITH LUXURY FEATURES AT
AN ECONOMY PRICE.

COLT NINE

The **Colt Nine** delivers the high-performance operation and feature-rich design of the Colt family of scooters. Feather-touch disassembly, exclusive low-profile tyres, a backlit battery gauge, and a wraparound delta tiller all combine to make the Colt Nine an exceptional value.



(Images are for illustration
purposes only)


Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —

(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pride-mobility.co.uk

Colt Nine

| | |
|---|---|
| Weight capacity: | 23 Stone 3 lbs. (147 kg) |
| Maximum speed ¹ : | Up to 4 mph (6.4 km) |
| Ground clearance ² : | 4.45 cm (1.75") at mid-frame |
| Overall length ² : | 114.3 cm (45") |
| Overall width ² : | 56 cm (22") |
| Turning radius ² : | 130.81 cm (51.5") |
| Tyre size: | Front: 22.86 cm (9") solid Rear: 22.86 cm (9") solid |
| Suspension: | N/A |
| Total weight without batteries ^{2,4} : | 61.5 kg (135.5 lbs.) |
| Heaviest piece when disassembled: | Front section: 28.1 kg (62 lbs.) |
| Standard seating: | Type: Foldable molded plastic Weight: 12.7 kg (28 lbs.) Material: Black vinyl Dimensions: 45.72 cm (18") overall width 43.18 cm (17") effective depth |
| Seat-to-deck range: | 44.5 cm - 48.3 cm (17.5" - 19") |
| Drive system: | Rear-wheel drive, sealed transaxle, 24-volt DC motor |
| Dual braking system: | Regenerative and electromechanical |
| Battery requirements ^{3,5} : | Two 12-volt, deep cycle Size: 30 AH (U-1) Weight: 9.5 kg (21 lbs.) each |
| Battery charger: | 5-amp, off-board |

Colour:



Features

- Feather-touch disassembly with only one hand, provides the ultimate quick transport system
- Exclusive, stylish, lightweight, low-profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Includes pinchless hinges & eight-position seat rotation
- 23 Stone (136 kg) weight capacity
- Suspension seat post

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Saddlebag

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition.

This specification can be subject to a variance of + 10%, -5%.

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Includes standard seat.

⁵ Battery weight may vary based on manufacturer.

⁶ Tested in accordance with ISO 7176-22, 2000 and ISO 7176-11, 1992 standards.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards. All specifications subject to change without notice.



32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.co.uk