

# Colt Plus

A FEATURE-RICH 10" SCOOTER,  
WITH LUXURY TOUCHES:  
THE INTELLIGENT CHOICE.

## COLT PLUS

The sleek, sporty **Colt Plus** scooter delivers high-performance operation, all-new features and feather-touch disassembly. With a full complement of advanced standard features including backlit battery gauge and a wraparound delta tiller, the Pride® Colt Plus offers simply the best value for your money.



(Images are for illustration purposes only)


**Pride**<sup>®</sup>  
Mobility Products Ltd.  
— LIVE YOUR BEST<sup>™</sup> —

(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)

# Colt Plus

|   |   |
|---|---|
| Weight capacity:                                | 25 Stone (158.75 kg)  |
| Maximum speed <sup>1</sup> :                    | Up to 4 mph (6.4 km/h)  |
| Ground clearance <sup>2</sup> :                 | 5.7 cm (2.25")  |
| Turning radius <sup>2</sup> :                   | 137 cm (54")  |
| Overall length <sup>2</sup> :                   | 119 cm (47")  |
| Overall width <sup>2</sup> :                    | 56 cm (22")   |
| Tyre size:                                      | Front: 9 cm x 25.4 cm (3.5" x 10") solid<br>Rear: 9 cm x 25.4 cm (3.5" x 10") solid   |
| Suspension:                                     | N/A   |
| Total weight without batteries <sup>2,4</sup> : | 62.5 kg (138 lbs.)  |
| Heaviest piece when disassembled:               | Front section: 29.7 kg (65.5 lbs.)  |
| Standard seating:                               | Type: CRS Junior Seat<br>Weight: 19 kg (42 lbs.)<br>Material: Black vinyl<br>Dimensions: 44.45 cm (17.5") overall width<br>43.18 cm (17") effective depth       |
| Seat-to-deck range:                             | 44.5 cm - 48.3 cm (17.5" - 19")   |
| Drive system:                                   | Rear-wheel drive, sealed transaxle, 24-volt DC motor  |
| Dual braking system:                            | Regenerative and electromechanical  |
| Battery requirements <sup>3</sup> :             | Two 12-volt, deep cycle<br>Size: 31-36 AH or 40-45 AH (optional)<br>Weight <sup>5</sup> : 31-36 AH: 10 kg (22.5 lbs.) each/ 40-45 AH: 14.74 kg (32.5 lbs.) each |
| Battery charger:                                | 5-amp, off-board  |

Colours:



## Features

- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- High-visibility automatic rear brake lighting system
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Exclusive stylish, lightweight, non-marking, low-profile wheels
- Backlit battery gauge
- Long-lasting, high-intensity LED headlight, taillight and indicators
- Lightweight, adjustable seat with viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual voltage off-board charger

## Options

- 40 amp batteries
- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of + 10%, -5%.

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Includes standard seat.

<sup>5</sup> Battery weight may vary based on manufacturer.

<sup>6</sup> Tested in accordance with ISO 7176-22, 2000 and ISO 7176-11, 1992 standards.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards. All specifications subject to change without notice.



**Pride**  
Mobility Products Ltd.  
— LIVE YOUR BEST —

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)

