

# Colt Pursuit

COMFORTABLE OUTDOOR  
PERFORMANCE WITH  
LUXURIOUS FEATURES

## COLT PURSUIT

The **Colt Pursuit** features a powerful drivetrain and large 33 cm (13") pneumatic wheels to handle rugged outdoor terrain. Luxury touches like a wraparound delta tiller and a new LCD console combine to give the Colt Pursuit a sleek, sporty style and a maximum level of comfort.



Image is for illustration  
purposes only.


**Pride**<sup>®</sup>  
Mobility Products Ltd.  
— LIVE YOUR BEST<sup>®</sup> —

(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)

# Colt Pursuit

Model number:	SC713
Weight capacity:	28 Stone 8 lbs. (181 kg)
Maximum speed <sup>1</sup> :	Up to 8 mph (12.87 km/h)
Ground clearance <sup>2</sup> :	9.2 cm (3.625") at motor
Turning radius:	167.5 cm (66")
Overall length <sup>2</sup> :	135 cm (53")
Overall width <sup>2</sup> :	63.5 cm (25")
Front tyres:	10 cm x 33 cm (4" x 13") pneumatic
Rear tyres:	10 cm x 33 cm (4" x 13") pneumatic
Suspension:	Front & rear
Total weight without batteries <sup>4</sup> :	97.5 kg (215 lbs.)
Heaviest piece when disassembled:	Front section: 43 kg (94.5 lbs.)
Standard seating:	Type: High-back, fold-flat seat with headrest and seat belt Weight: 20 kg (44 lbs.) Material: Black Vinyl Dimensions: Width - 48.25 cm (19") overall width Depth - 43.18 cm (17") effective depth
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	Two 12-volt, deep cycle Size: 50 AH or 70 AH (optional) Weight: 50 AH: 14.8 kg (32.5 lbs.) each/ 70 AH: 24 kg (54 lbs.)
Battery charger:	8-amp, off-board (standard)

Colours:



<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition.

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Includes standard seat.

<sup>5</sup> Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)

## Features

- LCD console with speed sensor
- Rear LED lights
- Front and rear suspension
- Large 33 cm (13") pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, LED headlight
- Delta tiller with wraparound handles
- Front basket standard
- Easily accessible tie down points (for transport of unoccupied scooter)
- Handbrake
- Rearview mirror standard

## Options

- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Rear basket
- Walker holder
- 70 AH batteries



**NEW LCD console with speed sensor (standard)**



**NEW rear LED lights**

