

# Colt Sport

## HIGH-PERFORMANCE OPERATION AT AN AFFORDABLE PRICE

### COLT SPORT

The all-new compact **Colt Sport** Scooter features a sleek, sporty design and delivers high-performance operation. Just some of the standard features include a handbrake, front and rear suspension, feather-touch disassembly, a wraparound delta tiller, high-visibility automatic brake lights and an LED kerb light. These features combine to make the Colt Sport a superb scooter at an exceptional price.



(Images are for illustration purposes only)


(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)

**Pride**<sup>®</sup>  
Mobility Products Ltd.  
— LIVE YOUR BEST —

## Colt Sport

Model number:	SC710DXW
Weight capacity:	25 Stone (158.75 kg)
Maximum speed <sup>1</sup> :	Up to 8 mph (12.9 km/h)
Ground clearance <sup>2</sup> :	5.7 cm (2.25")
Turning radius <sup>2</sup> :	133.35 cm (52.5")
Overall length <sup>2</sup> :	119.4 cm (47")
Overall width <sup>2</sup> :	65.7 cm (25.9")
Tyre size:	Front: 25.4 cm (10") solid Rear: 25.4 cm (10") solid
Suspension:	Front and rear
Total weight without batteries <sup>2,4</sup> :	68 kg (150 lbs.)
Heaviest piece when disassembled:	Front section: 29.48 kg (65 lbs.)
Standard seating:	Type: High-back, fold-flat seat with headrest and seat belt Weight: 20 kg (40 lbs.) Material: Black Vinyl Dimensions: Width - 43.18 cm (17") overall width Depth - 43.18 cm (17") effective depth
Seat-to-deck range:	47 cm - 49.5 cm (18.5" - 19.5")
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	Two 12-volt, deep cycle Size: 40 AH Weight <sup>5</sup> : 40 AH: 11 kg (24.3 lbs.) each
Battery charger:	5-amp, off-board

Colours:



### Features

- Front and rear suspension
- High-visibility automatic rear brake lighting system
- Exclusive feather-touch disassembly provides the ultimate quick transport system
- Long-lasting, high-intensity LED headlight, taillight and indicators
- Exclusive stylish, lightweight, non-marking, low-profile wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- LED battery meter for increased accuracy and visibility in any setting
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access, tiller-mounted charger port for standard dual voltage off-board charger
- Rearview mirror standard
- Handbrake

### Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Rear basket
- Walker holder

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of + 10%, -5%.

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Includes standard seat.

<sup>5</sup> Battery weight may vary based on manufacturer.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards. All specifications subject to change without notice.

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL

(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)



**Pride**  
Mobility Products Ltd.  
— LIVE YOUR BEST —

UK\_SC\_COLT SPORT\_SS\_4-1-14

DEALER INFO HERE

