

DUAL-FRONT WHEEL SCOOTER



HIGH-PERFORMANCE  
OPERATION AND  
ALL-NEW FEATURES

The sleek, sporty Colt Twin scooter delivers unique dual front wheels for the maximum stability of a four-wheel scooter with the superb manoeuvrability of a three-wheel scooter. The Colt Twin has a full complement of advanced standard features such as low-profile tyres and a wraparound delta tiller.



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

(Tel) 01869 324 600 • (Fax) 01869 323 070  
32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
[www.pridemobility.co.uk](http://www.pridemobility.co.uk)

# SPECIFICATIONS

## COLT TWIN

## dual front wheel Model # SC712DX

Weight capacity:	25 Stone (158.75 kg)
Maximum speed <sup>1</sup> :	Up to 4 mph (6.4 km/h)
Ground clearance <sup>2</sup> :	5.7 cm (2.25") at motor
Turning radius:	108.6 cm (42.75")
Overall length <sup>2</sup> :	116 cm (46")
Overall width <sup>2</sup> :	56 cm (22")
Front tyres:	Dual 23 cm (9") solid
Rear tyres:	25 cm (10") solid
Suspension:	N/A
Total weight without batteries <sup>5</sup> :	62.5 kg (138 lbs.)
Heaviest piece when disassembled:	Front section: 27.2 kg (60 lbs.)
Standard seating:	Type: High-back, fold-flat seat with headrest and seat belt Weight: 20 kg (44 lbs.) Material: Black Vinyl Dimensions: width - 48.25 cm (19") overall width depth - 43.18 cm (17") effective depth
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	(2) 12-volt, deep cycle Size: (2) 31-36 AH or (2) 40AH Weight: 31-36 AH: 9.8 kg (21.5 lbs.) 40 AH 14.74 kg (32.5 lbs.) each
Battery charger:	5-amp, off-board

Colour:



Red



Blue

- <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>3</sup> AGM or gel cell type required.
- <sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
- <sup>5</sup> Battery weight may vary based on manufacturer.

**NOTE: All specifications subject to change without notice.**

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

## FEATURES

- Dual front wheels for added stability
- Delta tiller with wraparound handles
- Exclusive stylish, lightweight, non-marking, low-profile wheels
- Backlit battery gauge
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual-voltage off-board charger

## OPTIONS

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Rear basket
- Walker holder



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

(Tel) 01869 324 600 • (Fax) 01869 323 070

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL

[www.pridemobility.co.uk](http://www.pridemobility.co.uk)