

GO-GO

ELITE TRAVELLER®

The Quickest Way to Travel.

The compact design of the **Go-Go Elite Traveller®** allows it to easily manoeuvre in tight spaces while providing stable outdoor performance. The Go-Go Elite Traveller disassembles into 5 pieces, enabling a safe and efficient transition from car boot to on the road. Go where you want to go, easily, with the Go-Go Elite Traveller.



Go-Go Elite Traveller®, 3-Wheel

Available as a 3 & 4-Wheel scooter



Feather-touch disassembly.



Go-Go Elite Traveller®, 4-Wheel



SPECIFICATIONS

Model numbers	SC40E (3-Wheel) SC44E (4-Wheel)
Max. weight capacity	21 Stone 6 lbs. (136 kg)
Maximum speed ²	Up to 4 mph (6.4 km/h)
Ground clearance ¹	5.7 cm (2.25") at midframe 3.5 cm (1.38") at motor
Turning radius ¹	3-Wheel: 83.8 cm (33") 4-Wheel: 112 cm (44.125")
Overall length ¹	3-Wheel: 94 cm (37") 4-Wheel: 100.3 cm (39.5")
Overall width ¹	49.5 cm (19.5")
TYRES:	
Front	3-W: 5 cm x 20.3 cm (2" x 8") solid 4-W: 5 cm x 17.8 cm (2" x 7") solid
Rear	6.35 cm x 20.3 cm (2.5" x 8") solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drivetrain	24-volt DC motor, sealed mini-transaxle
Dual braking system	Regenerative and electromechanical
Battery charger	Off-board, 2 amp
Battery requirements ³	(2) 12V deep-cycle (standard)
Compatible batteries	12 AH; 18 AH (optional)
WEIGHT:	
Batteries (2 required) ⁴	12 AH (8.8 kg (19.5 lbs.) battery pack) 18 AH (13.4 kg (29.5 lbs.) battery pack)
Weight w/o battery pack	3-Wheel: 34.2 kg (75.5 lbs.) 4-Wheel: 36 kg (79.5 lbs.)
Heaviest piece	3-Wheel: 13.4 kg (29.5 lbs.) (front section) 4-Wheel: 15.4 kg (34 lbs.) (front section)

FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black, non-scuffing tyres
- Auto-connecting front-to-rear cable
- Includes 3 sets of easily changeable coloured panels in Red, Blue & Silver
- Front frame-mounted seat post offers maximum stability
- Frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Standard front tiller-mounted basket included
- Removable deck basket for additional storage (4-Wheel only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board

¹ Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.
² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
³ AGM or gel-cell type recommended.
⁴ Battery weight may vary based on manufacturer.


NOTE: All specifications subject to change without notice.
 The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Standard Colours



3 sets of colours with every scooter!

Optional Colours available in colour pack (4-wheel version only):



(Tel) 01869 324 600 • (Fax) 01869 323 070 • 32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
www.pridemobility.co.uk



DEALER INFO HERE

