

# **Jazzy**<sup>®</sup> **SELECT**

The highly manoeuvrable Jazzy<sup>®</sup> Select delivers the advantages of in-line motor technology for enhanced efficiency, torque, range, and performance. Plus, it's loaded with performance and convenience features, like standard Active-Trac<sup>®</sup> Suspension, that make it as easy to enjoy as it is to use.



**Pride**<sup>®</sup>  
Mobility Products Ltd.  
— LIVE YOUR BEST —

## SPECIFICATIONS

<b>Weight capacity</b>	21 Stone 6 lbs. (136 kg)
<b>Maximum speed</b> <sup>1</sup>	Up to 6.44 km/h (4 mph)
<b>Ground clearance</b>	4.76 cm (1.875") at midframe
<b>Turning radius</b> <sup>4</sup>	50.8 cm (20")
<b>Overall length</b> <sup>4</sup>	90.80 cm (35.75") without front riggings
<b>Overall width</b> <sup>4</sup>	60.32 cm (23.75")
<b>Drive wheels</b>	Drive: 25.4 cm (10")
<b>Anti-tips</b>	12.7 cm (5") solid, front-mounted
<b>Castors</b>	15.24 cm (6") solid, rear-mounted
<b>Suspension</b>	Active-Trac®
<b>Drivetrain</b>	Two-motor, mid-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC 2 controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 5 amp
<b>Battery requirements</b>	(2) 12 volt deep cycle
<b>Battery size (largest)</b>	30 AH
<b>Battery weight</b> <sup>2</sup>	11.11 kg (24.5 lbs.) each
<b>Base weight</b>	40.82 kg (90 lbs.)
<b>Standard seat weight</b>	19 kg (42.5 lbs.) (high-back)
<b>Standard seating dimensions</b>	46 cm x 46 cm (18" x 18")
<b>Warranty</b>	Lifetime warranty on frame; 1-year limited warranty on electronics; 18-month warranty on PG Drive control system; 2-year limited warranty on drivetrain

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

1) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.

2) Battery weight may vary based on manufacturer.

3) These are optional (additional) for configuration 1.

4) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

## FEATURES

- 40 amp, PG GC 2 controller
- Built with ease of service in mind
- Easy side access to batteries
- Toolless front panel access to electronics to minimize service time
- Active-Trac® Suspension

## OPTIONS

- Cup holder
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## COLOUR CHOICES



Red



Blue



**PRIDE**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST® —

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL  
(Tel) 01869 324 600 • (Fax) 01869 323 070 • [www.pridemobility.co.uk](http://www.pridemobility.co.uk)

