

Jazzy® 600 XL
 Order Form / Price List
 Revised: November 2007



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:	Address:			Discount:
Phone Contact:	Marked For:			Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

JAZZY 600 XL-C CONFIGURATION #1

Jazzy 600 XL-C Red \$6,751 -20" W x 22" D Ltd. Recline Gray vinyl high back
 Group 2 / Heavy-Duty
K0825 w/ captain's seating
 400 lbs. weight capacity
 -Full length armrests
 -70 amp PG VR2 controller
 -14" Solid tires
 -4 amp on-board charger
 -60" lapbelt
 -Foot platform

NF-22 (pair) Installed AGM \$510 Gel \$590

JAZZY 600 XL-SS CONFIGURATION #2

Jazzy 600 XL-SS Red \$6,751 -20" W x 22" D Ltd. Recline Gray vinyl high back seat with solid seat pan
 Group 2 / Heavy-Duty
K0824 w/ sling/solid seating
 40 lbs. weight capacity
 *Select cushion for the solid seat pan in appropriate section or review specific cushion order form.
 -Full length armrests
 -70 amp PG VR2 controller
 -14" Solid tires
 -4 amp on-board charger
 -60" lapbelt
 -Foot platform

NF-22 (pair) Installed AGM \$510 Gel \$590

JAZZY 600 XL-C CONFIGURATION #3

Jazzy 600 XL-C Red \$6,751 -22" W x 20" D Gray vinyl high-back seat
 Group 2 / Heavy-Duty
K0825 w/ captain's seating
 40 lbs. weight capacity
 -Full length armrests
 -70 amp PG VR2 controller
 -14" Solid tires
 -4 amp on-board charger
 -60" lapbelt
 -Foot platform

NF-22 (pair) Installed AGM \$510 Gel \$590

JAZZY 600 XL-SS CONFIGURATION #4

Jazzy 600 XL-SS Red \$6,751 -22" W x 20" D Gray vinyl high-back seat with solid seat pan
 Group 2 / Heavy-Duty
K0824 w/ sling/solid seating
 400 lbs. weight capacity
 *Select cushion for the solid seat pan in appropriate section or review specific cushion order form.
 -Full length armrests
 -70 amp PG VR2 controller
 -14" Solid tires
 -4 amp on-board charger
 -60" lapbelt
 -Foot platform

NF-22 (pair) Installed AGM \$510 Gel \$590

SPECIFIC CONFIGURATION

Jazzy 600 XL \$6,751

COLORS

- Metallic Colors:**
- Champagne NC
 - Candy Apple Red \$52
 - Viper Blue \$52
 - Onyx Black \$52
 - Pewter \$52
 - Silver \$52
 - Emerald Green \$52
 - Graphite Blue \$52
 - Sapphire Blue \$52
- Solid Colors:**
- Kelly Green \$52
 - Tangerine \$52

BATTERIES (Installed)

- NF-22, 50 amp AGM E2361 \$255ea
- NF-22, 50 amp Gel Cell E2361 \$295ea

BATTERY CHARGER

- 4 amp On-board Charger Std.
- 4 amp On 5 amp Off-board Charger \$220
- Charger Plug Left Side Right Side

CONTROLLER

- PG VR2 70 amp⁹ Std.
- Pride Flight 75 amp NC
- (Flight only available as quick ship)

PROGRAMMER

- PG Programmer \$423
- Flight Programmer \$423

CONTROLLER MOUNT

- Fixed, In-Line Joystick Mount Std.
- Left Right
- Swing-away Joystick Mount¹⁶ E1028 \$196
- 4" L¹⁷ 7" L
- Left Right
- Vertical Adj. Clamp² K0108 \$95
- Provides vertical joystick adjustment w/o adjusting armrest height.

JOYSTICK HANDLE OPTIONS

- PG Controls**
- Std. Knob Std.
 - Mushroom Handle E2323 \$71
 - T-Bar E2323 \$107
 - Chin Cup E2323 \$71
 - Extended I-Handle E2323 \$71
 - Soft Ball Handle E2323 \$50
 - Medium Ball Handle E2323 \$50
 - Goal Post Handle E2323 \$66
 - 3" W 4" W

Flight Controls

- Flight Knob Std.
- Sponge Ball E2323 \$35
- Chin Cup E2323 \$30
- DL Knob E2323 \$20
- Carrot Knob E2323 \$20
- Straight Bar (Chrome) E2323 \$45
- Finger Tube (Chrome) E2323 \$45
- Round T-Bar (Chrome) E2323 \$45
- Flat T-Bar (Chrome) E2323 \$45
- Wrist Plate E2323 \$45

Jazzy® 600 XL
Order Form / Price List
 Revised: November 2007

Legrest measurements have changed. Lengths are now measured according to industry standard from the seat pan to the back of the footplate

NEW SWING-AWAY LEGRESTS

<input type="checkbox"/> 70° Swing-away footrests Style #8 ²⁴	\$162
<input type="checkbox"/> 70° Tapered S/A footrests Style #8 ²⁴	\$162
(N/C with power positioning systems)	
<input type="checkbox"/> Quad Release Lever	NC
□ S/A Standard Extensions:	
<input type="checkbox"/> 11"-13" Ext.	NC
<input type="checkbox"/> 13"-16.5" Ext	Std.
<input type="checkbox"/> 15.5"-19.5" Ext.	NC
(Lengths over 17" may cause interference)	
<input type="checkbox"/> Aluminum Footplate	Std.
<input type="checkbox"/> 7.25" x 6" ¹⁵ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 8.5 x 6" ¹⁵ (Tapered legs require 18"W or greater seat)	
□ Angle Adjustable Ext.: \$150	
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
(Lengths over 17" may cause interference)	
Angle Adjustable Footplates (pair)	
<input type="checkbox"/> 4"x6" ¹⁹ (Tapered legs require 14"W or greater seat)	
<input type="checkbox"/> 5"x6" ²⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x8" ²⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 6"x8" ¹⁵ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x11" ¹⁵ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8"x11" ⁵ (Tapered legs require 22"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25
□ Multi-Axis Ext. (ball & socket): \$495	
<input type="checkbox"/> 11"-16.5" Ext.	
<input type="checkbox"/> 15.5"-19.5" Ext.	
(Lengths over 17" may cause interference)	
Multi-Axis Footplates (pair)	
<input type="checkbox"/> 4"x6" ¹⁹ (Tapered legs require 14"W or greater seat)	
<input type="checkbox"/> 5"x6" ²⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x8" ²⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 6"x8" ¹⁵ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x11" ¹⁵ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8"x11" ⁵ (Tapered legs require 22"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

FRONT RIGGINGS

<input type="checkbox"/> Foot Platform ¹⁸	Std.
Rigging Type (pair)	
<input type="checkbox"/> Elevating Leg Rests E0990	\$325
Select footplate size for above:	
<input type="checkbox"/> 7 1/4 x 6	<input type="checkbox"/> 8 1/2 x 6
Adjustable Calf Pads w/ ELRS (1 hole std.)	
<input type="checkbox"/> 1 hole	<input type="checkbox"/> 3 holes* <input type="checkbox"/> 6 holes*
*Potential for interference w/ shrouds	
Footplates (pair)	
<input type="checkbox"/> Reinforced Footplate	Std.
<input type="checkbox"/> 7 1/4 x 6	<input type="checkbox"/> 8 1/2 x 6
<input type="checkbox"/> 5"x8" Angle-Adj. Footplate ²⁰	\$150
<input type="checkbox"/> 6"x8" Angle-Adj. Footplate ¹⁵	\$150
<input type="checkbox"/> 6"x11" Angle-Adj. Footplate ¹⁵	\$150
<input type="checkbox"/> 8"x11" Angle-Adj. Footplate ⁵	\$150
<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$150

FRONT RIGGINGS CONT.

Extensions¹²	
<input type="checkbox"/> 11"-14" Ext.	NC
<input type="checkbox"/> 14"-17" Ext	Std.
Rigging Options	
<input type="checkbox"/> Heel Loops (pair) ¹³ E0951	\$20
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ¹⁰ K0108	\$65
<input type="checkbox"/> ELR / ALR Adductor Buttons (pr) K0108	\$100
<input type="checkbox"/> POSAlinc Knee Adductor (pair) for S/A legrests K0108	\$464
<input type="checkbox"/> Residual limb support left E1020	\$300
<input type="checkbox"/> Residual limb support right E1020	\$300
□ HD Drop-Ins, 70° \$420	
N/A w/ contoured or Specialty seating	
Standard w/ Ang. dj. footplates	
Angle Extensions: NC	
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
<input type="checkbox"/> 17.5" - 21" Ext.	
(Lengths over 17" may cause interference)	

SEATING

<input type="checkbox"/> High-Back w/ Headrest	Std.
□ High-Back w/ Headrest & Solid Seat Pan⁸ Std.	
<input type="checkbox"/> 20" W x 22" D	
<input type="checkbox"/> 22" W x 20" D	
<input type="checkbox"/> 22" W x 22" D	
<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Beige Vinyl
Custom seat dimensions are available. Please contact the Custom Group for pricing and leadtime.	
22" W AND ABOVE WILL NOT BE LTD RECLINE.	

ARMRESTS

For High-Back Seats	
<input type="checkbox"/> Flip-up, removable	Std.
<input type="checkbox"/> Flip-up, height-adjustable removable ¹⁴ E0973	\$170
Armpads	
<input type="checkbox"/> Desk Length, 10"	NC
<input type="checkbox"/> Left NC	<input type="checkbox"/> Right NC
<input type="checkbox"/> Full Length, 14"	NC
<input type="checkbox"/> Left Std.	<input type="checkbox"/> Right Std.

MANUAL RECLINER^{1,3}

<input type="checkbox"/> Manual Recline High Back E1226	\$600
□ Manual Recline High Back⁸ w/ Solid Seat Pan E1226 \$600	
(180 degrees of Recline) (Dimensions 18"W x 18"D)	
<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Beige Vinyl
Width:	
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 20" NC
Depth:	
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 20" NC

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS¹¹

<input type="checkbox"/> Simplicity Cushion		
16 x 16 - 20 x 20	E2601	\$98
<input type="checkbox"/> Solution 1 Cushion		
14 x 14 - 20 x 18	E2603	\$380
20 x 20 - 20 x 24	E2603	\$430
22 x 16 - 24 x 24	E2604	\$430
<input type="checkbox"/> Solution Cushion		
10 x 10 - 20 x 18	E2607	\$380
20 x 20 - 20 x 24	E2607	\$430
22 x 16 - 24 x 24	E2608	\$430
<input type="checkbox"/> Spectrum Cushion		
<input type="checkbox"/> Gel		
10 x 10 - 20 x 18	E2605	\$480
20 x 20 - 20 x 24	E2605	\$530
22 x 16 - 24 x 24	E2606	\$530
<input type="checkbox"/> Visco		
10 x 10 - 20 x 18	E2607	\$480
20 x 20 - 20 x 24	E2607	\$530
22 x 16 - 24 x 24	E2608	\$530
<input type="checkbox"/> Structure 2 Cushion		
14 x 16 - 20 x 18	K0734	\$480
20 x 20 - 20 x 24	K0734	\$580
22 x 18 - 24 x 24	K0735	\$580

Footnotes:

1. Rear mount accessories are N/A.
2. N/A on contoured seating or waterfall armpads.
3. Special measurements are not available.
4. 18" x 18" std. Special order available. Call Customer Service.
5. Only avail. on 20 W" seats.
6. Not available w/ manual recline high-back.
7. Standard w/ heel loops and angle-adjustable footplates, check size below.
8. Includes black vinyl seat pan, back available in all listed colors.
9. Std. w/ one programmable drive profile.
10. Available w/ swing-away footrests only.
11. See Synergy cushion order form for sizes, pricing and coding.
12. Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and /or cushion.
13. Std. when angle-adjustable footplates are selected.
14. Not available on Manual Recliner.
15. Only avail on 16" W or wider seats.
16. Can only mount one armrest accessory in combination w/ S/A joystick.
17. N/A with waterfall armpads and dual post armrests.
18. Does not elevate.
19. Only avail. on 12" W seats.
20. Only avail. on 14" W seats.

Jazzy® 600 XL
Order Form / Price List
 Revised: November 2007

Synergy Seating Systems

Synergy Seating™ System \$434
 (Price with Base)
400 lbs weight capacity

WIDTH

Set Seat Width At:			
<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC		

Seat Width Package

<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"W-18"W
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W

DEPTH

Set Seat Depth At:			
<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC		

Seat Depth Package

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"D-20"D	<input type="checkbox"/> 20"W-24"W

BACK CANES

Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

BACK OPTIONS

- Solid Plate Back up to 20" W \$275
- Solid Plate Back 21"-24" W \$325
- Back Height:**¹²
- 18"H 19"H 20"H 21"H
- 22"H 23"H 24"H 25"H
- Sling Back¹³ \$98

*Back height will be the same as back cane height.
 Includes gap for cushion.
 Push Bar selection must match back width; center tube standard when selecting flip-up arms.*

- TRU-Comfort Back**¹⁶
- Back Height:** 22" 25"
- Back up to 20" W **E2620** \$595
- 22" W **E2621** \$695
- TRU-Comfort Cushion**
- 16 x 16 - 20 x 18 **E2601** \$380
- 20 x 20 - 20 x 22 **E2601** \$380
- 22 x 16 - 22 x 24 **E2602** \$480

See TRU-Comfort Back & Cushion order forms for specific sizes

ARMRESTS

- Flip-back, cane mounted ^{11,14} Std.
- 2-Post, flip-back, height ^{11,14} **E0973** \$150ea adjustable, removable
- (9"-13" of height adjustment in 1" increments)*
- Single Post height adj.^{10,15} **E0973** \$75ea quick height
- 11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.*
- Quad Release Lever NC
- Removable, height adjustable **E0973** \$75ea heavy duty armrests¹⁵
- 8"H-12"H 10"H-14"H
- Pediatric removable¹⁵ **E0973** \$75ea adjustable, quick height.
- Match Shroud Color Black NC
- Quad Release Lever NC
- (7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)*

ARMPADS

- Full length, Straight, 2" X 14" Std.
- Left Right
- Desk Length, Straight, 2" X 10" Std.
- Left Right
- Full Length, Waterfall, 2" X 14"⁹ NC
- Left Right
- Desk Length, Waterfall, 2" X 10"^{8,9} NC
- Left Right

HEAD REST ⁷

- Stealth** ^{5,6} **E0955**
- Large (11"wide x 6" high) \$210
- Stealth** ^{5,6} **E0955**
- Small (9 3/4" wide x 3" high) \$210
- Stealth**, swing-away (Flip-Down)
- Head Rest Bracket **E1028** \$150
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications** ^{5,6} **E0955**
- Large (11"wide x 6" high) \$275
- Medical Applications** ^{5,6} **E0955**
- Small (9 3/4" wide x 3" high) \$275
- Fixed, Offset Mounting Bracket NC
- TRU-Comfort** ^{11,12} **E0955**
- 8" Wide w/ Cobra hardware \$417
- 10" Wide w/ Cobra hardware \$417
- Bracket for after market headrest \$180

ACCESSORIES

- Rear Basket^{1,2} \$103
- Cane/Crutch Holder^{1,2} **E2207** \$98
- Oxygen Holder^{1,2} **E2208** \$103
- Walker Holder^{1,2,3} **K0108** \$92
- Flag Holder (SOF on rehab seating) \$17
- Weather Cover, Large, Black \$202
- Saddle bag(contoured seating only) \$28
- Cup holder \$17
- Left Right
- Lap Belt Std.
- 38" ⁴ 50" 60"
- 70" 80" 90"

Footnotes:

1. Rear mounted accessories are not available w/ manual recline high-back seating.
2. Includes mounting bracket and hardware.
3. Walker holder is available on contoured & semi-recline.
4. N/A w/ limited recline seating.
5. Comes with swing-away / removable hardware; HCPCS code is E1028.
6. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
7. Available w/ solid backs only.
8. N/A w/ S/A joystick mount.
9. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
10. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
11. S/A laterals interfere w/ flip-back arms.
12. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
13. N/A w/ Synergy manual tilt or recline.
14. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
15. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
16. N/A w/ flip-back arms.

Lowest seat-to-floor height
Synergy Seat: 17"
High-back Seat: 21"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Quantum Custom Delivery Service
 Includes: liftgate service and delivery of your product(s) inside your facility
 Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.