

Order Form

Discount: _____
 Terms: _____
 Freight: _____
 Taken By: _____



Account Number: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 E-Mail: _____
 Date: _____
 PO#: _____
 Marked For: _____
 Ship To Address: _____
 City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASE OPTIONS		
Description	Price	HCPCS Code
<input type="checkbox"/> Litestream XF <i>250lbs. Weight Capacity, Allows Seat-To-Floor Heights 15"-19.5"</i>	\$1,995.00	K0005

ADDITIONAL BASE OPTIONS		
Description	Price	HCPCS Code
<input type="checkbox"/> Ultra Hemi Package ¹ <i>250lbs. Weight Capacity, Allows Seat-To-Floor Heights 12.5"-16.5"</i>	\$120	K0056
<input type="checkbox"/> Heavy Duty Package ² <i>Increases Weight Capacity To 350lbs, Available 18-22"W</i>	\$495	K0108

1. Not available with HD Package or Occupied Transit Kit.
 2. 14"D-17"D Requires SOF. Not available with Occupied Transit Kit.

COLOR SELECTION		
Metallic Colors: N/C	Solid Colors: N/C	Custom Prints: \$795
<input type="checkbox"/> Cosmic Silver	<input type="checkbox"/> Green Mantis ¹	<input type="checkbox"/> Winter Camo
<input type="checkbox"/> Pacific Blue	<input type="checkbox"/> Pink Flamingo	<input type="checkbox"/> Digital Woodland Camo
<input type="checkbox"/> Rainforest Green	<input type="checkbox"/> Arctic White ¹	<input type="checkbox"/> Digital Marine Camo
<input type="checkbox"/> Plum Crazy ¹	<input type="checkbox"/> Orange Crush ¹	<input type="checkbox"/> American Flag
<input type="checkbox"/> Purple Haze	<input type="checkbox"/> Red Hot	<input type="checkbox"/> Woodland Camo
	<input type="checkbox"/> Blue Sky	<input type="checkbox"/> Pink Camo
	<input type="checkbox"/> Black Abyss	<input type="checkbox"/> Mossy Oak
	<input type="checkbox"/> Yellow Sunshine ¹	<input type="checkbox"/> Snake Skin
		<input type="checkbox"/> Digital Army Camo

1. Extended 12 day lead time.

SEATING DIMENSIONS		
Description	Price	HCPCS Code
Width:		
<input type="checkbox"/> 14"W	No Charge	
<input type="checkbox"/> 15"W	No Charge	
<input type="checkbox"/> 16"W	No Charge	
<input type="checkbox"/> 17"W	No Charge	
<input type="checkbox"/> 18"W	No Charge	
<input type="checkbox"/> 19"W	No Charge	
<input type="checkbox"/> 20"W	\$450	E2201
<input type="checkbox"/> 21"W		
<input type="checkbox"/> 22"W ¹		
Depth:		
<input type="checkbox"/> 14"D ²	No Charge	
<input type="checkbox"/> 15"D ²	No Charge	
<input type="checkbox"/> 16"D ²	No Charge	
<input type="checkbox"/> 17"D	No Charge	
<input type="checkbox"/> 18"D	No Charge	
<input type="checkbox"/> 19"D	No Charge	
<input type="checkbox"/> 20"D	\$500	E2203
Frame Length:		
<input type="checkbox"/> 16" (14"D-16"D)		
<input type="checkbox"/> 18"D (16"D-18"D)		
<input type="checkbox"/> 20"D (18"D-20"D)		

1. Requires HD package.
 2. 14"D-17"D with HD package requires SOF.

AXLE PLATE AND CAMBER		
Description	Price	HCPCS Code
<input type="checkbox"/> 3" Single Slot Axle Plate <i>Camber:</i> <input type="checkbox"/> No Camber <input type="checkbox"/> 1.5° Camber <input type="checkbox"/> 3° Camber	No Charge	
<input type="checkbox"/> 1-Arm Drive, 3" Single Slot Axle Plate	\$795	E0958
<input type="checkbox"/> Amputee Axle Plate ²	\$175	
<i>Camber:</i> <input type="checkbox"/> No Camber <input type="checkbox"/> 1.5° Camber <input type="checkbox"/> 3° Camber		

1. Dual Post, Flip Back Armrests are not available with Camber.
 2. Not available with HD Package or Occupied Transit Kit

SEAT-TO-FLOOR HEIGHT				
<i>Rear Seat-To-Floor Can Not Be Higher Than Front Seat-To-Floor Height Quick Release Fork Allows For A Maximum Dump Of 2" Bolt-On Fork Allows For A Maximum Dump Of 2.5"</i>				
Front:				
<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"
<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"
<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"
Rear:				
<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"
<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"
<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	

ARMRESTS		
Description	Price	HCPCS Code
<input type="checkbox"/> No Armrests		
Armrests:		
<input type="checkbox"/> Tubular, Swing-Away, 12"L, 10-11"H ¹	No Charge	
<input type="checkbox"/> Height Adjustable, 10"H-14"H, Removable, <i>Drop-in</i>	\$275	E0973
<input type="checkbox"/> Dual Post, Flip Back, Removable, Height- Adjustable, 10.5"-15.5" ^{1,2}	\$245	E0973
Armpads:		
<input type="checkbox"/> Straight Armpads	No Charge	
<i>Size:</i> <input type="checkbox"/> Desk Length <input type="checkbox"/> Full Length		
<input type="checkbox"/> Waterfall Armpads ³	No Charge	
<i>Size:</i> <input type="checkbox"/> Desk Length <input type="checkbox"/> Full Length		

1. Not available with Occupied Transit kit.
 2. Not available with Camber. Can not adjust to lower than 12.5"H when used with 24" or 26" rear wheels set in highest axle position.
 3. Not available w/ dual post armrests.

WHEEL LOCKS		
Description	Price	HCPCS Code
<input type="checkbox"/> Push-to-Lock ¹	No Charge	
<input type="checkbox"/> Push-to-Lock Extensions ²	\$75	E0961
<input type="checkbox"/> Pull-to-Lock	No Charge	
<input type="checkbox"/> Pull-to-Lock Extensions	\$75	E0961
<input type="checkbox"/> Hill Holder (<i>Push to Lock</i>) ³	\$295	E0974
<input type="checkbox"/> Scissor Locks ²	\$195	
<input type="checkbox"/> Unilateral Wheel Lock System ⁴	\$129	
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		

1. May interfere with ELR's.
 2. May interfere with Full Length Arms.
 3. Not available with polyurethane tires.
 4. Must Use Tubular or Flip-Up Height Adjustable Armrests. Not available with Side Guards. Not available with poly tires.

WHEEL LOCK CLAMP		
<input type="checkbox"/> Standard		
<input type="checkbox"/> Extended		

WHEEL TYPE	
<input type="checkbox"/> Pre-Configured Wheels	See Page 2
<input type="checkbox"/> Build-To-Order Wheels	See Pages 2 & 3

BACK CANES	
Description	Price
<input type="checkbox"/> 8° Back Canes ¹ w/ Push Handles	No Charge
Size: <input type="checkbox"/> 12"-16"H <input type="checkbox"/> 15"-19"H	
<input type="checkbox"/> Straight Back Canes ¹	No Charge
Size: <input type="checkbox"/> 12"-16"H <input type="checkbox"/> 15"-19"H	
Handles: <input type="checkbox"/> With Push Handles <input type="checkbox"/> Without Push Handles	
<input type="checkbox"/> Straight Reinforced Canes 15-19"H w/ Push Handles	No Charge
<input type="checkbox"/> Angle-Adjustable Back Canes ²	\$164
Size: <input type="checkbox"/> 12"-16"H <input type="checkbox"/> 15"-19"H	
Handles: <input type="checkbox"/> With Push Handles <input type="checkbox"/> Without Push Handles ¹	

1. Not available with HD package.
2. Not available with Occupied Transit Kit

BACK OPTIONS		
Description	Price	HCPCS Code
<input type="checkbox"/> Padded Sling Back	No Charge	
<input type="checkbox"/> Tension Adjustable Seat Back	\$350	E2611
<input type="checkbox"/> 22"W Tension Adjustable Seat Back	\$465	E2612
<input type="checkbox"/> Synergy Shape GP Back (19"H) ¹		
Width: <input type="checkbox"/> 14"W <input type="checkbox"/> 16"W <input type="checkbox"/> 18"W <input type="checkbox"/> 20"W	\$465	E2611
<input type="checkbox"/> 22"W	\$485	E2612
<input type="checkbox"/> Synergy Shape Back		
Hardware: <input type="checkbox"/> 2-Point ¹ <input type="checkbox"/> 4-Point		
Height: <input type="checkbox"/> Standard (19"H)		
Width: <input type="checkbox"/> 14"W <input type="checkbox"/> 16"W <input type="checkbox"/> 18"W <input type="checkbox"/> 20"W	\$495	E2613
<input type="checkbox"/> 22"W	\$595	E2614
<input type="checkbox"/> Synergy Shape Deep Contour ¹		
Height: <input type="checkbox"/> Medium <input type="checkbox"/> Tall		
Width: <input type="checkbox"/> 14"W <input type="checkbox"/> 16"W <input type="checkbox"/> 18"W <input type="checkbox"/> 20"W	\$595	E2620
<input type="checkbox"/> 22"W	\$695	E2621
<input type="checkbox"/> Omit Backrest ¹	(-\$15)	
<input type="checkbox"/> Comfort Company Back	Various	

[See Comfort Company Order Form](#)

1. Not available with Occupied Transit Kit.

MILITARY BACK PATCHES (FOR TENSION ADJUSTABLE BACK)

Branch Patches:

- U.S. Air Force No Charge U.S. Marines No Charge



- U.S. Army No Charge U.S. Navy No Charge



- U.S. Coast Guard No Charge



Purple Heart Patch:

- Purple Heart No Charge



HEADREST		
Description	Price	HCPCS Code
<i>Only compatible with Shape Back. Not configurable, must be line itemed. Refer to Positioning Components Order Form for details and part numbers.</i>		
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H)	\$205	E0955
<input type="checkbox"/> Large (11"W x 6"H)	\$205	E0955
<input type="checkbox"/> Small Comfort Plus (10"W)	\$205	E0955
<input type="checkbox"/> Large Comfort Plus (14"W)	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket	\$210	
<i>Interfaces with any headrest that uses 1/2" square down-tube.</i>		

REAR SEAT-TO FLOOR HEIGHT CHART

Drive Wheel	Rear STF
20"	12.5" - 16"
22"	13.5" - 17"
24"	14.5" - 18"
25" ¹	15" - 18.5"
26"	15.5" - 19"

1. When 25" rim is used with 26" mountain bike tire the STF height for 26" wheel applies.

PRECONFIGURED WHEEL ASSEMBLIES

Description	Price	HCPCS Code
<i>Includes Rim, Handrim and Tire</i>		
<i>For Customized Configurations See The Following Sections</i>		
<input type="checkbox"/> Composite Mag ¹ w/ Aluminum Handrim & 1 1/8" Poly-urethane Tire	No Charge	
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24"		
<input type="checkbox"/> Heavy Duty 24" Composite Wheels w/ Aluminum Handrim and Flat-Free Insert	No Charge	
<input type="checkbox"/> Composite Mag ¹ w/ Vinyl Coated Handrim & 1 1/8" Poly-urethane Tire	\$200	E2218
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24"		
<input type="checkbox"/> One-Arm Drive Wheels ^{1,2} w/ 24" Lightweight Quick-Release Wheel w/ Aluminum Handrim and Pneumatic Tire	No Charge	
<input type="checkbox"/> 24" Composite Mag ¹ w/ Vinyl Coated 45° Projection-Handrim & 1 1/8" Polyurethane Tire	\$295	E2218
<input type="checkbox"/> Spinergy Wire Spokes w/ Aluminum Handrim, & 1 1/8" Everyday Gray Tire		
Size: <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 26"	\$395	
	No Charge	
	\$395	
<input type="checkbox"/> 24" Magic Wheels ¹ 2-Speed Geared	\$5,495	E2227

1. Not available with HD package.
2. Requires 3" Single Slot Axle Plate for 1-Arm Drive. Not available with Occupied Transit Kit.

BUILD TO ORDER WHEEL CHARTS

Rear Seat-To-Floor Height:

Drive Wheel	Rear STF
20"	12.5" - 16"
22"	13.5" - 17"
24"	14.5" - 18"
25" ¹	15" - 18.5"
26"	15.5" - 19"

Tire/Rim Compatibility:

	Rim				
	Composite Mag	Spinergy Wire Spoke	Spinergy Flex Rim	Spinergy Lite Extreme	Sport Lite Extreme
Tires					
Everyday Gray	20", 22", 24"	20", 22", 24", 25"	24"	20", 22", 24", 25"	24", 25"
1.375" Knobby	24"	24", 26"	24"	24", 26"	24", 26"
Poly Touring	20", 22", 24"	N/A	N/A	N/A	N/A
145 psi Schwable	N/A	24", 25", 26"	24"	22", 25", 26"	24", 25", 26"
2.10" Knobby	N/A	25"	N/A	25"	25"
110 psi Primo	N/A	22", 24", 26"	24"	22", 24", 26"	24", 26"
Primo	20", 22", 24"	20", 22", 24", 26"	24"	20", 22", 24", 26"	24", 26"

1. When 25" rim is used with 26" mountain bike tire the STF height for 26" wheel applies.

RIMS WITH HANDRIMS

Description	Price
<i>Must Select Tire In Addition To Rim w/ Handrim.</i>	
<input type="checkbox"/> Composite Mag w/ Alum. Handrim, 24"	No Charge
<input type="checkbox"/> Spinergy Wire Spokes w/ Alum. Handrim, 24"	No Charge
<input type="checkbox"/> Spinergy Wire Spokes w/ Natural Fit Handrim, 24"	\$395
<input type="checkbox"/> Spinergy FlexRim, 24"	\$995

RIMS WITH TIRES

Description	Price	HCPCS Code
<i>Must Select Handrim In Addition To Rim w/ Tire.</i>		
<input type="checkbox"/> Composite Mags w/ 1 3/8" Gray Everyday Flat Free Tire	\$75	E2213
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24"		

RIMS

Description	Price
<i>Must Select Handrim & Tire In Addition To Rims.</i>	
<input type="checkbox"/> Composite Mags	No Charge
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24"	
<input type="checkbox"/> Spinergy Lite Extreme Carbon Fiber Spokes (12 Spoke)	\$895
Size: <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	
Color: <input type="checkbox"/> Blue <input type="checkbox"/> Black <input type="checkbox"/> White <input type="checkbox"/> Yellow <input type="checkbox"/> Red <input type="checkbox"/> Green <input type="checkbox"/> Pink <input type="checkbox"/> Orange	
<input type="checkbox"/> Spinergy Sport Lite Extreme Carbon Fiber Spokes (24 Spoke)	\$895
Size: <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	
Color: <input type="checkbox"/> Blue <input type="checkbox"/> Black <input type="checkbox"/> White <input type="checkbox"/> Yellow <input type="checkbox"/> Red <input type="checkbox"/> Green <input type="checkbox"/> Pink <input type="checkbox"/> Orange	
<input type="checkbox"/> Spinergy Wire Spokes	\$245
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	

TIRES

Description	Price	HCPCS Code
<i>Must Select Rim & Handrim In Addition To Tires. Consult Tire/Rim Compatibility Chart On Page 2.</i>		
<input type="checkbox"/> 1 3/8" Gray Everyday Use Tire		
Tube Type: <input type="checkbox"/> Pneumatic <input type="checkbox"/> Flat Free	No Charge \$145	E2213
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 26"		
<input type="checkbox"/> 1 3/8" Gray Knobby Tire		
Tube Type: <input type="checkbox"/> Pneumatic <input type="checkbox"/> Flat Free	No Charge \$145	E2213
Size: <input type="checkbox"/> 24" <input type="checkbox"/> 26"		
<input type="checkbox"/> 1 3/8" Polyurethane Full Profile Touring Tread ¹	\$180	E2218
Color: <input type="checkbox"/> Dark Charcoal Gray <input type="checkbox"/> Everyday Light Gray (24" Only)		
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25"		
<input type="checkbox"/> 1" Schwable 145 PSI Tube & Tire ²	\$180	
Color: <input type="checkbox"/> Black, Marathon Plus Evolution <input type="checkbox"/> Gray, Marathon Plus		
Size: <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"		
<input type="checkbox"/> 26" x 2.10" Knobby Mountain Bike Pneumatic Tire & Tube ³	\$145	
Color: <input type="checkbox"/> Black <input type="checkbox"/> Gray		
<input type="checkbox"/> 1 1/4" Gray Primo 110 PSI V-Trac Pneumatic Tire & Tube	\$95	
Size: <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"		
Color: <input type="checkbox"/> Gray <input type="checkbox"/> Red (24" & 26" Only) <input type="checkbox"/> Yellow (24" & 26" Only) <input type="checkbox"/> Blue (24" & 26" Only)		

1. Not available with HD package.
2. Must be used with spoke wheels.
3. Not available with HD package. Must use 25" spoke wheels.

HANDRIMS

Description	Price
<i>Must Select Tire and Rim In Addition To Handrim. Handrim Must Match Rim Size</i>	
For Any Rims:	
<input type="checkbox"/> Anodized Aluminum	No Charge
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	
<input type="checkbox"/> Vinyl Coated Aluminum	\$200
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	
<input type="checkbox"/> Projection Handrim ¹	\$195
Style: <input type="checkbox"/> 45° Oblique <input type="checkbox"/> Vertical	
Size: <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	
For Spoke Rims Only:	
<input type="checkbox"/> Natural Fit Handrim	
Grip: <input type="checkbox"/> Standard Grip <input type="checkbox"/> Super Grip	\$345 \$395
Size: <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26"	
<input type="checkbox"/> 24" Natural Fit LT Handrim	
Grip: <input type="checkbox"/> Standard Grip <input type="checkbox"/> Super Grip	\$345 \$395
For Mag Rims Only:	
<input type="checkbox"/> Chrome Handrim	No Charge

1. Not available with spoke protectors.

FRONT SEAT-TO-FLOOR HEIGHT/CASTER CHARTS

S = Short Fork L = Long Fork F = Frog Legs


		Standard Seat-to-Floor Heights									
		15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
Caster Size	4"	S3	S2	S1/F3	L3/F2	L2/F1	L1				
	5"			S2	S1/L4	L3/F2	L2/F1	L1			
	6"					S1/L4	L3	L2/F1	L1		
	8"									L2	L1

		Ultra Hemi Seat-to-Floor Heights								
		12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"
Caster Size	4"	S2	S1/F3	L3/F2	L2/F1	L1				
	5"		S2	S1/L4	L3/F2	L2/F1	L1			
	6"				S1/L4	L3	L2/F1	L1		
	8"								L2	L1

Hole Configuration

Short Fork

Long Fork



CASTER FORKS

Description Price HCPCS Code

Fork:			
<input type="checkbox"/> Short Fork		No Charge	
<input type="checkbox"/> Long Fork		No Charge	
<input type="checkbox"/> 2" W Long Fork, Quick Release ¹		No Charge	
<input type="checkbox"/> Frog Legs ²		\$395	E1015
<input type="checkbox"/> Frog Legs Unitine Fork ³		\$249	E1015
Accessories:			
<input type="checkbox"/> Extended Stem Bolt For Increased Front STF Height ⁴		No Charge	
<i>Allows For A Maximum Dump Of 2.5" With Bolt-On Fork Only</i>			
<input type="checkbox"/> .75" Higher			
<input type="checkbox"/> 1.5" Higher			
<input type="checkbox"/> Caster Pin Locks For Bolt-On Fork		\$195	K0073
<input type="checkbox"/> Short Fork			
<input type="checkbox"/> Long Fork			

1. 2"W Tall Fork is required with 2"W Casters.
 2. If an HD package is chosen with Frog Legs you must choose a 6" Caster and can only achieve an 18" STF Height.
 3. For Frog Legs Unitine Fork Seat-To-Floor Heights use Short Fork holes 2 & 3 (S2 & S3). Not available with HD package.
 4. Only available through SOF. Not available with HD package.

CASTER WHEELS

Description Price

4" Casters¹:	
<input type="checkbox"/> Lighted Caster 1"W, 10 Spoke, Clear Acrylic Smooth, Flat Free	\$68
<input type="checkbox"/> Black Soft Roll 1.5"W, 3 Spoke Alum. Rim, Smooth Polyurethane	\$230
<input type="checkbox"/> Black Urethane 1.25"W, 3 Spoke Flat Free Smooth	\$65
<input type="checkbox"/> Colored Casters 1"W, 10 Spoke Acrylic, Smooth, Flat Free	\$85
Color: <input type="checkbox"/> Orange <input type="checkbox"/> Yellow <input type="checkbox"/> Green <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Clear	
5" Casters²:	
<input type="checkbox"/> Black Urethane 1.5"W, 3 Spoke Flat Free Smooth	No Charge
<input type="checkbox"/> Black Soft Roll 1.5"W, 3 Spoke Alum. Rim, Smooth Polyurethane	\$230
<input type="checkbox"/> Gray Flat Free Aluminum Wheel, Smooth Flat Free Tire	\$65

1. Not available with HD package.
 2. 5" casters can not be used with 90 degree legrests.

CASTER WHEELS (CONT'D)

Description Price

6" Casters³:	
<input type="checkbox"/> Black Urethane 1.5"W, 3 Spoke Flat Free Smooth	No Charge
<input type="checkbox"/> Black Soft Roll 1.5"W, 3 Spoke Alum. Rim, Smooth Polyurethane	\$230
<input type="checkbox"/> Gray Pneumatic ¹ 1.25"W, 3 Spoke, Treaded Tire	\$85
8" Casters⁴:	
<input type="checkbox"/> Black Urethane 1.5"W, 3 Spoke Flat Free Smooth	No Charge
<input type="checkbox"/> Gray Pneumatic ^{1,5,2} W, 3 Spoke, Threaded Tire	\$85

3. 6" casters can not be used with 90 degree legrests. If an HD package is chosen with Frog Legs you must choose a 6" caster and can only achieve an 18" STF height.
 4. 8" casters can not be used with 90 degree or 80 degree legrests.
 5. 2"W caster wheels must be used with 2"W Tall Fork.

AXLE SELECTION

Description Price

<input type="checkbox"/> Quick Release	No Charge
<input type="checkbox"/> Quad Release ¹	\$90

1. Not available with HD package. Cosmetic cover will be removed w/ Mag wheel configuration.

BI-DIRECTIONAL FOOT RIGGINGS

Description Price HCPCS Code

<i>Refer To Charts On Page 5 For Maximum Available Leg Length Chart</i>		
<input type="checkbox"/> Aluminum Swing-Away Legrests	No Charge	
Style: <input type="checkbox"/> 60° <input type="checkbox"/> 70° <input type="checkbox"/> 70° Tapered ¹ <input type="checkbox"/> 80° ² <input type="checkbox"/> 90° ³		
<input type="checkbox"/> Elevating Legrests ⁴	\$295	E0990
Calf pad: <input type="checkbox"/> Standard <input type="checkbox"/> Offset ¹ <input type="checkbox"/> Depth Adjustable		
Extension: <input type="checkbox"/> 13.25"-17"		
Set Length at ⁵ : <input type="checkbox"/> 13.25" <input type="checkbox"/> 14" <input type="checkbox"/> 14.75" <input type="checkbox"/> 15.5" <input type="checkbox"/> 16.25" <input type="checkbox"/> 17"		
Extension: <input type="checkbox"/> 15.25"-19"		
Set Length at ⁵ : <input type="checkbox"/> 15.25" <input type="checkbox"/> 16" <input type="checkbox"/> 16.75" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18.25" <input type="checkbox"/> 19"		
Extension: <input type="checkbox"/> 18.25"-22"		
Set Length at ⁵ : <input type="checkbox"/> 18.25" <input type="checkbox"/> 19" <input type="checkbox"/> 19.75" <input type="checkbox"/> 20.5" <input type="checkbox"/> 21.25" <input type="checkbox"/> 22"		
Ultra Hemi:		
<input type="checkbox"/> Ultra Hemi Aluminum Swing-Aways	No Charge	
11.75"-14.75" Leg Length		
Style: <input type="checkbox"/> 60° <input type="checkbox"/> 70° ²		
Set Length at ⁵ : <input type="checkbox"/> 11.75" <input type="checkbox"/> 12.5" <input type="checkbox"/> 13.25" <input type="checkbox"/> 14" <input type="checkbox"/> 14.75"		
<input type="checkbox"/> High Pivot Elevating Legrests	\$295	E0990
11.75"-14.75" Leg Length		
Calf pad: <input type="checkbox"/> Standard <input type="checkbox"/> Offset ¹ <input type="checkbox"/> Depth Adjustable ¹		
Set Length at ⁵ : <input type="checkbox"/> 11.75" <input type="checkbox"/> 12.5" <input type="checkbox"/> 13.25" <input type="checkbox"/> 14" <input type="checkbox"/> 14.75"		

1. Only available on seats 16"W or greater. N/A in artic white.
 2. Not available with 8" casters.
 3. Can not be used with 5", 6", or 8" casters.
 4. Not available on 14"W seat. Not available with HD package. Not available with occupied transit kit.
 5. Measurement is taken from top of frame.

FOOTPLATES

Description Price HCPCS Code

<i>See Chart On Page 5 For Footplate Compatibility For Your Configuration</i>		
<input type="checkbox"/> Composite Footplate, Comfort Adjustable	No Charge	
Size: <input type="checkbox"/> 4.25"W <input type="checkbox"/> 4.25"W, Narrow ¹ <input type="checkbox"/> 5"W <input type="checkbox"/> 6"W		
<input type="checkbox"/> Angle Adjustable Aluminum Footplates	\$185	K0040
Size: <input type="checkbox"/> 5.5"W <input type="checkbox"/> 6.5"W <input type="checkbox"/> 7.5"W		

Accessories:		
<input type="checkbox"/> Heel Loops	\$38	E0951

1. Required on 14"W seat. Not available with HD package.

SLING SEAT

Description Price

<input type="checkbox"/> Adjustable Tension Sailcloth	No Charge
-------------------------------------------------------	-----------

SEAT CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required	No Charge	
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Simplicity Cushion <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$98	E2602
Solution@ 1 Cushion <i>Skin Protection</i>		
<input type="checkbox"/> 14" x 14" - 20" x 18"	\$380	E2603
<input type="checkbox"/> 20" x 20" - 20" x 20"	\$430	E2603
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$430	E2604
Solution@ Cushion Skin <i>Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 18"	\$380	E2607
<input type="checkbox"/> 20" x 20" - 20" x 20"	\$430	E2607
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$430	E2608
Spectrum@ Gel Cushion <i>Positioning</i>		
<input type="checkbox"/> 12" x 10" - 20" x 19"	\$480	E2605
<input type="checkbox"/> 20" x 20" - 20" x 20"	\$530	E2605
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$530	E2606
Spectrum@ Visco Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 18"	\$480	E2607
<input type="checkbox"/> 20" x 20" - 20" x 20"	\$530	E2607
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$530	E2608
Spectrum@ Air <i>Adjustable Skin Protection</i>		
<input type="checkbox"/> 14" x 12" - 20" x 20"	\$480	E2622
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$530	E2623
Spectrum@ Air Contour <i>Adjustable Skin Protection & Positioning</i>		
<input type="checkbox"/> 14" x 12" - 20" x 20"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 22" x 20"	\$549	E2625

FOOTPLATE COMPATABILITY CHART

*Maximum Footplate Size Is Measured To Provide 3/4" Minimum Between Footplates
Chart Shows Maximum Extension Length For Given Footplate Size And Seat Width*

Chair Width	Legrest	4.25" W (Narrow) Composite & 5.5" W Aluminum	4.25" W Composite	5" W Composite & 6.5" W Aluminum	6" W Composite & 7.5" W Aluminum
		14"	Straight	4.25" W Narrow	N/A
	70° Tapered	N/A	N/A	N/A	N/A
15"	Straight	ALL	ALL	N/A	N/A
	70° Tapered	N/A	N/A	N/A	N/A
16"	Straight	ALL	ALL	N/A	N/A
	70° Tapered	14"	N/A	N/A	N/A
17"	Straight	ALL	ALL	ALL	N/A
	70° Tapered	17.5" ext.	14.75" ext.	N/A	N/A
18"	Straight	ALL	ALL	ALL	N/A
	70° Tapered	ALL	19.75" ext.	N/A	N/A
19"	Straight	ALL	ALL	ALL	ALL
	70° Tapered	ALL	ALL	16" ext.	N/A
20"	Straight	ALL	ALL	ALL	ALL
	70° Tapered	ALL	ALL	ALL	14"
21"	Straight	ALL	ALL	ALL	ALL
	70° Tapered	ALL	ALL	ALL	16.75" ext.
22" (HD Only)	Straight	N/A	ALL	ALL	ALL
	70° Tapered	N/A	ALL	ALL	ALL

MAXIMUM LEG LENGTH CHARTS

Maximum Leg Length is measured with 3/4" minimum ground clearance with footplates set at 90° to legrest

Composite Footplates					
Front STF	Rear STF	60 Degree	70 Degree	80 Degree	90 Degree
15"	14.5"	16"	14.75"	14"	13.25"
	15"	15.5"	14.75"	14"	13.25"
15.5"	14.5"	16.25"	15.25"	14"	14"
	15"	16.25"	15.25"	14"	14"
	15.5"	16.25"	14.75"	14"	14"
16"	14.5"	17"	16"	15.25"	14.75"
	15"	16.75"	15.5"	14.75"	14"
	15.5"	16.75"	15.5"	14.75"	14"
	16"	16.75"	15.5"	14.75"	14"
16.5"	14.5"	18.25"	16.75"	16"	15.25"
	15"	17.5"	16.25"	15.5"	14.75"
	15.5"	17.5"	16.25"	15.5"	14.75"
	16"	17.5"	16.25"	15.5"	14.75"
	16.5"	16.75"	16"	15.25"	14.75"
17"	14.5"	19"	17.5"	16.25"	15.5"
	15"	19"	16.75"	16.25"	15.5"
	15.5"	18.25"	16.75"	16"	15.5"
	16"	18.25"	16.75"	16"	15.5"
	16.5"	17.5"	16.25"	16"	15.5"
	17"	17.5"	16.25"	15.5"	15.25"
17.5"	15"	19.75"	17.5"	16.75"	16"
	15.5"	19"	17.5"	16.75"	16"
	16"	19"	16.75"	16.25"	16"
	16.5"	18.25"	16"	16.25"	15.5"
	17"	18.25"	16"	16.25"	15.5"
	17.5"	18.25"	16"	16"	15.5"
18"	15.5"	19.75"	17.5"	17"	N/A
	16"	19.75"	17.5"	17"	N/A
	16.5"	19"	17.5"	17"	N/A
	17"	19"	17.5"	17"	N/A
	17.5"	18.25"	17.5"	16.75"	N/A
	18"	18.25"	16.75"	16.25"	N/A
18.5"	16"	20.5"	19"	17.5"	N/A
	16.5"	20.5"	19"	17.5"	N/A
	17"	19.75"	19"	17.5"	N/A
	17.5"	19.75"	17.5"	17.5"	N/A
	18"	19"	17.5"	17"	N/A
	18.5"	19"	17.5"	16.75"	N/A
19"	16.5"	21.25"	19"	N/A	N/A
	17"	20.5"	19"	N/A	N/A
	17.5"	19.75"	18.25"	N/A	N/A
	18"	19.75"	18.25"	N/A	N/A
	18.5"	19"	18.25"	N/A	N/A
	19"	19"	17.5"	N/A	N/A
19.5"	17"	21.25"	19.75"	N/A	N/A
	17.5"	20.5"	19.75"	N/A	N/A
	18"	20.5"	19"	N/A	N/A
	18.5"	20.5"	19"	N/A	N/A
	19"	19.75"	18.25"	N/A	N/A

Aluminum Footplates ²					
Front STF	Rear STF	60 Degree	70 Degree	80 Degree	90 Degree
15"	14.5	16.75	15.5	14	13.25
	15	16.25	15.5	14	13.25
15.5"	14.5	17.5	16.25	14	14
	15	17	16.25	14	14
	15.5	16.75	16	14	14
16"	14.5	18.25	16.75	15.25	14.75
	15	17.5	16.25	14.75	14
	15.5	17.5	16.25	14.75	14
	16	17.5	16	14.75	14
16.5"	14.5	19	17.5	16	15.25
	15	19	17.5	15.5	14.75
	15.5	18.25	16.75	15.5	14.75
	16	18.25	16.75	15.5	14.75
	16.5	17.5	16.75	15.25	14.75
17"	14.5	19.75	18.25	16.25	15.5
	15	19.75	18.25	16.25	15.5
	15.5	19.75	18.25	16	15.5
	16	19	17.5	16	15.5
	16.5	19	17.5	16	15.5
	17	18.25	17.5	15.5	15.25
17.5"	15	20.5	19	16.75	16
	15.5	20.5	19	16.75	16
	16	19.75	18.25	16.25	16
	16.5	19.75	18.25	16.25	15.5
	17	19	18.25	16.25	15.5
	17.5	19	17.5	16	15.5
18"	15.5	21.25	19.75	17	N/A
	16	20.5	19	17	N/A
	16.5	20.5	19	17	N/A
	17	20.5	19	17	N/A
	17.5	19.75	18.25	16.75	N/A
	18	19	18.25	16.25	N/A
18.5"	16	22	20.5	17.5	N/A
	16.5	21.25	19.75	17.5	N/A
	17	21.25	19.75	17.5	N/A
	17.5	20.5	19	17.5	N/A
	18	19.75	19	17	N/A
	18.5	19.75	18.25	16.75	N/A
19"	16.5	22	20.5	N/A	N/A
	17	22	20.5	N/A	N/A
	17.5	21.25	19.75	N/A	N/A
	18	20.5	19.75	N/A	N/A
	18.5	20.5	19	N/A	N/A
	19	19.75	19	N/A	N/A
19.5"	17	22	21.25	N/A	N/A
	17.5	22	20.5	N/A	N/A
	18	22	20.5	N/A	N/A
	18.5	21.25	20.5	N/A	N/A
	19	21.25	19.75	N/A	N/A

Ultra Hemi Frame					
Maximum Legrest Length					
Front STF	Rear STF	Composite Footplates		Angle Adjustable Footplates	
		60 degree	70 degree	60 degree	70 degree
12.5	12.5	12.5	N/A	12.5	N/A
13"	12.5	13.25	11.75	13.25	11.75
	13	13.25	11.75	13.25	11.75
13.5"	12.5	14	12.5	13.25	12.5
	13	14	12.5	13.25	12.5
	13.5	13.25	11.75	13.25	12.5
14"	12.5	14.75	13.25	14.75	13.25
	13	14.75	13.25	14	13.25
	13.5	14	13.25	14	13.25
	14	14	13.25	14	13.25
14.5"	12.5	14.75	14.75	14.75	14
	13	14.75	14.75	14.75	14
	13.5	14.75	14.75	14.75	14
	14	14.75	14.75	14.75	14
15"	12.5	14.75	14.75	14.75	14.75
	13	14.75	14.75	14.75	14.75
	13.5	14.75	14.75	14.75	14.75
	14	14.75	14.75	14.75	14.75
15.5"	14.5	14.75	14	14.75	14
	15	14.75	14	14.75	14
	15.5	14.75	14	14.75	14.75
	16	14.75	14	14.75	14.75
16"	13.5	14.75	N/A	14.75	N/A
	14	14.75	N/A	14.75	N/A
	14.5	14.75	N/A	14.75	N/A
	15	14.75	N/A	14.75	N/A
	16	14.75	N/A	14.75	N/A
16.5"	14	14.75	N/A	14.75	N/A
	14.5	14.75	N/A	14.75	N/A
	15	14.75	N/A	14.75	N/A
	15.5	14.75	N/A	14.75	N/A
	16	14.75	N/A	14.75	N/A
	16.5	14.75	N/A	14.75	N/A

1. Leg Length measured from top of seat pan to footplate.
2. Footplate in middle position.

ACCESSORIES

NOTES:

<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> Adjustable Rear Anti-Tip Wheels	\$105	E0971
<input type="checkbox"/> Removable Plastic Side Guards, Small ¹	\$90	
<input type="checkbox"/> Removable Aluminum Side Guards, Large ¹	\$120	
<input type="checkbox"/> Lap Belt ²	\$50	E0978
<input type="checkbox"/> Spinergy Spoke Protectors ³	\$90	
Size: <input type="checkbox"/> 24" <input type="checkbox"/> 26"		
<input type="checkbox"/> Gel Arm Covers For 12" Tubular Arms	\$99	
<input type="checkbox"/> Case Logic Under Seat Organizer	\$85	
<input type="checkbox"/> Touch-up Paint	\$25	
<input type="checkbox"/> 7-Piece Tool Kit	\$25	
<input type="checkbox"/> Combat Cane/Crutch Holder ¹	\$95	E2207
<input type="checkbox"/> Backpack w/ Quantum Logo	\$60	

- 1. Not available with camber.
- 2. Not available with occupied transit kit.
- 3. Not available with projection hand rims.

TRANSIT SECUREMENT PACKAGE

<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> Occupied ¹	\$315	
<input type="checkbox"/> with anti-tips	\$105	E0971
<input type="checkbox"/> without anti-tips	No Charge	
<input type="checkbox"/> Unoccupied	\$275	
<input type="checkbox"/> with anti-tips	\$105	E0971
<input type="checkbox"/> without anti-tips	No Charge	

- 1. Not available with HD or Ultra Hemi packages.

THIS SECTION MUST BE READ

The Transit Securement Package, when used in conjunction with a wheelchair was tested in the specific configuration of a 16"x16" seat size, with 15"-19" straight back canes, desk-length height-adjustable armrests, padded back sling upholstery, and 70-degree bi-directional legrests, conforms with the standard of ANSI/RESNA WC/ Vol. 4, Section 19, for the purpose of enabling wheelchairs to be secured in certain types of motor vehicles. The Transit Securement Package, including manufacturer-installed front and rear securement brackets and pelvic belt anchoring brackets, has been crash tested in accordance with ANSI/RESNA WC/Vol. 4, Section 19, Frontal Impact Test, with a 168-lbs. (76.2 kg) surrogate occupant corresponding to a user weight range of 115 lbs. (52.16 kg) to 210 lbs. (95.25 kg). While the Transit Securement Package is compatible with additional seat sizes, optional components, and varied configurations, and provides designated securement points, Pride Mobility makes no claims that any other components or configurations have been tested beyond the above described ANSI/RESNA conforming configuration. Please refer to the Transit Securement Package supplement (part number INFMANU3659) for additional information. **However, the transit kit can be purchased with any size chair.**