

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP. Options with a circular checkbox (○) are NOT available through Quantum West.

BASE MODEL

Description	Price	HCPCS Code
☐ Quantum Q6400Z 4S-SS	\$9,364	K0868
☐ Quantum Q6400Z 4S-C	\$9,364	K0869

COLORS**Metallic Colors:**

- ☐ Metallic Red ☐ Metallic Black
☐ Metallic Blue ☐ Metallic Silver
☐ Sunburst Orange

DRIVE WHEELS

Description	Price
☐ 14" Flat-Free (pair)	Standard
☐ 14" Pneumatic (pair)	No Charge

MOTOR OPTIONS

Description	Price	HCPCS Code
☐ High Speed Motor Package	\$1,574	
<i>Programmable Speed Up To 8 MPH, 300 lbs. Weight Capacity. HS Motors decrease per-charge range.</i>		

BASE MOUNTED FOOT PLATFORM

Description	Price
☐ Base Mounted Foot Platform Installed	Standard
☐ No Base Mounted Foot Platform	No Charge

BATTERY CHARGER

Description	Price
☐ 8 Amp Off-board Charger	Standard

TRANSIT KIT

Description	Price
<i>Refer to last page for Details Regarding Occupied and Unoccupied Transit Kits</i>	
☐ Occupied Transit Kit	No Charge
☐ Unoccupied Transit Kit ¹	No Charge
1. Unoccupied Transit Kit with a Vent Tray may not allow lowest seat to floor height, vent tray unavailable w/ Occupied Transit Kit.	

SUSPENSION OPTIONS

Description	Price
☐ Standard Suspension	Standard
☐ Pediatric Strut Suspension	No Charge
<i>For User Weights Less Than 120lbs</i>	

LED LIGHT PACKAGE

Description	Price
☐ Through Q-Logic EX or NE+ 8-Key w/ Lights	\$995

NON-EXPANDABLE CONTROLS

Description	Price
☐ 75 Amp NE	Standard
☐ 75 Amp NE+ 6-Key For Up To 2 Actuators Thru Joystick	\$200 <i>Standard on 3SP-C Package</i>
☐ 75 Amp NE+ 8-Key For Lights	\$400

EXPANDABLE CONTROLS

Description	Price	HCPCS Code
<i>Q-Logic EX Joystick or Enhanced Display must be ordered with a Q-Logic EX Controller and Harness For Expandable Controls</i>		
☐ 100 Amp Q-Logic EX Controller	\$695	E2377
☐ Expandable Controller Harness	\$350	E2313
☐ Q-Logic EX Joystick	\$700	
☐ Q-Logic EX Joystick for Lights	\$700	

PROGRAMMING

Description	Price
☐ Q-Logic Handheld Programmer	\$545
☐ Q-Logic PC Programming Cable with Disk ¹	\$750

1. Works with Q-Logic EX only.

JOYSTICK HANDLES

Description	Price	HCPCS Code
☐ Standard Carrot Handle	Standard	
☐ Mushroom Handle	\$71	E2323
☐ T-Bar	\$107	E2323
☐ Chin Cup	\$71	E2324 ¹
☐ Extended I-Handle	\$71	E2323
☐ Soft Ball Handle	\$50	E2323
☐ Medium Ball Handle	\$50	E2323
☐ Bodypoint 3" Goal Post Handle	\$69	E2323
☐ Bodypoint 4" Goal Post Handle	\$69	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

Options with a circular checkbox (○) are NOT available through Quantum West.

JOYSTICK MOUNTING BRACKETS

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
☐ Fixed Inline ¹	Standard	
Side: ☐ Left ☐ Right		
☐ Swing-away	\$210	E1028
Side: ☐ Left ☐ Right		
Style: ☐ Inline ☐ Offset		
Length: ☐ 4" Ext. ☐ 7" Ext. ²		
☐ Multi-axis Fixed	\$100	K0108
Side: ☐ Left ☐ Right		

1. Not available with Q-Logic EX Stand Alone Joystick
2. In iQuantum choose in addition to Swing-Away 4" mount.

CONTOURED SEATING

Description	Price	HCPCS Code
Depth Adjustable High Back:		
☐ Limited Recline Depth Adjustable	No Charge	
High Back w/ Headrest ¹		
Color: ☐ Black Vinyl		
Size: ☐ 18"W x 18-20"D ☐ 20"W x 18-20"D ☐ 22"W x 20-22"D		
High Back with Solid Seat Pan:		
☐ Limited Recline High-Back with Headrest & Solid Seat Pan	No Charge	
Color: ☐ Black Vinyl		
Size: ☐ 16"W x 16"D ☐ 16"W x 18"D ☐ 16"W x 20"D ☐ 16"W x 22"D		
☐ 18"W x 16"D ☐ 18"W x 18"D ☐ 18"W x 20"D ☐ 18"W x 22"D		
☐ 20"W x 16"D ☐ 20"W x 18"D ☐ 20"W x 20"D ☐ 20"W x 22"D		
☐ Extended Seat Mounts <i>Increase STF Height 3"</i>	No Charge	
Armrests:		
☐ Flip-Up, Removable	No Charge	
☐ Flip-Up, Height-Adjustable, Removable	\$85 ea.	E0973
Armrests:		
Left: ☐ Full Length ☐ Desk Length		
Right: ☐ Full Length ☐ Desk Length		

Special Order Seat

____ "W x ____ "D

Seat size can not exceed 24"W x 24"D

Contact your inside sales specialist for pricing

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part

MANUAL RECLINER

Description	Price	HCPCS Code
☐ Manual Recliner High Back with Headrest	\$600	E1226
☐ Manual Recliner High Back with Headrest and Solid Seat Pan ¹	\$600	E1226
Upholstery:		
☐ Gray Vinyl		
Seat Width:		
☐ 16" ☐ 18" ☐ 20"		
Seat Depth:		
☐ 16" ☐ 18" ☐ 20"		
Armrests:		
Left: ☐ Full Length ☐ Desk Length		
Right: ☐ Full Length ☐ Desk Length		

1. Includes black vinyl seatpan. Back available in listed colors.

SPECIALTY SEATING

Description	Price	HCPCS Code
☐ Specialty Seating Kit	No Charge	
Seat Width:		
☐ 16" ☐ 18" ☐ 20"		
Seat Depth:		
☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20"		
Back Cane Height:		
☐ 18" ☐ 20" ☐ 22"		
Back Rest Height:		
☐ 18" ☐ 20" ☐ 22"		
Armrests:		
☐ 2-Post, Fixed, Flip Back ¹	No Charge	
☐ Flip-Back, Cane Mounted ¹	\$95 ea.	
☐ Removable, Adjustable, Quick Height ²	\$275 ea.	E0973
☐ Removable, Single-Post, Fixed Height Adjustable Heavy Duty	\$150 ea.	E0973
☐ 8-12" Adjustment ☐ 10-14" Adjustment		
Armrests:		
Left: ☐ Full Length ☐ Desk Length		
Right: ☐ Full Length ☐ Desk Length		

1. Swing-away Laterals Interfere with flip-back arms. Call for configuration with Shape back.
2. Adjustable from 11" to 17"H in 1/2" increments.

SYNERGY® SEATING SYSTEMS

Description	Price	HCPCS Code
☐ Synergy Seating Frame	\$434	
☐ Standard Stem Mounts		
☐ 2" Extended Vertical Stem Mounts ¹		
☐ Synergy Manual Recline Package ²	\$700	E1226
☐ Synergy Manual Tilt Package ³	\$550	K0108
<i>Max. Seat Size 14"W x 14"D, 75 lbs Weight Capacity</i>		
Seat Width:		
☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"		
☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"		
☐ 20" ☐ 21" ☐ 22" ☐ 23" ☐ 24"		
Seat Width Package:		
☐ 10"W-14"W ☐ 14"W-18"W ☐ 16"W-20"W ☐ 20"W-24"W		
Seat Width Package (Power Elevate):		
☐ 14"W-18"W ☐ 18"W-22"W		
Seat Depth:		
☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"		
☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"		
☐ 20" ☐ 21" ☐ 22" ☐ 23" ☐ 24"		
Seat Depth Package:		
☐ 10"D-14"D ☐ 14"D-18"D ☐ 15"D-20"D ☐ 20"D-24"D		
Back Cane Height:		
☐ 14" ☐ 15" ☐ 16" ☐ 17"		
☐ 18" ☐ 19" ☐ 20" ☐ 21"		
☐ 22" ☐ 23" ☐ 24" ☐ 25"		

1. Increases seat to floor height by 2". Lowest seat-to-floor height not available.
2. Not available with Flip Back Arms, Cantilever Arms, Shape Back or Power Elevate. Requires 22"H canes. Pivot point 2.5" above seat pan.
3. In iQuantum this option is chosen under "Width Package Type"

SYNERGY® BACK OPTIONS

Description	Price	HCPCS Code
☐ Synergy® Shape Back <i>Positioning Back</i>	Various	E2613/E2614
<i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>		
☐ Sling Back	\$98	
☐ Omit Back	No Charge	
☐ Synergy® Back with 1.5" Thick Back Cushion	\$275	
Overall Back: ☐ 18" H ☐ 19" H ☐ 20" H ☐ 21" H		
Height: ☐ 22" H ☐ 23" H ☐ 24" H ☐ 25" H		
<i>Choose Lateral Extensions¹:</i>		
☐ 6" for 13"-16"H Back	\$25	
☐ 10" for 17"-18"H Back	\$25	
☐ 12" For 19"-22"H Back	\$25	
☐ 16" For 23"-25"H Back	\$25	

1. Lateral extensions are required to mount lateral thoracic supports.

Options with a circular checkbox (○) are NOT available through Quantum West.

SYNERGY® BACK OPTIONS

Description	Price	HCPCS Code
☐ TRU-Comfort™ Back ¹ 4-way Stretch Material Cover		
Width: ☐ 16"W ☐ 18" W ☐ 20"W	\$595	E2620
Overall Back Height ² : ☐ 20" H ☐ 21" H ☐ 22"H ☐ 23"H		
☐ 24"H ☐ 25"H ☐ 27"H ☐ 28"H		
Width: ☐ 17"W ☐ 19"W	\$595	E2620
Overall Back Height ² : ☐ 22"H ☐ 23"H ☐ 24"H ☐ 25"H		
☐ 27"H ☐ 28"H		
Width: ☐ 22"W ☐ 24"W ³	\$695	E2621
Overall Back Height ² : ☐ 24"H ☐ 25"H ☐ 27"H ☐ 28"H		

Tru Comfort Color Shroud Option:

☐ Black (Standard color option)

Metallic Colors:\$100

- ☐ Red
☐ Blue
☐ Sunburst Orange
☐ Silver

Synergy® Back Cane Options:

- ☐ Standard Back Canes Standard
☐ Omit Push Handles \$125
☐ Flip Forward Back Canes⁴ \$150
☐ Omit Push Handles \$125
☐ Offset (Meso)⁵ \$195
☐ 2" (1" Per Side) ☐ 4" (2" Per Side)
☐ Inset (Endo) \$195
☐ 2" (1" Per Side) ☐ 4" (2" Per Side)

1. Not available on Specialty Seating Kit.
 2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
 3. Available with offset (meso) brackets.
 4. Not available with Vent Tray, O2 Holder, cane mounted or cantilever armrests. Also not available w/ a Synergy manual recline or manual tilt. May interfere with Lateral Supports and Waterfall Armpads.
 5. Dual Post Flip-back Armrests interfere with offset canes.

SYNERGY® ARMRESTS

Description	Price	HCPCS Code
☐ Flip-back, Cane Mounted ¹	No Charge	
☐ 2-post, Flip-back, Height-adjustable, Removable ^{2,3}	\$150 ea.	E0973
☐ Single-post, Height-adjustable, Removable ^{3,4}	\$75 ea.	E0973
☐ Quad Release Lever		
☐ Removable, Single-post, Height Adjustable, Heavy Duty	\$75 ea.	E0973
☐ 8"H - 12"H ☐ 10"H-14"H		
☐ Pediatric, Removable ⁵	\$75 ea.	E0973
☐ Quad Release Lever		
☐ Bariatric, Removable, Height Adjustable ⁶	\$75 ea.	E0973
☐ 8"H - 12"H ☐ 10"H-14"H		
☐ Height Adjustable Cantilever Armrests ⁷	\$195	K0020
Angle-Adjustable: -30° through +30°		

1. Interferes with swing-away laterals. Not available on flip forward canes or manual recline.
 2. Interferes with swing-away laterals. Height adjustable from 9"-13" in 1" increments.
 3. Not available with a depth less than 16".
 4. Height adjustable from 11"-17" in 1/2" increments
 5. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armpads not available
 6. Bariatric Armrests require Bariatric Armpads
 7. Height adjustment from 11-17" depends on cane height. Not available with Manual Recline, Flip Forward back canes & Shape Back.

SYNERGY® ARMPADS

Description	Price
☐ Straight, Full Length (2"x14")	Standard
☐ Left ☐ Right	
☐ Straight, Desk Length (2"x10")	No Charge
☐ Left ☐ Right	
☐ Waterfall, Full Length (2"x14") ¹	No Charge
☐ Left ☐ Right	
☐ Waterfall, Desk Length (2"x10") ¹	No Charge
☐ Left ☐ Right	
☐ Falcon Gel, Full Length, (3.5"x14")*	\$160 pr.
☐ Left ☐ Right	
☐ Full Length Gel Armpad Mounting Plate (Required)*	\$56.69 ea.
☐ Falcon Gel, Desk Length, (3.5"x10")*	\$125 pr.
☐ Left ☐ Right	
☐ Desk Length Gel Armpad Mounting Plate (Required)*	\$53.52 ea.
☐ Flat, Bariatric ²	No Charge
☐ Left ☐ Right	

1. Not available w/ Cane Mounted. Flip Back Armrests combined with a Synergy Back.

2. Bariatric Armrests require Bariatric Armpads

HEADRESTS

Description	Price	HCPCS Code
☐ No Headrest Required		
TRU-Comfort™ Headrest¹:		
☐ 8" W Leatherette	\$205	E0955
☐ 8"W TRU-Comfort Plus	\$205	E0955
☐ 10" W Leatherette	\$205	E0955
☐ 10"W TRU-Comfort Plus	\$205	E0955
☐ Removable Headrest Bracket	\$210	E1028
Stealth:		
☐ Small (9.75"W x 3"H) (STL-OTSM250)	\$205	E0955
☐ Large (11"W x 6"H) (STL-OTLG150)	\$205	E0955
☐ Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
☐ Small Comfort Plus (10"W) (CP250)	\$142	E0955
☐ Large Comfort Plus (14"W) (CP150)	\$142	E0955
☐ Removable Headrest Bracket ² (HMO475P)	\$210	E1028
☐ Swing-away (Flip-down) Headrest Bracket (FDM350)	\$210	
Used in Addition to Removable Headrest Bracket		
Aftermarket Bracket:		
☐ Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028

1. Measurement is taken from the inside of the headrest.

NON-POWERED LEGRESTS

Description	Price	HCPCS Code
☐ 60-90° Swing-away Legrests¹	Standard	
(Adjustable in 5° increments)		
☐ Standard Extruded Footplates	No Charge	
Standard Extruded Extensions:		
☐ 13"-16.5" ☐ 15.5"-19.5"		
Standard Extruded Footplates:		
☐ 7.25" x 6" ☐ 8.5" x 6"		
☐ Angle-Adjustable Footplates ²	\$150	K0040
Angle Adjustable Extensions:		
☐ 13"-16.5" ☐ 15.5"-19.5"		
Angle-Adjustable Footplates:		
☐ 4" x 6" ☐ 5" x 6" ☐ 5" x 8"		
☐ 6" x 8" ³ ☐ 6" x 11" ³ ☐ 8" x 11" ⁴		
☐ Multi-Axis Footplates	\$495	K0040
Multi-Axis Extensions:		
☐ 15.5"-19.5"		
Multi-Axis Footplates:		
☐ 4" x 6" ☐ 5" x 6" ☐ 5" x 8"		
☐ 6" x 8" ³ ☐ 6" x 11" ³ ☐ 8" x 11" ⁴		

1. Possible caster wheel interference when adjusted to less than 70 degrees.

2. Includes Heel Loops

3. Causes interference on seat widths less than 18".

4. Causes interference on seat widths less than 20".

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

Options with a circular checkbox (☉) are NOT available through Quantum West.


NON-POWERED LEGRESTS (CONT'D)

Description	Price	HCPCS Code
☐ High Mount, Clamp-On Footrest*		
☐ Angle Adjustable		K0040
☐ Small, 4" x 7.75"	\$356	
☐ Large, 4.25" x 9.25"	\$356	
☐ Multi-Axis		K0040
☐ Small, 4" x 7.75"	\$456	
☐ Large, 4.25" x 9.25"	\$460	
☐ Heavy Duty Drop-In Legrests	\$420	
☐ Angle Adjustable Footplates ²	\$150	K0040
<i>Angle Adjustable Extension Length:</i>		
☐ 12"-14.5"	☐ 14.5"-18"	
<i>Angle Adjustable Footplate Size:</i>		
☐ 4" x 6"	☐ 5" x 6"	☐ 5" x 8"
☐ 6" x 8" ³	☐ 6" x 11" ³	☐ 8" x 11" ⁴
☐ Elevating Legrests¹	\$325	E0990
☐ Standard Footplates		
<i>Standard Extruded Extension Length:</i>		
☐ 11"-14"	☐ 14"-17"	
<i>Standard Extruded Footplate Size:</i>		
☐ 7.25" x 6"	☐ 8.5" x 6"	
☐ Angle Adjustable Footplates ²	\$150	K0040
<i>Angle Adjustable Extension Length:</i>		
☐ 12"-14.5"	☐ 14.5"-18"	
<i>Angle Adjustable Footplate Size:</i>		
☐ 4" x 6"	☐ 5" x 6"	☐ 5" x 8"
☐ 6" x 8" ³	☐ 6" x 11" ³	☐ 8" x 11" ⁴
☐ Calf Pad Adjustment		
☐ 1 Hole	☐ 3 Hole	☐ 6 Hole
☐ 90° Pediatric Swing-Away Legrests¹	\$162	
☐ Angle Adjustable Footplates ²		
<i>Angle Adjustable Extension Length:</i>		
☐ 7"-10"		
<i>Angle Adjustable Footplate Size:</i>		
☐ 4" x 6"	☐ 5" x 6"	

1. Extensions are measured from the seatpan to footplate.
 2. Includes heel loops.
 3. Causes interference on seat widths less than 18".
 4. Causes interference on seat widths less than 20".

ADJUSTABLE CENTER MOUNT FOOT PLATFORM


Description	Price
☐ Adjustable Center Mount Foot Platform	\$250
<i>Includes 10"W x 11"D Split Foot Platform</i>	
<i>Extensions:¹</i>	
☐ 8"-12" Ext.	☐ 10"-14" Ext.
☐ 10"-16" Ext.	☐ 13"-19" Ext. ²
<i>Calf Supports³:</i>	
☐ Calf Pads (pair)	
☐ Small (6"H x 6"W)	\$150
☐ Large (8"H x 8"W, Tapered to 6"W)	\$150
☐ Calf Panel (ea.)	
☐ Small (8"H x 18"W)	\$150
☐ Large (8"H x 24"W) ⁴	\$200
<i>Optional Modular Footplate Extension Components/Kits^{3,5}</i>	
☐ 4" Width Extension Kit	\$50
☐ 2" Bell Shaped, Length Extension Kit	\$50



1. Extension length is measured from top of the seat pan to footplate.
 2. Not available on Recline only.
 3. Choose in Rigging Accessory Section in iQuantum.
 4. Only available on 20"W or greater.
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference.

ADJUSTABLE CENTER MOUNT FOOT PLATFORM (CONT'D)

Description	Price
☐ 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75
☐ 6-Piece Extension Stock Kit (Includes all above pieces)	\$95
<i>Additional Optional Accessories⁶:</i>	
☐ Footplate Connector	\$49



6. Assembly required. May require seating configuration adjustments to avoid front caster interference.

ELECTRONICS FOR POWER FOOT RIGGINGS

Description	Price	HCPCS Code
☐ Single Actuator Through Joystick (AFP Only)	\$1,175	E2310 ¹
☐ Dual Actuator Through NE+ 6-Key (Power Legrests)	\$2,375	E2311 ²
☐ Multiple Actuator Through Q-Logic EX (Power Legrests)	\$2,375	E2311 ²

1. HCPCS used only with Power AFP.
 2. HCPCS used only with Power Legrests


POWER ELEVATING LEGRESTS

Description	Price	HCPCS Code
☐ 70° Articulating Swing-Away Legrests, 4" Articulation Range	\$2,395	E1010
☐ Quad Release Lever		
<i>Calf Pads¹:</i>		
☐ 6" x 5.35"	☐ 8" x 5.86"	
<i>Angle Adjustable Extension Length²:</i>		
☐ 11"-13"	☐ 14.5"-18"	
<i>Angle Adjustable Footplate Size:</i>		
☐ 4" x 6"	☐ 5" x 6"	☐ 5" x 8"
☐ 6" x 8" ³	☐ 6" x 11" ³	☐ 8" x 11" ⁴
<i>Operational Method:</i>		
☐ Independent	☐ Combined	
☐ 70° Articulating Drop-In Legrests	\$2,395	E1010
☐ 4" Range	☐ 6" Range	
<i>Calf Pads¹:</i>		
☐ 6" x 5.35"	☐ 8" x 5.86"	
<i>Angle Adjustable Extension Length²:</i>		
☐ 12"-14.5"	☐ 14.5"-18"	
<i>Angle Adjustable Footplate Size:</i>		
☐ 4" x 6"	☐ 5" x 6"	☐ 5" x 8"
☐ 6" x 8" ³	☐ 6" x 11" ³	☐ 8" x 11" ⁴
<i>Operational Method:</i>		
☐ Independent	☐ Combined	

1. Choose in Rigging Accessory Section in iQuantum.
 2. Extensions are measured from top of seat pan to footplate.
 3. Only available with seat width 18" or greater.
 4. Only available with seat width 20" or greater.

CENTER MOUNT ARTICULATING FOOT PLATFORM

Description	Price	HCPCS Code
☐ Center Mount Articulating Foot Platform	\$2895	K0108
<i>Includes 10"W x 11"D Split Foot Platform</i>		
<i>Extensions¹:</i>		
☐ 10"-16" Ext.	☐ 13"-19" Ext.	
<i>Calf Supports²:</i>		
☐ Calf Pads (pair)		
☐ Small (6"H x 6"W)		
☐ Large (8"H x 8"W, Tapered to 6"W)		
☐ Calf Panel (ea.)		
☐ Small (8"H x 18"W)	\$150	
☐ Large (8"H x 24"W) ³	\$200	
<i>Optional Modular Footplate Extension Components/Kits^{2,4}</i>		
☐ 4" Width Extension Kit	\$50	





1. Extension length is measured from top of the seat pan to footplate.
 2. Choose in Rigging Accessory Section in iQuantum.
 3. Only available on 20"W or greater.
 4. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

Options with a circular checkbox (○) are NOT available through Quantum West.

CENTER MOUNT ARTICULATING FOOT PLATFORM (CONT'D) LATERALS

Description	Price	Description	Price	HCPCS Code
☐ 2" Bell Shaped, Length Extension Kit	\$50	☐ Stealth Thoracic Laterals (pair)^{1*}		
		Bracket:		
☐ 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	☐ Swing-Away Rail Mounted	\$450	E1028
		For Use on Synergy Back and Tru-Comfort		
☐ 6-Piece Extension Stock Kit	\$95	Pad:		
(Includes all above pieces)		☐ Flat Pad, 3"x5", Leather-Like Cover	\$200	E0956
Additional Optional Accessories ¹ :		☐ Flat Pad, 4"x6", Leather-Like Cover	\$200	E0956
☐ Footplate Connector	\$49	☐ Curved Pad, 2"x3", 4-Way Stretch Cover	\$200	E0956
		☐ Curved Pad, 3"x5", 4-Way Stretch Cover	\$200	E0956
		☐ Curved Pad, 4"x6", 4-Way Stretch Cover	\$200	E0956
		☐ Stealth Pelvic/Thigh Laterals*		
		Side:		
		☐ Left Side		☐ Right Side
		Bracket:		
		☐ Removable Bracket (ea.)	\$250	E1028
		☐ Removable Bracket w/ 1" Offset (ea.)	\$325	E1028
		☐ Flip Down Bracket (ea.)	\$225	E1028
		Pad:		
		☐ Short Pad, 4"x4" (ea.)	\$50	E0956
		☐ Medium Pad, 4"x8" (ea.)	\$50	E0956
		☐ Long Pad, 4"x12" (ea.)	\$50	E0956
		Adapter:		
		☐ 2" Depth/Height Adapter Plate (ea.) ²	\$25	K0108

1. Assembly required. May require seating configuration adjustments to avoid front caster interference.

FRONT RIGGING ACCESSORIES

Description	Price	HCPCS Code
☐ Wheel Bumpers	\$25	
☐ Heel Loops (pair) ¹	\$38	E0951
☐ Gel Neoprene Footrest Sleeve (pair) ²	\$65	K0108
☐ ELR/ALR Adductor Buttons (pair)	\$125	K0108
☐ POSAlinc Knee Adductor for Swing-Away Legrests (pair)*	\$256	K0108
☐ POSAlinc Rail Mounted Knee Adductor (pair) *	\$287	K0108
☐ Residual Limb Support	\$300	E1020
☐ Left ☐ Right		
☐ Bodypoint Ankle Huggers (pair)*	\$134	K0108
☐ Small ☐ Medium		

1. Standard with Angle Adjustable or Multi-Axis footplates.

2. Only available with Swing-Away legrests.

SEAT CUSHIONS

Description	Price	HCPCS Code
☐ No Cushion Required		
Please refer to Synergy order forms for detailed pricing, sizes and codes.		
Simplicity Cushion General Use		
☐ 16" x 16" - 20" x 24"	\$98	E2601
☐ 22" x 16" - 24" x 24"	\$98	E2602
Solution® 1 Cushion Skin Protection:		
☐ 14" x 14" - 20" x 24"	\$380	E2603
☐ 22" x 16" - 24" x 24"	\$430	E2604
Solution® Cushion Skin Protection & Positioning		
☐ 10" x 10" - 20" x 24"	\$380	E2607
☐ 22" x 16" - 24" x 24"	\$430	E2608
Spectrum® Gel Cushion Positioning		
☐ 12" x 10" - 20" x 24"	\$480	E2605
☐ 22" x 16" - 24" x 24"	\$530	E2606
Spectrum® Visco Cushion Skin Protection & Positioning		
☐ 10" x 10" - 20" x 24"	\$480	E2607
☐ 22" x 16" - 24" x 24"	\$530	E2608
Spectrum® Air Adjustable Skin Protection		
☐ 14" x 12" - 20" x 24"	\$480	E2622
☐ 22" x 16" - 24" x 24"	\$530	E2623
Spectrum® Air Contour Adjustable Skin Protection & Positioning		
☐ 14" x 12" - 20" x 24"	\$499	E2624
☐ 22" x 16" - 24" x 24"	\$549	E2625
TRU-Comfort™ Seat Cushion General Use		
Includes Allure leather-like cover		
☐ 16" x 16" - 20" x 24"	\$380	E2601
☐ 22" x 16" - 22" x 24"	\$480	E2602
TRU-Comfort™ Plus Seat Cushion Skin Protection		
Includes 4-way stretch material cover		
☐ 16" x 16" - 20" x 24"	\$380	E2603
☐ 22" x 16" - 22" x 24"	\$480	E2604

1. Swing-away laterals interfere with flip-up armrests.

2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.

MISCELLANEOUS POSITIONING COMPONENTS

Description	Price	HCPCS Code
☐ Therafin Flip-down Knee Abductor Bracket*	\$172	E1028
Pad:		
☐ Small ☐ Medium ☐ Large	\$135	E0957
☐ Stealth Foot Boxes (pair)*	\$475	K0108
☐ Small ☐ Medium ☐ Large		
☐ Stealth Padded Shoe Holders (pair)*	\$166	K0108
☐ Small ☐ Medium ☐ X-Large		
☐ Otto Bock Arm Channels (ea.)*	\$104	E2209
Side:		
☐ Left ☐ Right		
Size:		
☐ Small ☐ Medium ☐ Large		
☐ Otto Bock Flat Hand Pad(ea.)*	\$52	
Side:		
☐ Left ☐ Right		
Size:		
☐ Medium ☐ Large		
Please see positioning components order form for additional types, sizes and details.		

ACCESSORIES

Description	Price	HCPCS Code
☐ Rear Basket ¹	\$103	
☐ Cane/Crutch Holder ¹	\$98	E2207
☐ Weather Cover, Medium, Black	\$158	
☐ Cup Holder	\$17 ²	
Side:		
☐ Left ☐ Right		
☐ Saddle Bag	\$28	
☐ Flag Holder	\$17	
☐ Quantum Backpack ³	\$60	
☐ Oxygen Holder, Contoured Seating ¹	\$103	E2208
☐ Oxygen Holder, Synergy Seat ⁴	\$150	E2208
☐ Walker Holder ⁴	\$92	K0108
☐ Tablet Holder ⁵	\$315	
☐ Deluxe Cell Phone/Media Holder ⁵	\$140	
☐ Cell Phone/Media Holder ⁵	No Charge	

1. Includes mounting hardware. Not available on Manual Recliner.

2. No charge on Synergy seating.

3. Backpack may interfere w/ other rear accessories. Not available w/ contoured seating.

4. Not available with Synergy manual recline or manual tilt. Includes mounting hardware.

5. Only available with Synergy or Specialty seating.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

Options with a circular checkbox (○) are NOT available through Quantum West.

LAPBELTS

Description	Price
☐ 38" Lap Belt ¹	Standard
☐ 50" Lap Belt ¹	Standard
☐ 60" Lap Belt	Standard
☐ 60" Non-Retractable (Sure-Lok) ²	Standard
☐ 70" Lap Belt	Standard
☐ 80" Lap Belt	Standard
☐ 90" Lap Belt	Standard
☐ Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width*	\$87
☐ Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width*	\$92
☐ Pelvic Positioning Belt, 4-Point Padded, Med, 9"-15" Hip Width*	\$112
☐ Pelvic Positioning Belt, 4-Point Padded, Lg, 13"-19" Hip Width*	\$116

1. Includes mounting hardware. Not available with Manual Recliner.
2. Only needed with the occupied transit kit.

BATTERY OPTIONS

Description	Price	HCPCS Code
☐ Group 24, 70 Amp, Gel (Universal)	\$352 ea.	E2363
☐ Group 24, 70 Amp, Gel (MK)	\$352 ea.	E2363
☐ NF-22, 50-Amp Batteries		
☐ Maintenance Free AGM (MK)	\$295 ea.	E2361
☐ Maintenance Free Gel Cell (Universal)	\$295 ea.	E2361
☐ Maintenance Free Gel Cell (MK)	\$295 ea.	E2361
Installation:		
☐ Batteries Installed ☐ Do Not Install		

COMPACT VENT TRAY

Description	Price
-------------	-------

The addition of a compact vent tray may increase the length & weight of your power chair.

☐ Compact Vent Tray	\$1,395
---	---------


Mounting Bracket: Bracket is pre-selected with TRU-Comfort Back

☐ Flush	☐ 3" Offset ¹	☐ 6" Offset
-----------------------------	--	---------------------------------

Battery Tray Selection:

☐ In Base Battery Tray Extension Kit		
☐ Include NF22 Battery		
☐ Maintenance Free AGM	\$255	E2361
☐ Maintenance Free Gel Cell	\$295	E2361

- Accommodates Pulmonetics LTV-style ventilators.
- Max Vent Dim. 3.25" H x 10.5" W x 13.5" / 20 lbs.
- Not Available with Shape Back, Sling Back, Manual Recline, Flip Down Back Canes or rear mounted accessories.



1. Not available with Flip-Up Armrests.

TRANSIT KIT INFO

Q6000Z

THIS SECTION MUST BE READ

The Transit Securement Package, when used in conjunction with a **wheelchair in the specific configuration of an 18"x 20" seat size with 25" high back canes, Full-length height-adjustable single post armrests, Synergy back 25"h, and style 8 swing away legrests and extruded footplates and heel-loops**, conforms with the standard of ANSI/RESNA WC/Vol. 4, Section 19, for the purpose of enabling wheelchairs to be secured in certain types of motor vehicles. The Transit Securement Package, including manufacturer-installed front and rear securement brackets and pelvic belt anchoring brackets, has been crash tested in accordance with ANSI/RESNA WC/Vol. 4, Section 19, Frontal Impact Test, with a 168-lbs. (76.2 kg) surrogate occupant corresponding to a user weight range of 115 lbs. (52.16 kg) to 210 lbs. (95.25 kg). While the Transit Securement Package is compatible with additional seat sizes, optional components, and varied configurations, and provides designated securement points, Pride Mobility makes no claims that any other components or configurations have been tested beyond the above described ANSI/RESNA conforming configuration. Please refer to the Transit Securement Package supplement (**part number INFMANU3658**) for additional information.

RULES:

- All options are available when unoccupied.

Notes when selecting occupied transit kit:

- It is recommended that usable back height be 2" lower than shoulder height.
- Available on Synergy Seat up to 20"W x 20"D
- Power AFP, ALRs & ELRs are available on Synergy Seating
- All legrests must have heel loops
- Synergy manual seating is N/A
- TRU-Balance Power Positioning seating, options and accessories are N/A
- TRU-Comfort seat cushion is N/A
- Contoured seating is N/A
- Specialty seating is N/A
- Vent option is N/A
- Omit Back option is N/A

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum