

# Tru Balance® Pediatric Power Tilt Order Form



Discount: \_\_\_\_\_  
Terms: \_\_\_\_\_  
Freight: \_\_\_\_\_  
Taken By: \_\_\_\_\_

Account Number: \_\_\_\_\_  
Provider Name: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
Date: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Ship To Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.*

## BASES: SYNERGY® ARMRESTS

- 250lb Weight Capacity. Max Dimensions 14"W x 18"D**  
 Quantum 610  
 Quantum 6000Z  
 Q6 Edge

## TRU BALANCE® PEDIATRIC POWER TILT

Description	Price	HCPCS Code
<input type="checkbox"/> Tru Balance® Pediatric Power Tilt List base serial # for off chair seating	\$5925	E1002

Description	Price
<input type="checkbox"/> Flip-back, Cane Mounted <sup>1</sup>	No Charge
<input type="checkbox"/> Pediatric, Removable <sup>2</sup> <input type="checkbox"/> Quad Release Lever	No Charge
<input type="checkbox"/> Height Adjustable Cantilever Armrests <sup>3</sup> Angle-Adjustable: -30° through +30°	\$195

1. Interferes with swing-away laterals. Not available on flip forward canes.  
 2. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armpads not available  
 3. Height adjustment from 11-17" depends on cane height. Not available with Flip Forward back canes or Shape Back.

S/N: \_\_\_\_\_

- Seat Width:**  
 11"     12"     13"     14"
- Seat Width Package:**  
 10"W-14"W
- Seat Depth:**  
 10"     11"     12"     13"     14"  
 15"     16"     17"     18"
- Seat Depth Package:**  
 10"D-14"D     14"D-18"D
- Back Cane Height:**  
 14"     15"     16"<sup>1</sup>     17"<sup>1</sup>  
 18"     19"     20"

1. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

## SYNERGY® BACK OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	No Charge	
Overall Back: <input type="checkbox"/> 12" H <input type="checkbox"/> 13" H <input type="checkbox"/> 14" H <input type="checkbox"/> 15" H		
Height: <input type="checkbox"/> 16" H <input type="checkbox"/> 17" H <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H		
Choose Lateral Extensions <sup>1</sup> :		
<input type="checkbox"/> 6" for 13"-16"H Back	\$26.25	
<input type="checkbox"/> 10" for 17"-18"H Back	\$26.25	
<input type="checkbox"/> 12" For 19"-22"H Back	\$26.25	
<input type="checkbox"/> 16" For 23"-25"H Back	\$26.25	

### Synergy® Back Cane Options:

- Standard Back Canes    Standard
- Omit Push Handles    \$125
- Flip Forward Back Canes<sup>2</sup>    \$150
- Omit Push Handles    \$125
- Offset (Meso)    \$195
- 2" (1" Per Side)     4" (2" Per Side)
- Inset (Endo)    \$195
- 2" (1" Per Side)     4" (2" Per Side)

1. Lateral extensions are required to mount lateral thoracic supports.  
 2. Not available with cane mounted or cantilever armrests. May interfere with Lateral Supports and Waterfall Armpads.

## SYNERGY® ARMPADS

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14") <input type="checkbox"/> Left <input type="checkbox"/> Right	Standard
<input type="checkbox"/> Straight, Desk Length (2"x10") <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Full Length (2"x14") <sup>1</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Desk Length (2"x10") <sup>1</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.

## JOYSTICK RECEIVER MOUNT

Description	Price
<input type="checkbox"/> Universal Joystick Receiver Mount <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Omit Universal Joystick Receiver Mount	No Charge

## ELECTRONICS FOR POWER OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Thru Toggle <input type="checkbox"/> Left Side <input type="checkbox"/> Right Side	Standard	
<input type="checkbox"/> Thru Joystick (Q-Logic 2 EX or Remote+ Only)	\$1,175	E2310
<input type="checkbox"/> Thru Single Pushbutton (VSI or Remote+ Only)	No Charge	
<input type="checkbox"/> Thru Dual Pushbutton (VSI or Remote+ Only)	No Charge	

*Available With NE, VSI, Remote+ or Q-Logic 2 Electronics*

<b>LAP BELTS</b>		
<i>Description</i>	<i>Price</i>	
<input type="checkbox"/> 38" Lap Belt <sup>1</sup>	Standard	
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width*	\$91	
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width	\$96	


1. Includes mounting hardware. Not available with Manual Recliner.

<b>HEADRESTS</b>		
<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> No Headrest Required		
<b>Stealth:</b>		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$64.20	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$151.94	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$151.94	E0955
<input type="checkbox"/> Removable Headrest Bracket <sup>2</sup> (HMO475P)	\$160.50	E1028
<b>Aftermarket Bracket:</b>		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$220.50	E1028

<b>NON-POWER LEGRESTS</b>		
<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> <b>90° Pediatric Swing-Away Legrests<sup>1</sup></b>	<b>Standard</b>	
<input type="checkbox"/> Angle Adjustable Footplates		
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 7"-10"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"	\$85 ea.	K0040
<input type="checkbox"/> <b>High Mount, Clamp-On Footplates*</b>		
<input type="checkbox"/> Angle Adjustable		K0040 <sup>2</sup>
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040 <sup>2</sup>
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	
<input type="checkbox"/> <b>60-90° Swing-away Legrests<sup>3</sup></b>	<b>\$162</b>	
<i>(Adjustable in 5° increments)</i>		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
<i>Standard Extruded Extensions:</i>		
<input type="checkbox"/> 13"-16.5"		
<i>Standard Extruded Footplates:</i>		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates <sup>4</sup>	\$85 ea.	K0040 <sup>2</sup>
<i>Angle Adjustable Extensions:</i>		
<input type="checkbox"/> 13"-16.5"		
<i>Angle-Adjustable Footplates:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> Multi-Axis Footplates	\$495	K0040 <sup>2</sup>
<i>Multi-Axis Extensions:</i>		
<input type="checkbox"/> 15.5"-19.5"		
<i>Multi-Axis Footplates:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		

1. Extensions are measured from the seatpan to footplate.  
 2. K0040 is not separately reimbursable on Group 2 units.  
 3. Possible caster wheel interference when adjusted to less than 70 degrees.  
 4. Includes Heel Loops

<b>ADJUSTABLE CENTER MOUNT FOOT PLATFORM</b>		
<i>Description</i>	<i>Price</i>	
<input type="checkbox"/> Adjustable Center Mount Foot Platform <sup>1</sup>	\$250	
<i>Footplates :</i>		
<input type="checkbox"/> 10"W x 11"D Split Footplates		
<input type="checkbox"/> 10"W x 8"D Split Footplates		
<i>Extensions:<sup>2</sup></i>		
<input type="checkbox"/> 8"-12" Ext. <input type="checkbox"/> 10"-14" Ext		
<i>Calf Supports<sup>3</sup>:</i>		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)	\$150	
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150	
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) <sup>4</sup>	\$200	
<i>Optional Modular Footplate Extension Components/Kits<sup>3,5</sup></i>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<i>Additional Optional Accessories<sup>5</sup>:</i>		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumpers, pair	\$12.50 ea.	



1. Not available on 614 or Select Series.  
 2. Extension length is measured from top of the seat pan to footplate.  
 3. Choose in Rigging Accessory Section in iQuantum.  
 4. Only available on 20"W or greater.  
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference. N/A with the 10"W x 8"D footplates.

<b>FOOT RIGGING ACCESSORIES</b>		
<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> Wheel Bumpers	\$12.50 ea	
<input type="checkbox"/> Heel Loops (pair) <sup>1</sup>	\$50	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) <sup>2</sup>	\$65	K0108

1. Standard with angle adjustable footplates.  
 2. Available only with swing-away legrests.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**LATERALS**

<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> <b>Stealth Thoracic Laterals (pair)*</b>		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away Rail Mounted	\$472.50	E1028
<i>For Use on Synergy Back</i>		
<i>Pad:</i>		
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> <b>Stealth Pelvic/Thigh Laterals*</b>		
<i>Side:</i>		
<input type="checkbox"/> Left Side		<input type="checkbox"/> Right Side
<i>Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$262.50	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$400	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$236.25	E1028
<i>Pad:</i>		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$75	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$75	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$75	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) <sup>2</sup>	\$26.25	K0108

1. Swing-away laterals interfere with flip-up armrests.
2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.

**SYNERGY® SEAT CUSHIONS**

<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
<b>Solution® Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 14" x 18"	\$430	E2607
<b>Spectrum® Visco Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 14" x 18"	\$530	E2607

**ACCESSORIES**

<i>Description</i>	<i>Price</i>
<input type="checkbox"/> Quantum Backpack <sup>1</sup>	\$75
<input type="checkbox"/> Weather Cover, Medium, Black	\$168
<input type="checkbox"/> Cup Holder	No Charge
<i>Side:</i>	
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Saddle Bag	\$35
<input type="checkbox"/> Tablet Holder	\$315
<input type="checkbox"/> Deluxe Cell Phone/Media Holder <sup>5</sup>	\$140

1. Backpack may interfere w/ other rear accessories.