



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

# Q6 Edge® 2.0 Order Form with Contoured and Synergy Seating - ODJFS Only

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_  
 Provider Name: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 Email: \_\_\_\_\_  
 PO Number: \_\_\_\_\_  
 Marked for: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_



## INTRODUCTION

This form is interactive when viewed as a PDF in most software. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Selections with HCPCS codes may have a part number noted with the selection. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email ([quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues.


If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum®!

## PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_  
 A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_  
 B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_  
 C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

## BASE MODEL

**Required.** Select a power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum® Q6 Edge® 2.0 3S-SS ..... \$7,595  
 Part: Q6 EDGE 2.0 3S-SS MFG Sugg HCPCS: K0848/K0011  
 Solid seat pan standard with no power option.

Quantum® Q6 Edge® 2.0 3S-C ..... \$7,595  
 Part: Q6 EDGE 2.0 3S-C MFG Sugg HCPCS: K0849/K0011  
 Contoured Seat (Captain) Standard.

Quantum® Q6 Edge® 2.0 3SP-SS ..... \$7,595  
 Part: Q6 EDGE 2.0 3SP-SS MFG Sugg HCPCS: K0856/K0011  
 Single power option.

## BASE COLOR

**Required.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



## DRIVE MOTORS

**Required.**

6 mph Drive Motors ..... Standard

Accu-Trac Motor Technology ..... \$1,295  
 HCPCS: K0108

Requires Q-Logic 2 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

## BASE MOUNTED FOOT PLATFORM

**Required.** May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform ..... Standard

Omit base mounted foot platform ..... No Charge

Bell Shaped, Base Mounted Foot Platform ..... \$235

## DRIVE WHEELS

**Required.** Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

3"x14" Flat Free Black Tires ..... Standard

3"x14" Pneumatic Black Tires ..... No Charge

4"x14" Knobby, Flat Free Black Tires ..... \$551

4"x14" Knobby, Pneumatic Black Tires ..... \$551

## Rim Color

Rim colors other than the standard silver bullet are available through a special order and extended lead time.

Silver Bullet ..... Standard

Match Rims to Base Color ..... \$85

Non-matching Rims ..... \$85

Select color below

Red Corvette

Back in Black

Pretty in Pink

Lime Light

Midnight Blue

Yellow Blaze

Electric Blue

White Satin

Deep Purple

Grey Street

Orange Crush

## POWER BASE CONTOURED/ SYNERGY OPTIONS

**Optional.** Select one option below. The power elevating seat actuator increases the seat to floor height. The elevate will be operated through the joystick and includes a 1/4 speed inhibit. Extended seat towers will increase seat to floor height.

Power Elevating Seat with 6" Stroke ..... \$1,895

Contoured w/ Solid Seat Pan Part: FRMASMB16309 HCPCS: E2300

Synergy Seat Part: FRMASMB16311 HCPCS: E2300

Power Elevating Seat with 7" Stroke ..... \$1,895

Contoured w/ Solid Seat Pan Part: FRMASMB16310 HCPCS: E2300

Synergy Seat Part: FRMASMB16312 HCPCS: E2300

Standard Seat Towers ..... Standard

Extended Seat Towers ..... No Charge

## CASTER WHEEL SIZE

**Required.**

6" Caster Wheels ..... Standard

8" Caster Wheels ..... \$150

Allows smoother transitions over thresholds and obstacles.

## BASE LIGHTING PACKAGES

**Optional.** Select one type of optional lighting packages below. Lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights ..... \$250  
Fender lights that operate only by pressing included single switch button.

Automatic Fender Lights..... \$275  
Fender lights that power on when the controller is activated and can be operated by the included single switch button.



Fender Lights

## ELECTRONICS

**Required.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Q-Logic Controls Order Form](#).

### Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick ..... Standard  
NE+ 6-Key Joystick (2 Actuators) ..... \$200  
Control up to 2 actuators through joystick.  
No Joystick, Setup for NE+ Joystick... No Charge

NE+ 8-Key Joystick (2 Actuators) ..... \$400  
Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

### Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

#### Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module ..... \$695  
HCPCS: E2377

Q-Logic 2 Expandable Controller..... \$425  
Harness  
Part: DWR1111H065 HCPCS: E2313



Q-Logic 2 EX Controller Module

#### Q-Logic 2 Joystick

Q-Logic 2 EX Joystick..... \$700  
Q-Logic 2 Ex Joystick for Seat..... \$700  
Mounted, LED Light Package  
No Joystick, Setup for Q-Logic

Stand Alone Joystick with ..... \$700  
External Charger Port  
Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.



Stand Alone Joystick

#### Q-Logic 2 Joystick Shroud

**Required with Q-Logic 2 EX Joystick selection above.** Reference colors shown on page 2.

Black, Soft Touch ..... Standard  
Non-matching Joystick Shroud ..... \$50  
Select color below

|               |                |               |               |
|---------------|----------------|---------------|---------------|
| Red Corvette  | White Satin    | Lime Light    | Midnight Blue |
| Electric Blue | Pretty in Pink | Grey Street   | Orange Crush  |
| Back in Black | Deep Purple    | Silver Bullet | Yellow Blaze  |



Q-Logic 2 EX Joystick with Lime Light

## JOYSTICK HANDLES

**Optional.** The NE Plus and Q-Logic 2 EX Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$71  
Part: ACCASMB1860 HCPCS: E2323  
Chin Cup (C) ..... \$71  
Part: ACCASMB1861 HCPCS: E2323  
(E2324 when chin control)  
T-Bar (D) ..... \$107  
Part: ACCASMB1858 HCPCS: E2323  
Medium Ball Handle (E)..... \$50  
Part: ACCASMB1859 HCPCS: E2323

Soft Ball Handle (F)..... \$50  
Part: ACCASMB1856 HCPCS: E2323  
Extended I-Handle (G)..... \$71  
Part: ACCASMB1857 HCPCS: E2323  
Bodypoint 3" Goal Post Handle..... \$69  
Part: PLSKNCP1057 HCPCS: E2323  
Bodypoint 4" Goal Post Handle..... \$69  
Part: PLSKNCP1058 HCPCS: E2323



Joystick Handles

## JOYSTICK MOUNTING BRACKETS

**Required.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket or See Positioning Positioning Components Page or [Order Form](#)

### Fixed, Inline Joystick Bracket

Fixed, Inline, Left Side Mount ..... Standard      Fixed, Inline, Right Side Mount ..... Standard

### Swing-Away (Retractable) Bracket


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount..... \$245  
NE/NE+ Part: ACC142526      HCPCS: E1028  
Q-Logic: ACC142524


Swing-away, 1" Offset, Left Side Mount .... \$245  
NE/NE+ Part: ACC142747      HCPCS: E1028  
Q-Logic: ACC142742

Swing-away, Inline, Right Side Mount ..... \$245  
NE/NE+ Part: ACC142527      HCPCS: E1028  
Q-Logic: ACC142525


Swing-away, 1" Offset, Right Side Mount.... \$245  
NE/NE+ Part: ACC142748      HCPCS: E1028  
Q-Logic: ACC142743



Swing Away Joystick  
Inline Mount (Right)




Swing Away Joystick  
Offset Mount (Right)



Retract4



Multi-Axis Joystick  
Mount



Swing-Away/Flip Down  
Bracket (Stand Alone  
Hardware Shown)

### Retract4 Swing-Away Arm Kit

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Retract4, Left Side Mount ..... \$250  
NE/NE+ Part: MEC138876      HCPCS: E1028  
Q-Logic: MEC138860

Retract4, Right Side Mount ..... \$250  
NE/NE+ Part: MEC138876      HCPCS: E1028  
Q-Logic: MEC138860

### Multi-Axis

Allows height adjustment of joystick.

Multi-axis, Left Side Mount .....\$100  
NE/NE+ Part: ACC142528      HCPCS: E1028  
Q-Logic: ACC131216

Multi-axis, Right Side Mount .....\$100  
NE/NE+ Part: ACC142528      HCPCS: E1028  
Q-Logic: ACC131216

### Stealth Swing Away/Flip Down Bracket

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket ..... \$275  
Left Side Mount  
Stand Alone Part: ST-TWBM480CJ-Q      HCPCS: E1028  
Q-Logic Part: ST-TWBM480QL-Q      HCPCS: E1028

Stealth Swing-Away/Flip Down Bracket ..... \$275  
Right Side Mount  
Stand Alone Part: ST-TWBM480CJ-Q      HCPCS: E1028  
Q-Logic Part: ST-TWBM480QL-Q      HCPCS: E1028

## ELECTRONICS ACCESSORIES

**Optional.**

USB Mobile Device Charger .....\$175  
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

Multiplier Harness .....\$221.21  
Q-Logic EX/NE Programmer ..... \$545  
Q-Logic PC Programming Disk with Cable ... \$750



USB Media  
Device Charger

## CONTOURED SEATING

**Optional.** Select only one seating system. Please see quote acknowledgment for seat assembly part numbers. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

### Black Vinyl High Back

High Back with Limited Recline ..... No Charge  
Select a size below.

|           |           |           |
|-----------|-----------|-----------|
| 16"Wx16"D | 16"Wx22"D | 20"Wx16"D |
| 16"Wx18"D | 18"Wx16"D | 20"Wx22"D |
| 16"Wx20"D | 18"Wx22"D |           |

High Back with Solid Seat Pan ..... No Charge  
and Limited Recline

Select a size below.

|             |      |      |      |      |
|-------------|------|------|------|------|
| Seat Width: | 16"W | 18"W | 20"W |      |
| Seat Depth: | 16"D | 18"D | 20"D | 22"D |

Depth-Adjustable High Back..... No Charge  
with Limited Recline

Seat is standard with 60" lap belt and full length armrests. Select a size below

|               |               |
|---------------|---------------|
| 18"Wx18"-20"D | 22"Wx20"-22"D |
| 20"Wx18"-20"D |               |

## CONTOURED SEATING CONTINUED

### Black Vinyl Medium Back

Medium Back ..... No Charge  
 Select a size below.  
 Seat Width: 16"W 18"W  
 Seat Depth: 16"D 18"D

Medium Back with Solid Seat Pan ..... No Charge  
 Select a size below.  
 Seat Width: 16"W 18"W  
 Seat Depth: 16"D 18"D

### Special Order Seat Size

Check the desired seat type above and complete width and depth below.

Width: \_\_\_\_\_ Depth: \_\_\_\_\_

### Contoured Seat Armrests

Flip-up, Height Adjustable, Removable ..... \$170  
 Part: See quote HCPCS: E0973

No Armrest ..... Standard

### Contoured Seat Armrest Size

14" Full Length   10" Desk Length  
 L R L R

### Contoured Seat Armpad

Black Bonded Armpad ..... No Charge

No Armpad ..... No Charge

## SYNERGY SEATING SYSTEM

**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges.

### Synergy Frame and Seat Selections

Synergy Seating Frame ..... \$434

### Synergy Seat Width

10"-14" Width 14"-18" Width 16"-20" Width 20"-24" Width

Set Width at: \_\_\_\_\_

### Synergy Seat Depth

10"-14" Depth 14"-18" Depth 16"-20" Depth 20"-24" Depth

Set Depth at: \_\_\_\_\_

### Synergy Transit Kit

Optional. The Q6 Edge® 2.0 is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option ..... \$200  
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://www.ercwts.org>.

### Synergy Back Options

**Required with Synergy Seat Frame.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with Special Order and extended lead time.

### Synergy Back Selection

Synergy Back ..... \$275  
 Select a back cane type and height and back height.

TRU-Comfort 4-way Stretch, 20" Width or less ..... \$595  
 Select back height. HCPCS: E2620

TRU-Comfort 4-way Stretch, 22" Width or greater ..... \$695  
 Select back height. HCPCS: E2621

Sling Back ..... \$98  
 Select a back height.

Available 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back ..... No Charge  
 Select a back cane type and height in following section.

|                |     | Width |     |     |      |      |      |      |      |     |      |
|----------------|-----|-------|-----|-----|------|------|------|------|------|-----|------|
|                |     | 12"   | 14" | 15" | 16"  | 17"  | 18"  | 19"  | 20"  | 21" | 22"  |
| Overall Height | 15" | S     | S   | S   | -    | -    | -    | -    | -    | -   | -    |
|                | 16" | S     | S   | S   | -    | -    | -    | -    | -    | -   | -    |
|                | 17" | -     | S   | S   | S    | S    | S    | S    | S    | S   | S    |
|                | 18" | -     | S   | S   | S    | S    | S    | S    | S    | S   | S    |
|                | 19" | -     | S   | S   | S    | S    | S    | S    | S    | S   | S    |
|                | 20" | -     | S   | S   | S TC | S    | S TC | S TC | S TC | S   | S    |
|                | 21" | -     | S   | S   | S TC | S    | S TC | S TC | S TC | S   | S    |
|                | 22" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S    |
|                | 23" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S    |
|                | 24" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S TC |
|                | 25" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S TC |

S Synergy Back S TC Synergy or TRU-Comfort Back

Synergy and TRU-Comfort Back Availability Chart

## SYNERGY SEATING SYSTEM CONTINUED

### Synergy Back Cane

Synergy Back Cane ..... No Charge  
Select a back cane height.

Flip Forward Back Cane ..... \$150  
Select a back cane height.

Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane..... \$195  
Select a cane height.

4" Mesomorph (2" per side) Back Cane ..... \$195  
Select a cane height.

Manual Recline Back ..... \$750  
Select an overall back height. HCPCS: E1226

Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, Shape Back, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.

2" Endomorph (2" per side) Back Cane ..... \$195  
Select a cane height.

4" Endomorph (2" per side) Back Cane ..... \$195  
Select a cane height.

### Synergy Push Handles

Push Handles ..... Standard

Omit Push Handles ..... \$125

### Synergy Accessory Bar

Accessory Bar ..... \$100

### Synergy Back Cane Height

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

Set Back Cane Height at: \_\_\_\_\_ (Range: 14"-25" High in increments of 1")

### Synergy Overall Back Height

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges.

Set Overall Back Height at: \_\_\_\_\_ (Range: 12"-25" High in increments of 1")

### Synergy Armrests

Required with Synergy Seat Frame.

Flip Back, Cane Mounted Armrests ..... No Charge  
May Interfere with swing-away laterals.

Two-Post, Flip-Back, Height Adjustable, ..... \$150 ea.  
Removable Armrests  
Part: See quote HCPCS: E0973  
Not available with a depth less than 16". Required with flip forward back and laptrays.

8"-12"H, Single Post, Height Adjustable, .....\$75 ea.  
Removable, Heavy Duty  
Part: See quote HCPCS: E0973

10"-14"H, Single Post, Height Adjustable, .....\$75 ea.  
Removable, Heavy Duty  
Part: See quote HCPCS: E0973

Pediatric, Removable, Height Adjustable .....\$73 ea.  
Part: See quote HCPCS: E0973  
Not available with full length armpads.

Height Adjustable, Cantilever Armrests (pair) ..... \$195  
Part: See quote  
Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, cane clamp pushbar, or Shape Back.

Single Post, Height Adjustable, Removable .....\$75 ea.  
Armrests  
Part: See quote HCPCS: E0973  
Height adjustment from 11"-17" in half inch increments.

No Armrests

### Synergy Armrests Size

Required with armrest selection above. Select a left (L) and right (R) armrest size. This also affects the armpad selection below.

L  R 14" Full Length

L  R 10" Desk Length

### Synergy Armpads

Required. Select to omit an armpad selection or select a left and right armpad length and a type of left and right armpad. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection and/or ..... No Charge  
See Positioning Components Page or Order Form

### Foam Armpads

Select a left (L) and right (R) foam type armpad (mix or match) OR choose the gel armpads that follow.


L  R Straight, Armpad ..... Standard

L  R Waterfall, Armpad..... No Charge

## JOYSTICK BRACKET RECEIVER

**Required.**

|   |              |
|---|--------------|
| Standard Fixed Receiver .....           | Standard     |
| Stealth Height Adjustable Receiver..... | \$200        |
| Tube Arm Part: ST-ARM250C-Q             | HCPCS: K0108 |
| Track Arm Part: ST-ARM250T-Q            | HCPCS: K0108 |



## Synergy Seat Mounted, LED Light Package

**Optional.** Available only with NE+ 8 Key Joystick or Q-Logic 2 Ex Joystick for LED Lights.

|                                      |       |
|--------------------------------------|-------|
| Seat Mounted, LED Light Package..... | \$995 |
|--------------------------------------|-------|



## Synergy Seat to Floor Height

**Required with Synergy Seat Frame.**

|         |         |         |
|---------|---------|---------|
| 15.875" | 16.875" | 17.875" |
|---------|---------|---------|

## LAP BELTS

**Required.** Select a lap belt below.

### Standard Lap Belts

|              |              |              |              |              |              |
|--------------|--------------|--------------|--------------|--------------|--------------|
| 38" Lap Belt | 50" Lap Belt | 60" Lap Belt | 70" Lap Belt | 80" Lap Belt | 90" Lap Belt |
|--------------|--------------|--------------|--------------|--------------|--------------|

### Optional Lap Belts


|  |      |  |       |
|--|------|--|-------|
| i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded, 1.5" Belt (IFPB5213) ..... | \$86 | i-Fit 60" Single Pull, Push Button, Extra Large Padded (11.5"), 2" Belt (IFPB5315) ..... | \$86  |
| i-Fit 60" Single Pull, Push Button, Large Padded (9.5"), 2" Belt (IFPB5314) .....    | \$86 | 60" Sure-Lok Lap Belt (Non-retractable).....   | \$120 |
|  |      | Required with ICO 7176-19 Compliant Transit Option.                                      |       |

## SYNERGY HEADRESTS

**Optional.** Select a headrest pad and bracket.

### Headrest Pads

|  |              |  |              |
|--|--------------|--|--------------|
| Stealth 8" Comfort Plus (CP450) .....      | \$142        | Stealth Small Contoured (9.75"Wx3"H).....    | \$60         |
| Part: POS148484                            | HCPCS: E0955 | Part: SETHEAD1052                            | HCPCS: E0955 |
| Stealth 10" Comfort Plus (CP250) (A) ..... | \$142        | Stealth Large Contoured (11"Wx6"H) (B) ..... | \$75         |
| Part: SETHEAD1057                          | HCPCS: E0955 | Part: SETHEAD1053                            | HCPCS: E0955 |
| Stealth 14" Comfort Plus (CP150) .....     | \$142        |  |              |
| Part: INDPART2327                          | HCPCS: E0955 |  |              |
| Pediatric Comfort Plus (6"Wx4.5"H) .....   | \$142        |  |              |
| (CPN350)                                   |              |  |              |
| Part: SETHEAD1069                          | HCPCS: E0955 |  |              |



### Headrest Bracket

|   |              |   |              |
|---|--------------|---|--------------|
| Stealth's Worlds Best Removable Headrest Bracket (TWB485-1) (C) ..... | \$210        | Pediatric Offset (HMO475P-17) (E) ..... | \$150        |
| Part: FRMASMB15999  | HCPCS: E1028 | Part: HDWASMB1166                       | HCPCS: E1028 |
| Stealth Multi-Axis Removable (D) .....                                | \$150        |   |              |
| Headrest Bracket (HMO475P)  |              |   |              |
| Part: HDWASMB1084   | HCPCS: E1028 |   |              |

## MANUFACTURER INSTALLED ACCESSORY

**Optional.**

|                         |              |
|-------------------------|--------------|
| Compact Vent Tray ..... | \$1,395      |
| Part: ACC133889         | HCPCS: E1029 |



## POWER LEGREST CONTROL

**Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

|   |  |
|---|--|
| Independent Legrests (Requires two actuator controls) | Articulating Foot Platform (Requires 1 actuator control) |
| Combined Legrests (Requires one actuator)             |  |





## LEGRESTS CONTINUED

### 70° Power Articulating, Swing Away Legrests

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper, lower extension, and footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

#### Configuration Verification

Left and Right Legrests will Match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

#### Power Articulating, Swing Away Legrest

L  R 70° Articulating, Swing-Away Legrest ..... \$1,197.50 ea.  
HCPCS: E1010

#### Lower Extension

L  R 11"-13" Length

L  R 14.5"-18" Length

Set Length at: \_\_\_\_\_

#### Angle Adjustable Footplate

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

L  R 4"Wx6"D

L  R 5"Wx6"D

L  R 5"Wx8"D

L  R 6"Wx8"D<sup>1</sup>

L  R 6"Wx11"D<sup>1</sup>

L  R 8"Wx11"D<sup>2</sup>

#### Accessories

Calf Pad (pair) 6"Wx6"H ..... No Charge

Gel Padded Knee Abductor Buttons (pair) ..... \$125  
INDPART2242 HCPCS: K0108

| 70° Power Articulating, Swing-Away Legrest Compatibility Matrix | Synergy Seat |        |      |  |
|---|--------------|--------|------|--|
|   | Low          | Medium | High |  |
| Angle Adjustable, 11" - 13"                                     | ✓            | ✓      | ✓    | ! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability. |
| Angle Adjustable, 14.5" - 18"                                   | !            | !      | !    | ✗ This rigging configuration does not work with this seat-to-floor height.   |
| Angle Adjustable, 17.5" - 21"                                   | ✗            | ✗      | ✗    |  |

### 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

#### Configuration Verification

Left and Right Legrests will match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

#### Swing-Away Legrest

L  R 60°-90° Swing-Away Legrest ..... \$81 ea.  
Lower extensions with minimum 11" length utilize a short style legrest upper.

#### Lower Extension

L  R 11"-13" Length  
N/A with Multi-axis

L  R 11"-16.5" Length  
Only with Multi-Axis

L  R 13"-16.5" Length  
N/A with Multi-Axis

L  R 15.5"-19.5" Length

Set Length at: \_\_\_\_\_

#### Extruded Footplate(s)

L  R 7.25"Wx6"D Extruded Footplate ..... Standard

L  R 8.5"Wx6"D Extruded Footplate ..... Standard  
May cause interference on seat widths less than 18."

#### Angle Adjustable or Multi Axis Selection

L  R Angle Adjustable with Heel Loops ..... \$75 ea.  
(Select footplate size) HCPCS: K0040

L  R Multi-Axis (Select Footplate Size) ..... \$247.50 ea  
HCPCS: K0040

#### Angle Adjustable or Multi-Axis Footplate Size

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

L  R 4"Wx6"D

L  R 5"Wx6"D

L  R 5"Wx8"D

L  R 6"Wx8"D<sup>1</sup>

L  R 6"Wx11"D<sup>1</sup>

L  R 8"Wx11"D<sup>2</sup>

# LEGRESTS CONTINUED

## 60°-90° Swing-Away Legrest Accessories

- Heel Loops (pair) .....\$38  
Part: FRMASMB7873 HCPCS: E0951  
Standard on angle adjustable and multi-axis footplates.
- Wheel Bumpers (pair) .....\$25
- Gel Neoprene Footrest Sleeve (pair) .....\$80  
Part: INDPART2243 HCPCS: K0108
- Therafin 31515 Multi-Axis Amputee Support, LEFT ..... \$230  
Part: INDPART2822 HCPCS: E1020

- Therafin 31515 Multi-Axis Amputee Support, RIGHT ..... \$230  
Part: INDPART2822 HCPCS: E1020
- Therafin 32387 Amputee Support Gel Cover, LEFT ..... \$150  
Part: INDPART2823 HCPCS: K0108
- Therafin 32387 Amputee Support Gel Cover, RIGHT .... \$150  
Part: INDPART2823 HCPCS: K0108

| 60°-90° Swing-Away Legrest Compatibility Matrix | Contoured/Synergy Seat |        |      |  |
|---|------------------------|--------|------|--|
|   | Low                    | Medium | High |  |
| Extruded, 13" - 16.5"                           | !                      | ✓      | ✓    | ✓ This rigging configuration causes no interference issues at this seat-to-floor height.   |
| Extruded, 15.5" - 19.5"                         | !                      | !      | !    | ! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability. |
| Angle Adjustable, 13" - 16.5"                   | !                      | ✓      | ✓    |  |
| Angle Adjustable, 15.5" - 19.5"                 | !                      | !      | !    |  |
| Multi-Axis, 15.5" - 19.5"                       | !                      | !      | !    | X This rigging configuration does not work with this seat-to-floor height.   |

## Elevating Legrests

Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

### Configuration Verification

Left and Right Legrests will match      Left and Right Legrests will not match      One Legrest (Left or Right as Marked)

### Elevating Legrest

Elevating Legrest Assembly ..... \$162.50 ea.  
L R Select lower extension and footplates HCPCS: E0990

### Lower Extension

11"-14" Length        12"-14.5" Length        14"-17" Length        14.5"-18" Length  
L R For reinforced only      L R For angle adjustable only      L R For reinforced only      L R For angle adjustable only

Set Length at: \_\_\_\_\_

### Reinforced Cast Footplate(s)

7.25"Wx6"D Reinforced Footplate..... Standard        8.5"Wx6"D Reinforced Footplate ..... Standard  
L R      L R May cause interference on seat widths less than 18."

### Angle Adjustable Footplate(s)

Angle Adjustable with Heel Loops ..... \$75 ea.  
L R (Select footplate size below) HCPCS: K0040

### Angle Adjustable Size

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D        5"Wx6"D        5"Wx8"D        6"Wx8"D<sup>1</sup>        6"Wx11"D<sup>1</sup>        8"Wx11"D<sup>2</sup>

### Accessories

Heel Loops (pair) .....\$38  
Part: FRMASMB7873 HCPCS: E0951  
Standard on angle adjustable footplates.

Wheel Bumpers (pair) .....\$25  
Falcon LPW Gel Padded Adductor Buttons (pair) ..... \$125  
Part: INDPART2242 HCPCS: K0108

| Elevating Legrest Compatibility Matrix | Contoured/Synergy Seat |        |      |  |
|--|------------------------|--------|------|--|
|  | Low                    | Medium | High |  |
| Standard, 11"-14"                      | ✓                      | ✓      | ✓    | ✓ This rigging configuration causes no interference issues at this seat-to-floor height.   |
| Standard, 14"-17"                      | !                      | !      | !    | ! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability. |
| Angle Adjustable, 12"-14.5"            | ✓                      | ✓      | ✓    |  |
| Angle Adjustable, 14.5"-18"            | !                      | !      | !    | X This rigging configuration does not work with this seat-to-floor height.   |

## CUSHION QUICK PICK

**Optional.** Available with solid seat pan configurations. General size selections are shown below. Refer to the [specific cushion order form](#) for detailed pricing, sizes, codes, and part numbers. Check the non-matching box and attached a cushion order form for a size other than the previously selected seat width and depth, or for a cushion not shown here.

Match seat selection width and depth

Non-matching cushion size, [see cushion order form](#)

### Simplicity Cushion (General Use)

16"Wx16"D to 20"Wx20"D ..... \$98  
*HCPCS: E2601*

22"Wx16"D to 24"Wx24"D ..... \$98  
*HCPCS: E2602*

### Solution® Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D ..... \$380  
*HCPCS: E2607*

22"Wx16"D to 24"Wx24"D ..... \$430  
*HCPCS: E2608*

### Solution® 1 Cushion (Skin Protection)

14"Wx14"D to 20"Wx24"D ..... \$380  
*HCPCS: E2603*

22"Wx16"D to 24"Wx24"D ..... \$430  
*HCPCS: E2604*

### Spectrum® Gel Cushion (Positioning)

12"Wx10"D to 20"Wx24"D ..... \$480  
*HCPCS: E2605*

22"Wx16"D to 24"Wx24"D ..... \$530  
*HCPCS: E2606*

### Spectrum® Foam Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D ..... \$480  
*HCPCS: E2607*

22"Wx16"D to 24"Wx24"D ..... \$530  
*HCPCS: E2608*

### TRU-Comfort™ 2 4-Way Stretch Cushion (Skin Protection)

16"Wx14"D to 20"Wx24"D ..... \$380  
*HCPCS: E2603*

22"Wx16"D to 22"Wx24"D ..... \$480  
*HCPCS: E2604*

### TRU-Comfort™ Sport Seat Cover

Available as an additional cover option for both TRU-Comfort seat cushions. Only available in sizes below.

16"Wx16"D ..... \$100

20"Wx20"D ..... \$100

18"Wx18"D ..... \$100

## PROVIDER INSTALLED ACCESSORIES

**Optional.**

Medium, Black Weather Cover ..... \$158

Rear Basket ..... \$103

Large, Black, Weather Cover ..... \$158

Walker Holder ..... \$92

Oxygen Holder ..... \$150  
 Not available with Synergy Manual Recline. *HCPCS: E2208*

Saddle Bag ..... \$28

Contoured Seat Cup Holder ..... \$17

Flag Holder ..... \$17

Synergy Seat Cup Holder ..... No Charge

Cane/Crutch Holder ..... \$98  
*HCPCS: E2207*

## BATTERIES

**Required.** Select a type of battery below, or to omit the batteries. These batteries are Group NF-22, maintenance free. The powerchair requires 2 batteries.

NF-22, 50 Amp AGM, UPG Battery ..... \$295 ea.  
Part: BATLIQ1018 HCPCS: E2361

NF-22, 50 Amp AGM, MK Battery..... \$295 ea.  
Part: BAT121346 HCPCS: E2361

NF-22, 50 Amp Gel Cell, UPG Battery ..... \$295 ea.  
Part: BATGEL1003 HCPCS: E2361

NF-22, 50 Amp Gel Cell, MK Battery..... \$295 ea.  
Part: BATGEL1008 HCPCS: E2361

No batteries. Provider will assume responsibility of batteries.

### Battery Installation

Please choose if batteries will be installed, if selected above.

Install Batteries ..... No charge

Do Not Install Batteries..... No Charge

## MANUFACTURER INSTALLED POSITIONING COMPONENTS

**Optional.** These components are for use with the Synergy Seating System. More options are available on the [Positioning Components Order Form](#).

### Stealth Thoracic Laterals

Swing-away laterals may interfere with TB3 and 2-post, flip up armrests. These laterals are not available with the back cane option.

#### Bracket

Swing-away, Rail Mounts, pair (TWBSTL)..... \$450  
Part: INDPART2896 HCPCS: E1028



#### Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$200  
Part: INDPART2663 HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$200  
Part: INDPART2664 HCPCS: E0956

#### Installation

Install Laterals

Do Not Install

### Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

#### Bracket

Flip Down Bracket (TWBTLADD-HW)..... \$225 ea.  
— Part: INDPART2804 HCPCS: E1028

Removable Bracket (LP430) ..... \$250 ea.  
— Part: INDPART2805 HCPCS: E1028

Flip Down with LEFT Full Surface Contact ..... \$225 ea.  
Feature (TWBADD-L/LPHW407)  
Part: POS140428 HCPCS: E1028

Removable Bracket with 1" Offset ..... \$325 ea.  
— (LP430 & LHW-121)  
Part: INDPART2854 HCPCS: E1028

Flip Down with RIGHT Full Surface Contact ..... \$225 ea.  
Feature (TWBADD-R/LPHW407)  
Part: POS140429 HCPCS: E1028

#### Adapter

2" Height Adapter Plate (LHW-122).....\$25  
Part: INDPART2821 HCPCS: K0108

Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) .....\$50  
Part: POS140430 HCPCS: K0108

Only available with the flip down full surface bracket.

#### Pad

Short Pad, 4"x4" (TWBADD-4X4) .....\$50 ea.  
— Part: INDPART2801 HCPCS: E0956

Long Pad, 4"x12" (TWBADD-4X12) .....\$50 ea.  
— Part: INDPART2803 HCPCS: E0956

Medium Pad, 4"x8" (TWBADD-4X8) .....\$50 ea.  
— Part: INDPART2802 HCPCS: E0956

#### Installation

Install on Left

Install on Right

Do Not Install

## NON-INSTALLED POSITIONING COMPONENTS

**Optional.** These options are for use with the Synergy Seating System. More options are available on the [Positioning Components Order Form](#).

### Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$475  
(SPLE107) (pair)  
Part: INDPART2786 HCPCS: K0108

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$475  
(SPLE109) (pair)  
Part: INDPART2788 HCPCS: K0108

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H..... \$475  
(SPLE108) (pair)  
Part: INDPART2787 HCPCS: K0108

### Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D..... \$266  
(SPLE111) (pair)  
Part: INDPART2783 HCPCS: K0108

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D..... \$266  
(SPLE114) (pair)  
Part: INDPART2785 HCPCS: K0108

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D..... \$266  
(SPLE112) (pair)  
Part: INDPART2784 HCPCS: K0108

### Stealth Elbow Support

Left, Swing-Away, Elbow Support Bracket..... \$175 ea.  
— Tube Style Arms (ARMS-ES-L)  
Part: POS124064 HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket..... \$175 ea.  
— Tube Style Arms (ARMS-ES-R)  
Part: POS124066 HCPCS: E1028

Left, Fixed Elbow Support Bracket..... \$175 ea.  
— TRU-Balance 3 Arms (ARMS-ESFTB3-L)  
Part: POS133465

Right, Fixed Elbow Support Bracket ..... \$175 ea.  
— TRU-Balance 3 Arms (ARMS-ESFTB3-R)  
Part: POS133466

### Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) .....\$50 ea.  
— Part: POS124065 HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) .....\$50 ea.  
— Part: POS124067 HCPCS: K0108

### Stealth i-Fit Dynamic Chest Straps

Upper Classic, Medium (IFPS1103) ..... \$142  
Part: ST-IFPS1103-Q HCPCS: E0960

Upper Contour, Large (IFPS2104) ..... \$142  
Part: ST-IFPS2104-Q HCPCS: E0960

Upper Classic, Large (IFPS1104) ..... \$142  
Part: ST-IFPS1104-Q HCPCS: E0960

Chest Strap, Medium, Belt Anchor (IFCS4113).....\$90

Upper Contour, Medium (IFPS2103) ..... \$142  
Part: ST-IFPS2103-Q HCPCS: E0960

Chest Strap, Large, Belt Anchor (IFCS4114).....\$90

### Stealth Swing Away Headrest Mounts

These Swing-Away headrest mounts are only available with the Stealth Comfort Series Headrests.

Facial Lateral, 4.5"x2" Pad, Right Side, ..... \$325  
for comfort Series (CPS9-R-HW & 900-Q)  
Part: INDPART2647 HCPCS: E1028

Facial Lateral, 4.5"x2" Pads, Left and Right Side ..... \$425  
for comfort Series (CPS9-B-HW & 900-Q)  
Part: INDPART2649 HCPCS: E1028

Facial Lateral, 4.5"x2" Pad, Left Side ..... \$325  
for comfort Series (CPS9-L-HW & 900-Q)  
Part: INDPART2648 HCPCS: E1028


### Other Accessories

Egg Switch Mount (ESM) (A) ..... \$50 ea.  
— Part: INDPART2815

Facial Pad, 6"x2" (920-Q) ..... \$50  
Part: ST-920-Q HCPCS: K0108

Figure 8 Bracket with 13" Rod and..... \$250  
Egg Switch Mount (B)  
Part: ST-F8678-13-Q HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) .... \$250  
Part: ST-SUS18-Q HCPCS: K0108  
Allows mounting of facial pads and attaches to the headrest



## SPECIALTY CONTROLS

**Optional.** Select a type of specialty controls below or to refer to the specialty controls order form.

### Enhanced Display

**Q-Logic 2 Enhanced Display Kit ..... \$1,884.93**

Select a color option below. Includes enhanced display, flexible goose-neck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

### Q-Logic 2 Enhanced Display Shroud

Required with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch ..... Standard

Match Joystick Shroud to Base Color ..... \$50

Part: See color below.

Non-matching Joystick Shroud ..... \$50

Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze

### Head Array

**Adult i-Drive Tri Array Series, ..... \$4,175**

3-Switch with Egg Switch for Mode (IDH200-1)

Part: CTL135940

HCPCS: E2330

Features 3 embedded proximity sensors for drive control and an egg switch for mode change.



Adult Tri-Array

### Other

See Specialty Controls Order Form

**Yellow Egg Switch (ESYEL) ..... \$77**

Multiplier Harness ..... \$221.21

**SCIM Module Kit ..... \$575**


Allows use of joystick display for specialty devices. Available with Q-Logic 2 EX Joystick selection above. Includes the TRU-Balance® 3 accessory bar. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

Red Egg Switch (ESRED) ..... \$77

Blue Egg Switch (ESBLU) ..... \$77

Black Egg Switch (ESBLK) ..... \$77

Green Egg Switch (ESGRN) ..... \$77



SCIM Module

### Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2016 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.