



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 Order Form with Pediatric TRU-Balance® 3 Positioning - ODJFS Only

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Marked for: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed as a PDF in most software. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Selections with HCPCS codes may have a part number noted with the selection. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues.


If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum®!

PATIENT INFORMATION

Weight: _____ Height: _____
 A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



BASE MODEL

Required. Select a power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum® Q6 Edge® 2.0 3S-SS \$7,595
 Part: Q6 EDGE 2.0 3S-SS MFG Sugg HCPCS: K0848/K0011
 Solid seat pan standard with no power option.

Quantum® Q6 Edge® 2.0 3MP-SS \$7,595
 Part: Q6 EDGE 2.0 3MP-SS MFG Sugg HCPCS: K0861/K0011
 Multiple actuator power option.

Quantum® Q6 Edge® 2.0 3SP-SS \$7,595
 Part: Q6 EDGE 2.0 3SP-SS MFG Sugg HCPCS: K0856/K0011
 Single power option.

BASE COLOR

Required. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



DRIVE MOTORS

Required.

6 mph Drive Motors Standard

Accu-Trac Motor Technology \$1,295
 HCPCS: K0108

Requires Q-Logic 2 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

BASE MOUNTED FOOT PLATFORM

Required. May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard

Omit base mounted foot platform No Charge

Bell Shaped, Base Mounted Foot Platform \$235

DRIVE WHEELS

Required. Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

3"x14" Flat Free Black Tires Standard

3"x14" Pneumatic Black Tires No Charge

4"x14" Knobby, Flat Free Black Tires \$551
 Not available with iLevel.

4"x14" Knobby, Pneumatic Black Tires \$551
 Not available with iLevel.

Rim Color

Rim colors other than the standard are available through a special order and extended lead time.

Standard Rims Standard

Match Rims to Base Color \$85

Non-matching Rims \$85
 Select color below

Red Corvette	Back in Black	Pretty in Pink	Lime Light	Midnight Blue	Yellow Blaze
Electric Blue	White Satin	Deep Purple	Grey Street	Orange Crush	

CASTER WHEEL SIZE

Required.

6" Caster Wheels Standard

8" Caster Wheels \$150
 Allows smoother transitions over thresholds and obstacles.

BASE LIGHTING PACKAGES

Optional. Select one type of optional lighting packages below. Lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights \$250
 Fender lights that operate only by pressing included single switch button.

Automatic Fender Lights \$275
 Fender lights that power on when the controller is activated and can be operated by the included single switch button.



Fender Lights

ELECTRONICS

Required. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Q-Logic Controls Order Form](#).

Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick Standard
Available with Tilt Only through Toggle.

NE+ 6-Key Joystick (1 Actuator)..... No Charge
Control of one actuator through Joystick on single power base model (3SP-5S) only.

No Joystick, Setup for NE+ Joystick... No Charge

NE+ 6-Key Joystick (2 Actuators) \$200
Control up to 2 actuators through joystick.

NE+ 8-Key Joystick (2 Actuators) \$400
Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module \$695
HCPCS: E2377

Q-Logic 2 Expandable Controller..... \$425
Harness
Part: DWR1111H065 HCPCS: E2313



Q-Logic 2 EX Controller Module

Q-Logic 2 Joystick

Q-Logic 2 EX Joystick..... \$700

Q-Logic 2 Ex Joystick for Seat Mounted, LED Light Package \$700

No Joystick, Setup for Q-Logic

Stand Alone Joystick with \$700
External Charger Port

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.



Stand Alone Joystick

Q-Logic 2 Joystick Shroud

Required with Q-Logic 2 EX Joystick selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Non-matching Joystick Shroud \$50

Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze



Q-Logic 2 EX Joystick with Lime Light

JOYSTICK HANDLES

Optional. The NE Plus and Q-Logic 2 EX Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$71
Part: ACCASMB1860 HCPCS: E2323

Chin Cup (C) \$71
Part: ACCASMB1861 HCPCS: E2323
(E2324 when chin control)

T-Bar (D) \$107
Part: ACCASMB1858 HCPCS: E2323

Medium Ball Handle (E)..... \$50
Part: ACCASMB1859 HCPCS: E2323

Soft Ball Handle (F)..... \$50
Part: ACCASMB1856 HCPCS: E2323

Extended I-Handle (G) \$71
Part: ACCASMB1857 HCPCS: E2323

Bodypoint 3" Goal Post Handle..... \$69
Part: PLSKNCP1057 HCPCS: E2323

Bodypoint 4" Goal Post Handle..... \$69
Part: PLSKNCP1058 HCPCS: E2323



Joystick Handles

JOYSTICK MOUNTING BRACKETS

Required. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket or See Positioning Positioning Components Page or [Order Form](#)

Fixed, Inline Joystick Bracket

Fixed, Inline, Left Side Mount Standard Fixed, Inline, Right Side Mount Standard

Swing-Away (Retractable) Bracket


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount..... \$245
NE/NE+ Part: ACC142526 HCPCS: E1028
Q-Logic: ACC142524


Swing-away, 1" Offset, Left Side Mount \$245
NE/NE+ Part: ACC142747 HCPCS: E1028
Q-Logic: ACC142742

Swing-away, Inline, Right Side Mount \$245
NE/NE+ Part: ACC142527 HCPCS: E1028
Q-Logic: ACC142525


Swing-away, 1" Offset, Right Side Mount.... \$245
NE/NE+ Part: ACC142748 HCPCS: E1028
Q-Logic: ACC142743



Swing Away Joystick
Inline Mount (Right)




Swing Away Joystick
Offset Mount (Right)



Retract4



Multi-Axis Joystick
Mount



Swing-Away/Flip Down
Bracket (Stand Alone
Hardware Shown)

Retract4 Swing-Away Arm Kit

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Retract4, Left Side Mount \$250
NE/NE+ Part: MEC138876 HCPCS: E1028
Q-Logic: MEC138860

Retract4, Right Side Mount \$250
NE/NE+ Part: MEC138876 HCPCS: E1028
Q-Logic: MEC138860

Multi-Axis

Allows height adjustment of joystick.

Multi-axis, Left Side Mount \$100
NE/NE+ Part: ACC142528 HCPCS: E1028
Q-Logic: ACC131216

Multi-axis, Right Side Mount \$100
NE/NE+ Part: ACC142528 HCPCS: E1028
Q-Logic: ACC131216

Stealth Swing Away/Flip Down Bracket

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket \$275
Left Side Mount
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028

Stealth Swing-Away/Flip Down Bracket \$275
Right Side Mount
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028

ELECTRONICS ACCESSORIES

Optional.

USB Mobile Device Charger \$175
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

Multiplier Harness \$221.21
Q-Logic EX/NE Programmer \$545
Q-Logic PC Programming Disk with Cable... \$750



USB Media
Device Charger

TRU-BALANCE® 3 SEATING

Required. Select one seating system. Be sure to complete the Seat Dimension, Seat to Floor, and back selections.

TRU-Balance® 3 Configuration

TRU-Balance® 3 Static Seat \$2,295
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance® 3 Power Recline \$7,595
Part: PTOINDV3444 HCPCS: E1004
175° of power recline.

TRU-Balance® 3 Power Tilt \$6,995
Part: PTOINDV3443 HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.

TRU-Balance® 3 Power Tilt & Recline \$11,895
Part: PTOINDV3446 HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

10" Power Lift or Manual Recline

Optional. Select one option below.

TRU-Balance® 3 10" Lift with iLevel® \$2,495
Part: PTOINDV3773 HCPCS: E2300
10" of scissor lift with iLevel® System. 22"W x 22"D maximum seat size.

Attendant Manual Recline \$500
Part: PTOINDV3501 HCPCS: E1226
Not available with NE+ joystick with Tilt through toggle. Available with Static Seat or Power Tilt. Not available with Power Lift.



TRU-BALANCE® 3 SEATING CONTINUED

Transit Kit

Optional. The Q6 Edge® 2.0 is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option..... \$200

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://www.ercwts.org>.

Seat Dimensions

Required with seat selection. Select width and depth package and enter desired set width in that package range, available in 1" increments. Configurations disallowed or not shown may be available through special order.

Seat Width

12"-16" Width Package

14"-16" Width Package

Set Width at: _____

Seat Depth

12"-16" Depth Package

14"-16" Depth Package

Set Depth at: _____

Seat to Floor Height

Required. Select a seat to floor height compatible with the seat functions previously selected.

16.5" Seat to Floor

17.5" Seat to Floor

18.5" Seat to Floor

19.5" Seat to Floor

Not available with Power Tilt or Power Lift configurations.

Not available with Static Seat or Power Recline Only.

Back Selection

Required. Match choice in this section with the options in the back cushion configuration section.

Static Back Standard
Select a back cushion configuration in the subsequent section.

Back Cane Option, Omit Back Pan/ Cushion(-\$125)
Available only with Static Back and Power Tilt (1 actuator) selections.
Select cane type and the back cane height below.

Flip Forward Back \$150
Available with Static Back and Power Tilt (1 actuator). Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width or with attendant manual recline. Select a back cushion type in the subsequent section.

Back Cane

Select with Back Cane Option selected above.

Angle Adjustable..... Standard

1" Endo \$35

Flip Forward \$25

2" Endo \$35

Back Cane Height

Select with Back Cane Option selected above.

16"

17"

20"

21"

24"

25"

Back Cushion Configuration

Required except when back cane option selected. Select a configuration type and applicable dimensions. Other size backs may be available with Special Order and extended lead time. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

Sport Back..... No Charge
Complete following width package and size..

TRU-Comfort 4-way Stretch, 20" Width or less..... \$595
Complete following width package and size. HCPCS: E2620

Omit Back Pan/Cushion..... No Charge
Complete following width package.

Aftermarket Back Interface \$200
Enter the preset back width below. Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount aftermarket backs. A Special Order is required on back widths less than 16".

Back Cushion Width

Available with +/-4" of seat width.

10"W-14"W Package

16"W-20"W Package

Set Width at: _____

	Tilt or Static Width			Recline Width		
	12"	14"	16"	12"	14"	16"
14"	S TC	S	-	-	-	-
15"	S TC	S	-	-	-	-
16"	S	S TC	S TC	S	S	S TC
17"	S	S TC	S TC	S	S	S TC
18"	-	-	S TC	-	-	S TC
19"	-	-	S TC	-	-	S TC
20"	-	-	S TC	-	-	S TC
21"	-	-	S TC	-	-	S TC
22"	-	-	S TC	-	-	S TC
23"	-	-	S TC	-	-	S TC
24"	-	-	S TC	-	-	S

S Sport Back S TC Sport or TRU-Comfort Back

Sport and TRU-Comfort Back Availability Chart

TRU-BALANCE® 3 SEATING CONTINUED

Back Cushion Overall Height

Static Seat and Power Tilt available height range: 14"-24". *Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.* **Power Recline available height range:** 16"-24". *Odd heights include 6" gap and even heights include 7" gap.* Please see the Q6Edge 2.0 TRU-Balance 3 order form for larger sizes. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

Set Overall Height at: _____

Aftermarket Back Interface Preset Width

Set Interface Width at: _____ (Range 16"-24" Wide)

TRU-Balance® 3 Back Shroud

Required. The back shroud is a feature of the TRU-Balance® 3 seating system and is available in matching shroud colors. The back shroud is not available with the back cane configuration. Reference colors shown on page 2.

Back in Black	Electric Blue	Pretty in Pink	Lime Light	Silver Bullet	Orange Crush
Red Corvette	White Satin	Deep Purple	Grey Street	Midnight Blue	Yellow Blaze

Additional TRU-Balance® 3 Back Options

Optional.

TRU-Comfort Lateral Positioning Wedge, pair\$30
Available with TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Privacy Flap, 9.75"\$35
Standard on recline seating systems.

Privacy Flap, 12.75"\$35
Standard on recline seating systems.

See Comfort Company Back order form

ARMRESTS

Required. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

TRU-Balance® 3 Flip-Back Armrests

Choose a left and right TRU-Balance® 3 armrest and set height and optional downpost if desired. Different Left (L) and Right (R) armrest height packages are allowed. Not available with Back Cane Option.

L R TRU-Balance® 3 Flip-Back, 8"-12"H..... Standard

L R TRU-Balance® 3 RIGHT Flip-Back, 12"-16"H Standard

Set LEFT height at: _____ Set RIGHT height at: _____

Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

Pediatric, Armrests\$75 ea.
HCPCS: E0973

Set LEFT height at: _____ Set RIGHT height at: _____

ARMREST SIZE

Required. Select a Left (L) and Right (R) armrest size. This also affects the armpad selection below.

L R 14" Full Length

L R 10" Desk Length

ARMPADS

Required. Select to omit an armpad selection or select a left and right armpad length and a type of left and right armpad. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection and/or No Charge
See Positioning Components Page or Order Form

Foam Armpads

Select a left and right standard type armpad (mix or match) OR choose the gel armpads that follow.

L R Straight, Armpad Standard

L R Waterfall, Armpad..... No Charge

ARMPADS CONTINUED

Falcon Gel Armpads

Sold only in pairs. Select an armpad set below. Armpads are 3" wide. The mounting plate allows the armpad to be offset from the armrest tube.

Falcon Gel Left and Right Armpad:
 Full Length\$92.50 ea
 Desk Length\$70.00 ea

For use with TB3 standard armrests.

Falcon Gel Left and Right Armpad with Mounting Plate:
 Full Length \$148.02 ea
 Desk Length \$126.20 ea

For use with 2-Post, Heavy Duty Removable, or Pediatric Armrests.

JOYSTICK BRACKET RECEIVER

Required.

Standard Fixed Receiver Standard

Stealth Height Adjustable Receiver \$200

Tube Arm Part: ST-ARM250C-Q

HCPCS: K0108

Track Arm Part: ST-ARM250T-Q

HCPCS: K0108



Height Adjustable Receiver (Tube Arm)



Height Adjustable Receiver (Track Arm)

HEADRESTS

Optional. Select a headrest pad and bracket.

Headrest Pads

Stealth 8" Comfort Plus (CP450)\$142
 Part: POS148484 HCPCS: E0955


Stealth 10" Comfort Plus (CP250) (A)\$142
 Part: SETHEAD1057 HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H)\$142
 (CPN350)
 Part: SETHEAD1069 HCPCS: E0955

Stealth Small Contoured (9.75"Wx3"H) \$60
 (STL-OTSM250)

Part: SETHEAD1052

HCPCS: E0955



Headrests

Headrest Bracket

Stealth's Worlds Best Removable.....\$210
 Headrest Bracket (TWB485-1) (C)
 Part: FRMASMB15999 HCPCS: E1028

Pediatric Offset (HMO475P-17) (E).....\$150
 Part: HDWASMB1166 HCPCS: E1028

Stealth Multi-axis Removable Headrest.....\$150
 Bracket (HMO475P) (D)
 Part: HDWASMB1084 HCPCS: E1028

LAP BELTS

Required. Select a lap belt below.

Standard Lap Belts

38" Lap Belt 50" Lap Belt 60" Lap Belt

Optional Lap Belts

i-Fit 60" Single Pull, Push Button, Medium\$86
 (7.5") Padded, 1.5" Belt (IFPB5213)

i-Fit 60" Single Pull, Push Button, Extra Large\$86
 Padded (11.5"), 2" Belt (IFPB5315)

i-Fit 60" Single Pull, Push Button, Large\$86
 Padded (9.5"), 2" Belt (IFPB5314)

60" Sure-Lok Lap Belt (Non-retractable)..... \$120
 Required with ICO 7176-19 Compliant Transit Option.

MANUFACTURER INSTALLED ACCESSORIES

Optional.

Trilogy Vent Tray \$400
 Part: ACC136791 HCPCS: E1029

Accessory Bar\$100

Standard with the SCIM module, Trilogy vent tray, and oxygen holder. May be added to mount additional accessory to back of unit. Not available with back cane option.

Compact Vent Tray \$1,395
 Part: ACC133889 HCPCS: E1029

Push Handles (pair)\$175

Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.



Trilogy Vent Tray



Compact Vent Tray

POWER LEGREST CONTROL

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)

Articulating Foot Platform (Requires 1 actuator control)

Combined Legrests (Requires 1 actuator control)

POWER POSITIONING ELECTRONICS

Required, except with Static Seat and no power options. Select the desired control type for the seating system. Please note the available actuator controls. The 4-way toggle and 4 function push button are mounted on the same side as the joystick.

Single Actuator Control, Through Toggle No Charge

Only available with Power Tilt only. Manual recline is not available with NE Tilt through toggle. Fender lights are not available with tilt through toggle. Select toggle side.

Left Side Toggle

Right Side Toggle

Through NE+ 6-Key Joystick with either:

Single Actuator Control..... \$1,895

Quote will select this price/code with 1 actuator. HCPCS: E2310

Two Actuator Control.....\$2,895

Quote will select this price/code with 2 actuators. HCPCS: E2311

Single Actuator Through Q-Logic 2 Controller \$1,895
and Interface (AM1) HCPCS: E2310

Multiple Actuator Through Q-Logic 2 Controller.....\$2,895
and Interface (AAM) HCPCS: E2311

Multiple Actuator Through 4-Way Toggle,\$2,895

Q-Logic 2 Controller, and Interface (AAM) HCPCS: E2311
Bracket mounts to outside of the armpad.

Multiple Actuators Through 4-Way Button Toggle,\$2,895

Q-Logic 2 Controller, and Interface (AAM) HCPCS: E2311
Bracket mounts to outside of the armpad.

Through 4 Function Push Button and Q-Logic 2 Controller
and Interface with either:

Single Actuator Control (AAM)\$2,895

Quote will select this price/code with 1 actuator. HCPCS: E2310

Multiple Actuator Control (AAM)\$2,895

Quote will select this with multiple actuators. HCPCS: E2311
Only available with TB3 armrests.

Attendant Control

Optional. Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$500

Part: ELE125264 HCPCS: E2331

Seat Mounted, LED Light Package

Optional. Available only with NE+ 8 Key Joystick or Q-Logic 2 Ex Joystick for LED Lights.

Seat Mounted, LED Light Package..... \$995



Seat Mounted LED Lights (One Shown)

LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

Center Mount Foot Platforms

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section.

Power Articulating Foot Platform \$2,995

Please verify the appropriate billing code to be used for this item with all third-party payors.

Manually Adjustable Center Mount \$250
Foot Platform



Power Articulating Foot Platform

Center Mount Lower Extension

8"-12" Extension

10"-16" Extension

13"-19" Extension

N/A with Power Articulation.

Set Length at: _____

Center Mount Footplates

5"x6" Footplates

5"x8" Footplates

5"x11" Footplates

Center Mount Foot Platform Compatibility Matrix	TRU-Balance® 3 Power Recline			TRU-Balance® 3 Power Tilt			TRU-Balance® 3 Power Lift			<input checked="" type="checkbox"/> No interference issues at this seat-to-floor height. <input type="checkbox"/> May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
	Low	Medium	High	Low	Medium	High	Low	Medium	High	
10"-16" Length Extension	!	!	✓	✓	✓	✓	✓	✓	✓	
13"-19" Length Extension	!	!	!	!	!	!	!	!	!	

LEGRESTS CONTINUED

Center Mount Calf Support

Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pad (pair) 6"Wx6"H.....	\$150	Small Calf Panel 18"Wx8"H	\$150
Large Calf Pad (pair) 8" taper to 6"Wx8"H	\$150	Large Calf Panel 24"Wx8"H	\$200

Only available on widths of 20" or greater.

Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension, 8"D and 11"D Footplates	\$50	4" Outside Edge	\$50
2" Bell-Shaped Length Extension	\$50	6" Outside Edge.....	\$50

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

Only available with 6"D footplates and 8"D footplate extension.

Available with 8"D and 11"D footplates and 11"D footplate extension.

Additional Accessories

Reversed Dartex Calf Pad or Panel Cover (pair)	\$120	Footplate Connector	\$15
Wheel Bumpers (pair)	\$25		

90° Pediatric Swing Away Legrests

Left (L) and right (R) configurations are selected by checking the appropriate boxes.

Configuration Verification

Left and Right Legrests will Match Left and Right Legrests will not match One Legrest (Left or Right as Marked)

Pediatric, Swing Away Legrest

90° Pediatric, Swing-Away Legrest..... \$81 ea.
L R

Lower Extension

7"-10" Length
L R

Set Length at: _____

Angle Adjustable Footplate with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D \$75 ea. 5"Wx6"D \$75 ea.
L R L R

Accessories

Wheel Bumpers (pair) \$25

Elevating Legrests

Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

Configuration Verification

Left and Right Legrests will match Left and Right Legrests will not match One Legrest (Left or Right as Marked)

Elevating Legrest

Elevating Legrest Assembly \$162.50 ea.
L R Select lower extension and footplates HCPCS: E0990

Lower Extension

11"-14" Length 12"-14.5" Length 14"-17" Length 14.5"-18" Length
L R For reinforced only L R For angle adjustable only L R For reinforced only L R For angle adjustable only

Set Length at: _____

Reinforced Cast Footplate(s)

7.25"Wx6"D Reinforced Footplate..... Standard 8.5"Wx6"D Reinforced Footplate Standard
L R L R May cause interference on seat widths less than 18."

LEGRESTS CONTINUED

Angle Adjustable Footplate(s)

L R Angle Adjustable with Heel Loops \$75 ea.
 (Select footplate size below) HCPCS: K0040

Angle Adjustable Size

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

L R 4"Wx6"D

L R 5"Wx6"D

L R 5"Wx8"D

L R 6"Wx8"D¹

L R 6"Wx11"D¹

L R 8"Wx11"D²

Accessories

Heel Loops (pair) \$38
 Part: FRMASMB7873 HCPCS: E0951
 Standard on angle adjustable footplates.

Wheel Bumpers (pair) \$25
 Gel Padded Adductor Buttons (pair) \$125
 Part: INDPART2242 HCPCS: K0108

CUSHION QUICK PICK

Optional. General size selections are shown below. Refer to the [specific cushion order form](#) for detailed pricing, sizes, codes, and part numbers. Check the non-matching box and attached a cushion order form for a size other than the previously selected seat width and depth, or for a cushion not shown here.

Match seat selection width and depth

Non-matching cushion size, [see cushion order form](#)

Simplicity Cushion (General Use)

16"Wx16"D to 20"Wx20"D \$98
 HCPCS: E2601

22"Wx16"D to 24"Wx24"D \$98
 HCPCS: E2602

Solution[®] Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D \$380
 HCPCS: E2607

22"Wx16"D to 24"Wx24"D \$430
 HCPCS: E2608

Solution[®] 1 Cushion (Skin Protection)

14"Wx14"D to 20"Wx24"D \$380
 HCPCS: E2603

22"Wx16"D to 24"Wx24"D \$430
 HCPCS: E2604

Spectrum[®] Gel Cushion (Positioning)

12"Wx10"D to 20"Wx24"D \$480
 HCPCS: E2605

22"Wx16"D to 24"Wx24"D \$530
 HCPCS: E2606

Spectrum[®] Foam Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D \$480
 HCPCS: E2607

22"Wx16"D to 24"Wx24"D \$530
 HCPCS: E2608

TRU-Comfort[™] 2 4-Way Stretch Cushion (Skin Protection)

16"Wx14"D to 20"Wx24"D \$380
 HCPCS: E2603

22"Wx16"D to 22"Wx24"D \$480
 HCPCS: E2604

TRU-Comfort[™] Sport Seat Cover

Available as an additional cover option for both TRU-Comfort seat cushions. Only available in sizes below.

16"Wx16"D \$100


20"Wx20"D \$100

18"Wx18"D \$100

PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU-Balance® 3 Cup Holder No Charge For use with TRU-Balance 3 (track style) armrests.	TRU-Balance® 3 Front Mounted..... \$50 Cell Phone Holder (E)
Cup Holder (Tube-style Armrest) No Charge For use with tube-style armrests.	Large, Black, Weather Cover.....\$158
TRU-Balance® 3 Backpack Holder (Pair) (A)... \$25 Not available with back cane option.	TRU-Balance® 3 Reflector Kit \$25
TRU-Balance® 3 Hydration System (B).....\$175	TRU-Balance® 3 Personal Item Hook \$25 Not available with back cane option.
TRU-Balance® 3 Glove Box (C)\$100	TRU-Balance® 3 Rear View Mirror..... \$25
TRU-Balance® 3 Transfer Bar,\$175 Bolt-On (Pair) HCPCS: E0705	TRU-Balance® 3 Oxygen Holder\$150 HCPCS: E2208
TRU-Balance® 3 Transfer Bar,\$175 Quick Release (Pair) (D) HCPCS: E0705	TRU-Balance® 3 Tablet Holder \$375
	TRU-Balance® 3 Clothing Guard (Pair).....\$175



BATTERIES

Required. Select a type of battery below, or to omit the batteries. These batteries are Group NF-22, maintenance free. The powerchair requires 2 batteries.

NF-22, 50 Amp AGM, UPG Battery \$295 ea. Part: BATLIQ1018 HCPCS: E2361	NF-22, 50 Amp AGM, MK Battery..... \$295 ea. Part: BAT121346 HCPCS: E2361
NF-22, 50 Amp Gel Cell, UPG Battery \$295 ea. Part: BATGEL1003 HCPCS: E2361	NF-22, 50 Amp Gel Cell, MK Battery..... \$295 ea. Part: BATGEL1008 HCPCS: E2361

No batteries. Provider will assume responsibility of batteries.

Battery Installation

Please choose if batteries will be installed, if selected above.

Install Batteries No charge	Do NOT install batteries..... No Charge
-----------------------------------	---

MANUFACTURER INSTALLED POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

Stealth Thoracic Laterals

Swing-away laterals may interfere with TB3 and 2-post, flip up armrests. These laterals are not available with the back cane option.

Bracket

Swing-away, Rail Mounts, pair (TWBSTL)..... \$450 Part: INDPART2896 HCPCS: E1028



Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$200 Part: INDPART2663 HCPCS: E0956	Curved Pads, 4"x6", pair (TWBL-4X6)..... \$200 Part: INDPART2664 HCPCS: E0956
--	--

Installation

Install Laterals	Do Not Install
------------------	----------------

Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

Bracket

Flip Down Bracket (TWBTLADD-HW)..... \$225 ea. — Part: INDPART2804 HCPCS: E1028	Removable Bracket (LP430)..... \$250 ea. — Part: INDPART2805 HCPCS: E1028
Flip Down with LEFT Full Surface Contact \$225 ea. Feature (TWBADD-L & LPHW407) Part: POS140428 HCPCS: E1028	Removable Bracket with 1" Offset \$325 ea. — (LP430 & LHW-121) Part: INDPART2854 HCPCS: E1028
Flip Down with RIGHT Full Surface Contact \$225 ea. Feature (TWBADD-R & LPHW407) Part: POS140429 HCPCS: E1028	

Adapter

2" Height Adapter Plate (LHW-122).....\$25 Part: INDPART2821 HCPCS: K0108 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.	2 1/4" Adapter Plate (LHW-129)\$50 Part: POS140430 HCPCS: K0108 Only available with the flip down full surface bracket.
--	---

MANUFACTURER INSTALLED POSITIONING COMPONENTS CONTINUED

Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4)\$50 ea.
— Part: INDPART2801 HCPCS: E0956

Medium Pad, 4"x8" (TWBADD-4X8)\$50 ea.
— Part: INDPART2802 HCPCS: E0956

Long Pad, 4"x12" (TWBADD-4X12)\$50 ea.
— Part: INDPART2803 HCPCS: E0956

Installation

Install on Left

Install on Right

Do Not Install

NON-INSTALLED POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$475
(SPLE107) (pair)
Part: INDPART2786 HCPCS: K0108

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$475
(SPLE109) (pair)
Part: INDPART2788 HCPCS: K0108

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$475
(SPLE108) (pair)
Part: INDPART2787 HCPCS: K0108

Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
(SPLE111) (pair)
Part: INDPART2783 HCPCS: K0108

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D..... \$266
(SPLE114) (pair)
Part: INDPART2785 HCPCS: K0108

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D..... \$266
(SPLE112) (pair)
Part: INDPART2784 HCPCS: K0108

Stealth Elbow Support

Left, Swing-Away, Elbow Support Bracket..... \$175 ea.
— Tube Style Arms (ARMS-ES-L)
Part: POS124064 HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket..... \$175 ea.
— Tube Style Arms (ARMS-ES-R)
Part: POS124066 HCPCS: E1028

Left, Fixed Elbow Support Bracket \$175 ea.
— TRU-Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465

Right, Fixed Elbow Support Bracket \$175 ea.
— TRU-Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466

Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$50 ea.
— Part: POS124065 HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$50 ea.
— Part: POS124067 HCPCS: K0108

Stealth i-Fit Dynamic Chest Straps

Upper Classic, Medium (IFPS1103) \$142
Part: ST-IFPS1103-Q HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142
Part: ST-IFPS2104-Q HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142
Part: ST-IFPS1104-Q HCPCS: E0960

Chest Strap, Medium, Belt Anchor (IFCS4113).....\$90

Upper Contour, Medium (IFPS2103) \$142
Part: ST-IFPS2103-Q HCPCS: E0960

Chest Strap, Large, Belt Anchor (IFCS4114).....\$90

Stealth Swing Away Headrest Mounts

These Swing-Away headrest mounts are only available with the Stealth Comfort Series Headrests.

Facial Lateral, 4.5"x2" Pad, Right Side, \$325
for comfort Series (CPS9-R-HW & 900-Q)
Part: INDPART2647 HCPCS: E1028

Facial Lateral, 4.5"x2" Pads, Left and Right Side \$425
for comfort Series (CPS9-B-HW & 900-Q)
Part: INDPART2649 HCPCS: E1028

Facial Lateral, 4.5"x2" Pad, Left Side \$325
for comfort Series (CPS9-L-HW & 900-Q)
Part: INDPART2648 HCPCS: E1028

Other Accessories


Egg Switch Mount (ESM) (A) \$50 ea.
— Part: INDPART2815

Facial Pad, 6"x2" (920-Q) \$50
Part: ST-920-Q HCPCS: K0108

Figure 8 Bracket with 13" Rod and..... \$250
Egg Switch Mount (B)
Part: ST-F8678-13-Q HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$250
Part: ST-SUS18-Q HCPCS: K0108

Allows mounting of facial pads and attaches to the headrest



SPECIALTY CONTROLS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

Enhanced Display

Q-Logic 2 Enhanced Display Kit \$1,884.93

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Q-Logic 2 Enhanced Display Shroud

Required with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Match Joystick Shroud to Base Color \$50
Part: See color below.

Non-matching Joystick Shroud \$50
Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze

Head Array

Pediatric i-Drive Tri Array Series, \$4,175
3-Switch with Egg Switch for Mode (IDH250-1)
Part: CTL135956 HCPCS: E2330

Features 3 embedded proximity sensors for drive control and an egg switch for mode change.



Other

See Specialty Controls Order Form

Yellow Egg Switch (ESYEL) \$77

Multiplier Harness \$221.21

SCIM Module Kit \$575

Red Egg Switch (ESRED) \$77

Allows use of joystick display for specialty devices. Available with Q-Logic 2 EX Joystick selection above. Includes the TRU-Balance® 3 accessory bar. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

Blue Egg Switch (ESBLU) \$77

Black Egg Switch (ESBLK) \$77

Green Egg Switch (ESGRN) \$77



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2016 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.