

J6 Series

Group 2 Standard
Order Form/ODJFS Only
300lb Weight Capacity

Pride[®]
Mobility Products Corp.

Pride Mobility Products Corporation

182 Susquehanna Ave., Exeter, PA 18643
Phone: (866) 800-2002 • Fax: (866) 707-3422

quantumorders@pridemobility.com

Discount: _____

Terms: _____

Freight: _____

Taken By: _____



Account Number: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

E-Mail: _____

Date: _____

PO#: _____

Marked For: _____

Ship To Address: _____

City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP. Items Designated with a Circular Checkbox (○) are not available through Quantum West

BASE MODEL

Description	Price	HCPCS Code
☐ J6 2S-C	\$5,971	K0823/K0011 MFG Sugg
☐ J6 2S-SS	\$5,971	K0822/K0011 MFG Sugg
☐ J6 2SP-SS	\$5,971	K0835/K0011 MFG Sugg

COLORS

- ☐ Candy Apple Red
☐ Viper Blue

BATTERY CHARGER

Description	Price
☐ 5 Amp Off-board Charger	Standard

DRIVE WHEELS

Description	Price
☐ 10" Pneumatic, pair	No Charge
☐ 10" Flat-Free, pair	Standard

POWER ELEVATING SEAT



Description	Price	HCPCS Code
☐ Power Elevating Seat Through VSI 6-Key ¹	\$975	E2300

1. Only available with VSI 6-Key. Not available with manual recliner, Synergy® Manual Recline or manual tilt. Max seat size is 20"W x 20"D. Increases STF ht. Includes 1/4 speed inhiit.

BASE MOUNTED FOOT PLATFORM

Description	Price
☐ Base Mounted Foot Platform	Standard
☐ No Foot Platform	No Charge
☐ Bell Shaped Foot Platform	\$235


NON-EXPANDABLE CONTROLS

Description	Price
☐ 75 Amp NE 	Standard
☐ 75 Amp NE+ 6-Key (For up to 2 Actuators Thru Joystick)	\$200
☐ 75 Amp NE+ 8-Key (For Lights) 	\$400
☐ 70 Amp PG VSI 6-Key ¹	Standard

1. Only available with a power elevating seat.

EXPANDABLE CONTROLS

Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls

Description	Price
☐ 100 Amp Q-Logic 2 EX Controller 	\$695
☐ Expandable Controller Harness	\$425
☐ Q-Logic 2 EX Joystick 	\$700
☐ Q-Logic 2 EX Joystick (For Lights) 	\$700
☐ Stand Alone Joystick w/ External Charger Port 	\$700

Armor Color Options:

- ☐ Black Soft Touch (Standard color option)
☐ Matching
☐ Non-Matching

Colors: \$50

- ☐ Candy Apple Red
☐ Viper Blue

Items Designated with a Circular Checkbox (○) are not available through Quantum West

JOYSTICK MOUNTING BRACKETS

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
☐ Fixed, Inline ¹	Standard	
☐ Left		
☐ Right		
☐ Swing-Away	\$245	E1028
SAMTIL7 ☐ Inline, Left		
SAMTIR7 ☐ Inline, Right		
SAMTOL7 ☐ Offset, Left		
SAMTOR7 ☐ Offset, Right		
☐ Retractable Swing-Away Joystick Arm Kit	\$250	E1028
☐ Left		
☐ Right		
☐ Multi-Axis	\$100	E1028
☐ Left		
☐ Right		
☐ No Joystick Mount		

1. Not available with Q-Logic 2 EX Stand Alone Joystick.

JOYSTICK HANDLES

Description	Price	HCPCS Code
☐ Standard Carrot Handle	Standard	
☐ Mushroom Handle	\$71	E2323
☐ T-Bar	\$107	E2323
☐ Chin Cup	\$71	E2324 ¹
☐ Extended I-Handle	\$71	E2323
☐ Soft Ball Handle	\$50	E2323
☐ Bodypoint 3" Goal Post Handle	\$69	E2323
☐ Bodypoint 4" Goal Post Handle	\$69	E2323
☐ Medium Ball Handle	\$50	E2323

1. Use HCPCS Code E2323 when not used in conjunction with a chin control.

PROGRAMMING

Description	Price
☐ Q-Logic EX/NE Programmer	\$545
	
☐ Q-Logic PC Programming Disk w/ Cable ¹	\$750
☐ P&G Programmer	\$423

1. Works with Q-Logic 2 EX only.

CONTOURED SEATING

Description	Price
Depth-Adjustable High Back: ¹	
☐ Black Vinyl, Depth-Adjustable High Back, Limited Recline	No Charge
☐ 18"W x 18-20"D	
☐ 20"W x 18-20"D	
Limited Recline High Back:	
☐ Black Vinyl, High Back with Limited Recline	No Charge
☐ 16"W x 16"D	☐ 18"W x 16"D
☐ 16"W x 18"D	☐ 20"W x 16"D
☐ 16"W x 20"D	
Limited Recline High Back w/ Solid Seat Pan:	
☐ Black Vinyl, High Back with Solid Seat pan, Limited Recline	No Charge
☐ 16"W x 16"D	☐ 18"W x 20"D
☐ 16"W x 18"D	☐ 20"W x 16"D
☐ 16"W x 20"D	☐ 20"W x 18"D
☐ 18"W x 16"D	☐ 20"W x 20"D
☐ 18"W x 18"D	
Medium Back:	
☐ Black Vinyl, Medium Back	No Charge
☐ 16"W x 16"D	☐ 18"W x 16"D
☐ 16"W x 18"D	☐ 18"W x 18"D
Medium Back w/ Solid Seat Pan:	
☐ Black Vinyl, Medium Back With Solid Seat Pan	No Charge
☐ 16"W x 16"D	☐ 18"W x 16"D
☐ 16"W x 18"D	☐ 18"W x 18"D
Armrest Type:	
☐ Flip up, Height-Adjustable, Removable	\$85 ea. E0973
Armpads:	
Left: ☐ Full Length	Right: ☐ Full Length
☐ Desk Length ²	☐ Desk Length ²

Special Order Seat

____ "W x ____ "D

Seat size can not exceed 20"W x 20"D

Contact your inside sales specialist for pricing

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part.
2. Selection of desk length armpads on a depth adjustable high back requires them to be added as an additional part.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

MANUAL RECLINER

Description	Price	HCPSC Code
☐ Manual Recliner High Back With Headrest ¹	\$600	E1226
☐ Manual Recliner High Back With Headrest and Solid Seat Pan ¹	\$600	E1226
Color: ☐ Gray Vinyl ☐ Beige Vinyl		
Sizes ☐ 16"W x 16"D ☐ 18"W x 20"D		
☐ 16"W x 18"D ☐ 20"W x 16"D		
☐ 16"W x 20"D ☐ 20"W x 18"D		
☐ 18"W x 16"D ☐ 20"W x 20"D		
☐ 18"W x 18"D		
Armrest Type:		
☐ Manual Recliner Armrests	No Charge	
Armpads:		
Left: ☐ Full Length ☐ Desk Length	Right: ☐ Full Length ☐ Desk Length	

1. Not available with a power elevating seat.

SYNERGY® SEATING SYSTEM

Description	Price	HCPSC Code
☐ Synergy® Seating Frame	\$434	
Towers:		
☐ Standard (15.75"-16.75")		
☐ Extended (16.5"-17.5") ¹		

SEATING SPECIFICATIONS

Seat Width:

- ☐ 10"W
☐ 11"W
☐ 12"W
☐ 13"W
☐ 14"W
☐ 15"W
☐ 16"W
☐ 17"W
☐ 18"W
☐ 19"W
☐ 20"W

Seat Width Package:

- ☐ 10"W-12"W
☐ 10"W-14"W
☐ 12"W-14"W
☐ 14"W-18"W
☐ 16"W-20"W

Seat Depth:

- ☐ 10"D
☐ 11"D
☐ 12"D
☐ 13"D
☐ 14"D
☐ 15"D
☐ 16"D
☐ 17"D
☐ 18"D
☐ 19"D
☐ 20"D

Seat Depth Package:

- ☐ 10"D-14"D
☐ 14"D-18"D
☐ 15"D-20"D

BACK CANE OPTIONS

Description	Price	HCPSC Code
Back Cane Height¹:		
☐ 14"H ² ☐ 15"H ² ☐ 16"H ³ ☐ 17"H ³ ☐ 18"H		
☐ 19"H ☐ 20"H ☐ 21"H ☐ 22"H ☐ 23"H		
☐ 24"H ☐ 25"H		
Cane Type¹:		
☐ Standard Back Canes	Standard	
☐ Omit Push Handles	\$125	
☐ Attendant Manual Recline ⁴	\$750	E1226
Cane Brackets:		
☐ Flip Forward Back Canes ⁵	\$150	
☐ Omit Push Handles	\$125	
☐ Offset (Meso) ⁶	\$195	
☐ 2" (1" per side) ☐ 4" (2" per side)		
☐ Inset (Endo) ⁶	\$195	
☐ 2" (1" per side) ☐ 4" (2" per side)		

1. No selection required with TRU-Comfort™ Back.
2. Selection of a Synergy® Back and cantilever armrests disqualifies 14" and 15"H back canes.
3. Selection of a Synergy® Back and cane mounted, flip up armrests disqualifies 16" and 17"H back canes.
4. Not available with Flip Back Arms, Cantilever Arms, Shape Back, TRU-Comfort back, AFP, rear accessories or power elevating seat. Requires 22"H canes. Pivot point 2.5" above seat pan.
5. Not available with vent tray, oxygen holder, cane mounted or cantilever armrests. May interfere with lateral supports and waterfall armpads.
6. Dual-Post, flip back armrests interfere with offset canes. Not available with attendant manual recline.

BACK OPTIONS

Description	Price	HCPSC Code
☐ Omit Back	No Charge	
Sling Backs:		
☐ Sling Back ¹	\$98	
Overall Height ² : ☐ 18"H ☐ 19"H ☐ 20"H ☐ 21"H ☐ 22"H		
☐ 23"H ☐ 24"H ☐ 25"H		
☐ Comfort Company Back	Various	

[See Comfort Company Order Form](#)**Synergy® Back:**

☐ Synergy® Back with 1.5" Thick Back Cushion	\$275	
Overall Height : ☐ 13"H ☐ 14"H ☐ 15"H ☐ 16"H ☐ 17"H		
☐ 18"H ☐ 19"H ☐ 20"H ☐ 21"H ☐ 22"H		
☐ 23"H ☐ 24"H ☐ 25"H		

Choose Lateral Extension³:

- ☐ 6" for 13-16"H Back \$25
☐ 10" for 17-18"H Back \$25
☐ 12" for 19-22"H Back \$25
☐ 16" for 23-25"H Back \$25

1. Not available with cantilever or flip, back, cane mounted armrests.
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
3. Lateral extensions are required to mount lateral thoracic supports.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

BACK OPTIONS (CONT'D)

Description	Price	HCPCS Code
TRU-Comfort™ Back: <i>Please refer to TRU-Comfort Back order form for detailed pricing, sizes, codes and part #'s</i>		
☐ TRU-Comfort™ 4-Way Stretch Back (20"W or less)	\$595	E2620
☐ TRU-Comfort™ 4-Way Stretch Back (22"W or greater)	\$695	E2621
Width: ☐ 16"W ☐ 18"W ☐ 20"W		
Height ¹ : ☐ 20"H ☐ 21"H ☐ 22"H ☐ 23"H ☐ 24"H		
☐ 25"H ☐ 27"H ☐ 28"H		
Width: ☐ 17"W ☐ 19"W		
Height ☐ 22"H ☐ 23"H ☐ 24"H ☐ 25"H ☐ 27"H		
☐ 28"H		
Width: ☐ 22"W ☐ 24"W		
Height : ☐ 24"H ☐ 25"H ☐ 27"H ☐ 28"H		

1. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".

ARMPADS

Option #	Description	Price
☐	Straight, Full Length (2"x14")	Standard
☐ Left ☐ Right		
☐	Straight, Desk Length (2"x10")	No Charge
☐ Left ☐ Right		
☐	Waterfall, Full Length (2"x14")	No Charge
☐ Left ☐ Right		
☐	Waterfall, Desk Length (2"x10")	No Charge
☐ Left ☐ Right		
☐	Flat, Bariatric ¹	No Charge
☐	Falcon, Gel, Full Length, Pair (3.5"x14")	\$92.50 ea.
☐	Gel Armpad Mounting Plate (Required)	\$56.69 ea.
☐	Falcon, Gel, Desk Length, Pair (3.5"x10")	\$70.00 ea.
☐	Gel Armpad Mounting Plate (Required)	\$53.52 ea.

1. Bariatric armrests are full length and require bariatric arm pads.

ARMRESTS

Option #	Description	Price	HCPCS Code
☐	Flip Back, Cane Mounted ¹	No Charge	
☐	2-Post, Flip-Back, Height Adjustable, Removable ²	\$150 ea.	E0973
☐	Single Post, Height Adjustable, Removable ³	\$75 ea.	E0973
☐	Quad Release Lever		
☐	Bariatric, Removable, Height Adjustable ⁴	\$75 ea.	E0973
Height: ☐ 8"H-12"H ☐ 10"H-14"H			
☐	Removable, Single Post, Height Adjustable, Heavy Duty	\$75 ea.	E0973
Height: ☐ 8"H-12"H ☐ 10"H-14"H			
☐	Pediatric, Removable Height Adjustable ⁵	\$75 ea.	E0973
☐	Quad Release Lever		
☐	Height Adjustable, Cantilever ⁶	\$195 pair	
	<i>Angle Adjustable: -30° through +30°</i>		

1. Interferes with Swing-Away laterals. Not available w/ manual recline, sling back or flip forward back canes
 2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16"
 3. Height adjustable from 11" - 17" in 1/2" increments.
 4. Bariatric armrests are full length and require bariatric arm pads.
 5. Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length arm pads.
 6. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, Shape Back or sling back.

LED LIGHT PACKAGE

Description	Price
☐ Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995

LAPBELTS

Option #	Description	Price
50BELT	☐ 50" Lap Belt	Standard
60BELT	☐ 60" Lap Belt	Standard
70BELT	☐ 70" Lap Belt	Standard
80BELT	☐ 80" Lap Belt	Standard
90BELT	☐ 90" Lap Belt	Standard
☐	60" Padded Push Button, Center Pull	\$115
☐	Large 13-19", Hip Padded, Push Button	\$92
☐	i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded 1.5" Belt	\$86
☐	i-Fit 60" Single Pull, Push Button, Large (9.5") Padded 2" Belt	\$86
☐	i-Fit 60" Single Pull, Push Button, Extra Large (11.5") Padded, 2" Belt	\$86

Items Designated with a Circular Checkbox (○) are not available through Quantum West

HEADRESTS

Option #	Description	Price	HCPSC Code
Headrest Pads:			
CP250	☐ Stealth 10" Comfort Plus	\$142	E0955
CP150	☐ Stealth 14" Comfort Plus	\$142	E0955
STL:OTSM250	☐ Stealth Small Contoured (9.75"W x 3"H)	\$60	E0955
STL:OTLD150	☐ Stealth Large Contoured (11"W x 6"H)	\$75	E0955
	☐ Omit Headrest		
Headrest Mounts:			
HMO475P	☐ Stealth Removable Headrest Bracket	\$150	E1028
TWB485-1	☐ Stealth's Worlds Best Removable Headrest Bracket	\$210	E1028
NOHR	☐ No Bracket		

ELECTRONICS FOR POWER FOOT RIGGINGS

Option #	Description	Price	HCPSC Code
THRUJOY	☐ Single Actuator Through Joystick (AFP Only)	\$1,895	E2310 ¹
THRUJOY	☐ Dual Actuator Through NE+ 6 Key	\$2,895	E2311 ¹
MTHRUJOY	☐ Multiple Actuator Through Q-Logic 2 EX	\$2,895	E2311 ¹

1. Applies to covered actuators only.

NON-POWERED LEGRESTS

Option #	Description	Price
----------	-------------	-------

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

☐ 60-90° Swing-Away Legrests \$162

(Adjustable in 5° Increments)

Rigging Selection:

- ☐ Standard Matching Legrests, pair
☐ Non-Matching Legrests, pair
☐ Individual Legrests (Left or Right Only)

Footplate Selection:

- ☐ Angle-Adjustable Footplates¹ \$75 ea.
☐ Standard, Extruded Footplates No Charge

Left Legrest:

Lower Extensions:

- ☐ 11"H-13"L
☐ 13"H-16.5"L²

Right Legrest:

Lower Extensions:

- ☐ 11"L-13"L
☐ 13"L-16.5"L

☐ Angle Adjustable/Multi-Axis Footplate Size

- | | |
|--|--|
| ☐ 4"W x 6"D | ☐ 4"W x 6"D |
| ☐ 5"W x 6"D | ☐ 5"W x 6"D |
| ☐ 5"W x 8"D | ☐ 5"W x 8"D |
| ☐ 6"W x 8"D ² | ☐ 6"W x 8"D ³ |
| ☐ 6"W x 11"D ³ | ☐ 6"W x 11"D ³ |
| ☐ 8"W x 11"D ³ | ☐ 8"W x 11"D ⁴ |

1. Includes Heel Loops.
 2. May cause interference on seat widths less than 18".
 3. May cause interference on seat widths less than 20".

NON-POWERED LEGRESTS (CONT'D)

Option #	Description	Price	HCPSC Code
☐ Standard, Extruded Footplate Size			
	☐ 7.25"W x 6"D ⁴	☐ 7.25"W x 6"D ⁵	
	☐ 8.5"W x 6"D ⁵	☐ 8.5"W x 6"D ⁶	
Accessories:			
	☐ Heel Loops, pair ⁶	\$38	E0951
BUMP	☐ Wheel Bumpers, pair	\$25	
LPW	☐ Gel Neoprene Footrest Sleeve	\$80 pr.	K0108
	☐ Therafin Multi-Axis Amputee Support ⁷	\$230	E1020
	☐ Amputee Sppt. Gel Cover	\$150	K0108

☐ Elevating Legrests \$325 K0195**Rigging Selection:**

- ☐ Standard Matching Legrests, pair
☐ Non-Matching Legrests, pair
☐ Individual Legrests (Left or Right Only)

Footplate Selection:

- ☐ Angle-Adjustable Footplates⁸ \$75 ea.
☐ Standard, Extruded Footplates No Charge

Left Legrest:

Lower Extensions:

- ☐ 11"-14"L⁹
☐ 12"H-14.5"L
☐ 14"-17"L¹⁰
☐ 14.5"-18"L

Right Legrest:

Lower Extensions:

- ☐ 11"-14"L¹⁰
☐ 12"L-14.5"L
☐ 14"-17"L¹⁰
☐ 14.5"-18"L

☐ Angle Adjustable Footplate Size

- | | |
|--|--|
| ☐ 4"W x 6"D | ☐ 4"W x 6"D |
| ☐ 5"W x 6"D | ☐ 5"W x 6"D |
| ☐ 5"W x 8"D | ☐ 5"W x 8"D |
| ☐ 6"W x 8"D ⁵ | ☐ 6"W x 8"D ⁵ |
| ☐ 6"W x 11"D ⁵ | ☐ 6"W x 11"D ⁵ |
| ☐ 8"W x 11"D ⁶ | ☐ 8"W x 11"D ⁶ |

☐ Standard, Extruded Footplate Size

- ☐ 7.25"W x 6"D⁵ ☐ 7.25"W x 6"D⁵
☐ 8.5"W x 6"D⁶ ☐ 8.5"W x 6"D⁶

Accessories:

- ☐ Heel Loops, pair⁷ \$38 E0951
 BUMP ☐ Wheel Bumpers, pair \$25
 LPW ☐ Gel Padded Adductor \$125 pr. K0108

4. May cause interference on seat widths less than 18".
 5. May cause interference on seat widths less than 20".
 6. Standard on Angle-Adjustable
 7. Swing-Away legrest upper assembly and receiver is required.
 8. Includes heel loops.
 9. Only available with standard, extruded footplates.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

NON-POWERED LEGRESTS (CONT'D)

Option #	Description	Price
☐ 90° Pediatric Swing-Away Legrests		\$162
Rigging Selection:		
☐	Standard Matching Legrests, pair	
☐	Non-Matching Legrests, pair	
☐	Individual Legrests (Left or Right Only)	
Footplate Selection:		
☐	Angle-Adjustable Footplates ¹	\$75 ea.
Left Legrest:		
Lower Extensions:		
☐	7"H-10"L	
Right Legrest:		
Lower Extensions:		
☐	7"L-10"L	
☐	Angle Adjustable Footplate Size	
☐	4"W x 6"D	☐ 4"W x 6"D
☐	5"W x 6"D	☐ 5"W x 6"D
Accessories:		
BUMP ☐	Wheel Bumpers, pair	\$25
☐ Heavy Duty Drop-In Legrests ¹		\$420
☐	70 Degree	
☐	80 degree	
Rigging Selection:		
☐	Standard Matching Legrests, pair	
☐	Non-Matching Legrests, pair	
☐	Individual Legrests (Left or Right Only)	
Footplate Selection:		
☐	Angle-Adjustable Footplates ²	\$75 ea.
Left Legrest:		
Lower Extensions:		
☐	12"H-14.5"L	
☐	14.5"-18"L	
Right Legrest:		
Lower Extensions:		
☐	12"L-14.5"L	
☐	14.5"-18"L	
☐	Angle Adjustable Footplate Size	
☐	4"W x 6"D	☐ 4"W x 6"D
☐	5"W x 6"D	☐ 5"W x 6"D
☐	5"W x 8"D	☐ 5"W x 8"D
☐	6"W x 8"D ³	☐ 6"W x 8"D ³
☐	6"W x 11"D ³	☐ 6"W x 11"D ³
☐	8"W x 11"D ⁴	☐ 8"W x 11"D ⁴
Accessories:		
BUMP ☐	Wheel Bumpers, pair	\$25

1. N/A with contoured seating.
2. Includes Heel Loops.
3. May cause interference on seat widths less than 18".
4. May cause interference on seat widths less than 20".

ADJUSTABLE CENTER MOUNT FOOT PLATFORM

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

☐ **Adjustable Center Mount Foot Platform** **\$250**

Footplates:

- ☐ 10"W x 11"D Split Footplates
- ☐ 10"W x 8"D Split Footplates
- ☐ 10"W x 6"D Split Footplates

Extensions:

CM1016 ☐ 10-16"L

Set Legrest Length:

- | | | |
|--|---|---|
| LL10 ☐ 10 | LL13 ☐ 13"H | LL16 ☐ 16"H |
| LL105 ☐ 10.5"H | LL135 ☐ 13.5"H | |
| LL11 ☐ 11"H | LL14 ☐ 14"H | |
| LL115 ☐ 11.5"H | LL145 ☐ 14.5"H | |
| LL12 ☐ 12"H | LL15 ☐ 15"H | |
| LL125 ☐ 12.5"H | LL15.5 ☐ 15.5"H | |

Calf Supports:

- ☐ **Calf Pads (pair)**
- ☐ Small (6"H x 6"W) \$150
- ☐ Large (8"H x 8"W, tapered to 6"W) \$150

☐ **Calf Pad Covers (pair)**

(1/2" Soft Visco Foam)

- ☐ Air Mesh \$120
- ☐ Reversed Dartex \$120

☐ **Calf Panel (ea.)**

- ☐ Small (8"H x 18"W) \$150
- ☐ Large (8"H x 24"W)¹ \$200

☐ **Calf Panel Covers**

(1/2" Soft Visco Foam)

- ☐ Air Mesh \$120
- ☐ Reversed Dartex \$120

Optional Modular Footplate Extension Components/Kits:²

- ☐ 4" Width Extension Kit \$50


- ☐ 2" Bell Shaped, Length Extension Kit \$50

- ☐ 4" Side Guard³ \$50

- ☐ 6" Side Guard⁴ \$50

Additional Optional Accessories:

- ☐ Wheel Bumpers, pair \$25
- ☐ Footplate Connector \$15



1. Only available on 20"W or greater.
2. Assembly required. May require seating configuration adjustments to avoid front caster interference.
3. Only available with 10"W x 6"D split footplates
4. Only available with 10"W x 8"D and 10"W x 11"D split footplates

Items Designated with a Circular Checkbox (○) are not available through Quantum West

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

SYNERGY® CUSHIONS

Option #	Description	Price	HCPCS Code
----------	-------------	-------	------------

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

☐ **Center Mount Power Articulating Foot Platform** **\$2995** See note¹

Footplates:

- ☐ 10"W x 11"D Split Footplates
☐ 10"W x 8"D Split Footplates
☐ 10"W x 6"D Split Footplates

Extensions:

CM1016 ☐ 10-16"L

Set Legrest Length:

- | | |
|---------------------------------------|---------------------------------------|
| LL10 ☐ 10"H | LL135 ☐ 13.5"H |
| LL105 ☐ 10.5"H | LL14 ☐ 14"H |
| LL11 ☐ 11"H | LL145 ☐ 14.5"H |
| LL115 ☐ 11.5"H | LL15 ☐ 15"H |
| LL12 ☐ 12"H | LL155 ☐ 15.5"H |
| LL125 ☐ 12.5"H | LL16 ☐ 16"H |
| LL13 ☐ 13"H | |

Calf Supports:

☐ **Calf Pads (pair)**

- ☐ Small (6"H x 6"W) No Charge
☐ Large (8"H x 8"W, tapered to 6"W) No Charge

☐ **Calf Pad Covers (pair)**

(1/2" Soft Visco Foam)

- ☐ Air Mesh \$120
☐ Reversed Dartex \$120

☐ **Calf Panel (ea.)**

- ☐ Small (8"H x 18"W) \$150
☐ Large (8"H x 24"W)² \$200

☐ **Calf Panel Covers**

(1/2" Soft Visco Foam)

- ☐ Air Mesh \$120
☐ Reversed Dartex \$120

Optional Modular Footplate Extension Components/Kits:³

- ☐ 4" Width Extension Kit \$50


- ☐ 2" Bell Shaped, Length Extension Kit \$50

- ☐ 4" Side Guard⁴ \$50

- ☐ 6" Side Guard⁵ \$50

Additional Optional Accessories:

- ☐ Footplate Connector \$15
☐ Wheel Bumpers, pair \$25



- ☐ Matching Seat Size
☐ Non-Matching Seat Size

Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s

Simplicity Cushion: General Use

- ☐ 16" x 16" - 20" x 20" \$98 E2601

Solution® Cushion: Skin Protection & Positioning

- ☐ 10" x 10" - 20" x 20" \$380 E2607

Solution®1 Cushion: Skin Protection

- ☐ 14" x 14" - 20" x 20" \$380 E2603

Spectrum® Gel Cushion: Positioning

- ☐ 12" x 10" - 20" x 20" \$480 E2605

Spectrum® Visco Cushion: Skin Protection & Positioning

- ☐ 10" x 10" - 20" x 20" \$480 E2607

Spectrum® Air Cushion: Adjustable Skin Protection

- ☐ 14" x 12" - 20" x 20" \$480 E2622

Air Tube:

- ☐ Left Side
☐ Right Side

Spectrum® Air Contour Cushion: Adjustable Skin Protection & Positioning

- ☐ 14" x 12" - 20" x 20" \$499 E2624

Air Tube:

- ☐ Left Side
☐ Right Side

TRU-Comfort™ 2 4-Way Stretch: Skin Protection:

- ☐ 16" x 14" - 20" x 20" \$380 E2603

CONTOURED SEATING/SYNERGY ACCESSORIES

Option #	Description	Price	HCPCS Code
----------	-------------	-------	------------

- ☐ Rear Basket¹ \$103
☐ Cane/Crutch Holder¹ \$98 E2207
☐ Weather Cover, Black \$158
☐ Cup Holder \$17²

Side:

- ☐ Left
☐ Right
☐ Saddle Bag \$28
☐ Flag Holder \$17
☐ Oxygen Holder
☐ Contoured Seating¹ \$103 E2208
☐ Synergy® Seat \$150 E2208
☐ Walker Holder \$92
☐ Tablet Holder³ \$315
☐ Deluxe Cell Phone/Media Holder³ \$140
☐ Transfer Bar, Pair³ \$175 E0705
 (Quick Release)

1. Please verify the appropriate billing code to be used for this item with all third-party payors.
 2. Only available on 20"W or greater.
 3. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 4. Only available with 10"W x 6"D split footplates
 5. Only available on 10"W x 8"D and 10"W x 11"D split footplates

1. Includes mounting hardware.
 2. No charge on Synergy seating
 3. Only available with Synergy seating.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

BATTERIES

Description	Price	HCPCS Code
☐ U-1, 32 Amp Batteries	\$133 ea.	E2365
☐ No Batteries		
Installation:		
☐ Batteries Installed		
☐ Do Not Install		
The following items can be chosen in the Positioning Components screen in iQuantum. Please refer to the Positioning Components order form for additional details, sizes and pricing		

LATERALS

Description	Price	HCPCS Code
☐ Stealth Thoracic Laterals¹		
<i>Bracket:</i>		
☐ Swing-Away, Rail Mounts	\$450 pr.	E1028
<i>Pads:</i>		
☐ Curved Pad, 3" x 5"	\$200 pr.	E0956
☐ Curved Pad, 4" x 6"	\$200 pr.	E0956
<i>Lateral Installation:</i>		
☐ Install		
☐ Do Not Install		
☐ Stealth Pelvic/Thigh Guides²		
<i>Side:</i>		
☐ Left		
☐ Right		
<i>Bracket:</i>		
☐ Flip Down	\$225 ea.	E1028
☐ Removable	\$250 ea.	E1028
☐ Flip Down w/ Full Surface Contact Feature	\$225 ea.	E1028
☐ Removable w/ 1" Offset	\$325 ea.	E1028
<i>Adapter:</i>		
☐ 2" Height Adapter Plate ³	\$25 ea.	K0108
☐ 2 1/4" Height Adapter Plate ⁴	\$25 ea.	K0108
<i>Pads:</i>		
☐ Short Pad, 4" x 4"	\$50 ea.	E0956
☐ Medium Pad, 4" x 8"	\$50 ea.	E0956
☐ Long Pad, 4" x 12"	\$50 ea.	E0956
<i>Lateral Installation:</i>		
☐ Install		
☐ Do Not Install		

1. Swing-away laterals may interfere with TB3 armrests. Not available with back cane option.

2. Pelvic/Thigh guides may interfere with the joystick receiver block.

3. Height adapter is needed w/ Tru Comfort seat cushions. Not available/necessary with the flip Down bracket.

4. Only available w/ the flip down brackets.

ADDITIONAL ACCESSORIES

Description	Price	HCPCS Code
☐ Stealth Foot Boxes (pair)	\$475	K0108
☐ Small		
☐ Medium		
☐ Large		
☐ Stealth Padded Shoe Holders (pair)	\$266	K0108
☐ Small		
☐ Medium		
☐ X-Large		
☐ Stealth Elbow Support	\$175 ea.	E1028
<i>Side:</i>		
☐ Bracket - Left		
☐ Bracket - Right		
<i>Pad:</i>		
☐ 3" x 3" (ARMS-ES-3x3)	\$50 ea.	K0108
☐ 4" x 4" (ARMS-ES-4x4)	\$50 ea.	K0108
☐ Stealth i-Fit Dynamic Chest Straps	\$142	E0960
☐ Upper Classic, Medium (IFPS1103)		
☐ Upper Classic, Large (IFPS1104)		
☐ Upper Contour, Medium (IFPS2103)		
☐ Upper Contour, Large (IFPS2104)		
☐ Stealth i-Fit Chest Strap Belt Anchor	\$90	
☐ Medium (IFPS4113)		
☐ Large (IFPS4114)		

Items Designated with a Circular Checkbox (○) are not available through Quantum West

ADDITIONAL HEADREST ACCESSORIES

SET-TO-SPEC*

Description Price HCPCS Code

- ☐ **Stealth Facial Laterals¹**
- ☐ Right Side, 4.5" x 2" Pad (CPS9-R-HW & 900-Q) \$325 E1028
- ☐ Left Side, 4.5" x 2" Pad (CPS9-L-HW & 900-Q) \$325 E1028
- ☐ Right and Left 4.5" x 2" Pad (CPS9-B-HW & 900-Q) \$425 E1028
- ☐ **Egg Switch Mount (ESM)** \$50 ea. K0108
- ☐ **Facial Pad, 6" x 2" (Q920-Q)** \$50 ea. K0108
- ☐ **Figure Eight Bracket with 13" Rod and Egg Switch Mount (F8678-13-Q)** \$250 K0108
- ☐ **Swing-away Joystick Mount, 18" Rod (SUS18-Q)** \$250 K0108

1. Only available with Stealth Comfort Series Headrests.

Please Supply Client Information:

Client Height: _____

Client Weight: _____

Mount Hip/Thigh Supports: ☐ Yes ☐ No

Mount Lateral Thoracic Supports: ☐ Yes ☐ No

Armrest Height from Pan: _____

Rigging Length (whole numbers only): _____

Seat to Floor Height (check one):

Seating System	Seat-To-Floor Heights	
Contoured Seating	20.5"	21.5"
Synergy Seat w/ Standard Towers	15.75"	16.75"
Synergy Seat w/ Extended Towers	16.5"	17.5"

Seat to floor height are +/- 1/8" with 18"W x 18"D with solid tires (add 3/8" for pneumatic tires).
Increased widths may increase lowest seat to floor heights.

* While every effort will be taken to ensure dimensions will be preset to your specifications prior to shipping variations in options may effect final settings.