

**Quantum® Q6 Edge with  
TB2 Power Positioning  
Group 3**



**QUANTUM REHAB®**

**Pride Mobility Products Corporation®**  
182 Susquehanna Ave., Exeter, PA 18643  
Phone: (866)800-2002 • Fax: (866) 707-3422  
Email: quantumorders@pridemobility.com

Discount: \_\_\_\_\_  
Terms: \_\_\_\_\_  
Freight: \_\_\_\_\_  
Taken By: \_\_\_\_\_

Account Number: \_\_\_\_\_  
Provider Name: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
Email: \_\_\_\_\_  
Date: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Ship to Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

**Note Section**

**For Quantum Use:**

**For Provider Use:**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Set-to-Spec**

*This section must be filled out, otherwise the selection will default to the mid-range measurements.*

Provide back height, armrest height, center mount height and legrest height.


**Set Back Height:** \_\_\_\_\_

**Set Armrest Height** Left: \_\_\_\_\_ Right: \_\_\_\_\_


**Set Center Mount Height:** \_\_\_\_\_

**Set Legrest Height** Left: \_\_\_\_\_ Right: \_\_\_\_\_

**Height:** \_\_\_\_\_ **Weight:** \_\_\_\_\_



- A. Shoulder Width \_\_\_\_\_
- B. Chest Width \_\_\_\_\_
- C. Pelvic Width \_\_\_\_\_
- D. Width at Knee \_\_\_\_\_



- E. Seat to Top of Head \_\_\_\_\_
- F. Seat to Top of Shoulder \_\_\_\_\_
- G. Seat to Axilla \_\_\_\_\_
- H. Seat to PSIS \_\_\_\_\_
- I. Chest Depth \_\_\_\_\_
- J. Back to ITS Anterior \_\_\_\_\_
- K. User Seat Depth \_\_\_\_\_
- L. Seat to Footplate \_\_\_\_\_

*DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.*



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Base Model

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Quantum Q6 Edge 3SP-SS</b> <i>Single power option. 300 lbs. wt capacity. Comes standard with 5 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24" W x 35.5" L.</i>	\$6,508	K0856
<input type="checkbox"/> <b>Quantum Q6 Edge 3MP-SS</b> <i>Multiple power option base. 300 lbs. wt capacity. Comes standard with 5 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24" W x 35.5" L.</i>	\$6,308	K0861

### Colors

**Metallic Colors:**

- ☐ ● Black  
☐ ● Orion Yellow  
☐ ● Silver  
☐ ● Viper Blue

- ☐ ● Candy Apple Red  
☐ ● Pearl Pink  
☐ ● Sunburst Orange

- ☐ ● Emerald Green  
☐ ● Sapphire Blue  
☐ ● True Blue

**Patterns: \$800** *(Available through SOF)*

- ☐ Mossy Oak

### Motor Upgrades

Description	Price	HCPCS Code
<input type="checkbox"/> Community Use Package <i>300 lbs. Weight Capacity; 4-Pole High Speed Motor Package; 4" wide flat free, or pneumatic, knobby, extended traction tires. Increases overall base width by 1.5"</i>	\$2,125	
<input type="checkbox"/> High Speed Motor Package <i>Programmable speed up to 6 Miles per Hour; High Speed Motors decrease per-charge range.</i>	\$1,574	
<input type="checkbox"/> Accu-Trac Motor Technology <i>Programmable speed up to 6 Miles per Hour; Delivers Enhanced Traction and Tracking, Especially with Switched Drive Controls. <b>Requires Q-Logic 2 EX.</b></i>	\$1,295	K0108

### Base Mounted Foot Platform

Description	Price	HCPCS Code
<input type="checkbox"/> Base Mounted Foot Platform	Standard	
<input type="checkbox"/> No Foot Platform	No Charge	

### Drive Wheels

Description	Price	HCPCS Code
<input type="checkbox"/> 3" x 14" Flat Free	Standard	
<input type="checkbox"/> 4" x 14" Flat Free, Knobby, Extended Traction Tires <i>Increases overall base width by 1.5"</i>	\$551	
<input type="checkbox"/> 3" x 14" Pneumatic	No Charge	
<input type="checkbox"/> 4" x 14" Pneumatic, Knobby, Extended Traction Tires <i>Increases overall base width by 1.5"</i>	\$551	

**Rim Color Options:**

- ☐ ● Silver *Standard*

- ☐ Matching

- ☐ Non-Matching

**Rim Color: \$50** *(Available through SOF)*

- ☐ ● Black  
☐ ● Orion Yellow  
☐ ● Sunburst Orange

- ☐ ● Candy Apple Red  
☐ ● Pearl Pink  
☐ ● True Blue

- ☐ ● Emerald Green  
☐ ● Sapphire Blue  
☐ ● Viper Blue

### Batteries

Description	Price	HCPCS Code
<input type="checkbox"/> NF-22, 50 Amp Batteries		
<input type="checkbox"/> Maintenance Free AGM (MK) (each)	\$295	E2361
<input type="checkbox"/> Maintenance Free Gel Cell (Universal) (each)	\$295	E2361
<input type="checkbox"/> Maintenance Free Gel Cell (MK) (each)	\$295	E2361
<input type="checkbox"/> No Batteries		

**Installation:**

- ☐ Batteries Installed  
☐ Do Not Install

### Rear Caster Beam

Description	Price	HCPCS Code
<input type="checkbox"/> Split Caster Beam	Standard	
<input type="checkbox"/> Extended Rear Beam <i>For Future Rear Accessory Upgrades; Automatically included with vent tray.</i>	No Charge	





### LED Light Package

Description	Price	HCPCS Code
<input type="checkbox"/> Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995	





# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Non-Expandable Controls

Description	Price	HCPCS Code
<input type="checkbox"/>  75 Amp NE	Standard	
<input type="checkbox"/>  75 Amp NE+ 6-Key on Tilt Only <i>For up to 2 actuators through joystick on single power base only</i>	No Charge	
<input type="checkbox"/>  75 Amp NE+ 6-Key <i>For up to 2 Actuators through joystick</i>	\$200	
<input type="checkbox"/>  75 Amp NE+ 8-Key <i>For Lights</i>	\$400	

## Expandable Controls

Description	Price	HCPCS Code
<input type="checkbox"/>  100 Amp Q-Logic 2 EX Controller <i>Required with Q-Logic Electronics</i>	\$695	E2377
<input type="checkbox"/>  Expandable Controller Harness <i>Required with Q-Logic Electronics</i>	\$350	E2313
<input type="checkbox"/>  Q-Logic 2 EX Joystick	\$700	
<input type="checkbox"/>  Q-Logic 2 EX Joystick <i>For Lights</i>	\$700	

### Color Options:

☐ ● Black Soft Touch *Standard*

☐ Matching



☐ Non-Matching

### Joystick Color Shroud: \$50

- ☐ ● Black  
☐ ● Orion Yellow  
☐ ● Silver  
☐ ● Viper Blue

- ☐ ● Candy Apple Red  
☐ ● Pearl Pink  
☐ ● Sunburst Orange

- ☐ ● Emerald Green  
☐ ● Sapphire Blue  
☐ ● True Blue

<input type="checkbox"/>  Stand Alone Joystick with External Charger Port	\$700	E2373
<b>Features:</b> <ul style="list-style-type: none"> <li>Up to 3 drive profiles and 1 seating profile when used alone.</li> <li>Use in conjunction w/ enhanced display for up to 5 profiles.</li> <li>Auxiliary functions including Bluetooth Mouse and IR.</li> </ul>		
<input type="checkbox"/>  Q-Logic 2 Enhanced Display Kit	\$700	
<b>Includes:</b> <ul style="list-style-type: none"> <li>Enhanced Display</li> <li>Flexible Gooseneck</li> <li>Gooseneck Mounting Bracket</li> </ul>		
<b>Features:</b> <ul style="list-style-type: none"> <li>Specialty controls interface w/ 3.5" LCD display</li> <li>Up to 5 profiles to control drive, seat and auxiliary functions</li> <li>Integrated single switch scanner, Bluetooth mouse emulator and IR controls</li> <li>Bluetooth connectivity for mobile devices</li> </ul>		

### Color Options:

☐ ● Black Soft Touch *Standard*

☐ Matching





☐ Non-Matching

### Display Color Shroud: \$50

- ☐ ● Black  
☐ ● Orion Yellow  
☐ ● Silver  
☐ ● Viper Blue

- ☐ ● Candy Apple Red  
☐ ● Pearl Pink  
☐ ● Sunburst Orange

- ☐ ● Emerald Green  
☐ ● Sapphire Blue  
☐ ● True Blue

<input type="checkbox"/>  Adult i-Drive Tri Array Series, 3-Switch w/ Egg Switch for Mode	\$4,175	E2330
<b>Stealth: IDH200-1</b> <b>Features:</b> 3 embedded proximity sensors for drive control, and an egg switch for mode change		
<input type="checkbox"/>  Egg Switches	\$77	
<input type="checkbox"/> ● Red <i>Stealth: ESRED</i> <input type="checkbox"/> ● Blue <i>Stealth: ESBLU</i> <input type="checkbox"/> ● Black <i>Stealth: ESBLK</i>	<input type="checkbox"/> ● Green <i>Stealth: ESGRN</i> <input type="checkbox"/> ● Yellow <i>Stealth: ESYEL</i>	
<input type="checkbox"/>  USB Mobile Device Charger	\$175	
Available with NE, NE+ and Q-Logic 2 EX controllers; Can only charge 1 device at a time. A multiplier harness is required with all Static and Tilt only seating systems. The selection of more than 3 modules requires an additional multiplier.		
<input type="checkbox"/>  <b>NEW!!! SCIM Module</b>	\$575	
Allows use of joystick display for specialty controls; Up to 5 profiles to control drive, seat and auxiliary functions; Excludes sip-n-puff and single switch functions; A switch is required to operate the input device attached to the SCIM module; Includes accessory bar on TB3; <i>Only available with Q-Logic 2.</i>		
<input type="checkbox"/> <a href="#">Click Here for Alternative Drive Devices</a>		



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Joystick Handles

	Standard	
<input type="checkbox"/> Standard Carrot Handle	\$71	E2323
<input type="checkbox"/> Mushroom Handle	\$107	E2323
<input type="checkbox"/> T-Bar	\$71	E2324
<input type="checkbox"/> Chin Cup <i>Use HCPCS Code E2323 when Chin Cup not used in conjunction with a chin control.</i>	\$71	E2323
<input type="checkbox"/> Extended I-Handle	\$50	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 3 Inch Goal Post Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 4 Inch Goal Post Handle		

### Joystick Mounting Brackets

	Standard	
<input type="checkbox"/> Fixed, Inline <i>Not available with Q-Logic 2 EX Stand Alone Joystick.</i>		
<input type="checkbox"/> Left		<input type="checkbox"/> Right
<input type="checkbox"/> Swing-Away	\$210	E1028
<input type="checkbox"/> Left		<input type="checkbox"/> Right
<input type="checkbox"/> <b>NEW!!!</b> Retractable Swing-Away Joystick Arm Kit <i>Full 180° Movement</i>	\$250	E1028
<input type="checkbox"/> Left		<input type="checkbox"/> Right
<input type="checkbox"/> Multi-Axis	\$100	
<input type="checkbox"/> Left		<input type="checkbox"/> Right
<input type="checkbox"/> No Joystick Mount		

### Programming

<input type="checkbox"/>  Q-Logic EX/NE Programmer	\$545
<input type="checkbox"/> Q-Logic PC Programming Disk w/ Cable <i>Works with Q-Logic 2 EX only.</i>	\$750

### TRU-Balance®2 Power Positioning

Description	Price	HCPCS Code
<input type="checkbox"/> Tru-Balance®2 Power Tilt <i>55 degree center of gravity (CG) power tilt seating system</i>	\$4,995	E1002
<input type="checkbox"/> Tru-Balance®2 Power Recline <i>168 degrees of power recline with full drive lockout at 20 degrees</i>	\$5,595	E1004
<input type="checkbox"/> Tru-Balance®2 Power Tilt & Recline <i>55 Degree center of gravity (CG) power tilt seating system with 168 degrees of power recline; Full drive lock-out at 20 degrees</i>	\$9,595	E1007

#### Power Seat Option:

<input type="checkbox"/> Power Elevating Seat <i>Maximum seat size 20"W x 20"D. Increases STF Height. Includes 1/4 speed inhibit.</i>	\$1,895
<input type="checkbox"/> 6" Stroke	
<input type="checkbox"/> 7" Stroke	
<input type="checkbox"/> Through NE+ 6-Key	
<input type="checkbox"/> Through Q-Logic 2 EX	

### TRU-Balance®2 Manual Positioning Option

Description	Price	HCPCS Code
<input type="checkbox"/> Attendant Manual Recline <i>Available w/ Tilt only. Not available w/ Flip Back armrests, Cantilever armrests, Shape Back, TRU-Comfort back, Rear accessories or power elevating seat. Requires 22"H Canes. Pivot point 2.5" above seat pan.</i>	\$750	E1226



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Transit Kit

Description	Price	HCPCS Code
<input type="checkbox"/> Unoccupied Transit Kit <i>Unoccupied transit kit w/ a vent tray may not allow the lowest seat to floor height.</i>	No Charge	

### Power Positioning Electronics

Description	Price	HCPCS Code
<b>Single Actuator Electronics:</b>		
<input type="checkbox"/> Through Toggle Box <i>Tilt Only. Manual Recline back is not available through NE tilt through toggle.</i>	No Charge	
<input type="checkbox"/> Left Mount		
<input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through NE+ 6 Key	\$1175	E2310
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$1175	E2310
<input type="checkbox"/> Through Single Switch <i>Tilt only w/ Q-Logic 2 EX. Single switch not included.</i>	\$500	
<b>Dual Actuator Electronics:</b>		
<input type="checkbox"/> Through NE+ 6-Key Controller	\$2,375	E2310/E2311 <sup>1</sup>
<b>Multiple Actuator Electronics:</b>		
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$2,375	E2311
<input type="checkbox"/> Q-Logic 2 EX Through 4-way Toggle <sup>2</sup>	\$2,375	
<input type="checkbox"/> Left Mount		
<input type="checkbox"/> Right Mount		
<input type="checkbox"/> Q-Logic 2 EX Through 4-way Button Toggle <sup>2</sup>	\$2,375	
<input type="checkbox"/> Left Mount		
<input type="checkbox"/> Right Mount		
<input type="checkbox"/> Q-Logic 2 EX Through Quad Push Button <sup>2</sup>	\$2,375	
<input type="checkbox"/> Left Mount		
<input type="checkbox"/> Right Mount		

1. Applies to covered actuators only
2. Bracket mounts to the outside of the armpad

### TRU-Balance®2 Seat Specifications

**Seat Width:**

- ☐ 14"W  
☐ 15"W  
☐ 16"W  
☐ 17"W  
☐ 18"W  
☐ 19"W  
☐ 20"W  
☐ 21"W  
☐ 22"W

**Seat Width Package:**

- ☐ 14"W-18"W  
☐ 16"W-20"W  
☐ 20"W-24"W

**Seat Depth:**

- ☐ 14"D  
☐ 15"D  
☐ 16"D  
☐ 17"D  
☐ 18"D  
☐ 19"D  
☐ 20"D  
☐ 21"D  
☐ 22"D

**Seat Depth Package:**

- ☐ 14"D-18"D  
☐ 16"D-20"D  
☐ 20"D-24"D

### TRU-Balance®2 Seat to Floor Height

**Tilt or Tilt/Recline:**

- ☐ 17.125" Seat to Floor  
☐ 18.125" Seat to Floor  
☐ 19.125" Seat to Floor

**Recline:**

- ☐ 15.875" Seat to Floor  
☐ 16.875" Seat to Floor  
☐ 17.875" Seat to Floor

**Tilt or Tilt/Recline w/ Power Elevating:**

- ☐ 18.5" -24.5" w/ 6" actuator  
☐ 19.5" -25.5" w/ 7" actuator

**Recline w/ Power Elevating Seat:**

- ☐ 15.5" -21.5" w/ 6" actuator  
☐ 16.5" -22.5" w/ 7" actuator



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### TRU-Balance® 2 TILT Back Cane Options

**Back Cane Height<sup>1</sup>:**

- |  |                               |                               |
|--|-------------------------------|-------------------------------|
| <input type="checkbox"/> 14"H <sup>2</sup> | <input type="checkbox"/> 18"H | <input type="checkbox"/> 22"H |
| <input type="checkbox"/> 15"H <sup>2</sup> | <input type="checkbox"/> 19"H | <input type="checkbox"/> 23"H |
| <input type="checkbox"/> 16"H <sup>3</sup> | <input type="checkbox"/> 20"H | <input type="checkbox"/> 24"H |
| <input type="checkbox"/> 17"H <sup>3</sup> | <input type="checkbox"/> 21"H | <input type="checkbox"/> 25"H |

**Cane Type<sup>1</sup>:**

- ☐ Standard Back Canes  
☐ Omit Push Handles

Standard  
\$125

**Cane Brackets:**

- ☐ Flip Forward Back Canes<sup>4</sup>  
☐ Omit Push Handles

\$150  
\$125

1. No selection required with TRU-Comfort™ Back.
2. Selection of a Synergy® Back and cantilever armrests disqualifies 14" and 15"H back canes.
3. Selection of a Synergy® Back and cane mounted, flip up armrests disqualifies 16" and 17"H back canes.
4. Not available with vent tray, oxygen holder, cane mounted or cantilever armrests. May interfere with lateral supports and waterfall arm pads.

### TRU-Balance® 2 TILT Back Specifications

+/-4" of Seat Width. Measurements listed is overall back height.



**TRU-Comfort Back** 4-way stretch fabric. \$595 E2620 for 20"W or less, \$695 E2621 for 22"W or greater

\$595/\$695 E2620/E2621



**Synergy® Back** 1.5" thick back cushion

No Charge

**Size:**

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 14"Wx15"H <sup>1</sup> | <input type="checkbox"/> 17"Wx22"H              | <input type="checkbox"/> 21"Wx19"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx16"H <sup>1</sup> | <input type="checkbox"/> 17"Wx23"H              | <input type="checkbox"/> 21"Wx20"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx17"H <sup>1</sup> | <input type="checkbox"/> 17"Wx24"H              | <input type="checkbox"/> 21"Wx21"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx18"H <sup>1</sup> | <input type="checkbox"/> 17"Wx25"H              | <input type="checkbox"/> 21"Wx22"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx19"H <sup>1</sup> | <input type="checkbox"/> 18"Wx17"H <sup>1</sup> | <input type="checkbox"/> 21"Wx23"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx20"H <sup>1</sup> | <input type="checkbox"/> 18"Wx18"H <sup>1</sup> | <input type="checkbox"/> 21"Wx24"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx21"H <sup>1</sup> | <input type="checkbox"/> 18"Wx19"H <sup>1</sup> | <input type="checkbox"/> 21"Wx25"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx22"H <sup>1</sup> | <input type="checkbox"/> 18"Wx20"H              | <input type="checkbox"/> 22"Wx17"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx23"H <sup>1</sup> | <input type="checkbox"/> 18"Wx21"H              | <input type="checkbox"/> 22"Wx18"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx24"H <sup>1</sup> | <input type="checkbox"/> 18"Wx22"H              | <input type="checkbox"/> 22"Wx19"H <sup>1</sup> |
| <input type="checkbox"/> 14"Wx25"H <sup>1</sup> | <input type="checkbox"/> 18"Wx23"H              | <input type="checkbox"/> 22"Wx20"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx17"H <sup>1</sup> | <input type="checkbox"/> 18"Wx24"H              | <input type="checkbox"/> 22"Wx21"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx18"H <sup>1</sup> | <input type="checkbox"/> 18"Wx25"H              | <input type="checkbox"/> 22"Wx22"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx19"H <sup>1</sup> | <input type="checkbox"/> 19"Wx17"H <sup>1</sup> | <input type="checkbox"/> 22"Wx23"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx20"H <sup>1</sup> | <input type="checkbox"/> 19"Wx18"H <sup>1</sup> | <input type="checkbox"/> 22"Wx24"H              |
| <input type="checkbox"/> 15"Wx21"H <sup>1</sup> | <input type="checkbox"/> 19"Wx19"H <sup>1</sup> | <input type="checkbox"/> 22"Wx25"H              |
| <input type="checkbox"/> 15"Wx22"H <sup>1</sup> | <input type="checkbox"/> 19"Wx20"H              | <input type="checkbox"/> 23"Wx17"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx23"H <sup>1</sup> | <input type="checkbox"/> 19"Wx21"H              | <input type="checkbox"/> 23"Wx18"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx24"H <sup>1</sup> | <input type="checkbox"/> 19"Wx22"H              | <input type="checkbox"/> 23"Wx19"H <sup>1</sup> |
| <input type="checkbox"/> 15"Wx25"H <sup>1</sup> | <input type="checkbox"/> 19"Wx23"H              | <input type="checkbox"/> 23"Wx20"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx17"H <sup>1</sup> | <input type="checkbox"/> 19"Wx24"H              | <input type="checkbox"/> 23"Wx21"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx18"H <sup>1</sup> | <input type="checkbox"/> 19"Wx25"H              | <input type="checkbox"/> 23"Wx22"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx19"H <sup>1</sup> | <input type="checkbox"/> 20"Wx17"H <sup>1</sup> | <input type="checkbox"/> 23"Wx23"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx20"H              | <input type="checkbox"/> 20"Wx18"H <sup>1</sup> | <input type="checkbox"/> 23"Wx24"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx21"H              | <input type="checkbox"/> 20"Wx19"H <sup>1</sup> | <input type="checkbox"/> 23"Wx25"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx22"H              | <input type="checkbox"/> 20"Wx20"H              | <input type="checkbox"/> 24"Wx17"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx23"H              | <input type="checkbox"/> 20"Wx21"H              | <input type="checkbox"/> 24"Wx18"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx24"H              | <input type="checkbox"/> 20"Wx22"H              | <input type="checkbox"/> 24"Wx19"H <sup>1</sup> |
| <input type="checkbox"/> 16"Wx25"H              | <input type="checkbox"/> 20"Wx23"H              | <input type="checkbox"/> 24"Wx20"H <sup>1</sup> |
| <input type="checkbox"/> 17"Wx17"H              | <input type="checkbox"/> 20"Wx24"H              | <input type="checkbox"/> 24"Wx21"H <sup>1</sup> |
| <input type="checkbox"/> 17"Wx18"H <sup>1</sup> | <input type="checkbox"/> 20"Wx25"H              | <input type="checkbox"/> 24"Wx22"H <sup>1</sup> |
| <input type="checkbox"/> 17"Wx19"H <sup>1</sup> | <input type="checkbox"/> 21"Wx17"H <sup>1</sup> | <input type="checkbox"/> 24"Wx23"H <sup>1</sup> |
| <input type="checkbox"/> 17"Wx20"H <sup>1</sup> | <input type="checkbox"/> 21"Wx18"H <sup>1</sup> | <input type="checkbox"/> 24"Wx24"H              |
| <input type="checkbox"/> 17"Wx21"H <sup>1</sup> |   |   |

1. Not available with TRU-Comfort Back

TILT BACK




Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

TRU-Balance® 2 RECLINE OR TILT/RECLINE Back Specifications

Measurements listed is overall back height.

☐



**TRU-Comfort Back** \$595 E2620 for 20"W or less, \$695 E2621 for 22"W or greater

☐ 4-way stretch

☐ Air mesh

\$595/\$695

E2620/E2621

☐



**Sport Back**

No Charge

Size:

- ☐
- 14"Wx17"H
- <sup>1</sup>
- ☐
- 14"Wx18"H
- <sup>1</sup>
- ☐
- 14"Wx19"H
- <sup>1</sup>
- ☐
- 14"Wx20"H
- <sup>1</sup>
- ☐
- 14"Wx21"H
- <sup>1</sup>
- ☐
- 14"Wx22"H
- <sup>1</sup>
- ☐
- 14"Wx23"H
- <sup>1</sup>
- ☐
- 14"Wx24"H
- <sup>1</sup>
- ☐
- 14"Wx25"H
- <sup>1</sup>
- ☐
- 14"Wx26"H
- <sup>1</sup>
- ☐
- 14"Wx27"H
- <sup>1</sup>
- ☐
- 14"Wx28"H
- <sup>1</sup>
- ☐
- 15"Wx19"H
- <sup>1</sup>
- ☐
- 15"Wx20"H
- <sup>1</sup>
- ☐
- 15"Wx21"H
- <sup>1</sup>
- ☐
- 15"Wx22"H
- <sup>1</sup>
- ☐
- 15"Wx23"H
- <sup>1</sup>
- ☐
- 15"Wx24"H
- <sup>1</sup>
- ☐
- 15"Wx25"H
- <sup>1</sup>
- ☐
- 15"Wx26"H
- <sup>1</sup>
- ☐
- 15"Wx27"H
- <sup>1</sup>
- ☐
- 15"Wx28"H
- <sup>1</sup>
- ☐
- 16"Wx19"H
- <sup>1</sup>
- ☐
- 16"Wx20"H
- <sup>1</sup>
- ☐
- 16"Wx21"H
- <sup>1</sup>
- ☐
- 16"Wx22"H
- <sup>1</sup>
- ☐
- 16"Wx23"H
- <sup>1</sup>
- ☐
- 16"Wx24"H
- <sup>1</sup>
- ☐
- 16"Wx25"H
- <sup>1</sup>
- ☐
- 16"Wx26"H
- <sup>1</sup>
- ☐
- 16"Wx27"H
- <sup>1</sup>
- ☐
- 16"Wx28"H
- <sup>1</sup>
- ☐
- 17"Wx19"H
- <sup>1</sup>
- ☐
- 17"Wx20"H
- <sup>1</sup>
- ☐
- 17"Wx21"H
- <sup>1</sup>
- ☐
- 17"Wx22"H
- <sup>1</sup>
- ☐
- 17"Wx23"H
- <sup>1</sup>
- ☐
- 17"Wx24"H
- <sup>1</sup>

1. Not available with TRU-Comfort Back

- ☐
- 17"Wx25"H
- ☐
- 17"Wx26"H
- ☐
- 17"Wx27"H
- ☐
- 17"Wx28"H
- ☐
- 18"Wx19"H
- <sup>1</sup>
- ☐
- 18"Wx20"H
- <sup>1</sup>
- ☐
- 18"Wx21"H
- <sup>1</sup>
- ☐
- 18"Wx22"H
- <sup>1</sup>
- ☐
- 18"Wx23"H
- <sup>1</sup>
- ☐
- 18"Wx24"H
- <sup>1</sup>
- ☐
- 18"Wx25"H
- <sup>1</sup>
- ☐
- 18"Wx26"H
- <sup>1</sup>
- ☐
- 18"Wx27"H
- <sup>1</sup>
- ☐
- 18"Wx28"H
- <sup>1</sup>
- ☐
- 19"Wx19"H
- <sup>1</sup>
- ☐
- 19"Wx20"H
- <sup>1</sup>
- ☐
- 19"Wx21"H
- <sup>1</sup>
- ☐
- 19"Wx22"H
- <sup>1</sup>
- ☐
- 19"Wx23"H
- <sup>1</sup>
- ☐
- 19"Wx24"H
- <sup>1</sup>
- ☐
- 19"Wx25"H
- <sup>1</sup>
- ☐
- 19"Wx26"H
- <sup>1</sup>
- ☐
- 19"Wx27"H
- <sup>1</sup>
- ☐
- 19"Wx28"H
- <sup>1</sup>
- ☐
- 20"Wx19"H
- <sup>1</sup>
- ☐
- 20"Wx20"H
- <sup>1</sup>
- ☐
- 20"Wx21"H
- <sup>1</sup>
- ☐
- 20"Wx22"H
- <sup>1</sup>
- ☐
- 20"Wx23"H
- <sup>1</sup>
- ☐
- 20"Wx24"H
- <sup>1</sup>
- ☐
- 20"Wx25"H
- <sup>1</sup>
- ☐
- 20"Wx26"H
- <sup>1</sup>
- ☐
- 20"Wx27"H
- <sup>1</sup>
- ☐
- 20"Wx28"H
- <sup>1</sup>
- ☐
- 21"Wx19"H
- <sup>1</sup>
- ☐
- 21"Wx20"H
- <sup>1</sup>
- ☐
- 21"Wx21"H
- <sup>1</sup>
- ☐
- 21"Wx22"H
- <sup>1</sup>

- ☐
- 21"Wx23"H
- <sup>1</sup>
- ☐
- 21"Wx24"H
- <sup>1</sup>
- ☐
- 21"Wx25"H
- <sup>1</sup>
- ☐
- 21"Wx26"H
- <sup>1</sup>
- ☐
- 21"Wx27"H
- <sup>1</sup>
- ☐
- 21"Wx28"H
- <sup>1</sup>
- ☐
- 22"Wx19"H
- <sup>1</sup>
- ☐
- 22"Wx20"H
- <sup>1</sup>
- ☐
- 22"Wx21"H
- <sup>1</sup>
- ☐
- 22"Wx22"H
- <sup>1</sup>
- ☐
- 22"Wx23"H
- <sup>1</sup>
- ☐
- 22"Wx24"H
- <sup>1</sup>
- ☐
- 22"Wx25"H
- <sup>1</sup>
- ☐
- 22"Wx26"H
- <sup>1</sup>
- ☐
- 22"Wx27"H
- <sup>1</sup>
- ☐
- 22"Wx28"H
- <sup>1</sup>
- ☐
- 23"Wx19"H
- <sup>1</sup>
- ☐
- 23"Wx20"H
- <sup>1</sup>
- ☐
- 23"Wx21"H
- <sup>1</sup>
- ☐
- 23"Wx22"H
- <sup>1</sup>
- ☐
- 23"Wx23"H
- <sup>1</sup>
- ☐
- 23"Wx24"H
- <sup>1</sup>
- ☐
- 23"Wx25"H
- <sup>1</sup>
- ☐
- 23"Wx26"H
- <sup>1</sup>
- ☐
- 23"Wx27"H
- <sup>1</sup>
- ☐
- 23"Wx28"H
- <sup>1</sup>
- ☐
- 24"Wx19"H
- <sup>1</sup>
- ☐
- 24"Wx20"H
- <sup>1</sup>
- ☐
- 24"Wx21"H
- <sup>1</sup>
- ☐
- 24"Wx22"H
- <sup>1</sup>
- ☐
- 24"Wx23"H
- <sup>1</sup>
- ☐
- 24"Wx24"H
- <sup>1</sup>
- ☐
- 24"Wx25"H
- <sup>1</sup>
- ☐
- 24"Wx26"H
- <sup>1</sup>
- ☐
- 24"Wx27"H
- <sup>1</sup>
- ☐
- 24"Wx28"H
- <sup>1</sup>

TILT/RECLINE BACK  
RECLINE OR

TRU-Balance® 2 Back Specifications

Color Shroud Option:

- ☐ ● Black *Standard*
- Color Options: \$100
- ☐ ● Black
- ☐ ● Orion Yellow
- ☐ ● Silver
- ☐ ● Viper Blue

- ☐ ● Candy Apple Red
- ☐ ● Pearl Pink
- ☐ ● Sunburst Orange

- ☐ ● Emerald Green
- ☐ ● Sapphire Blue
- ☐ ● True Blue



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Additional TRU-Balance® 2 Back Options

<input type="checkbox"/>	Aftermarket Back Interface	\$200
	<i>Overall seat depth may need to be adjusted to accommodate certain backs. Contact your Quantum specialist for details. Does not include brackets to mount aftermarket back. SOF is required on back widths less than 16".</i>	
<input type="checkbox"/>	Omit Back	(-\$100)
<input type="checkbox"/>	Comfort Company Back <i>See Comfort Company order form</i>	Varies
<input type="checkbox"/>	Offset (Meso) <i>Dual-Post, flip back armrests interfere with offset canes. Not available with attendant manual recline.</i>	\$195
<input type="checkbox"/>	<input type="checkbox"/> 2" (1" per side) <input type="checkbox"/> 4" (2" per side)	
<input type="checkbox"/>	Inset (Endo) <i>Dual-Post, flip back armrests interfere with offset canes. Not available with attendant manual recline.</i>	\$195
<input type="checkbox"/>	<input type="checkbox"/> 2" (1" per side) <input type="checkbox"/> 4" (2" per side)	

### Armrests

#### Tilt Armrests:

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Flip Back, Cane Mounted</b> <i>Interferes with Swing-Away laterals. Not available w/ manual recline, flip forward back canes, or TRU-Comfort back.</i>	No Charge	
<input type="checkbox"/> <b>2-Post, Flip-Back, Height Adjustable, Removable</b> <i>Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16"</i>	No Charge	
<input type="checkbox"/> <b>Single Post, Height Adjustable, Removable</b> <i>Height adjustment from 11" - 17" in 1/2" increments</i>	No Charge	
<input type="checkbox"/> <input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> <b>Removable, Single Post, Height Adjustable, Heavy Duty</b>	No Charge	
<b>Armrest Height Package:</b>		
<input type="checkbox"/> 8"-12" Height		
<input type="checkbox"/> 10"-14" Height		
<input type="checkbox"/> <b>Pediatric, Removable, Height Adjustable</b> <i>Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length armpads.</i>	No Charge	
<input type="checkbox"/> <input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> <b>Height Adjustable, Cantilever</b> <i>Angle adjustable: -30° through +30°</i> <i>Height adjustment from 11" - 17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, or Shape Back.</i>	\$195	

#### Recline Armrests:

<input type="checkbox"/> <b>Reclining, Flip-Back</b>	Standard
<b>Armrest Height Package:</b>	
<input type="checkbox"/> 8"-12" Height	
<input type="checkbox"/> 10"-14" Height	
<input type="checkbox"/> <b>Height Adjustable, Cantilever</b> <i>Angle adjustable: -30° through +30°</i> <i>Height adjustment from 11" - 17" depends on cane height.</i>	\$195

#### Armpad:

##### Left Armpad Length:

- ☐ 14"L Full Length  
☐ 10"L Desk Length

##### Left Armpad Style:

- ☐ Straight  
☐ Waterfall

<input type="checkbox"/> Falcon, Gel, Full Length 3"W x 14"L (each)	\$92.50
<input type="checkbox"/> Mounting Plate, Full Length <sup>1</sup> (each)	\$56.52
<input type="checkbox"/> Falcon, Gel, Desk Length 3"W x 10"L (each)	\$70.00
<input type="checkbox"/> Mounting Plate, Desk Length <sup>1</sup> (each)	\$56.20

##### Right Armpad Length:

- ☐ 14"L Full Length  
☐ 10"L Desk Length

##### Right Armpad Style:

- ☐ Straight  
☐ Waterfall





<input type="checkbox"/> Falcon, Gel, Full Length 3"W x 14"L (each)	\$92.50
<input type="checkbox"/> Mounting Plate, Full Length <sup>1</sup> (each)	\$56.52
<input type="checkbox"/> Falcon, Gel, Desk Length 3"W x 10"L (each)	\$70.00
<input type="checkbox"/> Mounting Plate, Desk Length <sup>1</sup> (each)	\$56.20

1. Required only with 2-post, flip back, height adjustable, removable armrests.





## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Vent Trays

	Description	Price	HCPCS Code
<input type="checkbox"/>	 <b>Compact Vent Tray</b> The addition of a compact vent tray will increase the length & weight of your power chair. Not available with Flip-Up Armrests or power elevating seat <ul style="list-style-type: none"> <li>Accommodates Pulmonetics LTV-style ventilators</li> <li>Max Vent. Dim. 3.25"H x 10.5"W x 13.5"D / 20 lbs.</li> <li>Not available with Shape Back, Sling Back, Manual Recline, Flip Down Back Canes or rear mounted accessories.</li> </ul> <b>Mounting Bracket:</b> Bracket is pre-selected with TRU-Comfort Back <ul style="list-style-type: none"> <li><input type="checkbox"/> Flush</li> <li><input type="checkbox"/> 3" Offset</li> <li><input type="checkbox"/> 6" Offset</li> </ul>	\$1,395	E1029
<input type="checkbox"/>	 <b>Articulating Vent Tray</b> 14.5"W x 10.5-15.5"D Not available with TRU-Balance power elevating seat, Manual Recline or with rear mounted accessories. Must choose the extended rear beam on page 2 of base order form.	\$1,450	E1030
<input type="checkbox"/>	 <b>Articulating Bi-Pap Tray</b> 9.5"W x 9.5-13.5"D Not available with TRU-Balance power elevating seat, Manual Recline or with rear mounted accessories. Must choose the extended rear beam on page 2 of base order form.	\$1,450	E1030
<b>Battery Tray Selection:</b>			
<input type="checkbox"/>	 <b>Battery Tray</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Include U1 Battery Box</li> <li><input type="checkbox"/> Include U1 Battery</li> </ul>	\$450 \$50 \$433	K0108 K0108 E2365

### Laterals




<input type="checkbox"/>	 <b>Stealth Thoracic Lateral</b> Swing-Away laterals may interfere with cantilever and 2-post, flip-up armrest. <i>Not available w/ Back Cane Option</i>		
	<b>Bracket:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Swing-Away, Rail Mounts (pair)</li> <li><input type="checkbox"/> 125° Swing-Away Bracket For use on TRU-Comfort Recline Back (pair)</li> </ul>	\$450 \$450	E1028 E1028
	<b>Pads:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Curved Pad, 2" x 3" (pair)</li> <li><input type="checkbox"/> Curved Pad, 3" x 5" (pair)</li> <li><input type="checkbox"/> Curved Pad, 4" x 6" (pair)</li> </ul>	\$200 \$200 \$200	E0956 E0956 E0956
	<b>Lateral Installation:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Install</li> <li><input type="checkbox"/> Do not Install</li> </ul>		
<input type="checkbox"/>	 <b>Stealth Pelvic/Thigh Guides</b> Pelvic/Thigh guides may interfere with the joystick receiver block.		
	<b>Side:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Left</li> <li><input type="checkbox"/> Right</li> </ul>		
	<b>Bracket:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Flip Down (each)</li> <li><input type="checkbox"/> <b>NEW!!!</b> Flip Down with Full Surface Contact Feature (each)</li> <li><input type="checkbox"/> Removable (each)</li> <li><input type="checkbox"/> Removable w/ 1" Offset (each)</li> </ul>	\$225 \$225 \$250 \$325	E1028 E1028 E1028 E1028
	<b>Adapter:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> 2"Height Adapter Plate (each)</li> <li><input type="checkbox"/> <b>NEW!!!</b> 2 1/4" Adapter Plate (each) <i>Only available/necessary with the flip down bracket</i></li> </ul>	\$25 \$50	K0108 K0108
	<b>Pads:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Short Pad, 4" x 4" (each)</li> <li><input type="checkbox"/> Medium Pad, 4" x 8" (each)</li> <li><input type="checkbox"/> Long Pad, 4" x 12" (each)</li> </ul>	\$50 \$50 \$50	E0956 E0956 E0956
	<b>Lateral Installation:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Install</li> <li><input type="checkbox"/> Do not Install</li> </ul>		



# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Center Mount Foot Platforms

All efforts will be made to accommodate legrest length selection, but minor changes may be needed to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of unit.  
Additional foot riggings can be found on page 11

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Center Mount Power Articulating Foot Platform</b> Includes 10"W x 11"D Split Footplates	\$2,895	K0108
<input type="checkbox"/> <b>Adjustable Center Mount Foot Platform</b> Includes 10"W x 11"D Split Footplates	\$250	
<b>Extensions:</b>		
<input type="checkbox"/> 10-16"L		
<input type="checkbox"/> 13-19"L		
<b>Calf Supports:</b>		
<input type="checkbox"/> <b>Calf Pads (pair)</b>		
<input type="checkbox"/> Small <sup>1</sup> (6"H x 6"W)	\$150	
<input type="checkbox"/> Large <sup>1</sup> (8"H x 8"W, tapered to 6"W)	\$150	
<input type="checkbox"/> <b>NEW!!! Calf Pad Covers (pair) (1/2" Soft Visco Foam)</b>		
<input type="checkbox"/> Air Mesh	\$120	
<input type="checkbox"/> Reversed Dartex	\$120	
<input type="checkbox"/> <b>Calf Panel (each)</b>		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) Only available on 20"W or greater.	\$200	
<input type="checkbox"/> <b>NEW!!! Calf Panel Covers (1/2" Soft Visco Foam)</b>		
<input type="checkbox"/> Air Mesh	\$120	
<input type="checkbox"/> Reversed Dartex	\$120	
<b>Optional Modular Footplate Extension Components/Kits:</b> Assembly required. May require seating configuration adjustments to avoid front caster interference.		
<input type="checkbox"/>  4" Width Extension Kit	\$50	
<input type="checkbox"/>  2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/>  4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<input type="checkbox"/>  6-Piece Extension Stock Kit (Includes all above pieces)	\$95	
<b>Additional Optional Accessories:</b>		
<input type="checkbox"/> Wheel Bumpers (pair)	\$25	
<input type="checkbox"/> Footplate Connector	\$49	
<input type="checkbox"/> <b>NEW!!! Footplate Cover, Reversed Dartex (1/2" Soft Visco Foam)</b> Not available with modular footplate extension components/kits.	\$120	

1. Comes standard with articulating foot platform.

Adjustable Center Mount Foot Platform Compatibility Matrix	TRU-Balance® 2 Recline			TRU-Balance® 2 Tilt			TRU-Balance® 2 Recline w/PES		TRU-Balance® 2 Tilt w/PES	
	Low	Medium	High	Low	Medium	High	Low (6")	High(7")	Low (6")	High(7")
10"-16" Extension	!	✓	✓	✓	✓	✓	!	✓	✓	✓
13"-19" Extension	X	X	X	X	!	✓	X	X	!	✓



This rigging Configuration causes no interference issues at this seat-to-floor height.



This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will



This rigging configuration does not work with this seat-to-floor height.



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Headrests

#### Headrest Pads

Description	Price	HCPCS Code
<input type="checkbox"/> 8"W TRU-Comfort™ Air Mesh <sup>1</sup>	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort™ <sup>1</sup>	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort™ Air Mesh <sup>1</sup>	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort™ <sup>1</sup>	\$205	E0955
<input type="checkbox"/> Stealth 10" Comfort Plus	\$142	E0955
<input type="checkbox"/> Stealth 14" Comfort Plus	\$142	E0955
<input type="checkbox"/> Stealth Small Contoured (9.75"W x 3"H)	\$205	E0955
<input type="checkbox"/> Stealth Large Contoured (11"W x 6"H)	\$205	E0955





#### Headrest Mounts

Description	Price	HCPCS Code
<input type="checkbox"/> TRU-Comfort™ Removable Headrest Bracket <sup>2</sup>	\$210	E1028
<input type="checkbox"/> Stealth Removable Headrest Bracket	\$210	E1028
<input type="checkbox"/> Stealth's Worlds Best Removable Headrest Bracket	\$210	E1028
<input type="checkbox"/> No Bracket		

1. Measurement is taken from the inside of the headrest.

2. Only available w/ TRU-Comfort™ headrest pads.

### Swing-Away Headrest Mounts

Description	Price	HCPCS Code
<input type="checkbox"/>  Facial Lateral, 4.5" x 2", Right Side, for Comfort Series <sup>1</sup> (Stealth: CPS9-R)	\$325.00	E2018
<input type="checkbox"/>  Facial Lateral, 4.5" x 2", Left Side, for Comfort Series <sup>1</sup> (Stealth: CPS9-L)	\$325.00	E2018
<input type="checkbox"/>  Facial Lateral, 4.5" x 2", Left and Right Side, for Comfort Series <sup>1</sup> (Stealth: CPS9-B)	\$425.00	E2018
<input type="checkbox"/>  Egg Switch Mount (each) (Stealth: ESM)	\$50.00	

1. Only available with the Stealth Products Comfort series headrests.

### Lap Belts

Description	Price	HCPCS Code
<input type="checkbox"/> 38" Lap Belt	Standard	
<input type="checkbox"/> 50" Lap Belt	Standard	
<input type="checkbox"/> 60" Lap Belt	Standard	
<input type="checkbox"/> 70" Lap Belt	Standard	
<input type="checkbox"/> 80" Lap Belt	Standard	
<input type="checkbox"/> 90" Lap Belt	Standard	
<input type="checkbox"/> 60" Padded Push Button, center Pull	\$115	
<input type="checkbox"/> Large 13-19", Hip Padded, Push Button	\$92	



## Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

### Synergy® Cushions







Description	Price	HCPCS Code
<input type="checkbox"/> Matching Seat Size		
<input type="checkbox"/> Non-Matching Seat Size		
<i>Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s</i>		
<b>Simplicity Cushion: General Use</b>		
<input type="checkbox"/> 16" x 16" - 20" x 24"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602
<b>Solution® Cushion: Skin Protection &amp; Positioning</b>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608
<b>Solution® 1 Cushion: Skin Protection</b>		
<input type="checkbox"/> 14" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604
<b>Spectrum® Gel Cushion: Positioning</b>		
<input type="checkbox"/> 12" x 10" - 20" x 24"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2606
<b>Spectrum® Foam Cushion: Skin Protection &amp; Positioning</b>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2608
<b>Spectrum® Air Contour Cushion: Adjustable Skin Protection</b>		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$500	E2622
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$550	E2623
<b>Air Tube:</b>		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
<b>Spectrum® Air Flat Cushion: Adjustable Skin Protection &amp; Positioning</b>		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$549	E2625
<b>Air Tube:</b>		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
<b>TRU-Comfort™ 2 4-Way Stretch: Skin Protection</b>		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$480	E2604
<b>TRU-Comfort™ 2 Air Mesh: Skin Protection</b>		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$480	E2604
<b>TRU-Comfort™ Accessory:</b>		
<input type="checkbox"/> Sport Seat Cushion Cover	\$100	
<i>Available as an additional cover option for both TRU-Comfort seat cushion styles. Only available in sizes listed.</i>		
<input type="checkbox"/> 16"W x 16"D		
<input type="checkbox"/> 18"W x 18"D		
<input type="checkbox"/> 20"W x 20"D		

### Additional Accessories




Description	Price	HCPCS Code
<input type="checkbox"/> Quantum Backpack <i>May interfere with other rear accessories</i>	\$60	
<input type="checkbox"/> Oxygen Holder <i>Not available with a Power Elevating Seat. Oxygen holder on Tilt with TRU-Comfort back requires SOF. N/A on Power Recline.</i>	\$150	E2208
<input type="checkbox"/> Recline Accessory Mounting Bracket <i>1" adapter tube allows cane mounted accessories to mount to power recline back or side rail.</i>	\$50	
<input type="checkbox"/> Cup Holder	No Charge	
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> Tablet Holder	\$315	
<input type="checkbox"/> Deluxe Cellphone/Media Holder	\$140	
<input type="checkbox"/> Transfer Bar (pair) <i>Quick Release</i>	\$175	

# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Additional Foot Rigging Accessories

Description	Price	HCPCS Code
<input type="checkbox"/>  <b>Therafin Flip-down Knee Abductor bracket (each)</b> <b>Pads:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Small (3.5" x 3.25" x 4")</li> <li><input type="checkbox"/> Small, Gel (3.5" x 3.25" x 4")</li> <li><input type="checkbox"/> Medium (4" x 3.25" x 6")</li> <li><input type="checkbox"/> Large (4.5" x 2.75" x 6")</li> </ul>	\$172  \$135 \$145 \$135 \$135	E1028  E0957 E0957 E0957 E0957
<input type="checkbox"/>  <b>POSAlinc Knee Adductor for Swing-Away Legrests (pair)</b> <input type="checkbox"/> <b>POSAlinc Rail Mounted Knee Adductors (pair)</b>	\$512 \$574	K0108 K0108
<input type="checkbox"/>  <b>Stealth Foot Boxes (pair)</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Small</li> <li><input type="checkbox"/> Medium</li> <li><input type="checkbox"/> Large</li> </ul>	\$475	K0108
<input type="checkbox"/>  <b>Stealth Padded Shoe Holders (pair)</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Small</li> <li><input type="checkbox"/> Medium</li> <li><input type="checkbox"/> Large</li> </ul>	\$166	K0108
<input type="checkbox"/>  <b>Ankle Huggers (pair)</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Medium</li> <li><input type="checkbox"/> Large</li> </ul>	\$156	K0108
<input type="checkbox"/>  <b>Bodypoint High Mount, Clamp-On Footplates (pair)</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Angle Adjustable: <ul style="list-style-type: none"> <li><input type="checkbox"/> Small, 4"x 7.75"</li> <li><input type="checkbox"/> Large, 4.25" x 9.25"</li> </ul> </li> <li><input type="checkbox"/> Multi-Axis: <ul style="list-style-type: none"> <li><input type="checkbox"/> Small, 4"x 7.75"</li> <li><input type="checkbox"/> Large, 4.25" x 9.25"</li> </ul> </li> </ul>	 \$356 \$356  \$456 \$460	 K0040 K0040  K0040 K0040

## Additional Armrest Accessories

<input type="checkbox"/>  <b>Ottobock Arm Channels (each)</b> <i>A 2" narrower back width is recommended to avoid interference w/ the back cushion.</i> <ul style="list-style-type: none"> <li><input type="checkbox"/> Small (4.75" x 8")</li> <li><input type="checkbox"/> Medium (4.75" x 11.5")</li> <li><input type="checkbox"/> Large (4.75" x 14")</li> <li><input type="checkbox"/> One Piece (4.75" x 18")</li> </ul> <b>Side:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Left</li> <li><input type="checkbox"/> Right</li> </ul>	\$104	E2209
<input type="checkbox"/>  <b>Ottobock Hand Pads (each)</b> <i>Not available on same side as controller</i> <b>Side:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Left</li> <li><input type="checkbox"/> Right</li> </ul> <b>Size:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Medium</li> <li><input type="checkbox"/> Large</li> </ul>	\$52	
<input type="checkbox"/>  <b>Stealth Elbow Support (each)</b> <b>Side:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> Swing-Away Bracket-Left (each)</li> <li><input type="checkbox"/> Swing-Away Bracket-Right (each)</li> </ul> <b>Pad Size:</b> <ul style="list-style-type: none"> <li><input type="checkbox"/> 3"W x 3"H (each)</li> <li><input type="checkbox"/> 4"W x 4"H (each)</li> </ul>	 \$175 \$175  \$50 \$50	 E1028 E1028  K0108 K0108



# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Powered Legrests

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

Description	Price	HCPCS Code
<input type="checkbox"/> <b>70° Power Articulating Swing-Away Legrests</b> 4" articulation range <i>Use code K0108 for an individual legrest.</i>	\$2,395	E1010

### Rigging Selection:

- ☐ Standard Matching Legrests (pair)
- ☐ Non-Matching Legrests (pair)
- ☐ Individual Legrests (Left or Right Only)

### Operation Method:

- ☐ Independent
- ☐ Combined

### Footplate Selection:

- ☐ Angle Adjustable Footplates (pair) *Includes heel loops.*

\$150

K0040

### Left Legrest

#### Lower Extensions:

- ☐ 11"-13"L
- ☐ 14.5"-18"L
- ☐ 17.5"-21"L

### Right Legrest

#### Lower Extensions:

- ☐ 11"-13"L
- ☐ 14.5"-18"L
- ☐ 17.5"-21"L

### ☐ Angle Adjustable Footplate Size:

- ☐ 4"W x 6"D
- ☐ 5"W x 6"D
- ☐ 5"W x 8"D
- ☐ 6"W x 8"D<sup>1</sup>
- ☐ 6"W x 11"D<sup>1</sup>
- ☐ 8"W x 11"D<sup>2</sup>

- ☐ 4"W x 6"D
- ☐ 5"W x 6"D
- ☐ 5"W x 8"D
- ☐ 6"W x 8"D<sup>1</sup>
- ☐ 6"W x 11"D<sup>1</sup>
- ☐ 8"W x 11"D<sup>2</sup>

### Accessories:

Description	Price	HCPCS Code
<input type="checkbox"/> Small Calf Pad (pair) 6"W - 5.35"H	No Charge	
<input type="checkbox"/> Large Calf Pad (pair) 8"W x 5.86"H	No Charge	
<input type="checkbox"/> Gel Padded Knee Adductor (pair)	\$100	K0108

1. May cause interference on seat widths less than 18"

2. May cause interference on seat widths less than 20"

Power Elevating Legrest Compatibility Matrix	TRU-Balance® 2 Recline			TRU-Balance® 2 Tilt			TRU-Balance® 2 Recline w/PES		TRU-Balance® 2 Tilt w/PES	
	Low	Medium	High	Low	Medium	High	Low (6")	High(7")	Low (6")	High(7")
Swing-Away, 11" - 13"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Swing-Away, 14.5" - 18"	!	!	!	!	!	!	X	!	!	!



This rigging Configuration causes no interference issues at this seat-to-floor height.



This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will



This rigging configuration does not work with this seat-to-floor height.



# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Non-Powered Legrests

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

Description	Price	HCPCS Code
<input type="checkbox"/> <b>60-90° Swing-Away Legrests</b> adjustable in 5° Increments	Standard	
<b>Rigging Selection:</b>		
<input type="checkbox"/> Standard Matching Legrests (pair)		
<input type="checkbox"/> Non-Matching Legrests (pair)		
<input type="checkbox"/> Individual Legrests Left or Right Only		
<b>Footplate Selection:</b>		
<input type="checkbox"/> Angle-Adjustable Footplates (pair) Includes Heel Loops.	\$150	K0040
<input type="checkbox"/> Multi-Axis Footplates (pair) Includes Heel Loops.	\$495	K0040
<input type="checkbox"/> Standard, Extruded Footplates (pair)	No Charge	
<b>Left Legrest:</b>		
<b>Lower Extensions:</b>		
<input type="checkbox"/> 11"-13"L <sup>1</sup>		
<input type="checkbox"/> 11"-16.5"L <sup>2</sup>		
<input type="checkbox"/> 13"-16.5"L <sup>1</sup>		
<input type="checkbox"/> 15.5"-19.5"L <sup>2</sup>		
<input type="checkbox"/> <b>Angle Adjustable/Multi-Axis Footplate Size:</b>		
<input type="checkbox"/> 4"W x 6"D		
<input type="checkbox"/> 5"W x 6"D		
<input type="checkbox"/> 5"W x 8"D		
<input type="checkbox"/> 6"W x 8"D <sup>3</sup>		
<input type="checkbox"/> 6"W x 11"D <sup>3</sup>		
<input type="checkbox"/> 8"W x 11"D <sup>4</sup>		
<input type="checkbox"/> <b>Standard, Extruded Footplate Size:</b>		
<input type="checkbox"/> 7.25"W x 6"D <sup>3</sup>		
<input type="checkbox"/> 8.5"W x 6"D <sup>4</sup>		
<b>Right Legrest:</b>		
<b>Lower Extensions:</b>		
<input type="checkbox"/> 11"-13"L <sup>1</sup>		
<input type="checkbox"/> 11"-16.5"L <sup>2</sup>		
<input type="checkbox"/> 13"-16.5"L <sup>1</sup>		
<input type="checkbox"/> 15.5"-19.5"L <sup>2</sup>		
<input type="checkbox"/> 4"W x 6"D		
<input type="checkbox"/> 5"W x 6"D		
<input type="checkbox"/> 5"W x 8"D		
<input type="checkbox"/> 6"W x 8"D <sup>3</sup>		
<input type="checkbox"/> 6"W x 11"D <sup>3</sup>		
<input type="checkbox"/> 8"W x 11"D <sup>4</sup>		
<input type="checkbox"/> 7.25"W x 6"D <sup>3</sup>		
<input type="checkbox"/> 8.5"W x 6"D <sup>4</sup>		

### Accessories:

Description	Price	HCPCS Code
<input type="checkbox"/> Heel Loops (pair) Standard on Angle Adjustable and Multi-Axis footplates	\$38	E0951
<input type="checkbox"/> Wheel Bumpers (pair)	\$25	
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair)	\$100	K0108
<input type="checkbox"/> Therafin Multi-Axis Amputee Support (each) <i>Swing-Away legrest upper assembly and receiver is required</i>	\$230	E1020
<input type="checkbox"/> Amputee Support Gel Cover (each)	\$150	K0108

1. Not available with multi-axis footplates.
2. Not available with angle adjustable footplates.
3. May cause interference on seat widths less than 18"
4. May cause interference on seat widths less than 20"

60-90° Swing-Away Legrest Compatibility Matrix	TRU-Balance® 2 Recline			TRU-Balance® 2 Tilt			TRU-Balance® 2 Recline w/PES		TRU-Balance® 2 Tilt w/PES	
	Low	Medium	High	Low	Medium	High	Low (6")	High(7")	Low (6")	High(7")
Standard, 13" - 16.5"	!	!	✓	!	✓	✓	!	!	✓	✓
Standard, 15.5" - 19.5"	X	X	!	!	!	!	X	X	!	!
Angle-Adjustable, 13" - 16.5"	!	!	✓	!	✓	✓	!	!	✓	✓
Angle-Adjustable, 15.5" - 19.5"	X	X	!	!	!	!	X	X	!	!
Multi-Axis, 11" - 16.5"	!	!	✓	!	✓	✓	!	!	✓	✓
Multi-Axis, 15.5" - 19.5"	X	X	!	!	!	!	X	X	!	!



This rigging Configuration causes no interference issues at this seat-to-floor height.



This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.



This rigging configuration does not work with this seat-to-floor height.



# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Non-Powered Legrests (cont.)

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

Description	Price	HCPCS Code
<input type="checkbox"/> <b>90° Pediatric Swing-Away Legrests</b>	Standard	
<b>Rigging Selection:</b>		
<input type="checkbox"/> Standard Matching Legrests (pair)		
<input type="checkbox"/> Non-Matching Legrests (pair)		
<input type="checkbox"/> Individual Legrests, left or right only		
<b>Footplate Selection:</b>		
<input type="checkbox"/> Angle-Adjustable Footplates (pair) <i>Includes heel loops.</i>	\$150	K0040
<b>Left Legrest:</b>		<b>Right Legrest:</b>
<b>Lower Extensions:</b>		<b>Lower Extensions:</b>
<input type="checkbox"/> 7"-10"L		<input type="checkbox"/> 7"-10"L
<input type="checkbox"/> <b>Angle Adjustable Footplate Size:</b>		
<input type="checkbox"/> 4"W x 6"D		<input type="checkbox"/> 4"W x 6"D
<input type="checkbox"/> 5"W x 6"D		<input type="checkbox"/> 5"W x 6"D
<b>Accessories:</b>		
<input type="checkbox"/> Wheel Bumpers (pair)	\$25	

Pediatric Swing-Away Compatibility Matrix	TRU-Balance® 2 Recline			TRU-Balance® 2 Tilt		
	Low	Medium	High	Low	Medium	High
10"-16" Extension	✓	✓	✓	✓	✓	✓



This rigging Configuration causes no interference issues at this seat-to-floor height.

<input type="checkbox"/> <b>Elevating Legrests</b>	\$325	E0990
<b>Rigging Selection:</b>		
<input type="checkbox"/> Standard Matching Legrests (pair)		
<input type="checkbox"/> Non-Matching Legrests (pair)		
<input type="checkbox"/> Individual Legrests <i>Left or Right Only</i>		
<b>Footplate Selection:</b>		
<input type="checkbox"/> Angle-Adjustable Footplates (pair) <i>Includes Heel Loops.</i>	\$150	K0040
<input type="checkbox"/> Standard, Extruded Footplates (pair)	No Charge	
<b>Left Legrest Lower Extensions:</b>		<b>Right Legrest Lower Extensions:</b>
<input type="checkbox"/> 11"-14" <sup>1</sup>		<input type="checkbox"/> 11"-14" <sup>1</sup>
<input type="checkbox"/> 12"-14.5" <sup>1</sup>		<input type="checkbox"/> 12"-14.5" <sup>1</sup>
<input type="checkbox"/> 14"-17" <sup>1</sup>		<input type="checkbox"/> 14"-17" <sup>1</sup>
<input type="checkbox"/> 14.5"-18" <sup>1</sup>		<input type="checkbox"/> 14.5"-18" <sup>1</sup>
<input type="checkbox"/> <b>Angle Adjustable Footplate Size:</b>		
<b>Left Legrest:</b>		<b>Right Legrest:</b>
<input type="checkbox"/> 4"W x 6"D		<input type="checkbox"/> 4"W x 6"D
<input type="checkbox"/> 5"W x 6"D		<input type="checkbox"/> 5"W x 6"D
<input type="checkbox"/> 5"W x 8"D		<input type="checkbox"/> 5"W x 8"D
<input type="checkbox"/> 6"W x 8"D <sup>2</sup>		<input type="checkbox"/> 6"W x 8"D <sup>2</sup>
<input type="checkbox"/> 6"W x 11"D <sup>2</sup>		<input type="checkbox"/> 6"W x 11"D <sup>2</sup>
<input type="checkbox"/> 8"W x 11"D <sup>3</sup>		<input type="checkbox"/> 8"W x 11"D <sup>3</sup>
<input type="checkbox"/> <b>Standard, Extruded Footplate Size:</b>		
<b>Left Legrest:</b>		<b>Left Legrest:</b>
<input type="checkbox"/> 7.5"W x 6"D <sup>2</sup>		<input type="checkbox"/> 7.5"W x 6"D <sup>2</sup>
<input type="checkbox"/> 8.5"W x 6"D <sup>3</sup>		<input type="checkbox"/> 8.5"W x 6"D <sup>3</sup>
<b>Accessories:</b>		
<b>Description</b>	<b>Price</b>	<b>HCPCS Code</b>
<input type="checkbox"/> Heel Loops (pair) <i>Standard on Angle Adjustable and Multi-Axis footplates</i>	\$38	E0951
<input type="checkbox"/> Wheel Bumpers (pair)	\$25	
<input type="checkbox"/> Gel Padded Adductor (pair)	\$100	K0108

1. Only available with standard, extruded footplates.

2. May cause interference on seat widths less than 18"

3. May cause interference on seat widths less than 20"





# Quantum Q6 Edge with TRU-Balance® 2 Power Positioning

## Non-Powered Legrests (cont.)

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

Power Elevating Legrest Compatibility Matrix	TRU-Balance® 2 Recline			TRU-Balance® 2 Tilt			TRU-Balance® 2 Recline w/PES		TRU-Balance® 2 Tilt w/PES	
	Low	Medium	High	Low	Medium	High	Low (6")	High(7")	Low (6")	High(7")
Standard, 11" - 14"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Standard, 14" - 17"	!	!	!	!	!	✓	!	!	!	✓
Angle-Adjustable, 12" - 14.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle-Adjustable, 14.5" - 18"	!	!	!	!	!	!	X	!	!	!



This rigging Configuration causes no interference issues at this seat-to-floor height.



This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will



This rigging configuration does not work with this seat-to-floor height.



### Description

#### Heavy Duty, Drop-In Legrests

Price

\$325

HCPCS Code

E0990

#### Rigging Selection:

☐ Standard Matching Legrests (pair)☐ Non-Matching Legrests (pair)☐ Individual Legrests *Left or Right Only*

#### Footplate Selection:

☐ Angle-Adjustable Footplates (pair) *Includes Heel Loops.*

\$150

K0040

#### Left Legrest:

##### Lower Extensions:

☐ 12"-14.5"L☐ 14.5"-18"L

#### Right Legrest:

##### Lower Extensions:

☐ 12"-14.5"L☐ 14.5"-18"L

#### Angle Adjustable Footplate Size:

#### Left Legrest

☐ 4"W x 6"D☐ 5"W x 6"D☐ 5"W x 8"D☐ 6"W x 8"D<sup>1</sup>☐ 6"W x 11"D<sup>1</sup>☐ 8"W x 11"D<sup>2</sup>

#### Right Legrest

☐ 4"W x 6"D☐ 5"W x 6"D☐ 5"W x 8"D☐ 6"W x 8"D<sup>1</sup>☐ 6"W x 11"D<sup>1</sup>☐ 8"W x 11"D<sup>2</sup>

### Accessories:

#### Description

☐ Wheel Bumpers (pair)

Price

\$25

HCPCS Code

1. May cause interference on seat widths less than 18"

3. May cause interference on seat widths less than 20"

Heavy Duty Drop-In Compatibility Matrix	TRU-Balance® 2 Recline			TRU-Balance® 2 Tilt			TRU-Balance® 2 Recline w/PES		TRU-Balance® 2 Tilt w/PES	
	Low	Medium	High	Low	Medium	High	Low (6")	High(7")	Low (6")	High(7")
12" - 14.5" Extension	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
14.5" - 18" Extension	!	!	!	!	!	!	!	!	!	!



This rigging Configuration causes no interference issues at this seat-to-floor height.



This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.