

J6 Series

Group 2 Standard
Order Form
300lb Weight Capacity



Pride Mobility Products Corporation
182 Susquehanna Ave., Exeter, PA 18643
Phone: (866) 800-2002 • Fax: (866) 707-3422
quantumorders@pridemobility.com

Discount: _____
Terms: _____
Freight: _____
Taken By: _____



Account Number: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
E-Mail: _____
Date: _____
PO#: _____
Marked For: _____
Ship To Address: _____
City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP. Items Designated with a Circular Checkbox (○) are not available through Quantum West

BASE MODEL

Description	Price	HCPCS Code
<input type="checkbox"/> J6 2S-C	\$5,971	K0823
<input type="checkbox"/> J6 2S-SS	\$5,971	K0822

COLORS

- Candy Apple Red
- Viper Blue

BATTERY CHARGER

Description	Price
<input type="checkbox"/> 5 Amp Off-board Charger	Standard

DRIVE WHEELS

Description	Price
<input type="checkbox"/> 10" Pneumatic, pair	No Charge
<input type="checkbox"/> 10" Flat-Free, pair	Standard

POWER ELEVATING SEAT



Description	Price
<input type="checkbox"/> Power Elevating Seat Through VSI 6-Key ¹	\$975

1. Only available with VSI 6-Key. Not available with manual recliner, Synergy® Manual Recline or manual tilt. Max seat size is 20"W x 20"D. Increases STF ht. Includes 1/4 speed inhiit.

BASE MOUNTED FOOT PLATFORM

Description	Price
<input type="checkbox"/> Base Mounted Foot Platform	Standard
<input type="checkbox"/> No Foot Platform	No Charge





NON-EXPANDABLE CONTROLS

Description	Price
<input type="checkbox"/> 75 Amp NE 	Standard
<input type="checkbox"/> 75 Amp NE+ 6-Key (For up to 2 Actuators Thru Joystick)	\$200
<input type="checkbox"/> 75 Amp NE+ 8-Key (For Lights) 	\$400
<input type="radio"/> 70 Amp PG VSI 6-Key ¹	Standard

1. Only available with a power elevating seat.

EXPANDABLE CONTROLS

Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls

Description	Price
<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller 	\$695
<input type="checkbox"/> Expandable Controller Harness	\$350
<input type="checkbox"/> Q-Logic 2 EX Joystick 	\$700
<input type="radio"/> Q-Logic 2 EX Joystick (For Lights) 	\$700
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port 	\$700

Armor Color Options:

- Black Soft Touch (Standard color option)
- Matching
- Non-Matching

Colors: \$50

- Candy Apple Red
- Viper Blue

Items Designated with a Circular Checkbox (☐) are not available through Quantum West

JOYSTICK MOUNTING BRACKETS

Description Price HCPCS Code

Comes Standard with a Quick Release, Adjustable-Height Joystick Block

- Fixed, Inline¹ Standard
 - Left
 - Right
- Swing-Away \$210 E1028
 - Left
 - Right
- New!!!** Retract4 Swing-Away Joystick Arm Kit \$250 E1028
 - Left
 - Right
- Multi-Axis, Fixed \$100
 - Left
 - Right
- No Joystick Mount

1. Not available with Q-Logic 2 EX Stand Alone Joystick.

JOYSTICK HANDLES

Description Price HCPCS Code

- Standard Carrot Handle Standard
- Mushroom Handle \$71 E2323
- T-Bar \$107 E2323
- Chin Cup \$71 E2324¹
- Extended I-Handle \$71 E2323
- Soft Ball Handle \$50 E2323
- Bodypoint 3" Goal Post Handle \$69 E2323
- Bodypoint 4" Goal Post Handle \$69 E2323
- Medium Ball Handle \$50 E2323

1. Use HCPCS Code E2323 when not used in conjunction with a chin control.

PROGRAMMING

Description Price

- Q-Logic EX/NE Programmer \$545
- Q-Logic PC Programming Disk w/ Cable¹ \$750
- P&G Programmer \$423



1. Works with Q-Logic 2 EX only.

CONTOURED SEATING

Description Price

Depth-Adjustable High Back:¹

- Black Vinyl, Depth-Adjustable High Back, Limited Recline No Charge
 - 18"W x 18-20"D
 - 20"W x 18-20"D

Limited Recline High Back:

- Black Vinyl, High Back with Limited Recline No Charge
 - 16"W x 16"D 18"W x 16"D
 - 16"W x 18"D 20"W x 16"D
 - 16"W x 20"D

Limited Recline High Back w/ Solid Seat Pan:

- Black Vinyl, High Back with Solid Seat pan, Limited Recline No Charge
 - 16"W x 16"D 18"W x 20"D
 - 16"W x 18"D 20"W x 16"D
 - 16"W x 20"D 20"W x 18"D
 - 18"W x 16"D 20"W x 20"D
 - 18"W x 18"D

Medium Back:

- Black Vinyl, Medium Back No Charge
 - 16"W x 16"D 18"W x 16"D
 - 16"W x 18"D 18"W x 18"D

Medium Back w/ Solid Seat Pan:

- Black Vinyl, Medium Back With Solid Seat Pan No Charge
 - 16"W x 16"D 18"W x 16"D
 - 16"W x 18"D 18"W x 18"D

Armrest Type:

- Flip up, Removable No Charge
- Flip up, Height-Adjustable, Removable \$85 ea. E0973

Armpads:

- | | |
|--|---|
| <p><i>Left:</i></p> <ul style="list-style-type: none"> <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length² | <p><i>Right:</i></p> <ul style="list-style-type: none"> <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length² |
|--|---|

Special Order Seat

_____ "W x _____ "D

Seat size can not exceed 20"W x 20"D

Contact your inside sales specialist for pricing

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part.
 2. N/A with depth-adjustable high back seat.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

MANUAL RECLINER

Description	Price	HCPCS Code
<input type="checkbox"/> Manual Recliner High Back With Headrest ¹	\$600	E1226
<input type="checkbox"/> Manual Recliner High Back With Headrest and Solid Seat Pan ¹	\$600	E1226
Color: <input type="checkbox"/> Gray Vinyl <input type="checkbox"/> Beige Vinyl		
Sizes <input type="checkbox"/> 16"W x 16"D <input type="checkbox"/> 18"W x 20"D		
<input type="checkbox"/> 16"W x 18"D <input type="checkbox"/> 20"W x 16"D		
<input type="checkbox"/> 16"W x 20"D <input type="checkbox"/> 20"W x 18"D		
<input type="checkbox"/> 18"W x 16"D <input type="checkbox"/> 20"W x 20"D		
<input type="checkbox"/> 18"W x 18"D		
Armrest Type:		
<input type="checkbox"/> Manual Recliner Armrests	No Charge	
Armpads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length	Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length	

1. Not available with a power elevating seat.

BACK CANE OPTIONS

Description	Price	HCPCS Code
Back Cane Height¹:		
<input type="checkbox"/> 14"H ² <input type="checkbox"/> 15"H ² <input type="checkbox"/> 16"H ³ <input type="checkbox"/> 17"H ³ <input type="checkbox"/> 18"H		
<input type="checkbox"/> 19"H <input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H		
<input type="checkbox"/> 24"H <input type="checkbox"/> 25"H		
Cane Type¹:		
<input type="checkbox"/> Standard Back Canes	Standard	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Attendant Manual Recline ⁴	\$750	E1226
Cane Brackets:		
<input type="checkbox"/> Flip Forward Back Canes ⁵	\$150	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Offset (Meso) ⁶	\$195	
<input type="checkbox"/> 2" (1" per side) <input type="checkbox"/> 4" (2" per side)		
<input type="checkbox"/> Inset (Endo) ⁶	\$195	
<input type="checkbox"/> 2" (1" per side) <input type="checkbox"/> 4" (2" per side)		

1. No selection required with TRU-Comfort™ Back.
 2. Selection of a Synergy® Back and cantilever armrests disqualifies 14" and 15"H back canes.
 3. Selection of a Synergy® Back and cane mounted, flip up armrests disqualifies 16" and 17"H back canes.
 4. Not available with Flip Back Arms, Cantilever Arms, Shape Back, TRU-Comfort back, rear accessories or power elevating seat. Requires 22"H canes. Pivot point 2.5" above seat pan.
 5. Not available with vent tray, oxygen holder, cane mounted or cantilever armrests. May interfere with lateral supports and waterfall armpads.
 6. Dual-Post, flip back armrests interfere with offset canes. Not available with attendant manual recline.

SYNERGY® SEATING SYSTEM

Description	Price
<input type="checkbox"/> Synergy® Seating Frame	\$434
<input type="checkbox"/> Synergy® Manual Tilt Package ¹	\$434
<i>Max seat size 14"W x 14"D, 75 lbs weight capacity</i>	
Towers:	
<input type="checkbox"/> Standard (15.75"-16.75")	
<input type="checkbox"/> Extended (16.5"-17.5") ²	

1. In iQuantum this option is chosen under "width package type".
 2. Not available with Synergy® manual tilt or recline.

BACK OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Omit Back	No Charge	
Sling Backs:		
<input type="checkbox"/> Sling Back ¹	\$98	
Overall Height ² : <input type="checkbox"/> 18"H <input type="checkbox"/> 19"H <input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H		
<input type="checkbox"/> 23"H <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H		
<input type="checkbox"/> Synergy Shape Back	Various	E2613/E2614
<i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>		
<input type="checkbox"/> Comfort Company Back	Various	
See Comfort Company Order Form		
Synergy® Back:		
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	\$275	
Overall Height : <input type="checkbox"/> 13"H <input type="checkbox"/> 14"H <input type="checkbox"/> 15"H <input type="checkbox"/> 16"H <input type="checkbox"/> 17"H		
<input type="checkbox"/> 18"H <input type="checkbox"/> 19"H <input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H		
<input type="checkbox"/> 23"H <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H		
Choose Lateral Extension ³ :		
<input type="checkbox"/> 6" for 13-16"H Back	\$25	
<input type="checkbox"/> 10" for 17-18"H Back	\$25	
<input type="checkbox"/> 12" for 19-22"H Back	\$25	
<input type="checkbox"/> 16" for 23-25"H Back	\$25	

SEATING SPECIFICATIONS

Seat Width:	Seat Depth:
<input type="checkbox"/> 10"W	<input type="checkbox"/> 10"D
<input type="checkbox"/> 11"W	<input type="checkbox"/> 11"D
<input type="checkbox"/> 12"W	<input type="checkbox"/> 12"D
<input type="checkbox"/> 13"W	<input type="checkbox"/> 13"D
<input type="checkbox"/> 14"W	<input type="checkbox"/> 14"D
<input type="checkbox"/> 15"W	<input type="checkbox"/> 15"D
<input type="checkbox"/> 16"W	<input type="checkbox"/> 16"D
<input type="checkbox"/> 17"W	<input type="checkbox"/> 17"D
<input type="checkbox"/> 18"W	<input type="checkbox"/> 18"D
<input type="checkbox"/> 19"W	<input type="checkbox"/> 19"D
<input type="checkbox"/> 20"W	<input type="checkbox"/> 20"D
Seat Width Package:	Seat Depth Package:
<input type="checkbox"/> 10"W-12"W	<input type="checkbox"/> 10"D-14"D
<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 12"W-14"W	<input type="checkbox"/> 15"D-20"D
<input type="checkbox"/> 14"W-18"W	
<input type="checkbox"/> 16"W-20"W	

1. Not available with cantilever or flip, back, cane mounted armrests.
 2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
 3. Lateral extensions are required to mount lateral thoracic supports.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

BACK OPTIONS (CONT'D)

ARMPADS

Description	Price	HCPSC Code
TRU-Comfort™ Back:		
<i>Please refer to TRU-Comfort Back order form for detailed pricing, sizes, codes and part #'s</i>		
<input type="checkbox"/> TRU-Comfort™ 4-Way Stretch Back (20"W or less)	\$595	E2620
<input type="checkbox"/> TRU-Comfort™ 4-Way Stretch Back (22"W or greater)	\$695	E2621
Width: <input type="checkbox"/> 16"W <input type="checkbox"/> 18"W <input type="checkbox"/> 20"W		
Height ¹ : <input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H <input type="checkbox"/> 24"H		
<input type="checkbox"/> 25"H <input type="checkbox"/> 27"H <input type="checkbox"/> 28"H		
Width: <input type="checkbox"/> 17"W <input type="checkbox"/> 19"W		
Height <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H <input type="checkbox"/> 27"H		
<input type="checkbox"/> 28"H		
Width: <input type="checkbox"/> 22"W <input type="checkbox"/> 24"W		
Height : <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H <input type="checkbox"/> 27"H <input type="checkbox"/> 28"H		
Color Shroud Option²:		
<input type="checkbox"/> Standard Black	No Charge	
Color Options: \$100		
<input type="radio"/> Candy Apple Red		
<input type="radio"/> Viper Blue		

1. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
2. Only available with TRU-Comfort back.

Option #	Description	Price
<input type="checkbox"/>	Straight, Full Length (2"x14")	Standard
<input type="checkbox"/>	Left <input type="checkbox"/> Right	
<input type="checkbox"/>	Straight, Desk Length (2"x10")	No Charge
<input type="checkbox"/>	Left <input type="checkbox"/> Right	
<input type="checkbox"/>	Waterfall, Full Length (2"x14")	No Charge
<input type="checkbox"/>	Left <input type="checkbox"/> Right	
<input type="checkbox"/>	Waterfall, Desk Length (2"x10")	No Charge
<input type="checkbox"/>	Left <input type="checkbox"/> Right	
<input type="checkbox"/>	Flat, Bariatric ¹	No Charge
<input type="checkbox"/>	Falcon, Gel, Full Length, Pair (3.5"x14")	\$92.50 ea.
<input type="checkbox"/>	Gel Armpad Mounting Plate (Required)	\$56.69 ea.
<input type="checkbox"/>	Falcon, Gel, Desk Length, Pair (3.5"x10")	\$70.00 ea.
<input type="checkbox"/>	Gel Armpad Mounting Plate (Required)	\$53.52 ea.

1. Bariatric armrests are full length and require bariatric arm pads.

LED LIGHT PACKAGE

Description	Price
<input type="radio"/> Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995

ARMRESTS

LAPBELTS

Option #	Description	Price	HCPSC Code
<input type="checkbox"/>	Flip Back, Cane Mounted ¹	No Charge	
<input type="checkbox"/>	2-Post, Flip-Back, Height Adjustable, Removable ²	\$150 ea.	E0973
<input type="checkbox"/>	Single Post, Height Adjustable, Removable ³	\$75 ea.	E0973
<input type="checkbox"/>	Quad Release Lever		
<input type="checkbox"/>	Bariatric, Removable, Height Adjustable ⁴	\$75 ea.	E0973
Height: <input type="checkbox"/>	8"H-12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/>	Removable, Single Post, Height Adjustable, Heavy Duty	\$75 ea.	E0973
Height: <input type="checkbox"/>	8"H-12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/>	Pediatric, Removable Height Adjustable ⁵	\$75 ea.	E0973
<input type="checkbox"/>	Quad Release Lever		
<input type="checkbox"/>	Height Adjustable, Cantilever ⁶	\$195 pair	
	<i>Angle Adjustable: -30° through +30°</i>		

Option #	Description	Price
<input type="checkbox"/>	50" Lap Belt	Standard
<input type="checkbox"/>	60" Lap Belt	Standard
<input type="checkbox"/>	70" Lap Belt	Standard
<input type="checkbox"/>	80" Lap Belt	Standard
<input type="checkbox"/>	90" Lap Belt	Standard
<input type="checkbox"/>	60" Padded Push Button, Center Pull	\$115
<input type="checkbox"/>	Large 13-19", Hip Padded, Push Button	\$92

HEADRESTS

Option #	Description	Price	HCPSC Code
Headrest Pads:			
<input type="checkbox"/>	8"W TRU-Comfort™ Air Mesh ¹	\$205	E0955
<input type="checkbox"/>	8"W TRU-Comfort™ ¹	\$205	E0955
<input type="checkbox"/>	10" TRU-Comfort™ Air Mesh ¹	\$205	E0955
<input type="checkbox"/>	10" TRU-Comfort™ ¹	\$205	E0955
<input type="checkbox"/>	Stealth 10" Comfort Plus	\$142	E0955
<input type="checkbox"/>	Stealth 14" Comfort Plus	\$142	E0955
<input type="checkbox"/>	Stealth Small Contoured (9.75"W x 3"H)	\$205	E0955
<input type="checkbox"/>	Stealth Large Contoured (11"W x 6"H)	\$205	E0955
<input type="checkbox"/>	Omit Headrest		

1. Measurement is taken from the inside of the headrest.

1. Interferes with Swing-Away laterals. Not available w/ manual recline, sling back or flip forward back canes.
2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16".
3. Height adjustable from 11" - 17" in 1/2" increments.
4. Bariatric armrests are full length and require bariatric arm pads.
5. Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length arm pads.
6. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, Shape Back or sling back.

Items Designated with a Circular Checkbox (☐) are not available through Quantum West

HEADRESTS (CONT'D)

Option #	Description	Price	HCPCS Code
Headrest Mounts:			
<i>TCBRKT</i>	<input type="checkbox"/> TRU-Comfort™ Removable Headrest Bracket ¹	\$210	E1028
<i>HMO475P</i>	<input type="checkbox"/> Stealth Removable Headrest Bracket	\$210	E1028
<i>TWB485-1</i>	<input type="checkbox"/> Stealth's Worlds Best Removable Headrest Bracket	\$210	E1028
<i>NOHR</i>	<input type="checkbox"/> No Bracket		

1. Only available w/ TRU-Comfort™ headrest pads.

ELECTRONICS FOR POWER FOOT RIGGINGS

Option #	Description	Price	HCPCS Code
<i>THRUJOY</i>	<input type="checkbox"/> Single Actuator Through Joystick (AFP Only)	\$1,175	E2310 ¹
<i>THRUJOY</i>	<input type="checkbox"/> Dual Actuator Through NE+ 6 Key	\$2,375	E2311 ¹
<i>MTHRUJOY</i>	<input type="checkbox"/> Multiple Actuator Through Q-Logic 2 EX	\$2,375	E2311 ¹

1. Applies to covered actuators only.

NON-POWERED LEGRESTS

Option #	Description	Price
----------	-------------	-------

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

60-90° Swing-Away Legrests **\$162**
(Adjustable in 5° Increments)

Rigging Selection:

- Standard Matching Legrests, pair
- Non-Matching Legrests, pair
- Individual Legrests (Left or Right Only)

Footplate Selection:

- Angle-Adjustable Footplates¹ \$150
- Multi-Axis Footplates¹ \$495
- Standard, Extruded Footplates No Charge

Left Legrest:

Lower Extensions:

- 11"H-13"L²
- 11"-16.5"L
- 13"H-16.5"L²
- 15.5"-19.5"L

Right Legrest:

Lower Extensions:

- 11"L-13"L²
- 11"-16.5"L
- 13"L-16.5"L²
- 15.5"-19.5"L

Angle Adjustable/Multi-Axis Footplate Size

- | | |
|--|--|
| <input type="checkbox"/> 4"W x 6"D | <input type="checkbox"/> 4"W x 6"D |
| <input type="checkbox"/> 5"W x 6"D | <input type="checkbox"/> 5"W x 6"D |
| <input type="checkbox"/> 5"W x 8"D | <input type="checkbox"/> 5"W x 8"D |
| <input type="checkbox"/> 6"W x 8"D ³ | <input type="checkbox"/> 6"W x 8"D ³ |
| <input type="checkbox"/> 6"W x 11"D ³ | <input type="checkbox"/> 6"W x 11"D ³ |
| <input type="checkbox"/> 8"W x 11"D ⁴ | <input type="checkbox"/> 8"W x 11"D ⁴ |

1. Includes Heel Loops.
2. Not available with multi-axis footplates.
3. May cause interference on seat widths less than 18".
4. May cause interference on seat widths less than 20".

NON-POWERED LEGRESTS (CONT'D)

Option #	Description	Price	HCPCS Code
<input type="checkbox"/> Standard, Extruded Footplate Size			
	<input type="checkbox"/> 7.25"W x 6"D ⁵	<input type="checkbox"/> 7.25"W x 6"D ⁵	
	<input type="checkbox"/> 8.5"W x 6"D ⁶	<input type="checkbox"/> 8.5"W x 6"D ⁶	
<i>Accessories:</i>			
	<input type="checkbox"/> Heel Loops, pair ⁷	\$38	E0951
<i>BUMP</i>	<input type="checkbox"/> Wheel Bumpers, pair	\$25	
<i>LPW</i>	<input type="checkbox"/> Gel Neoprene Footrest Sleeve	\$100 pr.	K0108
	<input type="checkbox"/> Therafin Multi-Axis Amputee Support ⁸	\$230	E1020
	<input type="checkbox"/> Amputee Sppt. Gel Cover	\$150	K0108

Elevating Legrests **\$325** **K0195**

Rigging Selection:

- Standard Matching Legrests, pair
- Non-Matching Legrests, pair
- Individual Legrests (Left or Right Only)

Footplate Selection:

- Angle-Adjustable Footplates⁹ \$150
- Standard, Extruded Footplates No Charge

Left Legrest:

Lower Extensions:

- 11"-14"L¹⁰
- 12"H-14.5"L
- 14"-17"L¹⁰
- 14.5"-18"L

Right Legrest:

Lower Extensions:

- 11"-14"L¹⁰
- 12"L-14.5"L
- 14"-17"L¹⁰
- 14.5"-18"L

Angle Adjustable Footplate Size

- | | |
|--|--|
| <input type="checkbox"/> 4"W x 6"D | <input type="checkbox"/> 4"W x 6"D |
| <input type="checkbox"/> 5"W x 6"D | <input type="checkbox"/> 5"W x 6"D |
| <input type="checkbox"/> 5"W x 8"D | <input type="checkbox"/> 5"W x 8"D |
| <input type="checkbox"/> 6"W x 8"D ⁵ | <input type="checkbox"/> 6"W x 8"D ⁵ |
| <input type="checkbox"/> 6"W x 11"D ⁵ | <input type="checkbox"/> 6"W x 11"D ⁵ |
| <input type="checkbox"/> 8"W x 11"D ⁶ | <input type="checkbox"/> 8"W x 11"D ⁶ |

Standard, Extruded Footplate Size

- | | |
|--|--|
| <input type="checkbox"/> 7.25"W x 6"D ⁵ | <input type="checkbox"/> 7.25"W x 6"D ⁵ |
| <input type="checkbox"/> 8.5"W x 6"D ⁶ | <input type="checkbox"/> 8.5"W x 6"D ⁶ |

Accessories:

- Heel Loops, pair⁷ \$38 E0951
- BUMP* Wheel Bumpers, pair \$25
- LPW* Gel Padded Adductor \$100 pr. K0108

5. May cause interference on seat widths less than 18".
6. May cause interference on seat widths less than 20".
7. Standard on Angle-Adjustable and Multi-Axis footplates.
8. Swing-Away legrest upper assembly and receiver is required.
9. Includes heel loops.
10. Only available with standard, extruded footplates.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

NON-POWERED LEGRESTS (CONT'D)

ADJUSTABLE CENTER MOUNT FOOT PLATFORM




Option #	Description	Price
<input type="checkbox"/>	90° Pediatric Swing-Away Legrests	\$162
Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair	
<input type="checkbox"/>	Non-Matching Legrests, pair	
<input type="checkbox"/>	Individual Legrests (Left or Right Only)	
Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ¹	\$150
Left Legrest:		
Lower Extensions:		
<input type="checkbox"/>	7"H-10"L	
Right Legrest:		
Lower Extensions:		
<input type="checkbox"/>	7"L-10"L	
<input type="checkbox"/>	Angle Adjustable Footplate Size	
<input type="checkbox"/>	4"W x 6"D	<input type="checkbox"/> 4"W x 6"D
<input type="checkbox"/>	5"W x 6"D	<input type="checkbox"/> 5"W x 6"D
Accessories:		
BUMP	<input type="checkbox"/> Wheel Bumpers, pair	\$25

<input type="checkbox"/>	Heavy Duty Drop-In Legrests¹	\$420
<input type="checkbox"/>	70 Degree	
<input type="checkbox"/>	80 degree	
Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair	
<input type="checkbox"/>	Non-Matching Legrests, pair	
<input type="checkbox"/>	Individual Legrests (Left or Right Only)	
Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ²	\$150
Left Legrest:		
Lower Extensions:		
<input type="checkbox"/>	12"H-14.5"L	
<input type="checkbox"/>	14.5"-18"L	
Right Legrest:		
Lower Extensions:		
<input type="checkbox"/>	12"L-14.5"L	
<input type="checkbox"/>	14.5"-18"L	
<input type="checkbox"/>	Angle Adjustable Footplate Size	
<input type="checkbox"/>	4"W x 6"D	<input type="checkbox"/> 4"W x 6"D
<input type="checkbox"/>	5"W x 6"D	<input type="checkbox"/> 5"W x 6"D
<input type="checkbox"/>	5"W x 8"D	<input type="checkbox"/> 5"W x 8"D
<input type="checkbox"/>	6"W x 8"D ³	<input type="checkbox"/> 6"W x 8"D ³
<input type="checkbox"/>	6"W x 11"D ³	<input type="checkbox"/> 6"W x 11"D ³
<input type="checkbox"/>	8"W x 11"D ⁴	<input type="checkbox"/> 8"W x 11"D ⁴
Accessories:		
BUMP	<input type="checkbox"/> Wheel Bumpers, pair	\$25

1. N/A with contoured seating.
 2. Includes Heel Loops.
 3. May cause interference on seat widths less than 18".
 4. May cause interference on seat widths less than 20".

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

<input type="checkbox"/>	Adjustable Center Mount Foot Platform¹ (Includes 10"W x 11"D Split Footplates)	\$250
Extensions:		
CM1016	<input type="checkbox"/> 10-16"L	
Set Legrest Length:		
LL10	<input type="checkbox"/> 10"H	LL13 <input type="checkbox"/> 13"H
LL105	<input type="checkbox"/> 10.5"H	LL135 <input type="checkbox"/> 13.5"H
LL11	<input type="checkbox"/> 11"H	LL14 <input type="checkbox"/> 14"H
LL115	<input type="checkbox"/> 11.5"H	LL145 <input type="checkbox"/> 14.5"H
LL12	<input type="checkbox"/> 12"H	LL15 <input type="checkbox"/> 15"H
LL125	<input type="checkbox"/> 12.5"H	LL15.5 <input type="checkbox"/> 15.5"H
		LL16 <input type="checkbox"/> 16"H
		LL165 <input type="checkbox"/> 16.5"H
		LL17 <input type="checkbox"/> 17"H
		LL175 <input type="checkbox"/> 17.5"H
		LL18 <input type="checkbox"/> 18.5"H
		LL1 <input type="checkbox"/> 19"H
Calf Supports:		
<input type="checkbox"/>	Calf Pads (pair)	
<input type="checkbox"/>	Small (6"H x 6"W)	\$150
<input type="checkbox"/>	Large (8"H x 8"W, tapered to 6"W)	\$150
New!!!	<input type="checkbox"/> Calf Pad Covers (pair) (1/2" Soft Visco Foam)	
<input type="checkbox"/>	Air Mesh	\$120
<input type="checkbox"/>	Reversed Dartex	\$120
<input type="checkbox"/>	Calf Panel (ea.)	
<input type="checkbox"/>	Small (8"H x 18"W)	\$150
<input type="checkbox"/>	Large (8"H x 24"W) ²	\$200
New!!!	<input type="checkbox"/> Calf Panel Covers (1/2" Soft Visco Foam)	
<input type="checkbox"/>	Air Mesh	\$120
<input type="checkbox"/>	Reversed Dartex	\$120
Optional Modular Footplate Extension Components/Kits: ³		
<input type="checkbox"/>	4" Width Extension Kit	\$50
<input type="checkbox"/>	2" Bell Shaped, Length Extension Kit	\$50
<input type="checkbox"/>	4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75
<input type="checkbox"/>	6-Piece Extension Stock Kit (Includes all above pieces)	\$95
Additional Optional Accessories:		
<input type="checkbox"/>	Wheel Bumpers, pair	\$25
<input type="checkbox"/>	Footplate Connector	\$49
New!!!	<input type="checkbox"/> Footplate Cover, Reversed Dartex ⁴ (1/2" Soft Visco Foam)	\$120



1. N/A with contoured seating.
 2. Only available on 20"W or greater.
 3. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 4. Not available with modular footplate extension components/kits.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

SYNERGY® CUSHIONS

Option #	Description	Price	HCPCS Code
----------	-------------	-------	------------

Description	Price	HCPCS Code
-------------	-------	------------

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

Center Mount Power Articulating Foot Platform¹ **\$2895** **K0108**
(Includes 10"W x 11"D Split Footplates)

Extensions:

10-16"L

Set Legrest Length:

LL10 <input type="checkbox"/> 10"H	LL135 <input type="checkbox"/> 13.5"H	LL165 <input type="checkbox"/> 16.5"H
LL105 <input type="checkbox"/> 10.5"H	LL14 <input type="checkbox"/> 14"H	LL17 <input type="checkbox"/> 17"H
LL11 <input type="checkbox"/> 11"H	LL145 <input type="checkbox"/> 14.5"H	LL175 <input type="checkbox"/> 17.5"H
LL115 <input type="checkbox"/> 11.5"H	LL15 <input type="checkbox"/> 15"H	LL18 <input type="checkbox"/> 18"H
LL12 <input type="checkbox"/> 12"H	LL155 <input type="checkbox"/> 15.5"H	LL185 <input type="checkbox"/> 18.5"H
LL125 <input type="checkbox"/> 12.5"H	LL16 <input type="checkbox"/> 16"H	LL19 <input type="checkbox"/> 19"H
LL13 <input type="checkbox"/> 13"H		

Calf Supports:

Calf Pads (pair)

- Small (6"H x 6"W) No Charge
- Large (8"H x 8"W, tapered to 6"W) No Charge

New!!!

Calf Pad Covers (pair)
(1/2" Soft Visco Foam)

- Air Mesh \$120
- Reversed Dartex \$120

Calf Panel (ea.)

- Small (8"H x 18"W) \$150
- Large (8"H x 24"W)² \$200

New!!!

Calf Panel Covers
(1/2" Soft Visco Foam)

- Air Mesh \$120
- Reversed Dartex \$120


Optional Modular Footplate Extension Components/Kits:³

4" Width Extension Kit \$50

2" Bell Shaped, Length Extension Kit \$50

4" Width Extension with 2" Bell Shaped, Length Extension Kit \$75

6-Piece Extension Stock Kit \$95
(Includes all above pieces)



Additional Optional Accessories:

- Footplate Connector \$49
- Wheel Bumpers, pair \$25
- Footplate Cover, Reversed Dartex⁴ (1/2" Soft Visco Foam) \$120

New!!!

1. N/A with contoured seating.
 2. Only available on 20"W or greater.
 3. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 4. Not available with modular footplate extension components/kits.

- Matching Seat Size
- Non-Matching Seat Size

Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s

Simplicity Cushion: General Use

16" x 16" - 20" x 20" \$98 E2601

Solution® Cushion: Skin Protection & Positioning

10" x 10" - 20" x 20" \$380 E2607

Solution®1 Cushion: Skin Protection

14" x 14" - 20" x 20" \$380 E2603

Spectrum® Gel Cushion: Positioning

12" x 10" - 20" x 20" \$480 E2605

Spectrum® Visco Cushion: Skin Protection & Positioning

10" x 10" - 20" x 20" \$480 E2607

Spectrum® Air Cushion: Adjustable Skin Protection

14" x 12" - 20" x 20" \$500 E2622

Air Tube:

- Left Side
- Right Side

Spectrum® Air Contour Cushion: Adjustable Skin Protection & Positioning

14" x 12" - 20" x 20" \$499 E2624

Air Tube:

- Left Side
- Right Side

TRU-Comfort™ 2 4-Way Stretch: Skin Protection:

16" x 14" - 20" x 20" \$380 E2603

TRU-Comfort™ 2 Air Mesh: Skin Protection:

16" x 14" - 20" x 20" \$380 E2603

CONTOURED SEATING/SYNERGY ACCESSORIES

Option #	Description	Price	HCPCS Code
----------	-------------	-------	------------

Rear Basket¹ \$103

Cane/Crutch Holder¹ \$98 E2207

Weather Cover, Black \$158

Cup Holder \$17²

Side:

- Left
- Right

Saddle Bag \$28

Flag Holder \$17

Oxygen Holder

Contoured Seating¹ \$103 E2208

Synergy® Seat \$150 E2208

Walker Holder \$92

Quantum Back Pack³ \$60

Tablet Holder⁴ \$315

Deluxe Cell Phone/Media Holder⁴ \$140

New!!!

Transfer Bar, Pair⁴ \$175
(Quick Release)

1. Includes mounting hardware.
 2. No charge on Synergy seating
 3. Backpack may interfere w/ other rear accessories. Not available w/ contoured seating.
 4. Only available with Synergy seating.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

BATTERIES **ADDITIONAL FRONT RIGGING ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> U-1, 32 Amp Batteries	\$133.ea.	E2365
<input type="checkbox"/> No Batteries		
Installation:		
<input type="checkbox"/> Batteries Installed		
<input type="checkbox"/> Do Not Install		

The following items can be chosen in the Positioning Components screen in iQuantum. Please refer to the Positioning Components order form for additional details, sizes and pricing

LATERALS

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals¹		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away, Rail Mounts	\$450 pr.	E1028
<i>Pads:</i>		
<input type="checkbox"/> Curved Pad, 3" x 5"	\$200 pr.	E0956
<input type="checkbox"/> Curved Pad, 4" x 6"	\$200 pr.	E0956
<i>Lateral Installation:</i>		
<input type="checkbox"/> Install		
<input type="checkbox"/> Do Not Install		
<input type="checkbox"/> Stealth Pelvic/Thigh Guides²		
<i>Side:</i>		
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<i>Bracket:</i>		
<input type="checkbox"/> Flip Down	\$225 ea.	E1028
<input type="checkbox"/> Removable	\$250 ea.	E1028
<input type="checkbox"/> Flip Down w/ Full Surface Contact Feature	\$225 ea.	E1028
<input type="checkbox"/> Removable w/ 1" Offset	\$325 ea.	E1028
<i>Adapter:</i>		
<input type="checkbox"/> 2" Height Adapter Plate ³	\$25 ea.	K0108
<input type="checkbox"/> 2 1/4" Height Adapter Plate ⁴	\$25 ea.	K0108
<i>Pads:</i>		
<input type="checkbox"/> Short Pad, 4" x 4"	\$50 ea.	E0956
<input type="checkbox"/> Medium Pad, 4" x 8"	\$50 ea.	E0956
<input type="checkbox"/> Long Pad, 4" x 12"	\$50 ea.	E0956
<i>Lateral Installation:</i>		
<input type="checkbox"/> Install		
<input type="checkbox"/> Do Not Install		

New!!!

1. Swing-away laterals may interfere with TB3 armrests. Not available with back cane option.
 2. Pelvic/Thigh guides may interfere with the joystick receiver block.
 3. Height adapter is needed w/ Tru Comfort seat cushions. Not available/necessary with the flip Down bracket.
 4. Only available w/ the flip down bracket

Description	Price	HCPCS Code
<input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket	\$172	E1028
<i>Pads:</i>		
<input type="checkbox"/> Small (3.5" x 3.25" x 4")	\$135	E0957
<input type="checkbox"/> Small, Gel (3.5" x 3.25" x 4")	\$145	E0957
<input type="checkbox"/> Medium (4" x 3.25" x 6")	\$135	E0957
<input type="checkbox"/> Large (4.5" x 2.75" x 6")	\$135	E0957
<input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair)	\$512	K0108
<input type="checkbox"/> POSAlinc Rail Mounted Knee Adductors (pair)	\$574	K0108
<input type="checkbox"/> Stealth Foot Boxes (pair)	\$475	K0108
<i>Small</i>		
<input type="checkbox"/> Medium		
<input type="checkbox"/> Large		
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)	\$166	K0108
<i>Small</i>		
<input type="checkbox"/> Medium		
<input type="checkbox"/> X-Large		
<input type="checkbox"/> Ankle Huggers	\$156	K0108
<i>Medium</i>		
<input type="checkbox"/> Large		
<input type="checkbox"/> Bodypoint High Mount, Clamp-on Footplates		
<i>Angle-Adjustable, pair</i>		
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<i>Multi-Axis, pair</i>		
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	

Items Designated with a Circular Checkbox (○) are not available through Quantum West

ADDITIONAL ARMREST ACCESSORIES

SET-TO-SPEC*

Description Price HCPCS Code

- Otto Bock Arm Channels (ea.)¹ \$104** E2209
 - Small (4.75" x 8")
 - Medium (4.75" x 11.5")
 - Large (4.75" x 14")
 - One Piece (4.75" x 18")

Side:

- Left
- Right

- Otto Bock Hand Pads (ea.)² \$52**

Side:

- Left
- Right

Size:

- Medium
- Large

- Stealth Elbow Support**

Bracket Side:

- Bracket-Left \$170 ea. E1028
- Bracket-Right \$170 ea. E1028

Pad Size:

- 3"W x 3"H \$50 ea. K0108
- 4"W x 4"H \$50 ea. K0108

1. A 2" narrower back width is recommended to avoid interference w/ the back cushion.
 2. Not available on same side as controller.

Please Supply Client Information:

Client Height: _____

Client Weight: _____

Mount Hip/Thigh Supports: Yes No

Mount Lateral Thoracic Supports: Yes No

Armrest Height from Pan: _____

Rigging Length (whole numbers only): _____

Seat to Floor Height (check one):

Seating System	Seat-To-Floor Heights	
Contoured Seating	20.5"	21.5"
Synergy Seat w/ Standard Towers	15.75"	16.75"
Synergy Seat w/ Extended Towers	16.5"	17.5"

*Seat to floor height are +/- 1/8" with 18"W x 18"D with solid tires (add 3/8" for pneumatic tires).
 Increased widths may increase lowest seat to floor heights.*

* While every effort will be taken to ensure dimensions will be preset to your specifications prior to shipping variations in options may effect final settings.