

Jazzy® 614 Series Group 2

Order Form / Price List
Revised: June 2011



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at www.mypride.com

Account #: _____ Date: _____ Order No. _____ Fax No. _____ Phone No. _____

Bill To: _____ Ship To: _____ P.O. #: _____

Address: _____ Marked For: _____ Phone Contact: _____ Discount: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.

All pricing is MSRP.

Terms: _____
Freight: _____
Taken By: _____

BASE MODEL

- Jazzy 614 2S-C K0823 \$5,363
- Jazzy 614 2S -SS K0822 \$5,363

COLORS

- Jazzy Red NC
- Jazzy Blue NC

BATTERIES (2 required)

- NF-22, 50 AGM (Installed) E2361 \$255
- NF-22, 50 Gel Cell (Installed) E2361 \$295

BATTERY CHARGER

- 5 amp Off-board Charger
- 50 amp controller¹ Std.

NON-EXPANDABLE CONTROLS

- 50 amp controller¹ Std.

JOYSTICK MOUNT

- Fixed, In-Line Joystick Mount Std.
 - Left Right
- Swing-away Joystick Mount¹⁶ E1028 \$210
(available thru SOF) 4" ¹⁸ 7"
 Left Right
- Quick Disconnect, Joystick Mount^{3,20} NC
 Left Right
- Multi-axis Joystick Mount K0108 \$100
 Left Right
(Fixed, does not swing-away)
- Vertical Adj. Clamp^{6,8} K0108 \$95

PROGRAMMER

- Dynamic Programmer \$423

DRIVE WHEELS

- 14" Pneumatic (pair) NC
- 14" Flat-Free (pair) \$213

JOYSTICK HANDLE OPTIONS

Dynamic Controls

- Standard Knob Std.
- Sponge Ball E2323 \$50
- DL Knob E2323 \$20
- Straight Bar (Chrome) E2323 \$45
- Finger Tube (Chrome) E2323 \$45
- Round T-Bar (Chrome) E2323 \$45
- Flat T-Bar (Chrome) E2323 \$45
- Wrist Plate E2323 \$45

SWING-AWAY LEGRESTS

- 70° Swing-away footrests Style #8¹⁹ \$162
- 70° Tapered S/A footrests Style #8¹⁹ \$162
(N/C with power positioning systems)
- Quad Release Lever NC
- S/A Standard Extensions: NC
 - 11"-13" Ext. NC
 - 13"-16.5" Ext. Std.
 - 15.5"-19.5" Ext. NC

- Extruded Footplate Std.
- 7.25" x 6" ¹⁰ (Tapered legs require 18"W or greater seat)
- 8.5 x 6" ⁹ (Tapered legs require 20"W or greater seat)

Angle Adjustable Footplates (pair) \$150

- 4"x6" ⁵ (Tapered legs require 16"W or greater seat)
- 5"x6" ¹⁰ (Tapered legs require 18"W or greater seat)
- 5"x8" ¹⁰ (Tapered legs require 18"W or greater seat)
- 6"x8" ⁹ (Tapered legs require 20"W or greater seat)
- 6"x11" ⁹ (Tapered legs require 20"W or greater seat)
- 8"x11" ¹¹ (Tapered legs require 24"W or greater seat)
- Wheel Bumpers \$25

Multi-Axis Ext. (ball & socket):

- 11"-16.5" Ext.
- 15.5"-19.5" Ext.
(Lengths over 17" may cause interference)
- Multi-Axis Footplates (pair) \$495
 - 4"x6" ⁵ (Tapered legs require 16"W or greater seat)
 - 5"x6" ¹⁰ (Tapered legs require 18"W or greater seat)
 - 5"x8" ¹⁰ (Tapered legs require 18"W or greater seat)
 - 6"x8" ⁹ (Tapered legs require 20"W or greater seat)
 - 6"x11" ⁹ (Tapered legs require 20"W or greater seat)
 - 8"x11" ¹¹ (Tapered legs require 24"W or greater seat)
 - Wheel Bumpers \$25

FRONT RIGGINGS

- Foot Platform¹⁷ Std.

Rigging Type (pair)

- Elevating Leg Rests²⁶ E0990/K0195 \$325
Select footplate size for above:
 - 71/4 x 6 81/2 x 6
 Adjustable Calf Pads w/ ELRS (1 hole std.)
 1 hole 3 holes* 6 holes*
**Potential for interference w/ shrouds*
- Extensions:**
 - 11" - 14" Ext.
 - 14" - 17" Ext.
 (Lengths over 17" may cause interference)

- HD Drop-Ins, 70^{07,8,13} \$420

- N/A w/ contoured or Specialty seating*
- Angle Extensions:** NC
 - 12" - 14.5" Ext.
 - 14.5" - 18" Ext.
 (Lengths over 17" may cause interference)

- 90° Pediatric Swing-Away footrests Style #8¹⁹ \$162
Extension Range 7" - 10"

- Footplates, Ang. Adj.: 4"x6" ⁵ 5"x6" ¹⁰

- Footplates (pair)**
 - Reinforced Footplate Std.
 - 71/4 x 6¹⁰ 81/2 x 6⁹
 - 5"x8" Angle-Adj. Footplate¹⁰ \$150
 - 6"x8" Angle-Adj. Footplate⁹ \$150
 - 6"x11" Angle-Adj. Footplate⁹ \$150
 - 8"x11" Angle-Adj. Footplate¹¹ \$150

- Wheel Bumpers (Ang-adj. only) \$150

Extensions¹²

- 11"-14" Ext. NC
- 14"-17" Ext. Std.
- High mount Clamp-on; Ang. Adj. \$326
 - Small Large
- High mount Clamp-on; Multi-axis \$412
 - Small Large

- High mount footplatform¹⁵ K0108 \$450

- Mid-wheel 6 foot platform Lower Extensions
 - 10"-14" L
 - 14"-17" L

Rigging Options

- Heel Loops (pair)¹³ E0951 \$38
- Gel Neoprene Footrest K0108 \$65
Sleeve (pair)²
- ELR / ALR Adductor Buttons (pair) K0108 \$100
- POSAlinc Knee Adductor (pair) for S/A legrests K0108 \$464
- Residual limb support left E1020 \$300
- Residual limb support right E1020 \$300

Jazzy® 614 Series Group 2

Order Form / Price List
Revised: June 2011

SEATING

Med-Back Gray Vinyl Std.

Med-Back, Gray Vinyl
& Solid Seat Pan ²¹ NC
 18" W x 18" D

Ltd. Recline, Gray Vinyl, High-Back
w/ Headrest, **Depth Adjustable** \$275

18" W x 18" - 20" D \$275
 20" W x 20" - 18" D \$275

SOLID SEAT PAN SELECTION

Ltd. Recline, Gray Vinyl, High-Back
w/ Headrest & Solid Seat Pan ²¹ \$275
 18" W x 18" D
 18" W x 20" D
 20" W x 18" D
 20" W x 20" D

Custom seat dimensions are available.
Please contact the Custom Group for pricing and leadtime.

22" W AND ABOVE WILL NOT BE LTD RECLINE.

Special Order Seat

W x D

Seat sizes can not exceed 22 W x 22 D
\$450 MSRP

ARMRESTS

For Medium and High-Back Seats

Flip-up, removable Std.
 Flip-up, height-adjustable E0973 \$170
removable²⁴

Armpads

Desk Length, 10" NC
 Left NC Right NC
 Full Length, 14" NC
 Left Std. Right Std.

MANUAL RECLINER ^{22,25}

Manual Recline High Back E1226 \$600

Manual Recline High Back²¹ E1226 \$600
w/ Solid Seat Pan
(180 degrees of Recline)
(Dimensions 18"W x 18"D)

Gray Vinyl Beige Vinyl

Width: Depth:
 16" NC 16" NC
 18" Std. 18" Std.
 20" NC 20" NC

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS²³

Simplicity Cushion
16 x 16 - 20 x 20 E2601 \$98
22 x 16 - 24 x 24 E2602 \$98
 Solution 1 Cushion
14 x 14 - 20 x 18 E2603 \$380
20 x 20 - 20 x 24 E2603 \$430
22 x 16 - 24 x 24 E2604 \$430

Solution Cushion
10 x 10 - 20 x 18 E2607 \$380
20 x 20 - 20 x 20 E2607 \$430
22 x 16 - 22 x 20 E2608 \$430

Spectrum Cushion
 Gel
10 x 10 - 20 x 18 E2605 \$480
20 x 20 - 20 x 20 E2605 \$530
22 x 16 - 22 x 20 E2606 \$530

Visco
10 x 10 - 20 x 18 E2607 \$480
20 x 20 - 20 x 20 E2607 \$530
22 x 16 - 22 x 20 E2608 \$530

Spectrum Air
14 x 12 - 20 x 24 E2622 \$480
22 x 16 - 24 x 24 E2623 \$530

Spectrum Air Contour
14 x 12 - 20 x 24 E2624 \$499
22 x 16 - 24 x 24 E2625 \$549

No Cushion Required

Footnotes for pages 1 & 2:

- Std. w/ one programmable drive profile.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.
- Only avail. on 12" W or greater seats.
- Provides vertical joystick adjustment w/o adjusting armrest height.
- May cause interference with front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The over all length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Please note swing away footrests are standard when ordered with Power Positioning.
- Std. calf pad is not depth adjustable.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- Fixed bracket, does not elevate.
- N/A with waterfall armpads and dual post armrests.
- Please note swing away footrests are standard when ordered with Power Positioning.
- N/A w/ 2-post, desk length arms w/ a toggle.
- Includes black vinyl seat pan, back available in all listed colors.
- 18" x 18" std. Special measurements are not available.
- See Synergy Cushion order forms for detailed pricing, sizes & codes.
- N/A on manual recliner.
- Rear mount accessories are N/A.
- E0990 Code used for purchase in CBA. K0195 code used for capped rental.

Jazzy® 614 Series Group 2

Order Form / Price List
Revised: June 2011

SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) \$275

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

- 16" NC
- 18" Std.
- 20" NC

Specialty Seating Kit Depth

- 16" Std. 17" NC
- 18" NC 19" NC
- 20" NC

Specialty Seating Kit Back Cane Height

- 18" Std. 20" NC
- 22" NC

Back Selection

See Synergy Shape Back Order Form^{1,2,3}

Synergy Back up to 20"W¹⁸ \$275

See Positioning Components Order Form for additional backs

Back Height Solid Pan⁴

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

Sling Back

- 16" wide \$98
- 18" wide \$98
- 20" wide \$98

Sling Back Height

(Hgt must match cane hgt.)

- 18" height NC
- 20" height NC
- 22" height NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

2-Post, fixed, flip-back^{3,5} Std.

Flip-back, cane mounted^{3,5} \$95ea

Removable, adjustable, quick height^{6,19} E0973 \$275ea
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.

Removable, height adjustable¹⁹ E0973 \$150ea heavy duty armrests

- 8"H-12"H 10"H-14"H

Armrest Pads

Full Length, straight, 2"x14" NC

- Left Right

Desk Length, straight, 2"x10" NC

- Left Right

Full Length, Waterfall, 2"x14"⁷ NC

- Left Right

Desk Length, Waterfall, 2"x10"⁷ NC

- Left Right

HEAD REST⁸

Stealth E0955

- Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205

Stealth Comfort Plus E0955

- Large (14"W) Small (10" W) \$205

Removable Head Rest Bracket E1028 \$210

Stealth, swing-away (Flip-Down) Head Rest Bracket E1028 \$210

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications E0955

Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$230

Removable Head Rest Bracket E1028 \$235

Fixed, Offset Mounting Bracket NC

Synergy Seating Systems¹¹

Synergy SeatingTM Frame \$434
(On chair MSRP)

Synergy Manual Recline Pkg.^{12,13,14,15} \$700
Maximum weight 300 lbs. E1226

Synergy Manual Tilt Pkg.^{13,16} K0108 \$550
(Max Dimensions 14"W x 14"D)
Maximum weight 75 lbs.

Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 10" NC | <input type="checkbox"/> 11" NC |
| <input type="checkbox"/> 12" NC | <input type="checkbox"/> 13" NC |
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | |

Seat Width Package

- | | |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"W-14"W | <input type="checkbox"/> 14"W-18"W |
| <input type="checkbox"/> 16"W-20"W | <input type="checkbox"/> 20"W-22"W |

Seat Width Package Manual Tilt

- | | |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"W-12"W | <input type="checkbox"/> 12"W-14"W |
|------------------------------------|------------------------------------|

DEPTH

Set Seat Depth At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 10" NC | <input type="checkbox"/> 11" NC |
| <input type="checkbox"/> 12" NC | <input type="checkbox"/> 13" NC |
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | |

Seat Depth Package

- | | |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"D-14"D | <input type="checkbox"/> 14"D-18"D |
| <input type="checkbox"/> 16"D-20"D | <input type="checkbox"/> 20"D-22"D |

Seat Depth Package Manual Tilt

- 10"D-14"D

BACK OPTIONS

See Synergy Shape Back Order Form^{1,3}

Synergy Back up to 20" W \$275

Omit Back

Back Height:⁴

18"H 19"H 20"H 21"H

22"H 23"H 24"H 25"H

Sling Back² \$98

Back height will be the same as back cane height.

Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

Please consult your Quantum Specialist or contact the Custom Group before ordering

BACK CANES²²

Back Canes w/ push handles Std.

Back Canes w/out push handles \$125

Flip Forward²⁴ \$150

Offset(Meso) \$195 Inset(Endo) \$195

2" (1" per side) 4" (2" per side)

N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests .

N/A w/ any reclining systems.

Actual Back Cane Height is reflected in measurement

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

Footnotes for pages 5:

1. 14" & 22" width are N/A.
2. N/A w/ Synergy Manual Systems.
3. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
4. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
5. S/A laterals interfere w/ flip-back arms.
6. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
7. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
8. Available w/ solid backs only.
9. Comes with swing-away / removable hardware; HCPCS code is E1028.
10. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
11. Manual systems do not include a headrest.
12. N/A w/ flip-back, cane mounted arms.
13. Synergy Back not included.
14. Includes removable, adjustable, quick-height armrests.
15. N/A w/ Shape Back.
16. Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
17. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
18. 20"D is N/A.
19. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
21. N/A w/ manual or PES systems.
22. Height is measured from seat pan to top of cane.
23. N/A w/ TRU-Comfort.
24. N/A w/ TB power recline, Endo/Meso canes, Synergy Manual Systems, vent trays, bariatric armpads, occupied transit kit and TRU-Comfort or sling backs. May also interfere w/ lateral supports and

Jazzy® 614 Series Group 2

Order Form / Price List
Revised: June 2011

BACK OPTIONS

- See Synergy Shape Back Order Form¹
 Synergy Back up to 20" W \$275
 Omit Back

Back Height:²

- 18"H 19"H 20"H 21"H
 22"H 23"H 24"H 25"H
 Sling Back³ \$98

*Back height will be the same as back cane height.
Includes gap for cushion.*

*Push Bar selection must match back width; center
tube standard when selecting flip-up arms.*

TRU-Comfort Back¹⁴

(Even sizes only 16" W - 20"W)

- Std. New Plus
Back Height: 22" 25"
 Back up to 20" W **E2620** \$595
 22" W **E2621** \$695
 TRU-Comfort Cushion
16 x 16 - 20 x 18 **E2601** \$380
20 x 20 - 20 x 22 **E2601** \$380
22 x 16 - 22 x 24 **E2602** \$480
 TRU-Comfort Plus Cushion
16 x 16 - 20 x 18 **E2603** \$380
20 x 20 - 20 x 22 **E2603** \$380
22 x 16 - 22 x 24 **E2604** \$480
 Omit Back -(\$100)

See TRU-Comfort Back & Cushion
order forms for specific sizes

*TRU-Comfort Plus-utilizes a 4-way stretch
material cover on cushions & backs*

ARMRESTS

- Flip-back, cane mounted^{4,5,6} Std.
 2-Post, flip-back, height^{4,5} **E0973** \$150ea
adjustable, removable armrests
(9"-13" of height adjustment in 1" increments)
 Single Post height adj.^{7,12} **E0973** \$75ea
quick height,removable
11"-17" of height adj. 3" on top in 1/2"
increments & 3" on bottom in 1" increments.
 Quad Release Lever NC
 Removable, height adjustable **E0973** \$75ea
heavy duty armrests¹²
 8"H-12"H 10"H-14"H
 Pediatric removable¹² **E0973** \$75ea
adjustable, quick height.
 Match Shroud Color Black
 Quad Release Lever NC
(7 3/4" - 9 3/4" height adjustment in 1/2"
increments from seat pan to top of armrest)
 Cantilever Armrests²⁴ SOF
Angle-Adjustable: - 30° through + 30°

Armrest Pads

- Full Length, straight, 2"x14" NC
 Left Right
 Desk Length, straight, 2"x10" NC
 Left Right
 Full Length, Waterfall, 2"x14"¹⁷ NC
 Left Right
 Desk Length, Waterfall, 2"x10"^{17,18} NC
 Left Right

HEAD REST¹⁰

- Stealth **E0955**
 Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205
 Stealth Comfort Plus **E0955**
 Large (14"W) Small (10" W) \$205
 Removable Head Rest Bracket **E1028** \$210
 Stealth, swing-away (Flip-Down)
Head Rest Bracket **E1028** \$210
Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications **E0955**

- Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$230
 Removable Head Rest Bracket **E1028** \$235
 Fixed, Offset Mounting Bracket NC

TRU-Comfort **E0955**

- 8" W Leatherette \$205
 10" W Plus \$205
 Removable Head Rest Bracket **E1028** \$210
 Removable Bracket for
after market headrest **E1028** \$210
 No Headrest Required

ACCESSORIES

- Rear Basket^{20,21} \$103
 Cane/Crutch Holder^{20,21,22} **E2207** \$98
 Oxygen Holder^{20,21,22} **E2208** \$103
 Oxygen Holder^{18,19,20,21} **E2208** \$150
(Synergy Seat only)
 Stroller handle \$85
(Synergy Manual Tilt only)
 Walker Holder^{20,21,22} **K0108** \$92
 Flag Holder (SOFO Power Positioning) \$17
 Weather Cover, Medium, Black \$158
 Saddle bag(contoured seating only) \$28
 Cup holder \$17
 Left Right
N/A on desk length arm with S/A
mount or toggle box.
 Lap Belt Std.
 Pediatric 38" N/A w/ limited recline seating
 50"²¹ 60" 70" 80" 90"

Special Instructions

Footnotes for page 4:

- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- N/A w/ Synergy manual tilt or recline.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- S/A laterals interfere w/ flip-back arms.
- N/A w/ Synergy manual recline.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- N/A w/ S/A joystick mount.
- Available w/ solid backs only.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- N/A w/ flip-back, cane mounted arms.
- N/A w/ Shape Back.
- Includes price for removable, adjustable, quick-height arms & 2-post arms.
- Height is measured from seat pan to top of cane.
- Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
- N/A w/ flip-back arms.
- Includes mounting bracket and hardware.
- N/A on manual recliner high-back
- N/A on Synergy Manual Tilt.
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting after market arm troughs or Universal Swivel w/ desk length arms.

Lowest seat-to-floor height

Specialty Seat:	16.25"
Synergy Seat:	17.75"
High-back Seat:	21.25"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility. Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.