

Quantum® Rival

Group 3 Standard
Order Form
300 lbs. weight capacity



ACCOUNT INFO

Account Number: _____
Provider Name: _____
Contact Name: _____
Phone #: _____
Fax #: _____
E-mail Address: _____

Date: _____
PO#: _____
Marked For: _____
3 Column Pricing on Quote/Order: Yes No
Quote #: _____

Ship To: _____
Address: _____
City: _____ State: _____ Zip: _____
Lift Gate Needed: Yes No

Taken By: _____

SET-TO-SPEC*

Please Supply Client Information:

Client Height: _____
Client Weight: _____
Mount Hip/Thigh Supports: Yes No
Mount Lateral Thoracic Supports: Yes No

Armrest Height from Pan: _____
Rigging Length (whole numbers only): _____

Seat to Floor Height (check one):

Seating System	Seat-To-Floor Heights		
	20.75"	21.5"	22.25"
Contoured Seating	20.75"	21.5"	22.25"
Synergy Seat	16.5"	17.25"	18"
TRU-Balance®3 Seat	16.5"	17.3"	18.1"
TRU-Balance®3 Seat with 10" Lift	18.5"	19.3"	20.1"

*Seat to floor height are +/- 1/8" with 18"W x 18"D with solid tires (add 3/8" for pneumatic tires).
Increased widths may increase lowest seat to floor heights.*

* While every effort will be taken to ensure dimensions will be preset to your specifications prior to shipping variations in options may effect final settings.

NOTES:

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASE MODEL **NON-EXPANDABLE CONTROLS**

Description	Price	HCPCS Code
<input type="checkbox"/> Quantum Rival 3S-SS	\$6,608	K0848
<input type="checkbox"/> Quantum Rival 3S-C	\$6,608	K0849

COLORS

Metallic Colors:

<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Pearl Pink
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Emerald Green
<input type="checkbox"/> Black	<input type="checkbox"/> Sapphire Blue
<input type="checkbox"/> Sunburst Orange	<input type="checkbox"/> Orion Yellow
<input type="checkbox"/> Silver	<input type="checkbox"/> True Blue

BATTERY CHARGER

Description	Price
<input type="checkbox"/> 8 Amp Off-board Charger	Standard

MOTOR OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Community Use Package ¹ <u>Package Includes:</u> - 300 lbs. Weight Capacity - High Speed Motor Package - 4" Wide Flat Free, or Pneumatic, Knobby, Extended Traction Tires	\$2,125	
<input type="checkbox"/> 6 MPH Motor Package	Standard	
<input type="checkbox"/> 6 MPH Accu-Trac Motor Package	\$1,295	K0108
<input type="checkbox"/> 8 MPH Song Motor Package	\$1,574	

1. Increases overall base width by 1.5".

DRIVE WHEELS

Description	Price
<input type="checkbox"/> 14" Flat Free	Standard
<input type="checkbox"/> 4"W Flat Free, Knobby, Extended Traction Tires ¹	\$551
<input type="checkbox"/> 14" Pneumatic	No Charge
<input type="checkbox"/> 4"W Pneumatic, Knobby, Extended Traction Tires ¹	\$551

1. Increases overall base width by 1.5"

FRONT CASTER WHEELS

Description	Price
<i>*For Optimal ride quality, pneumatic casters are recommended*</i>	
<input type="checkbox"/> 9" Flat Free	\$105
<input type="checkbox"/> 9" Pneumatic	No Charge

SUSPENSION OPTIONS



Description	Price
<input type="checkbox"/> Standard Suspension	Standard
<input type="checkbox"/> Heavy Suspension ¹	No Charge

1. Required with a TB3 Lift or vent tray.

LED LIGHT PACKAGE







Description	Price
<input type="checkbox"/> Through Q-logic 2 EX or NE+ 8-Key w/ Lights	\$995
<input type="checkbox"/> Seat Mounted	
<input type="checkbox"/> Base Mounted ¹	

1. Only available w/ the Q-logic 2 EX joystick.

Description	Price
<input type="checkbox"/> 75 Amp NE 	Standard
<input type="checkbox"/> 75 Amp NE+ 6-Key (For up to 2 Actuators Thru Joystick) 	\$200 <i>Standard on 3SP-SS Package</i>
<input type="checkbox"/> 75 Amp NE+ 8-Key For Lights	\$400

EXPANDABLE CONTROLS

Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls

Description	Price
<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller 	\$695
<input type="checkbox"/> Expandable Controller Harness 	\$350
<input type="checkbox"/> Q-Logic 2 EX Joystick 	\$700
<input type="checkbox"/> Q-Logic 2 EX Joystick (For Lights) 	\$700
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port 	\$700
New!!! <input type="checkbox"/> SCIM Module ¹ Specialty Control Interface Module 	\$575
<i>Features:</i> - Allows use of joystick display for specialty controls - Up to 5 profiles to control drive, seat and auxiliary functions - Excludes sip-n-puff and single switch functions - A switch is required to operate the input device attached to the SCIM module - Includes Accessory bar	
<input type="checkbox"/> USB Mobile Device Charger ²	\$175

Armor Color Options:

<input type="checkbox"/> Black Soft Touch (Standard color option)	
<input type="checkbox"/> Matching	
<input type="checkbox"/> Non-Matching	

Colors: \$50

<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Pearl Pink
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Emerald Green
<input type="checkbox"/> Black	<input type="checkbox"/> Sapphire Blue
<input type="checkbox"/> Sunburst Orange	<input type="checkbox"/> Orion Yellow
<input type="checkbox"/> Silver	<input type="checkbox"/> True Blue

1. Only available with Q-logic 2.
2. Available with NE, NE+ and Q-Logic 2 EX controllers. Can only charge 1 device at a time. A multiplier harness is required with all Static and Tilt only seating systems. The selection of more than 3 modules required an additional multiplier.

JOYSTICK MOUNTING BRACKETS

CONTOURED SEATING

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed, Inline ¹	Standard	
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> Swing-Away	\$210	E1028
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
New!!! <input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit	\$250	E1028
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> Multi-Axis, Fixed	\$100	
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> No Joystick Mount		


1. Not available with Q-Logic 2 EX Stand Alone Joystick.

JOYSTICK HANDLES

Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$71	E2323
<input type="checkbox"/> T-Bar	\$107	E2323
<input type="checkbox"/> Chin Cup	\$71	E2324 ¹
<input type="checkbox"/> Extended I-Handle	\$71	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$50	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$69	E2323

1. Use HCPCS Code E2323 when not used in conjunction with a chin control.

PROGRAMMING

Description	Price
<input type="checkbox"/> Q-Logic EX/NE Programmer	\$545
	
<input type="checkbox"/> Q-Logic PC Programming Disk w/ Cable ¹	\$750

1. Works with Q-Logic 2 EX only.

Description	Price	HCPCS Code
Depth-Adjustable High Back: ¹		
<input type="checkbox"/> Black Vinyl, Depth-Adjustable High Back, Limited Recline	No Charge	
<input type="checkbox"/> 18"W x 18-20"D		
<input type="checkbox"/> 20"W x 18-20"D		
<input type="checkbox"/> 22"W x 20-22"D		
Limited Recline High Back:		
<input type="checkbox"/> Black Vinyl, High Back with Limited Recline	No Charge	
<input type="checkbox"/> 16"W x 16"D		<input type="checkbox"/> 18"W x 16"D
<input type="checkbox"/> 16"W x 18"D		<input type="checkbox"/> 20"W x 16"D
<input type="checkbox"/> 16"W x 20"D		
Limited Recline High Back w/ Solid Seat Pan:		
<input type="checkbox"/> Black Vinyl, High Back with Solid Seat pan, Limited Recline	No Charge	
<input type="checkbox"/> 16"W x 16"D		<input type="checkbox"/> 18"W x 20"D
<input type="checkbox"/> 16"W x 18"D		<input type="checkbox"/> 20"W x 16"D
<input type="checkbox"/> 16"W x 20"D		<input type="checkbox"/> 20"W x 18"D
<input type="checkbox"/> 18"W x 16"D		<input type="checkbox"/> 20"W x 20"D
<input type="checkbox"/> 18"W x 18"D		
Medium Back:		
<input type="checkbox"/> Black Vinyl, Medium Back	No Charge	
<input type="checkbox"/> 16"W x 16"D		<input type="checkbox"/> 18"W x 16"D
<input type="checkbox"/> 16"W x 18"D		<input type="checkbox"/> 18"W x 18"D
Medium Back w/ Solid Seat Pan:		
<input type="checkbox"/> Black Vinyl, Medium Back With Solid Seat Pan	No Charge	
<input type="checkbox"/> 16"W x 16"D		<input type="checkbox"/> 18"W x 16"D
<input type="checkbox"/> 16"W x 18"D		<input type="checkbox"/> 18"W x 18"D
Armrest Type:		
<input type="checkbox"/> Flip up, Removable	No Charge	
<input type="checkbox"/> Flip up, Height-Adjustable, Removable	\$85 ea.	E0973
Armpads:		
<i>Left:</i>		<i>Right:</i>
<input type="checkbox"/> Full Length		<input type="checkbox"/> Full Length
<input type="checkbox"/> Desk Length ²		<input type="checkbox"/> Desk Length ²

Special Order Seat

_____ "W x _____ "D

Seat size can not exceed 22"W x 22"D

****Contact your inside sales specialist for pricing****

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part
 2. N/A with depth-adjustable high back seat.

MANUAL RECLINER

SYNERGY® SEAT (CONT'D)

Description	Price	HCPCS Code
<input type="checkbox"/> Manual Recliner High Back With Headrest	\$600	E1226
<input type="checkbox"/> Manual Recliner High Back With Headrest and Solid Seat Pan	\$600	E1226
Color: <input type="checkbox"/> Gray Vinyl <input type="checkbox"/> Beige Vinyl		
Sizes <input type="checkbox"/> 16"W x 16"D <input type="checkbox"/> 18"W x 20"D		
<input type="checkbox"/> 16"W x 18"D <input type="checkbox"/> 20"W x 16"D		
<input type="checkbox"/> 16"W x 20"D <input type="checkbox"/> 20"W x 18"D		
<input type="checkbox"/> 18"W x 16"D <input type="checkbox"/> 20"W x 20"D		
<input type="checkbox"/> 18"W x 18"D		
Armrest Type:		
<input type="checkbox"/> Manual Recliner Armrests	No Charge	
Armpads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length

Description	Price	HCPCS Code
Back Cane Options:		
<input type="checkbox"/> Standard Back Canes	Standard	
<input type="checkbox"/> Flip Forward Back Canes ³	\$150	
Armrests:		
<input type="checkbox"/> Height Adjustable Cantilever Armrests ⁴	Standard	
<input type="checkbox"/> Bariatric, Removable, Height Adjustable ⁵	\$75 ea.	E0973
<input type="checkbox"/> 8"H-12"H <input type="checkbox"/> 10"H-14"H		
Armpads:		
<input type="checkbox"/> Straight, Full Length	Standard	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Falcon, Gel, Full Length	\$80.00 ea.	
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)	\$56.69 ea.	
<input type="checkbox"/> Flat, Bariatric ⁵	No Charge	

3. Not available w/ cantilever armrests or O2 holder. May interfere with lateral thoracic supports.
 4. Not available w/ flip forward back canes. May interfere with lateral thoracic supports.
 5. Bariatric armrests require bariatric armpads.

For additional Synergy® Seat Options, please see the Synergy® Seat Order Form

SYNERGY® SEAT

Description	Price	HCPCS Code
-------------	-------	------------

Choose Synergy® Seat options below along with additional items starting on page 6

Synergy® Seating Frame:

<input type="checkbox"/> Synergy Seating Frame	\$434	
Seat Width:		
<input type="checkbox"/> 16"W <input type="checkbox"/> 17"W <input type="checkbox"/> 18"W <input type="checkbox"/> 19"W <input type="checkbox"/> 20"W		
<input type="checkbox"/> 21"W <input type="checkbox"/> 22"W		
Seat Width Package:		
<input type="checkbox"/> 16"W-20"W <input type="checkbox"/> 20"W-24"W		
Seat Depth:		
<input type="checkbox"/> 16"D <input type="checkbox"/> 17"D <input type="checkbox"/> 18"D <input type="checkbox"/> 19"D <input type="checkbox"/> 20"D		
<input type="checkbox"/> 21"D <input type="checkbox"/> 22"D		
Seat Depth Package:		
<input type="checkbox"/> 15"D-20"D <input type="checkbox"/> 20"D-24"D		
Back Cane Height:		
<input type="checkbox"/> 18"H <input type="checkbox"/> 19"H <input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H		
<input type="checkbox"/> 24"H <input type="checkbox"/> 25"H		
Back Type:		

See Synergy® Shape Back & TRU-Comfort™ Back Order Forms for Detailed Pricing and Codes

<input type="checkbox"/> Synergy® Shape Back	Various	E2613/E2614
<input type="checkbox"/> Sling Back ¹	\$98	
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> TRU-Comfort™ Back (20"W or less)	\$595	E2620
<input type="checkbox"/> TRU-Comfort™ Back (22"W)	\$695	E2621
Width: <input type="checkbox"/> 16"W <input type="checkbox"/> 18"W <input type="checkbox"/> 20"W		
Height ² : <input type="checkbox"/> 20"H <input type="checkbox"/> 22"H <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H		
<input type="checkbox"/> 28"H		
Width: <input type="checkbox"/> 22"W		
Height ² : <input type="checkbox"/> 24"H <input type="checkbox"/> 28"H		

1. Overall back height will be the same height as the back canes. N/A w/ cantilever armrests.
 2. Overall back height measurement is measured from the seat pan to the top of the back, +/- 1/2".

TRU-BALANCE® 3 SEATING

Option #	Description	Price
TB3STATIC	<input type="checkbox"/> TRU-Balance®3 Static Seat	\$2295
TB3LIFT	<input type="checkbox"/> TRU-Balance®3 Lift ¹	\$2295

Must be picked in conjunction with the TRU-Balance® 3 Static Seat

1. 10" of scissor lift, 1/4 speed inhibit.

TRU-BALANCE® 3 MANUAL POSITIONING OPTION

Option #	Description	Price	HCPCS Code
TB3MRECL	<input type="checkbox"/> Attendant Manual Recline ¹	\$500	E1226

1. Available w/ tilt only. N/A w/ NE tilt through toggle or back cane option. Push handles are required.

TRU-BALANCE® 3 SEATING SPECIFICATIONS

Seat Width:		Seat Depth:	
12W	<input type="checkbox"/> 12"W	12D	<input type="checkbox"/> 12"D
13W	<input type="checkbox"/> 13"W	13D	<input type="checkbox"/> 13"D
14W	<input type="checkbox"/> 14"W	14D	<input type="checkbox"/> 14"D
15W	<input type="checkbox"/> 15"W	15D	<input type="checkbox"/> 15"D
16W	<input type="checkbox"/> 16"W	16D	<input type="checkbox"/> 16"D
17W	<input type="checkbox"/> 17"W	17D	<input type="checkbox"/> 17"D
18W	<input type="checkbox"/> 18"W	18D	<input type="checkbox"/> 18"D
19W	<input type="checkbox"/> 19"W	19D	<input type="checkbox"/> 19"D
20W	<input type="checkbox"/> 20"W	20D	<input type="checkbox"/> 20"D
21W	<input type="checkbox"/> 21"W	21D	<input type="checkbox"/> 21"D
22W	<input type="checkbox"/> 22"W	22D	<input type="checkbox"/> 22"D

Choose a width and depth package and combine both option #'s for the full option #

Seat Width Package:		Seat Depth Package:	
1216W	<input type="checkbox"/> 12"W-16"W	1216D	<input type="checkbox"/> 12"D-16"D
1620W	<input type="checkbox"/> 16"W-20"W	1620D	<input type="checkbox"/> 16"D-20"D
2024W	<input type="checkbox"/> 20"W-24"W	2024D	<input type="checkbox"/> 20"D-24"D

TRU-BALANCE® 3 SEAT TO FLOOR HEIGHT

Please refer to page 1 for seat to floor heights available on each setting system

- STF165 16.5" STF185 18.5"
- STF173 17.3" STF193 19.3"
- STF181 18.1" STF201 20.1"

TRU-BALANCE® 3 BACK TYPE

Description Price

- Static Back Standard
- Flip Forward Back¹ \$150
- Back Cane Option² (-\$125)
- Back Cane Type:**
 - Angle Adjustable Standard
 - Flip Forward \$35
 - 1" Endo \$25
 - 2" Endo \$25
- Back Cane Height:**
 - 16"H 21"H
 - 17"H 24"H
 - 20"H 25"H

New!!!

1. Only available w/ 2-post, flip back armrests. Not available on 12 or 14" widths. N/A w/ attendant manual recline.
 2. Only available with "Omit Back" and 2-post, flip back armrests.

TRU-BALANCE® 3 BACK OPTIONS

Option # Description Price HCPCS Code

- Sport Back¹ No Charge
- TRU-Comfort™ 2 4-Way Stretch Back^{1,2} (20"W or less) \$595 E2620
- TRU-Comfort™ 2 4-Way Stretch Back^{1,2} (22"W or greater) \$695 E2621
- LGWG** Large Lateral Positioning Wedge, pair^{1,3} \$30
- TRU-Comfort™ 2 Air Mesh Back^{1,2} (20"W or less) \$595 E2620
- TRU-Comfort™ 2 Air Mesh Back^{1,2} (22"W or greater) \$695 E2621
- LGWG** Large Lateral Positioning Wedge, pair^{1,2} \$30
- ABI** Aftermarket Back Interface⁴ \$200
Please Indicate which aftermarket back is being used
- NOBK** Omit Back (-\$100)
- Comfort Company Back Various

See Comfort Company Order Form

1. Not available with back cane option.
 2. Standard with bolster kit and small lateral positioning wedge.
 3. Not available on Back heights 16-17"H w/ limited adjustability on 18"H.
 4. Overall seat depth may need to be adjusted to accommodate certain backs. Contact your Quantum Specialist for details. Does not include brackets to mount aftermarket back. SOF is required on back width less than 16".

TRU-BALANCE® 3 BACK SPECIFICATIONS

Seat Width to Back Width Measurements:

- MSIZE** Matching
- NMSIZE** Non-Matching¹

Sport Back:

- Width: 12"W 14"W
- Height²: 14"H 15"H 16"H 17"H
- Width: 16"W 17"W 18"W 19"W 20"W
- Height²: 21"W 22"W 24"W
- Height²: 16"H 17"H 18"H 19"H 20"H
- 21"H 22"H 23"H 24"H 25"H

Back Width Package:

- 10"W"-14"W 20"W"-24"W
- 16"W"-20"W 24"W"-28"W

TRU-Comfort™ Back:

Please refer to TRU-Comfort Back order form for detailed pricing, sizes, codes and part #'s

- Width: 12"W³
- Height²: 14"H 15"H
- Width: 14"W
- Height²: 16"H 17"H
- Width: 16"W 17"W 19"W 22"W
- Height²: 16"H 17"H 18"H 19"H 20"H
- 21"H 22"H 23"H
- Width: 18"W 20"W
- Height²: 16"H 17"H 18"H 19"H 20"H
- 21"H 22"H 23"H 24"H 25"H
- Width: 22"W 24"W
- Height²: 16"H 17"H 18"H 19"H 20"H
- 21"H 22"H 23"H

Back Width Package:

- 10"W"-14"W 20"W"-24"W
- 16"W"-20"W 24"W"-28"W

1. +/- 4" of seat width.
 2. Odd heights = 4" gap, Even heights = 3" gap. Measurement listed is overall back height.
 3. Not available in Air Mesh cover.

ADDITIONAL TRU-BALANCE® 3 BACK OPTIONS

ARMRESTS (CONT'D)

Option #	Description	Price
Color Shroud Option:		
	<input type="checkbox"/> Standard Black	No Charge
MSH	<input type="checkbox"/> Matching	
NMSH	<input type="checkbox"/> Non-Matching	
Color Options:		
BBLACK	<input type="checkbox"/> Black	
BRED	<input type="checkbox"/> Candy Apple Red	
BVBLUE	<input type="checkbox"/> Viper Blue	
BSPINK	<input type="checkbox"/> Pearl Pink	
BSSBLUE	<input type="checkbox"/> Sapphire Blue	
BSEGREEN	<input type="checkbox"/> Emerald Green	
BSYELLOW	<input type="checkbox"/> Orion Yellow	
BORANGE	<input type="checkbox"/> Sunburst Orange	
BSSILVER	<input type="checkbox"/> Silver	
BSTBLUE	<input type="checkbox"/> True Blue	
Privacy Flap:		
PFLAP	<input type="checkbox"/> 9.75"W ¹	\$35
PFLAP	<input type="checkbox"/> 12.75"W	\$35

1. Recommended on back widths 14"W or less.

ARMRESTS

Option #	Description	Price	
<input type="checkbox"/> TB3 Armrests ¹		No Charge	
<input type="checkbox"/> Down Post		\$175	
<i>+/- up to 2" difference between seat and back width</i>			
Left Armrest:		Right Armrest:	
<i>Height Package:</i>		<i>Height Package:</i>	
ARM812	<input type="checkbox"/> 8"H-12"H	ARM812	<input type="checkbox"/> 8"H-12"H
ARM1216	<input type="checkbox"/> 12"-16"H	ARM1216	<input type="checkbox"/> 12"-16"H
<i>Set Armrest Height:</i>		<i>Set Armrest Height:</i>	
AH8	<input type="checkbox"/> 8"H ²	AH8	<input type="checkbox"/> 8"H ²
AH9	<input type="checkbox"/> 9"H ²	AH9	<input type="checkbox"/> 9"H ²
AH10	<input type="checkbox"/> 10"H	AH10	<input type="checkbox"/> 10"H
AH11	<input type="checkbox"/> 11"H	AH11	<input type="checkbox"/> 11"H
AH12	<input type="checkbox"/> 12"H	AH12	<input type="checkbox"/> 12"H
AH13	<input type="checkbox"/> 13"H	AH13	<input type="checkbox"/> 13"H
AH14	<input type="checkbox"/> 14"H	AH14	<input type="checkbox"/> 14"H
AH15	<input type="checkbox"/> 15"H	AH15	<input type="checkbox"/> 15"H
AH16	<input type="checkbox"/> 16"H	AH16	<input type="checkbox"/> 16"H
Left Armpad:		Right Armpad:	
<i>Choose an armpad length and style, then combine the Abbreviated option #'s for the full option #</i>			
<i>Length:</i>		<i>Length:</i>	
FL	<input type="checkbox"/> Full Length	FR	<input type="checkbox"/> Full Length
DL	<input type="checkbox"/> Desk Length ³	DR	<input type="checkbox"/> Desk Length ³
<i>Style:</i>		<i>Style:</i>	
s	<input type="checkbox"/> Straight	s	<input type="checkbox"/> Straight
w	<input type="checkbox"/> Waterfall	w	<input type="checkbox"/> Waterfall
TB3314G	<input type="checkbox"/> Falcon, Gel, Full Length	\$92.50 ea.	
TB3314G	<input type="checkbox"/> Falcon, Gel, Desk Length	\$70.00 ea.	

1. May interfere with swing away laterals. Not available with back cane option.
 2. Not available with the 28-32" width package.
 3. Not available on depths over 22"D.

Option #	Description	Price	HCPCS Code
<input type="checkbox"/> 2-Post, Flip Back, Height Adjustable, Removable ⁴		\$150 ea.	E0973
Left Armrest:		Right Armrest:	
<i>Set Armrest Height:</i>		<i>Set Armrest Height:</i>	
AH9	<input type="checkbox"/> 9"H	AH9	<input type="checkbox"/> 9"H
AH10	<input type="checkbox"/> 10"H	AH10	<input type="checkbox"/> 10"H
AH11	<input type="checkbox"/> 11"H	AH11	<input type="checkbox"/> 11"H
AH12	<input type="checkbox"/> 12"H	AH12	<input type="checkbox"/> 12"H
AH13	<input type="checkbox"/> 13"H	AH13	<input type="checkbox"/> 13"H
Left Armpad:		Right Armpad:	
<i>Choose an armpad style and length, then combine the Abbreviated option #'s for the full option #</i>			
<i>Length:</i>		<i>Length:</i>	
FL	<input type="checkbox"/> Full Length	FR	<input type="checkbox"/> Full Length
DL	<input type="checkbox"/> Desk Length	DR	<input type="checkbox"/> Desk Length
<i>Style:</i>		<i>Style:</i>	
s	<input type="checkbox"/> Straight	s	<input type="checkbox"/> Straight
w	<input type="checkbox"/> Waterfall	w	<input type="checkbox"/> Waterfall
TB3314G	<input type="checkbox"/> Falcon, Gel, Full Length	\$92.50 ea.	
	<input type="checkbox"/> Armpad Mounting Plate, Full Length ⁵	\$56.52 ea.	
TB3310G	<input type="checkbox"/> Falcon, Gel, Desk Length	\$70.00 ea.	
	<input type="checkbox"/> Armpad Mounting Plate, Desk Length ⁵	\$56.20 ea.	

4. Not available with a depth less than 16". Required with flip forward back and laptrays.
 5. Required only with 2-post, flip back, height adjustable, removable armrests.

HEADRESTS

Option #	Description	Price	HCPCS Code
Headrest Pads:			
	<input type="checkbox"/> 8"W TRU-Comfort™ Air Mesh ¹	\$205	E0955
8TCPAD	<input type="checkbox"/> 8"W TRU-Comfort™ ¹	\$205	E0955
	<input type="checkbox"/> 10" TRU-Comfort™ Air Mesh ¹	\$205	E0955
10TCPAD	<input type="checkbox"/> 10" TRU-Comfort™ ¹	\$205	E0955
CP250	<input type="checkbox"/> Stealth 10" Comfort Plus	\$142	E0955
CP150	<input type="checkbox"/> Stealth 14" Comfort Plus	\$142	E0955
STL:OTSM250	<input type="checkbox"/> Stealth Small Contoured (9.75"W x 3"H)	\$205	E0955
STL:OTLD150	<input type="checkbox"/> Stealth Large Contoured (11"W x 6"H)	\$205	E0955
	<input type="checkbox"/> Omit Headrest		
Headrest Mounts:			
TCBRKT	<input type="checkbox"/> TRU-Comfort™ Removable Headrest Bracket ²	\$210	E1028
HMO475P	<input type="checkbox"/> Stealth Removable Headrest Bracket	\$210	E1028
TWB485-1	<input type="checkbox"/> Stealth's Worlds Best Removable Headrest Bracket	\$210	E1028
NOHR	<input type="checkbox"/> No Bracket		

1. Measurement is taken from the inside of the headrest.
 2. Only available w/ TRU-Comfort™ headrest pads.

LAPBELTS **NON-POWERED LEGRESTS**

Option #	Description	Price
38BELT	<input type="checkbox"/> 38" Lap Belt	Standard
50BELT	<input type="checkbox"/> 50" Lap Belt	Standard
60BELT	<input type="checkbox"/> 60" Lap Belt	Standard
70BELT	<input type="checkbox"/> 70" Lap Belt	Standard
80BELT	<input type="checkbox"/> 80" Lap Belt	Standard
90BELT	<input type="checkbox"/> 90" Lap Belt	Standard
	<input type="checkbox"/> 60" Padded Push Button, Center Pull	\$115
	<input type="checkbox"/> Large 13-19", Hip Padded, Push Button	\$92

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

Option #	Description	Price	HCPCS Code
<input type="checkbox"/>	60-90° Swing-Away Legrests <i>(Adjustable in 5" Increments)</i>	\$162	
Rigging Selection:			
<input type="checkbox"/>	Standard Matching Legrests, pair		
<input type="checkbox"/>	Non-Matching Legrests, pair		
<input type="checkbox"/>	Individual Legrests (Left or Right Only)		
Footplate Selection:			
<input type="checkbox"/>	Angle-Adjustable Footplates ¹	\$150	K0040
<input type="checkbox"/>	Multi-Axis Footplates ¹	\$495	K0040
<input type="checkbox"/>	Standard, Extruded Footplates	No Charge	
Left Legrest:		Right Legrest:	
<i>Lower Extensions:</i>		<i>Lower Extensions:</i>	
<input type="checkbox"/>	11"H-13"L ²	<input type="checkbox"/>	11"L-13"L ²
<input type="checkbox"/>	11"-16.5"L	<input type="checkbox"/>	11"-16.5"L
<input type="checkbox"/>	13"H-16.5"L ²	<input type="checkbox"/>	13"L-16.5"L ²
<input type="checkbox"/>	15.5"-19.5"L	<input type="checkbox"/>	15.5"-19.5"L

TRU-BALANCE® 3 PREMIUM INSTALLED ACCESSORIES

Option #	Description	Price
<i>Available on TB3 Seating Only</i>		
TB3PUSH	<input type="checkbox"/> Push Handles (pair) ¹	\$175
	<input type="checkbox"/> Accessory Bar ²	\$100



1. Not available with O2 holder, back cane option or vent tray. Required with attendant manual recline.
2. Comes standard with SCIM Module, Trilogy vent tray and oxygen holder. May be added to mount additional accessories to back of unit. Not available with back cane option.

ELECTRONICS FOR POWER FOOT RIGGINGS

Option #	Description	Price	HCPCS Code
THRUJOY	<input type="checkbox"/> Single Actuator Through Joystick (AFP or Lift Only)	\$1175	E2310 ¹
THRUJOY	<input type="checkbox"/> Dual Actuator Through NE+ 6 Key ³	\$2375	E2311 ²
MTHRUJOY	<input type="checkbox"/> Multiple Actuator Through Q-Logic 2 EX	\$2375	E2311 ²
	<input type="checkbox"/> Through Q-Logic 2 EX Controller and Single Function with 4-Button Actuator <i>(Mounts to same side as controller)</i> ⁴	\$2375	E2310

New!!!

1. HCPCS used only with power AFP.
2. HCPCS used only with power legrests.
3. If a lift is chosen an AFP is available, not power ALR's.
4. Includes AAM and 4-button actuator modular. Only available with TB3 armrests.

<input type="checkbox"/>	Angle Adjustable/Multi-Axis Footplate Size		
<input type="checkbox"/>	4"W x 6"D	<input type="checkbox"/>	4"W x 6"D
<input type="checkbox"/>	5"W x 6"D	<input type="checkbox"/>	5"W x 6"D
<input type="checkbox"/>	5"W x 8"D	<input type="checkbox"/>	5"W x 8"D
<input type="checkbox"/>	6"W x 8"D ³	<input type="checkbox"/>	6"W x 8"D ³
<input type="checkbox"/>	6"W x 11"D ³	<input type="checkbox"/>	6"W x 11"D ³
<input type="checkbox"/>	8"W x 11"D ⁴	<input type="checkbox"/>	8"W x 11"D ⁴
<input type="checkbox"/>	Standard, Extruded Footplate Size		
<input type="checkbox"/>	7.25"W x 6"D ³	<input type="checkbox"/>	7.25"W x 6"D ³
<input type="checkbox"/>	8.5"W x 6"D ⁴	<input type="checkbox"/>	8.5"W x 6"D ⁴
<i>Accessories:</i>			
<input type="checkbox"/>	Heel Loops, pair ⁵	\$38	E0951
BUMP	<input type="checkbox"/> Wheel Bumpers, pair	\$25	
LPW	<input type="checkbox"/> Gel Neoprene Footrest Sleeve	\$100 pr.	K0108
	<input type="checkbox"/> Therafin Multi-Axis Amputee Support ⁶	\$230	E1020
	<input type="checkbox"/> Amputee Sppt. Gel Cover	\$150	K0108

60-90° Swing-Away Legrest Compatability Matrix

	Contoured Seat			Synergy Seat			TRU-Balance®3 Static Seat			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High	Low	Medium	High
Standard, 13" - 16.5"	!	!	!	!	!	!	!	!	!	!	!	!
Standard, 15.5" - 19.5"	X	X	!	!	!	!	X	!	!	!	!	!
Angle-Adjustable, 13" - 16.5"	!	!	!	!	!	!	!	!	!	!	!	!
Angle-Adjustable, 15.5" - 19.5"	X	X	!	!	!	!	X	!	!	!	!	!
Multi-Axis, 15.5" - 19.5"	X	X	!	!	!	!	X	!	!	!	!	!

1. Includes Heel Loops.
2. Not available with multi-axis footplates.
3. May cause interference on seat widths less than 18".
4. May cause interference on seat widths less than 20".
5. Standard on Angle-Adjustable and Multi-Axis footplates.
6. Swing-Away legrest upper assembly and receiver is required.

NON-POWERED LEGRESTS (CONT'D)

NON-POWERED LEGRESTS (CONT'D)

Option #	Description	Price	HCPCS Code
<input type="checkbox"/>	Elevating Legrests	\$325	E0990
	Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair		
<input type="checkbox"/>	Non-Matching Legrests, pair		
<input type="checkbox"/>	Individual Legrests (Left or Right Only)		
	Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ¹	\$150	K0040
<input type="checkbox"/>	Standard, Extruded Footplates	No Charge	
	Left Legrest:		
	<i>Lower Extensions:</i>		
<input type="checkbox"/>	11"-14" ^{L2}		
<input type="checkbox"/>	12"H-14.5"L		
<input type="checkbox"/>	14"-17" ^{L2}		
<input type="checkbox"/>	14.5"-18"L		
	Right Legrest:		
	<i>Lower Extensions:</i>		
<input type="checkbox"/>	11"-14" ^{L2}		
<input type="checkbox"/>	12"L-14.5"L		
<input type="checkbox"/>	14"-17" ^{L2}		
<input type="checkbox"/>	14.5"-18"L		
<input type="checkbox"/>	Angle Adjustable Footplate Size		
<input type="checkbox"/>	4"W x 6"D		
<input type="checkbox"/>	5"W x 6"D		
<input type="checkbox"/>	5"W x 8"D		
<input type="checkbox"/>	6"W x 8"D ³		
<input type="checkbox"/>	6"W x 11"D ³		
<input type="checkbox"/>	8"W x 11"D ⁴		
<input type="checkbox"/>	Standard, Extruded Footplate Size		
<input type="checkbox"/>	7.25"W x 6"D ³		
<input type="checkbox"/>	8.5"W x 6"D ⁴		
	<i>Accessories:</i>		
<input type="checkbox"/>	Heel Loops, pair ⁵	\$38	E0951
BUMP <input type="checkbox"/>	Wheel Bumpers, pair	\$25	
LPW <input type="checkbox"/>	Gel Padded Adductor	\$100 pr.	K0108

Option #	Description	Price	HCPCS Code
<input type="checkbox"/>	90° Pediatric Swing-Away Legrests	\$162	
	Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair		
<input type="checkbox"/>	Non-Matching Legrests, pair		
<input type="checkbox"/>	Individual Legrests (Left or Right Only)		
	Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ¹	\$150	K0040
	Left Legrest:		
	<i>Lower Extensions:</i>		
<input type="checkbox"/>	7"H-10"L		
	Right Legrest:		
	<i>Lower Extensions:</i>		
<input type="checkbox"/>	7"L-10"L		
<input type="checkbox"/>	Angle Adjustable Footplate Size		
<input type="checkbox"/>	4"W x 6"D		
<input type="checkbox"/>	5"W x 6"D		
<input type="checkbox"/>	4"W x 6"D		
<input type="checkbox"/>	5"W x 6"D		
	<i>Accessories:</i>		
BUMP <input type="checkbox"/>	Wheel Bumpers, pair	\$25	

Pediatric Swing-Away Compatibility Matrix	Contoured Seat			Synergy Seat			TRU-Balance® Static Seat			TRU-Balance® Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High	Low	Medium	High
	7"-10" Extension	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

Elevating Legrest Compatibility Matrix	Contoured Seat			Synergy Seat			TRU-Balance®3 Static Seat			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High	Low	Medium	High
	Standard, 11" - 14"	!	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Standard, 14" - 17"	✗	!	!	!	!	!	!	!	!	!	!	!
Angle-Adjustable, 12" - 14.5"	!	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle-Adjustable, 14.5" - 18"	✗	!	!	!	!	!	✗	!	!	✗	!	✓

1. Includes heel loops.
 2. Only available with standard, extruded footplates.
 3. May cause interference on seat widths less than 18".
 4. May cause interference on seat widths less than 20".
 5. Standard on Angle-Adjustable footplates.

1. Includes Heel Loops.

- This rigging configuration causes no interference issues at this seat-to-floor height.
- This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.
- This rigging configuration does not work with this seat-to-floor height.

POWERED LEGRESTS

ADJUSTABLE CENTER MOUNT FOOT PLATFORM

Option # Description Price HCPCS Code

Description Price

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

May require seating configuration adjustments to avoid front caster interference. All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

70° Power Articulating Swing-Away Legrests (4" Articulation Range) **\$2395 E1010¹**

Rigging Selection:

- Standard Matching Legrests, pair
- Non-Matching Legrests, pair
- Individual Legrests (Left or Right Only)

Operational Method:

- Independent
- Combined

Footplate Selection:

- Angle-Adjustable Footplates²

Left Legrest:

Lower Extensions:

- 11"H-13"L
- 14.5"H-18"L
- 17.5"H-21"L

Right Legrest:

Lower Extensions:

- 11"L-13"L
- 14.5"L-18"L
- 17.5"L-21"L

Angle Adjustable Footplate Size

- | | |
|--|--|
| <input type="checkbox"/> 4"W x 6"D | <input type="checkbox"/> 4"W x 6"D |
| <input type="checkbox"/> 5"W x 6"D | <input type="checkbox"/> 5"W x 6"D |
| <input type="checkbox"/> 5"W x 8"D | <input type="checkbox"/> 5"W x 8"D |
| <input type="checkbox"/> 6"W x 8"D ³ | <input type="checkbox"/> 6"W x 8"D ³ |
| <input type="checkbox"/> 6"W x 11"D ³ | <input type="checkbox"/> 6"W x 11"D ³ |
| <input type="checkbox"/> 8"W x 11"D ⁴ | <input type="checkbox"/> 8"W x 11"D ⁴ |

Accessories:

- Small Calf Pad (6"W-5.35"H) No Charge
- Large Calf Pad (8"W x 5.86"H) No Charge

KBEALR Gel Padded Knee Adductor \$100 pr. **K0108**

Power Elevating Legrest Compatibility Matrix	Synergy Seat			TRU-Balance® Static Seat			TRU-Balance® Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High
Angle Adjustable, 11" - 13"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5" - 18"	!	!	!	!	!	!	!	!	!
Angle Adjustable, 17.5" - 21"	✗	✗	✗	✗	✗	✗	!	✗	!

1. Use code K0108 for an individual legrest.
 2. Includes Heel Loops.
 3. May cause interference on seat widths less than 18".
 4. May cause interference on seat widths less than 20".

Adjustable Center Mount Foot Platform (Includes 10"W x 11"D Split Footplates) **\$250**

Extensions:

- CM812** 8-12"L
- CM1014** 10-14"L
- CM1016** 10-16"L
- CM1319** 13-19"L

Set Legrest Length:

- | | | |
|--|---|--|
| LL10 <input type="checkbox"/> 10"H | LL13 <input type="checkbox"/> 13"H | LL16 <input type="checkbox"/> 16"H |
| LL105 <input type="checkbox"/> 10.5"H | LL135 <input type="checkbox"/> 13.5"H | LL165 <input type="checkbox"/> 16.5"H |
| LL11 <input type="checkbox"/> 11"H | LL14 <input type="checkbox"/> 14"H | LL17 <input type="checkbox"/> 17"H |
| LL115 <input type="checkbox"/> 11.5"H | LL145 <input type="checkbox"/> 14.5"H | LL175 <input type="checkbox"/> 17.5"H |
| LL12 <input type="checkbox"/> 12"H | LL15 <input type="checkbox"/> 15"H | LL18 <input type="checkbox"/> 18.5"H |
| LL125 <input type="checkbox"/> 12.5"H | LL15.5 <input type="checkbox"/> 15.5"H | LL1 <input type="checkbox"/> 19"H |

Calf Supports:

- Calf Pads (pair)**
 - Small (6"H x 6"W) \$150
 - Large (8"H x 8"W, tapered to 6"W) \$150


- New!!!** **Calf Pad Covers (pair)** (1/2" Soft Visco Foam)
- Air Mesh \$120
 - Reversed Dartex \$120

- Calf Panel (ea.)**
 - Small (8"H x 18"W) \$150
 - Large (8"H x 24"W)¹ \$200


- New!!!** **Calf Panel Covers** (1/2" Soft Visco Foam)
- Air Mesh \$120
 - Reversed Dartex \$120

Optional Modular Footplate Extension Components/Kits:²


- 4" Width Extension Kit \$50



- 2" Bell Shaped, Length Extension Kit \$50



- 4" Width Extension with 2" Bell Shaped, Length Extension Kit \$75



- 6-Piece Extension Stock Kit \$95 (Includes all above pieces)

Additional Optional Accessories:

- Wheel Bumpers, pair \$25
- Footplate Connector \$49
- New!!!** Footplate Cover, Reversed Dartex³ (1/2" Soft Visco Foam) \$120

1. Only available on 20"W or greater.
 2. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 3. Not available with modular footplate extension components/kits.

ADJUSTABLE CENTER MOUNT FOOT PLATFORM (CONT'D) CENTER MOUNT ARTICULATING FOOT PLATFORM (CONT'D)

Adjustable Center Mount Foot Platform Compatibility Matrix

	Synergy Seat			TRU-Balance®3 Static Seat			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High
10"-16" Extension	!	✓	✓	!	!	✓	✓	✓	✓
13"-19" Extension	!	!	!	!	!	!	!	!	!

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

Description Price HCPCS Code

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

Important Provider Note: The addition of the Power Center Mount Articulating Foot Platform to the Quantum Rival can effect the performance of the power chair due to added weight to the front end. This option will require Rehab Professionals to make the proper programming and seating configuration adjustments to meet the specific client needs, avoid front caster interference and enhance performance.

Center Mount Power Articulating Foot Platform **\$2895 K0108** (Includes 10"W x 11"D Split Footplates)

Extensions:

- 10-16"L
- 13-19"L

Set Legrest Length:

LL10 <input type="checkbox"/> 10"H	LL135 <input type="checkbox"/> 13.5"H	LL165 <input type="checkbox"/> 16.5"H
LL105 <input type="checkbox"/> 10.5"H	LL14 <input type="checkbox"/> 14"H	LL17 <input type="checkbox"/> 17"H
LL11 <input type="checkbox"/> 11"H	LL145 <input type="checkbox"/> 14.5"H	LL175 <input type="checkbox"/> 17.5"H
LL115 <input type="checkbox"/> 11.5"H	LL15 <input type="checkbox"/> 15"H	LL18 <input type="checkbox"/> 18"H
LL12 <input type="checkbox"/> 12"H	LL155 <input type="checkbox"/> 15.5"H	LL185 <input type="checkbox"/> 18.5"H
LL125 <input type="checkbox"/> 12.5"H	LL16 <input type="checkbox"/> 16"H	LL19 <input type="checkbox"/> 19"H
LL13 <input type="checkbox"/> 13"H		

Calf Supports:

- Calf Pads (pair)**
 - Small (6"H x 6"W) No Charge
 - Large (8"H x 8"W, tapered to 6"W) No Charge

- New!!!** **Calf Pad Covers (pair)** (1/2" Soft Visco Foam)
- Air Mesh \$120
 - Reversed Dartex \$120

- Calf Panel (ea.)**
 - Small (8"H x 18"W) \$150
 - Large (8"H x 24"W)¹ \$200

- New!!!** **Calf Panel Covers** (1/2" Soft Visco Foam)
- Air Mesh \$120
 - Reversed Dartex \$120

1. Only available on 20"W or greater.

Description Price

Optional Modular Footplate Extension Components/Kits:²

4" Width Extension Kit \$50

2" Bell Shaped, Length Extension Kit \$50

4" Width Extension with 2" Bell Shaped, Length Extension Kit \$75

6-Piece Extension Stock Kit \$95 (Includes all above pieces)

Additional Optional Accessories:

Footplate Connector \$49

Wheel Bumpers, pair \$25

Footplate Cover, Reversed Dartex³ (1/2" Soft Visco Foam) \$120

New!!!

	Synergy Seat			TRU-Balance®3 Static Seat			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High
10"-16" Extension	!	✓	✓	!	!	✓	✓	✓	✓

2. Assembly required. May require seating configuration adjustments to avoid front caster interference.
3. Not available with modular footplate extension components/kits.

SYNERGY® CUSHIONS

Description Price HCPCS Code

- Matching Seat Size
- Non-Matching Seat Size

Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s

Simplicity Cushion: General Use

<input type="checkbox"/> 16" x 16" - 20" x 24"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602

Solution® Cushion: Skin Protection & Positioning

<input type="checkbox"/> 10" x 10" - 20" x 24"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608

Solution@1 Cushion: Skin Protection

<input type="checkbox"/> 14" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604

Spectrum® Gel Cushion: Positioning

<input type="checkbox"/> 12" x 10" - 20" x 24"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2606

SYNERGY® CUSHIONS (CONT'D)

TRU-BALANCE® 3 COMPLIMENTARY ACCESSORIES

Description	Price	HCPCS Code
Spectrum® Visco Cushion: Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2608
Spectrum® Air Cushion: Adjustable Skin Protection		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$500	E2622
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$550	E2623
Air Tube:		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
Spectrum® Air Contour Cushion: Adjustable Skin Protection & Positioning		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$549	E2625
Air Tube:		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
TRU-Comfort™ 2 4-Way Stretch: Skin Protection:		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™ 2 Air Mesh: Skin Protection:		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™ Accessory:		
<input type="checkbox"/> Sport Seat Cushion Cover ¹	\$100	
<input type="checkbox"/> 16"W x 16"D		
<input type="checkbox"/> 18"W x 18"D		
<input type="checkbox"/> 20"W x 20"D		

1. Available as an additional cover option for both TRU-Comfort seat cushion styles. Only available in sizes listed.






CONTOURED SEATING/SYNERGY ACCESSORIES

Option #	Description	Price	HCPCS Code
<input type="checkbox"/>	Rear Basket ¹	\$103	
<input type="checkbox"/>	Cane/Crutch Holder ¹	\$98	E2207
<input type="checkbox"/>	Cup Holder	\$17 ²	
Side:			
<input type="checkbox"/>	Left		
<input type="checkbox"/>	Right		
<input type="checkbox"/>	Saddle Bag	\$28	
<input type="checkbox"/>	Flag Holder	\$17	
<input type="checkbox"/>	Oxygen Holder		
<input type="checkbox"/>	Contoured Seating ¹	\$103	E2208
<input type="checkbox"/>	Synergy® Seat	\$150	E2208
<input type="checkbox"/>	Walker Holder	\$92	
<input type="checkbox"/>	Quantum Back Pack ³	\$60	
<input type="checkbox"/>	Tablet Holder ⁴	\$315	
<input type="checkbox"/>	Deluxe Cell Phone/Media Holder ⁴	\$140	
<input type="checkbox"/>	Transfer Bar, Pair ⁴ (Quick Release)	\$175	

New!!!

1. Includes mounting hardware.
 2. No charge on Synergy seating
 3. Backpack may interfere w/ other rear accessories. Not available w/ contoured seating.
 4. Only available with Synergy seating.








Choose two complimentary items in this section. Each additional item will be \$25 MSRP. Available on TB3 Seating Only

Option #	Description
<input type="checkbox"/> CUP	Cup Holder ¹ 
<input type="checkbox"/> PIH	Personal Item Hook ² 
<input type="checkbox"/> RKIT	Reflector Kit 
<input type="checkbox"/> CELL	Cell Phone/Media Holder ¹ 
<input type="checkbox"/> BPH New!!!	Backpack Holder, Pair ² 

1. Not available w/ 2-post, flip back armrests or L...
 2. Not available with back cane option.

TRU-BALANCE® 3 PREMIUM ACCESSORIES

Premium accessories will be shipped with the unit, not installed. Multiple accessories can be mounted along the various tracks, some interference may occur

Option #	Description	Price	HCPCS Code
<input type="checkbox"/> CGUAR	Clothing Guard, Pair 	\$175	
<input type="checkbox"/> GBOX	Glove Box 	\$100	
<input type="checkbox"/> H2O	Hydration System ¹ 	\$175	
<input type="checkbox"/> OXYGEN	Oxygen Holder ¹ 	\$150	E2208
<input type="checkbox"/> OXYGEN	Oxygen Holder (For use w/ back cane option) 	\$150	E2208
<input type="checkbox"/>	Transfer Bar, Pair (Quick Release) 	\$175	
<input type="checkbox"/> QRTB	Transfer Bar, Pair (Bolt-On) 	\$175	
<input type="checkbox"/>	iPad Holder	\$375	

1. Not available with back cane option.

BATTERIES

Description	Price	HCPCS Code
<input type="checkbox"/> Group 24, 70 Amp Gel (Universal)	\$352 ea.	E2363
<input type="checkbox"/> Group 24, 70 Amp Gel (MK)	\$352 ea.	E2363
<input type="checkbox"/> No Batteries		

Installation:

Batteries Installed
 Do Not Install

The following items can be chosen in the Positioning Components screen in iQuantum
Please refer to the Positioning Components order form for additional details, sizes and pricing

LATERALS **ADDITIONAL FRONT RIGGING ACCESSORIES (CONT'D)**

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals¹		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away, Rail Mounts	\$450 pr.	E1028
<i>Pads:</i>		
<input type="checkbox"/> Curved Pad, 3" x 5"	\$200 pr.	E0956
<input type="checkbox"/> Curved Pad, 4" x 6"	\$200 pr.	E0956
<i>Lateral Installation:</i>		
<input type="checkbox"/> Install		
<input type="checkbox"/> Do Not Install		
<input type="checkbox"/> Stealth Pelvic/Thigh Guides²		
<i>Side:</i>		
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<i>Bracket:</i>		
<input type="checkbox"/> Flip Down	\$225 ea.	E1028
<input type="checkbox"/> Removable	\$250 ea.	E1028
<input type="checkbox"/> Flip Down w/ Full Surface Contact Feature	\$225 ea.	E1028
<input type="checkbox"/> Removable w/ 1" Offset	\$325 ea.	E1028
<i>Adapter:</i>		
<input type="checkbox"/> 2" Height Adapter Plate ³	\$25 ea.	K0108
<input type="checkbox"/> 2 1/4" Height Adapter Plate ⁴	\$25 ea.	K0108
<i>Pads:</i>		
<input type="checkbox"/> Short Pad, 4" x 4"	\$50 ea.	E0956
<input type="checkbox"/> Medium Pad, 4" x 8"	\$50 ea.	E0956
<input type="checkbox"/> Long Pad, 4" x 12"	\$50 ea.	E0956
<i>Lateral Installation:</i>		
<input type="checkbox"/> Install		
<input type="checkbox"/> Do Not Install		

New!!!

1. Swing-away laterals may interfere with TB3 armrests. Not available with back cane option.
2. Pelvic/Thigh guides may interfere with the joystick receiver block.
3. Height adapter is needed w/ Tru Comfort seat cushions. Not available/necessary with the flip Down bracket.
4. Only available w/ the flip down bracket

ADDITIONAL FRONT RIGGING ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket	\$172	E1028
<i>Pads:</i>		
<input type="checkbox"/> Small (3.5" x 3.25" x 4")	\$135	E0957
<input type="checkbox"/> Small, Gel (3.5" x 3.25" x 4")	\$145	E0957
<input type="checkbox"/> Medium (4" x 3.25" x 6")	\$135	E0957
<input type="checkbox"/> Large (4.5" x 2.75" x 6")	\$135	E0957
<input type="checkbox"/> POSAInc Knee Adductor for Swing-Away Legrests (pair)	\$512	K0108
<input type="checkbox"/> POSAInc Rail Mounted Knee Adductors (pair)	\$574	K0108
<input type="checkbox"/> Stealth Foot Boxes (pair)	\$475	K0108
<input type="checkbox"/> Small		
<input type="checkbox"/> Medium		
<input type="checkbox"/> Large		

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)	\$166	K0108
<input type="checkbox"/> Small		
<input type="checkbox"/> Medium		
<input type="checkbox"/> X-Large		
<input type="checkbox"/> Ankle Huggers	\$156	K0108
<input type="checkbox"/> Medium		
<input type="checkbox"/> Large		
<input type="checkbox"/> Bodypoint High Mount, Clamp-on Footplates		
<input type="checkbox"/> Angle-Adjustable, pair		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis, pair		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	

ADDITIONAL ARMREST ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Otto Bock Arm Channels (ea.)¹	\$104	E2209
<input type="checkbox"/> Small (4.75" x 8")		
<input type="checkbox"/> Medium (4.75" x 11.5")		
<input type="checkbox"/> Large (4.75" x 14")		
<input type="checkbox"/> One Piece (4.75" x 18")		
<i>Side:</i>		
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> Otto Bock Hand Pads (ea.)²	\$52	
<i>Side:</i>		
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Medium		
<input type="checkbox"/> Large		
<input type="checkbox"/> Stealth Elbow Support Bracket		
<i>Side:</i>		
<input type="checkbox"/> Bracket-Left	\$170 ea.	E1028
<input type="checkbox"/> Bracket-Right	\$170 ea.	E1028
<i>Pad Size:</i>		
<input type="checkbox"/> 3"W x 3"H	\$50 ea.	K0108
<input type="checkbox"/> 4"W x 4"H	\$50 ea.	K0108

1. A 2" narrower back width is recommended to avoid interference w/ the back cushion.
2. Not available on same side as controller.