

Tru Balance® Pediatric Power Tilt Order Form

Discount: _____
 Terms: _____
 Freight: _____
 Taken By: _____



Account Number: _____
 Provider Name: _____
 Contact: _____
 Phone: _____
 Fax: _____
 Date: _____
 PO#: _____
 Marked For: _____
 Ship To Address: _____
 City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASES: SYNERGY® ARMRESTS

250lb Weight Capacity. Max Dimensions 14"W x 18"D
 Quantum 610
 Quantum 6000Z
 Q6 Edge

Description	Price
<input type="checkbox"/> Flip-back, Cane Mounted ¹	No Charge
<input type="checkbox"/> Pediatric, Removable ² <input type="checkbox"/> Quad Release Lever	No Charge
<input type="checkbox"/> Height Adjustable Cantilever Armrests ³ Angle-Adjustable: -30° through +30°	\$195

1. Interferes with swing-away laterals. Not available on flip forward canes.
 2. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armpads not available
 3. Height adjustment from 11-17" depends on cane height. Not available with Flip Forward back canes or Shape Back.

TRU BALANCE® PEDIATRIC POWER TILT

Description	Price	HCPCS Code
<input type="checkbox"/> Tru Balance® Pediatric Power Tilt List base serial # for off chair seating	\$5925	E1002

S/N: _____

Seat Width:
 11" 12" 13" 14"

Seat Width Package:
 10"W-14"W

Seat Depth:
 10" 11" 12" 13" 14"
 15" 16" 17" 18"

Seat Depth Package:
 10"D-14"D 14"D-18"D

Back Cane Height:
 14" 15" 16"¹ 17"¹
 18" 19" 20"

1. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

SYNERGY® ARMPADS

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14") <input type="checkbox"/> Left <input type="checkbox"/> Right	Standard
<input type="checkbox"/> Straight, Desk Length (2"x10") <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Full Length (2"x14") ¹ <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Desk Length (2"x10") ¹ <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.

SYNERGY® BACK OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy® Shape Back <i>Positioning Back</i> <i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>	Various	E2613/E2614
<input type="checkbox"/> Omit Back	(-\$100)	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	No Charge	
<i>Overall Back</i> <input type="checkbox"/> 12" H <input type="checkbox"/> 13" H <input type="checkbox"/> 14" H <input type="checkbox"/> 15" H		
<i>Height:</i> <input type="checkbox"/> 16" H <input type="checkbox"/> 17" H <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H		
<i>Choose Lateral Extensions¹:</i>		
<input type="checkbox"/> 6" for 13"-16"H Back	\$25	
<input type="checkbox"/> 10" for 17"-18"H Back	\$25	
<input type="checkbox"/> 12" For 19"-22"H Back	\$25	
<input type="checkbox"/> 16" For 23"-25"H Back	\$25	

Synergy® Back Cane Options:

Standard Back Canes Standard
 Omit Push Handles \$125

Flip Forward Back Canes² \$150
 Omit Push Handles \$125

Offset (Meso) \$195
 2" (1" Per Side) 4" (2" Per Side)

Inset (Endo) \$195
 2" (1" Per Side) 4" (2" Per Side)

1. Lateral extensions are required to mount lateral thoracic supports.
 2. Not available with cane mounted or cantilever armrests. May interfere with Lateral Supports and Waterfall Armpads.

JOYSTICK RECEIVER MOUNT

Description	Price
<input type="checkbox"/> Universal Joystick Receiver Mount <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Omit Universal Joystick Receiver Mount	No Charge

ELECTRONICS FOR POWER OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Thru Toggle <input type="checkbox"/> Left Side <input type="checkbox"/> Right Side	Standard	
<input type="checkbox"/> Thru Joystick (Q-Logic 2 EX or Remote+ Only)	\$1,175	E2310
<input type="checkbox"/> Thru Single Pushbutton (VSI or Remote+ Only)	No Charge	
<input type="checkbox"/> Thru Dual Pushbutton (VSI or Remote+ Only)	No Charge	

Available With NE, VSI, Remote+ or Q-Logic 2 Electronics

LAP BELTS		
Description	Price	
<input type="checkbox"/> 38" Lap Belt ¹	Standard	
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width*	\$87	
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width	\$92	
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Med, 9"-15" Hip Width*	\$112	
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Lg, 13"-19" Hip Width*	\$116	

1. Includes mounting hardware. Not available with Manual Recliner.

HEADRESTS		
Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
TRU-Comfort™ Headrest¹:		
<input type="checkbox"/> 8" W Leatherette	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> 10" W Leatherette	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTS250)	\$205	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket ² (HMO475P)	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket (FDM350) Used in Addition to Removable Headrest Bracket	\$210	
Aftermarket Bracket:		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028


1. Measurement is taken from the inside of the headrest.

NON-POWER LEGRESTS		
Description	Price	HCPCS Code
<input type="checkbox"/> 90° Pediatric Swing-Away Legrests¹	Standard	
<input type="checkbox"/> Angle Adjustable Footplates		
Angle Adjustable Extension Length:		
<input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		
<input type="checkbox"/> High Mount, Clamp-On Footplates*		
<input type="checkbox"/> Angle Adjustable		K0040 ²
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040 ²
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	
<input type="checkbox"/> 60-90° Swing-away Legrests³	\$162	
(Adjustable in 5° increments)		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
Standard Extruded Extensions:		
<input type="checkbox"/> 13"-16.5"		
Standard Extruded Footplates:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates ⁴	\$150	K0040 ²
Angle Adjustable Extensions:		
<input type="checkbox"/> 13"-16.5"		
Angle-Adjustable Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> Multi-Axis Footplates	\$495	K0040 ²
Multi-Axis Extensions:		
<input type="checkbox"/> 15.5"-19.5"		
Multi-Axis Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		

1. Extensions are measured from the seatpan to footplate.
 2. K0040 is not separately reimbursable on Group 2 units.
 3. Possible caster wheel interference when adjusted to less than 70 degrees.
 4. Includes Heel Loops.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

ADJUSTABLE CENTER MOUNT FOOT PLATFORM		
Description	Price	
<input type="checkbox"/> Adjustable Center Mount Foot Platform ¹	\$250	
Includes 10"W x 11"D Split Foot Platform		
Extensions: ² <input type="checkbox"/> 8"-12" Ext. <input type="checkbox"/> 10"-14" Ext		
Calf Supports ³ :		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)	\$150	
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150	
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) ⁴	\$200	
Optional Modular Footplate Extension Components/Kits ^{3,5}		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<input type="checkbox"/> 6-Piece Extension Stock Kit (Includes all above pieces)	\$95	
Additional Optional Accessories ⁵ :		
<input type="checkbox"/> Footplate Connector	\$49	
<input type="checkbox"/> Wheel Bumpers, pair	\$25	



1. Not available on 614 or Select Series.
 2. Extension length is measured from top of the seat pan to footplate.
 3. Choose in Rigging Accessory Section in iQuantum.
 4. Only available on 20"W or greater.
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference.

FOOT RIGGING ACCESSORIES		
Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) ¹	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	\$65	K0108

1. Standard with angle adjustable footplates.
 2. Available only with swing-away legrests.

LATERALS		
Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals (pair)^{1*}		
Bracket:		
<input type="checkbox"/> Swing-Away Rail Mounted	\$450	E1028
<i>For Use on Synergy Back</i>		
Pad:		
<input type="checkbox"/> Flat Pad, 3"x5", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Flat Pad, 4"x6", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 2"x3", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Stealth Pelvic/Thigh Laterals*		
Side:		
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		
Bracket:		
<input type="checkbox"/> Removable Bracket (ea.)	\$250	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$325	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$225	E1028
Pad:		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$50	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$50	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$50	E0956
Adapter:		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) ²	\$25	K0108

1. Swing-away laterals interfere with flip-up armrests.
 2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.

SYNERGY® SEAT CUSHIONS

<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> No Cushion Required <i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Solution® Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 14" x 18"	\$380	E2607
Spectrum® Visco Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 14" x 18"	\$480	E2607

ACCESSORIES

<i>Description</i>	<i>Price</i>
<input type="checkbox"/> Quantum Backpack ¹	\$158
<input type="checkbox"/> Weather Cover, Medium, Black	\$158
<input type="checkbox"/> Cup Holder	No Charge
<i>Side:</i>	
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Saddle Bag	\$28
<input type="checkbox"/> Tablet Holder	\$315
<input type="checkbox"/> Deluxe Cell Phone/Media Holder ⁵	\$140

1. Backpack may interfere w/ other rear accessories.