

J6 Series

Group 2 Single Power
Order Form
250lb Weight Capacity



Pride Mobility Products Corporation
182 Susquehanna Ave., Exeter, PA 18643
Phone: (866) 800-2002 • Fax: (866) 707-3422
quantumorders@pridemobility.com



Discount: _____
Terms: _____
Freight: _____
Taken By: _____

Account Number: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
E-Mail: _____
Date: _____
PO#: _____
Marked For: _____
Ship To Address: _____
City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP. Items Designated with a Circular Checkbox (○) are not available through Quantum West

BASE MODEL	NON-EXPANDABLE CONTROLS
------------	-------------------------

Description	Price	HCPCS Code
-------------	-------	------------

<input type="checkbox"/> J6 2SP-SS	\$6,171	K0835
------------------------------------	---------	-------

Description	Price
-------------	-------

<input type="checkbox"/> 75 Amp NE	Standard
<input type="checkbox"/> 75 Amp NE+ 6-Key <small>(For up to 2 Actuators Thru Joystick)</small>	Standard
<input type="checkbox"/> 75 Amp NE+ 8-Key <small>(For Lights)</small>	\$400

COLORS

Candy Apple Red
 Viper Blue

BATTERY CHARGER

Description	Price
-------------	-------

<input type="checkbox"/> 5 Amp Off-board Charger	Standard
--	----------

EXPANDABLE CONTROLS

Description	Price
-------------	-------

Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls

<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller	\$695
<input type="checkbox"/> Expandable Controller Harness	\$425
<input type="checkbox"/> Q-Logic 2 EX Joystick	\$700
<input type="radio"/> Q-Logic 2 EX Joystick <small>(For Lights)</small>	\$700
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port	\$700

DRIVE WHEELS

Description	Price
-------------	-------

<input type="checkbox"/> 10" Pneumatic, pair	No Charge
<input type="checkbox"/> 10" Flat-Free, pair	Standard

BASE MOUNTED FOOT PLATFORM

Description	Price
-------------	-------

<input type="checkbox"/> Base Mounted Foot Platform	Standard
<input type="checkbox"/> No Foot Platform	No Charge
<input type="checkbox"/> Bell Shaped Foot Platform	\$235

Armor Color Options:

Black Soft Touch *(Standard color option)*
 Matching
 Non-Matching

Colors: \$50

Candy Apple Red
 Viper Blue

Items Designated with a Circular Checkbox (☐) are not available through Quantum West

JOYSTICK MOUNTING BRACKETS

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed, Inline ¹	Standard	
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> Swing-Away	\$245	E1028
SAMTIL7 <input type="checkbox"/> Inline, Left		
SAMTIR7 <input type="checkbox"/> Inline, Right		
SAMTOL7 <input type="checkbox"/> Offset, Left		
SAMTOR7 <input type="checkbox"/> Offset, Right		
<input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit	\$350	E1028
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> Multi-Axis	\$100	E1028
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<input type="checkbox"/> No Joystick Mount		

1. Not available with Q-Logic 2 EX Stand Alone Joystick.

TRU-BALANCE® 2 SEATING SYSTEM

Description	Price	HCPCS Code
<input type="checkbox"/> TRU-Balance® 2 Power Tilt (55°)	\$4,995	E1002
<i>55° Center of Gravity Power Tilt Seating System.</i>		

SEATING SPECIFICATIONS

Seat Width:	Seat Depth:
<input type="checkbox"/> 14"W	<input type="checkbox"/> 14"D
<input type="checkbox"/> 15"W	<input type="checkbox"/> 15"D
<input type="checkbox"/> 16"W	<input type="checkbox"/> 16"D
<input type="checkbox"/> 17"W	<input type="checkbox"/> 17"D
<input type="checkbox"/> 18"W	<input type="checkbox"/> 18"D
<input type="checkbox"/> 19"W	<input type="checkbox"/> 19"D
<input type="checkbox"/> 20"W	<input type="checkbox"/> 20"D
Seat Width Package:	Seat Depth Package:
<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 15"D-20"D

JOYSTICK HANDLES

Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$85	E2323
<input type="checkbox"/> T-Bar	\$120	E2323
<input type="checkbox"/> Chin Cup	\$85	E2324 ¹
<input type="checkbox"/> Extended I-Handle	\$85	E2323
<input type="checkbox"/> Soft Ball Handle	\$65	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$99	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$99	E2323
<input type="checkbox"/> Medium Ball Handle	\$65	E2323

1. Use HCPCS Code E2323 when not used in conjunction with a chin control.

BACK CANE OPTIONS

Description	Price	HCPCS Code
Back Cane Height¹:		
<input type="checkbox"/> 14"H ²	<input type="checkbox"/> 15"H ²	<input type="checkbox"/> 16"H ³
<input type="checkbox"/> 17"H ³	<input type="checkbox"/> 18"H	<input type="checkbox"/> 19"H
<input type="checkbox"/> 20"H	<input type="checkbox"/> 21"H	<input type="checkbox"/> 22"H
<input type="checkbox"/> 23"H	<input type="checkbox"/> 24"H	<input type="checkbox"/> 25"H
Cane Type¹:		
<input type="checkbox"/> Standard Back Canes	Standard	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="radio"/> Attendant Manual Recline ⁴	\$750	E1226
Cane Brackets:		
<input type="checkbox"/> Flip Forward Back Canes ⁵	\$150	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Offset (Meso) ⁶	\$195	
<input type="checkbox"/> 2" (1" per side)	<input type="checkbox"/> 4" (2" per side)	
<input type="checkbox"/> Inset (Endo) ⁶	\$195	
<input type="checkbox"/> 2" (1" per side)	<input type="checkbox"/> 4" (2" per side)	

1. No selection required with TRU-Comfort™ Back.
 2. Selection of a Synergy® Back and cantilever armrests disqualifies 14" and 15"H back canes.
 3. Selection of a Synergy® Back and cane mounted, flip up armrests disqualifies 16" and 17"H back canes.
 4. Not available with Flip Back Arms, Cantilever Arms, Shape Back, TRU-Comfort back, AFP or rear accessories. Requires 22"H canes. Pivot point 2.5" above seat pan.
 5. Not available with vent tray, oxygen holder, cane mounted or cantilever armrests. May interfere with lateral supports and waterfall armpads.
 6. Dual-Post, flip back armrests interfere with offset canes. Not available with attendant manual recline.

PROGRAMMING

Description	Price
<input type="checkbox"/> Q-Logic EX/NE Programmer	\$545
<input type="checkbox"/> Q-Logic PC Programming Disk w/ Cable ¹	\$750



1. Works with Q-Logic 2 EX only.

Items Designated with a Circular Checkbox (☐) are not available through Quantum West

BACK OPTIONS			ARMRESTS		
Description	Price	HCPCS Code	Option #	Description	Price
<input type="checkbox"/> Omit Back	No Charge		<input type="checkbox"/> Flip Back, Cane Mounted ¹		No Charge
<input type="checkbox"/> Comfort Company Back	Various		<input type="checkbox"/> 2-Post, Flip-Back, Height Adjustable, Removable ²		No Charge
See Comfort Company Order Form			<input type="checkbox"/> Single Post, Height Adjustable, No Charge Removable ³		No Charge
Synergy® Back:			<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	No Charge		<input type="checkbox"/> Bariatric, Removable, Height Adjustable ⁴		No Charge
Overall Height : <input type="checkbox"/> 15"H <input type="checkbox"/> 16"H <input type="checkbox"/> 17"H <input type="checkbox"/> 18"H <input type="checkbox"/> 19"H			Height: <input type="checkbox"/> 8"H-12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H <input type="checkbox"/> 24"H			<input type="checkbox"/> Removable, Single Post, Height Adjustable, Heavy Duty		No Charge
<input type="checkbox"/> 25"H			Height: <input type="checkbox"/> 8"H-12"H <input type="checkbox"/> 10"H-14"H		
Choose Lateral Extension ¹ :			<input type="checkbox"/> Pediatric, Removable ⁵ Height Adjustable		No Charge
<input type="checkbox"/> 6" for 13-16"H Back	\$26.25		<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> 10" for 17-18"H Back	\$26.25		<input type="checkbox"/> Height Adjustable, Cantilever ⁶		\$195 pair
<input type="checkbox"/> 12" for 19-22"H Back	\$26.25		<i>Angle Adjustable: -30° through +30°</i>		
<input type="checkbox"/> 16" for 23-25"H Back	\$26.25		1. Interferes with Swing-Away laterals. Not available w/ manual recline, sling back or flip forward back canes 2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16". 3. Height adjustable from 11" - 17" in 1/2" increments. 4. Bariatric armrests are full length and require bariatric armpads. 5. Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length armpads. 6. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, Shape Back or sling back.		
TRU-Comfort™ Back:			ARMPADS		
<i>Please refer to TRU-Comfort Back order form for detailed pricing, sizes, codes and part #'s</i>			Option #	Description	Price
<input type="checkbox"/> TRU-Comfort™	\$675	E2620	<input type="checkbox"/> Straight, Full Length (2"x14")	Standard	
4-Way Stretch Back (20"W or less)			<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> TRU-Comfort™	\$775	E2621	<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge	
4-Way Stretch Back (22"W or greater)			<input type="checkbox"/> Left <input type="checkbox"/> Right		
Width: <input type="checkbox"/> 16"W <input type="checkbox"/> 18"W <input type="checkbox"/> 20"W			<input type="checkbox"/> Waterfall, Full Length (2"x14")	No Charge	
Height ² : <input type="checkbox"/> 20"H <input type="checkbox"/> 21"H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H <input type="checkbox"/> 24"H			<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> 25"H <input type="checkbox"/> 27"H <input type="checkbox"/> 28"H			<input type="checkbox"/> Waterfall, Desk Length (2"x10")	No Charge	
Width: <input type="checkbox"/> 17"W <input type="checkbox"/> 19"W			<input type="checkbox"/> Left <input type="checkbox"/> Right		
Height: <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H <input type="checkbox"/> 27"H			<input type="checkbox"/> Flat, Bariatric ¹	No Charge	
<input type="checkbox"/> 28"H			<input type="checkbox"/> Falcon, Gel, Full Length, Pair (3.5"x14")	\$92.50 ea.	
Width: <input type="checkbox"/> 22"W <input type="checkbox"/> 24"W			<input type="checkbox"/> Gel Armpad Mounting Plate (Required)	\$56.69 ea.	
Height: <input type="checkbox"/> 24"H <input type="checkbox"/> 25"H <input type="checkbox"/> 27"H <input type="checkbox"/> 28"H			<input type="checkbox"/> Falcon, Gel, Desk Length, Pair (3.5"x10")	\$70.00 ea.	
			<input type="checkbox"/> Gel Armpad Mounting Plate (Required)	\$53.52 ea.	
1. Lateral extensions are required to mount lateral thoracic supports. 2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".			1. Bariatric armrests are full length and require bariatric armpads.		

Items Designated with a Circular Checkbox (○) are not available through Quantum West

POWER POSITIONING ELECTRONICS

Option #	Description	Price	HCPCS Code
Single Actuator Electronics:			
THRUTOG	<input type="checkbox"/> Through Toggle Box ¹ (Tilt Only)	No Charge	
TOGL	<input type="checkbox"/> Left Mount		
TOGR	<input type="checkbox"/> Right Mount		
THRUJOY	<input type="checkbox"/> Through NE+ 6 Key	\$1895	E2310
THRUJOY	<input type="checkbox"/> Through Single Switch ² (Tilt Only with Q-Logic 2 EX)	\$500	
THRUJOY	<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$1895	E2310
Multiple Actuator Electronics:			
MTHRUJOY	<input type="checkbox"/> Through NE+ 6-Key Controller	\$2895	E2310/E2311 ³
MTHRUJOY	<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$2895	E2311
MTHRUTOG	<input type="checkbox"/> Q-Logic 2 EX Through 4-way Toggle	\$2895	
TOGL	<input type="checkbox"/> Left Mount		
TOGR	<input type="checkbox"/> Right Mount		
	<input type="checkbox"/> Q-Logic 2 EX Through 4-way Button Toggle	\$2895	
TOGL	<input type="checkbox"/> Left Mount		
TOGR	<input type="checkbox"/> Right Mount		
	<input type="checkbox"/> Q-Logic 2 EX Through Quad Push Button	\$2895	
TOGL	<input type="checkbox"/> Left Mount		
TOGR	<input type="checkbox"/> Right Mount		

1. Manual Recline back is not available through NE tilt through toggle.
 2. Single switch not included. See specialty controls order form for a list of available switches.
 3. Applies to covered actuators only.

I-FIT LAPBELTS

Option #	Description	Price	HCPCS Code
<input type="checkbox"/>	i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded 1.5" Belt	\$86	
<input type="checkbox"/>	i-Fit 60" Single Pull, Push Button, Large (9.5") Padded 2" Belt	\$86	
<input type="checkbox"/>	i-Fit 60" Single Pull, Push Button, Extra Large (11.5") Padded, 2" Belt	\$86	

HEADRESTS

Option #	Description	Price	HCPCS Code
Headrest Pads:			
CP250	<input type="checkbox"/> Stealth 10" Comfort Plus	\$151.94	E0955
CP150	<input type="checkbox"/> Stealth 14" Comfort Plus	\$151.94	E0955
STL:OTS250	<input type="checkbox"/> Stealth Small Contoured (9.75"W x 3"H)	\$64.20	E0955
STL:OTLD150	<input type="checkbox"/> Stealth Large Contoured (11"W x 6"H)	\$80.25	E0955
	<input type="checkbox"/> Omit Headrest		
Headrest Mounts:			
HMO475P	<input type="checkbox"/> Stealth Removable Headrest Bracket	\$160.50	E1028
TWB485-1	<input type="checkbox"/> Stealth's Worlds Best Removable Headrest Bracket	\$224.70	E1028
NOHR	<input type="checkbox"/> No Bracket		

LED LIGHT PACKAGE

Description	Price
<input type="checkbox"/> Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995

LAPBELTS

Option #	Description	Price
50BELT	<input type="checkbox"/> 50" Lap Belt	Standard
60BELT	<input type="checkbox"/> 60" Lap Belt	Standard
70BELT	<input type="checkbox"/> 70" Lap Belt	Standard
80BELT	<input type="checkbox"/> 80" Lap Belt	Standard
90BELT	<input type="checkbox"/> 90" Lap Belt	Standard
	<input type="checkbox"/> 60" Padded Push Button, Center Pull	\$125
	<input type="checkbox"/> Large 13-19", Hip Padded, Push Button	\$92

Items Designated with a Circular Checkbox (⊖) are not available through Quantum West

NON-POWERED LEGRESTS

NON-POWERED LEGRESTS (CONT'D)

Option #	Description	Price	
<input type="checkbox"/>	60-90° Swing-Away Legrests <i>(Adjustable in 5° Increments)</i>	Standard	
	Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair		
<input type="checkbox"/>	Non-Matching Legrests, pair		
<input type="checkbox"/>	Individual Legrests (Left or Right Only)		
	Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ¹	\$85 ea.	
<input type="checkbox"/>	Standard, Extruded Footplates	No Charge	
	Left Legrest:		Right Legrest:
	<i>Lower Extensions:</i>		<i>Lower Extensions:</i>
<input type="checkbox"/>	11"H-13"L		<input type="checkbox"/>
<input type="checkbox"/>	13"H-16.5"L		<input type="checkbox"/>
	Angle Adjustable/Multi-Axis Footplate Size		
<input type="checkbox"/>	4"W x 6"D		<input type="checkbox"/>
<input type="checkbox"/>	5"W x 6"D		<input type="checkbox"/>
<input type="checkbox"/>	5"W x 8"D		<input type="checkbox"/>
<input type="checkbox"/>	6"W x 8"D ²		<input type="checkbox"/>
<input type="checkbox"/>	6"W x 11"D ³		<input type="checkbox"/>
<input type="checkbox"/>	8"W x 11"D ³		<input type="checkbox"/>
	Standard, Extruded Footplate Size		
<input type="checkbox"/>	7.25"W x 6"D ³		<input type="checkbox"/>
<input type="checkbox"/>	8.5"W x 6"D ⁴		<input type="checkbox"/>
	<i>Accessories:</i>		
<input type="checkbox"/>	Heel Loops, pair ⁴	\$50	E0951
BUMP <input type="checkbox"/>	Wheel Bumpers, pair	\$25	
LPW <input type="checkbox"/>	Gel Neoprene Footrest Sleeve	\$80 pr.	K0108
<input type="checkbox"/>	Therafin Multi-Axis Amputee Support ⁵	\$230	E1020
<input type="checkbox"/>	Amputee Sppt. Gel Cover	\$155	K0108

1. Includes Heel Loops.
 2. May cause interference on seat widths less than 18".
 3. May cause interference on seat widths less than 20".
 4. Standard on Angle-Adjustable.
 5. Swing-Away legrest upper assembly and receiver is required.

Option #	Description	Price	HCPCS Code
<input type="checkbox"/>	Elevating Legrests	\$325	E0990
	Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair		
<input type="checkbox"/>	Non-Matching Legrests, pair		
<input type="checkbox"/>	Individual Legrests (Left or Right Only)		
	Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ⁶	\$85 ea.	
<input type="checkbox"/>	Standard, Extruded Footplates	No Charge	
	Left Legrest:		Right Legrest:
	<i>Lower Extensions:</i>		<i>Lower Extensions:</i>
<input type="checkbox"/>	11"-14"L ⁷		<input type="checkbox"/>
<input type="checkbox"/>	12"H-14.5"L		<input type="checkbox"/>
<input type="checkbox"/>	14"-17"L ⁸		<input type="checkbox"/>
<input type="checkbox"/>	14.5"-18"L		<input type="checkbox"/>
	Angle Adjustable Footplate Size		
<input type="checkbox"/>	4"W x 6"D		<input type="checkbox"/>
<input type="checkbox"/>	5"W x 6"D		<input type="checkbox"/>
<input type="checkbox"/>	5"W x 8"D		<input type="checkbox"/>
<input type="checkbox"/>	6"W x 8"D ¹		<input type="checkbox"/>
<input type="checkbox"/>	6"W x 11"D ¹		<input type="checkbox"/>
<input type="checkbox"/>	8"W x 11"D ²		<input type="checkbox"/>
	Standard, Extruded Footplate Size		
<input type="checkbox"/>	7.25"W x 6"D ³		<input type="checkbox"/>
<input type="checkbox"/>	8.5"W x 6"D ⁴		<input type="checkbox"/>
	<i>Accessories:</i>		
<input type="checkbox"/>	Heel Loops, pair ⁵	\$50	E0951
BUMP <input type="checkbox"/>	Wheel Bumpers, pair	\$25	
LPW <input type="checkbox"/>	Gel Padded Adductor	\$125 pr.	K0108
<input type="checkbox"/>	90° Pediatric Swing-Away Legrests	\$162	
	Rigging Selection:		
<input type="checkbox"/>	Standard Matching Legrests, pair		
<input type="checkbox"/>	Non-Matching Legrests, pair		
<input type="checkbox"/>	Individual Legrests (Left or Right Only)		
	Footplate Selection:		
<input type="checkbox"/>	Angle-Adjustable Footplates ⁷	\$85 ea.	
	Left Legrest:		Right Legrest:
	<i>Lower Extensions:</i>		<i>Lower Extensions:</i>
<input type="checkbox"/>	7"H-10"L		<input type="checkbox"/>
	Angle Adjustable Footplate Size		
<input type="checkbox"/>	4"W x 6"D		<input type="checkbox"/>
<input type="checkbox"/>	5"W x 6"D		<input type="checkbox"/>
	<i>Accessories:</i>		
BUMP <input type="checkbox"/>	Wheel Bumpers, pair	\$25	

6. Includes heel loops.
 7. Only available with standard, extruded footplates.

Items Designated with a Circular Checkbox (○) are not available through Quantum West

NON-POWERED LEGRESTS (CONT'D)


Option #	Description	Price
<input type="checkbox"/>	Heavy Duty Drop-In Legrests	\$420
<input type="checkbox"/>	70 Degree	
<input type="checkbox"/>	80 degree	
	Rigging Selection:	
<input type="checkbox"/>	Standard Matching Legrests, pair	
<input type="checkbox"/>	Non-Matching Legrests, pair	
<input type="checkbox"/>	Individual Legrests (Left or Right Only)	
	Footplate Selection:	
<input type="checkbox"/>	Angle-Adjustable Footplates ¹	\$85 ea.
	Left Legrest:	Right Legrest:
<i>Lower Extensions:</i>		<i>Lower Extensions:</i>
<input type="checkbox"/> 12"H-14.5"L		<input type="checkbox"/> 12"L-14.5"L
<input type="checkbox"/> 14.5"-18"L		<input type="checkbox"/> 14.5"-18"L
<input type="checkbox"/>	Angle Adjustable Footplate Size	
<input type="checkbox"/>	4"W x 6"D	<input type="checkbox"/> 4"W x 6"D
<input type="checkbox"/>	5"W x 6"D	<input type="checkbox"/> 5"W x 6"D
<input type="checkbox"/>	5"W x 8"D	<input type="checkbox"/> 5"W x 8"D
<input type="checkbox"/>	6"W x 8"D ²	<input type="checkbox"/> 6"W x 8"D ²
<input type="checkbox"/>	6"W x 11"D ²	<input type="checkbox"/> 6"W x 11"D ²
<input type="checkbox"/>	8"W x 11"D ³	<input type="checkbox"/> 8"W x 11"D ³
<i>Accessories:</i>		
BUMP	<input type="checkbox"/> Wheel Bumpers, pair	\$25

1. Includes Heel Loops.
 2. May cause interference on seat widths less than 18".
 3. May cause interference on seat widths less than 20".

ADJUSTABLE CENTER MOUNT FOOT PLATFORM

Option #	Description	Price
<i>All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.</i>		
<input type="checkbox"/>	Adjustable Center Mount Foot Platform	\$250

<i>Footplates:</i>		
<input type="checkbox"/>	10"W x 11"D Split Footplates	
<input type="checkbox"/>	10"W x 8"D Split Footplates	
<input type="checkbox"/>	10"W x 6"D Split Footplates	
<i>Extensions:</i>		
CM1016	<input type="checkbox"/> 10-16"L	
<i>Set Legrest Length:</i>		
LL10	<input type="checkbox"/> 10	LL13 <input type="checkbox"/> 13"H
LL105	<input type="checkbox"/> 10.5"H	LL135 <input type="checkbox"/> 13.5"H
LL11	<input type="checkbox"/> 11"H	LL14 <input type="checkbox"/> 14"H
LL115	<input type="checkbox"/> 11.5"H	LL145 <input type="checkbox"/> 14.5"H
LL12	<input type="checkbox"/> 12"H	LL15 <input type="checkbox"/> 15"H
LL125	<input type="checkbox"/> 12.5"H	LL15.5 <input type="checkbox"/> 15.5"H
<i>Calf Supports:</i>		
<input type="checkbox"/>	Calf Pads (pair)	
<input type="checkbox"/>	Small (6"H x 6"W)	\$170
<input type="checkbox"/>	Large (8"H x 8"W, tapered to 6"W)	\$170
<input type="checkbox"/>	Calf Pad Covers (pair) (1/2" Soft Visco Foam)	
<input type="checkbox"/>	Air Mesh	\$150
<input type="checkbox"/>	Reversed Dartex	\$150
<input type="checkbox"/>	Calf Panel (ea.)	
<input type="checkbox"/>	Small (8"H x 18"W)	\$150
<input type="checkbox"/>	Large (8"H x 24"W) ¹	\$200
<input type="checkbox"/>	Calf Panel Covers (1/2" Soft Visco Foam)	
<input type="checkbox"/>	Air Mesh	\$150
<input type="checkbox"/>	Reversed Dartex	\$150
<i>Optional Modular Footplate Extension Components/Kits:²</i>		
<input type="checkbox"/>	4" Width Extension Kit	\$50
<input type="checkbox"/>	2" Bell Shaped, Length Extension Kit	\$50
<input type="checkbox"/>	4" Side Guard ³	\$50
<input type="checkbox"/>	6" Side Guard ⁴	\$50
<i>Additional Optional Accessories:</i>		
<input type="checkbox"/>	Wheel Bumpers, pair	\$25
<input type="checkbox"/>	Footplate Connector	\$15



1. Only available on 20"W or greater.
 2. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 3. Only available with 10"W x 6"D split footplates
 4. Only available with 10"W x 8"D and 10"W x 11"D split footplates

Items Designated with a Circular Checkbox (○) are not available through Quantum West

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

SYNERGY® CUSHIONS

Option #	Description	Price	HCPCS Code
----------	-------------	-------	------------

Description	Price	HCPCS Code
-------------	-------	------------

All efforts will be made to accommodate legrest length selection but minor changes may need to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of the unit.

See cushion order form or details below

Center Mount Power Articulating Foot Platform **\$2995** See note¹

Footplates:

- 10"W x 11"D Split Footplates
- 10"W x 8"D Split Footplates
- 10"W x 6"D Split Footplates

Extensions:

CM1016 10-16"L

Set Legrest Length:

- | | |
|---------------------------------------|---------------------------------------|
| LL10 <input type="checkbox"/> 10"H | LL135 <input type="checkbox"/> 13.5"H |
| LL105 <input type="checkbox"/> 10.5"H | LL14 <input type="checkbox"/> 14"H |
| LL11 <input type="checkbox"/> 11"H | LL145 <input type="checkbox"/> 14.5"H |
| LL115 <input type="checkbox"/> 11.5"H | LL15 <input type="checkbox"/> 15"H |
| LL12 <input type="checkbox"/> 12"H | LL155 <input type="checkbox"/> 15.5"H |
| LL125 <input type="checkbox"/> 12.5"H | LL16 <input type="checkbox"/> 16"H |
| LL13 <input type="checkbox"/> 13"H | |

Calf Supports:

- Calf Pads (pair)
 - Small (6"H x 6"W) No Charge
 - Large (8"H x 8"W, tapered to 6"W) No Charge


- Calf Pad Covers (pair)
(1/2" Soft Visco Foam)
 - Air Mesh \$150
 - Reversed Dartex \$150

- Calf Panel (ea.)
 - Small (8"H x 18"W) \$150
 - Large (8"H x 24"W)² \$200

- Calf Panel Covers
(1/2" Soft Visco Foam)
 - Air Mesh \$150
 - Reversed Dartex \$150

Optional Modular Footplate Extension Components/Kits:³

- 4" Width Extension Kit \$50
- 2" Bell Shaped, Length Extension Kit \$50
- 4" Side Guard⁴ \$50
- 6" Side Guard⁵ \$50



Additional Optional Accessories:

- Footplate Connector \$15
- Wheel Bumpers, pair \$25

ADDITIONAL ACCESSORIES

Description	Price	HCPCS Code
-------------	-------	------------

- Oxygen Holder¹ \$200 E2208
- Weather Cover, Medium, Black \$168
- Cup Holder No Charge
 - Left
 - Right
- Tablet Holder \$315
- Deluxe Cell Phone/Media Holder \$140
- Transfer Bar, Pair \$200 E0705
(Quick Release)

1. Oxygen holder on Tilt with TRU-Comfort back requires SOF.

1. Please verify the appropriate billing code to be used for this item with all third-party payors.
 2. Only available on 20"W or greater.
 3. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 4. Only available on 10"W x 6"D split footplates
 5. Only available on 10"W x 8"D and 10"W x 11"D split footplates

Items Designated with a Circular Checkbox (○) are not available through Quantum West

BATTERIES **ADDITIONAL ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> U-1, 32 Amp Batteries	\$133 ea.	E2365
<input type="checkbox"/> No Batteries		
Installation:		
<input type="checkbox"/> Batteries Installed		
<input type="checkbox"/> Do Not Install		

The following items can be chosen in the Positioning Components screen in iQuantum. Please refer to the Positioning Components order form for additional details, sizes and pricing

LATERALS

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals¹		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away, Rail Mounts	\$472.50 pr.	E1028
<i>Pads:</i>		
<input type="checkbox"/> Curved Pad, 3" x 5"	\$210 pr.	E0956
<input type="checkbox"/> Curved Pad, 4" x 6"	\$210 pr.	E0956
<i>Lateral Installation:</i>		
<input type="checkbox"/> Install		
<input type="checkbox"/> Do Not Install		
<input type="checkbox"/> Stealth Pelvic/Thigh Guides²		
<i>Side:</i>		
<input type="checkbox"/> Left		
<input type="checkbox"/> Right		
<i>Bracket:</i>		
<input type="checkbox"/> Flip Down	\$262.50 ea.	E1028
<input type="checkbox"/> Removable	\$262.50 ea.	E1028
<input type="checkbox"/> Flip Down w/ Full Surface Contact Feature	\$236.25 ea.	E1028
<input type="checkbox"/> Removable w/ 1" Offset	\$400 ea.	E1028
<i>Adapter:</i>		
<input type="checkbox"/> 2" Height Adapter Plate ³	\$26.25 ea.	K0108
<input type="checkbox"/> 2 1/4" Height Adapter Plate ⁴	\$53.50 ea.	K0108
<i>Pads:</i>		
<input type="checkbox"/> Short Pad, 4" x 4"	\$75 ea.	E0956
<input type="checkbox"/> Medium Pad, 4" x 8"	\$75 ea.	E0956
<input type="checkbox"/> Long Pad, 4" x 12"	\$75 ea.	E0956
<i>Lateral Installation:</i>		
<input type="checkbox"/> Install		
<input type="checkbox"/> Do Not Install		

1. Swing-away laterals may interfere with flip back armrests.
 2. Pelvic/Thigh guides may interfere with the joystick receiver block.
 3. Height adapter is needed w/ Tru Comfort seat cushions. Not available/necessary with the flip Down bracket.
 4. Only available w/ the flip down brackets

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Foot Boxes (pair)	\$498.75	K0108
<input type="checkbox"/> Small		
<input type="checkbox"/> Medium		
<input type="checkbox"/> Large		
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)	\$266	K0108
<input type="checkbox"/> Small		
<input type="checkbox"/> Medium		
<input type="checkbox"/> X-Large		
<input type="checkbox"/> Stealth Elbow Support	\$183.75 ea.	E1028
<i>Side:</i>		
<input type="checkbox"/> Bracket - Left		
<input type="checkbox"/> Bracket - Right		
<i>Pad:</i>		
<input type="checkbox"/> 3" x 3" (ARMS-ES-3x3)	\$52.50 ea.	K0108
<input type="checkbox"/> 4" x 4" (ARMS-ES-4x4)	\$52.50 ea.	K0108
<input type="checkbox"/> Stealth i-Fit Dynamic Chest Straps	\$142	E0960
<input type="checkbox"/> Upper Classic, Medium (IFCS1103)		
<input type="checkbox"/> Upper Classic, Large (IFCS1104)		
<input type="checkbox"/> Upper Contour, Medium (IFCS2103)		
<input type="checkbox"/> Upper Contour, Large (IFCS2104)		
<input type="checkbox"/> Stealth i-Fit Chest Strap Belt Anchor	\$90	
<input type="checkbox"/> Medium (IFPS4113)		
<input type="checkbox"/> Large (IFPS4114)		

Items Designated with a Circular Checkbox (○) are not available through Quantum West

ADDITIONAL HEADREST ACCESSORIES

<i>Description</i>	<i>Price</i>	<i>HCPCS Code</i>
<input type="checkbox"/> Stealth Facial Laterals¹		
<input type="checkbox"/> Right Side, 4.5" x 2" Pad (CPS9-R-HW & 900-Q)	\$347.75	
<input type="checkbox"/> Left Side, 4.5" x 2" Pad (CPS9-L-HW & 900-Q)	\$347.75	
<input type="checkbox"/> Right and Left 4.5" x 2" Pad (CPS9-B-HW & 900-Q)	\$453.75	
<input type="checkbox"/> Egg Switch Mount (ESM)	\$53.50 ea.	
<input type="checkbox"/> Facial Pad, 6" x 2" (Q920-Q)	\$53.50 ea.	
<input type="checkbox"/> Figure Eight Bracket with 13" Rod and Egg Switch Mount (F8678-13-Q)	\$267.50	
<input type="checkbox"/> Swing-away Joystick Mount, 18" Rod (SUS18-Q)	\$250	

1. Only available with Stealth Comfort Series Headrests.

SET-TO-SPEC*

Please Supply Client Information:

Client Height: _____

Client Weight: _____

Mount Hip/Thigh Supports: Yes No

Mount Lateral Thoracic Supports: Yes No

Armrest Height from Pan: _____

Rigging Length (whole numbers only): _____

Seat to Floor Height (check one):

Seating System	Seat-To-Floor Heights
TRU-Balance® 2 Tilt	16.875"

Seat to floor height are +/- 1/8" with 18"W x 18"D with solid tires (add 3/8" for pneumatic tires). Increased widths may increase lowest seat to floor heights.

* While every effort will be taken to ensure dimensions will be preset to your specifications prior to shipping variations in options may effect final settings.