

PRODUCT SPECIFICATIONS SHEET**Model:** Pride Boom 250: Manual Base/Power Base**Manufacturer:** Pride Mobility Products Corporation**Internet Address:** www.pridemobility.com**Document Information:** INFINFB2739/Rev B/April 2012

Straight Arm Weight:	Manual Base: 58 lbs. (26 kg) Power Base: 98 lbs. (44.45 kg)
Offset Arm Weight:	Manual Base: 65 lbs. (29 kg) Power Base: 105 lbs. (47.6 kg)
Lifting Capacity:	250 lbs. (113 kg)
Fuse Rating:	15 amp
Circuit Breaker Rating:	Manual Base: 10 amp Power Base: 12 amp
Minimum Height:	30.5 in. (77 cm)
Maximum Height:	38 in. (97 cm)
Arm Length:¹	36 in. (91 cm)
Angle Adjustment:	0°, 5°, 10°, 12°, 15°
Strap Adjustment:	21 in. (53 cm), 22 in. (56 cm), 23 in. (58 cm), 24 in. (61 cm), 25 in. (64 cm), 26 in. (66 cm), 27 in. (69 cm), 28 in. (71 cm), 29 in. (74 cm), 30 in. (76 cm), 31 in. (79 cm), 32 in. (81 cm), 33 in. (84 cm), 34 in. (86 cm), 35 in. (89 cm)

¹ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

NOTE: All specifications subject to change without notice.