

PRODUCT SPECIFICATIONS SHEET

Model: Jazzy Select Series – Jazzy Select

Manufacturer: Pride Mobility Products Australia Pty. Ltd.

Internet Address: www.pridemobility.com

Document Information: INFINF2767/Rev A/February 2012

| | |
|------------------------------------|---|
| Class of Use: | A |
| Weight Capacity: | 136 kg (300 lbs.) |
| Maximum Speed: ¹ | Up to 6.44 km/h (4 mph) |
| Ground Clearance: ² | 0.794 cm (0.313 in.) at motor mount skid plate 3.334 cm (1.313 in.) at gusset 4.604 cm (1.813 in.) at rear frame skid plate |
| Turning Radius: ² | 55.25 cm (21.75 in.) |
| Overall Size: ^{2, 6} | Length: 107 cm (42.125 in.) with front riggings 91.44 cm (36 in.) without front riggings Width: 60.33 cm (23.75 in.) |
| Drive Wheels: | 25.4 cm (10 in.), solid |
| Caster Wheels: | 15.24 cm (6 in.), solid, rear-mounted |
| Anti-tip Wheels: | 12.7 cm (5 in.), solid, front-mounted |
| Suspension: | Active-Trac |
| Range: ^{1, 4} | Up to 27 km (17 miles) |
| Component Weights: ⁵ | Base: 46.5 kg (102.5 lbs.) Batteries: 9.75 kg (21.5 lbs.) each |
| Drivetrain: | Two motor, mid-wheel drive |
| Brakes: | “Intelligent Braking” electronic regenerative, disc park brake |
| Batteries: ³ | Two 12-volt, U-1 batteries |
| Battery Charger: | 3.5 amp, off-board |
| Maximum Allowable Incline: | 10.5% (6°) |
| Maximum Angle Capacity: | 17.6% (10°) |
| Maximum Obstacle Climbing Ability: | 6.68 cm (2.63 in.) |

¹ Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tyre type, and tyre condition. This specification is subject to a variance of +10%, -5%.

² Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.

³ AGM or Gel-Cell type required. See VI. “Batteries and Charging.”


⁴ Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 Series EN12184 standards. All specifications subject to change without notice.

*** I N F I N F B 2 7 6 7 ***



Jazzy Select Dimensions and Ground Clearance

PRODUCT SPECIFICATIONS SHEET

Model: Jazzy Select Series – Jazzy Select 6

Manufacturer: Pride Mobility Products Australia Pty. Ltd.

Internet Address: www.pridemobility.com

Document Information: INFNFB2767/Rev A/February 2012

| | |
|------------------------------------|---|
| Class of Use: | A |
| Weight Capacity: | 136 kg (300 lbs.) |
| Maximum Speed: ¹ | Up to 6.8 km/h (4.2 mph) |
| Ground Clearance: ² | 3.175 cm (1.25 in.) at caster brackets 4.445 cm (1.75 in.) at motors and at center frame |
| Turning Radius: ² | 54.6 cm (21.5 in.) |
| Overall Size: ^{2, 6} | Length: 106.7 cm (42 in.) with front riggings 90.2 cm (35.5 in.) without front riggings Width: 59.7 cm (23.5 in.) |
| Drive Wheels: | 25.4 cm (10 in.), solid |
| Caster Wheels: | 15.24 cm (6 in.), solid, front-mounted and rear-mounted |
| Suspension: | Active-Trac with Mid-Wheel 6 Technology |
| Range: ^{1, 4} | Up to 20.8 km (12.92 miles) |
| Component Weights: ⁵ | Base: 45.4 kg (100 lbs.) Batteries: 9.75 kg (21.5 lbs.) each |
| Drivetrain: | Two motor, mid-wheel-6 drive |
| Brakes: | “Intelligent Braking” electronic regenerative, disc park brake |
| Batteries: ³ | Two 12-volt, U-1 batteries |
| Battery Charger: | 3.5 amp, off-board |
| Maximum Allowable Incline: | 10.5% (6°) |
| Maximum Angle Capacity: | 17.6% (10°) |
| Maximum Obstacle Climbing Ability: | 6.5 cm (2.55 in.) |

¹ Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tyre type, and tyre condition. This specification is subject to a variance of +10%, -5%.

² Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.



³ AGM or Gel-Cell type required. See VI. “Batteries and Charging.”

⁴ Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 Series EN12184 standards. All specifications subject to change without notice.



Jazzy Select 6 Dimensions and Ground Clearance