



SAFETY

THANK YOU for choosing your new lift chair to assist you for years in ease and comfort. The product you have purchased combines state-of-the-art components with **safety**, comfort, and styling in mind. **Read and follow** all instructions, warnings, and notes in this manual and all other accompanying literature before attempting to operate this product for the first time. If there is any information in this manual which you do not understand, or if you require additional assistance for setup or operation, please contact your authorized provider. **Failure to follow the instructions, warnings, and notes in this manual and those located on your product can result in personal injury or product damage and will void the product warranty.**

PURCHASER'S AGREEMENT

By accepting delivery of this product, you promise that you will not change, alter, or modify this product or remove or render inoperable or unsafe any guards, shields, or other safety features of this product; fail, refuse, or neglect to install any retrofit kits from time to time provided to enhance or preserve the safe use of this product.

SHIPPING AND DELIVERY

Before using your lift chair, make sure your delivery is complete as some components may be individually packaged. If you do not receive a complete delivery, please contact your authorized provider immediately. Where damage has occurred during transport, either to the packaging or content, please contact the delivery company responsible.

LANGUAGE USAGE

Many of the English terms in the owner's manual are written in language and/or spelling common to the USA. Please refer to the chart below to clarify some of the terms if you use a different word variation.

USA	Variation		USA	Variation		USA	Variation
asphalt	tarmac		backward(s)	rearward(s)		center	centre
sidewalk	pavement		authorized	authorised		curb	kerb
path	footpath		counterclockwise	anticlockwise		provider	dealer; agent
path	bridleway		tire	tyre		trunk	boot
yard	grounds		cord	lead		wrench	spanner
turn signal	turn indicator		pocketbook	handbag		color	colour
meter	metre		elevator	lift		labor	labour

Pride Mobility Products Australia Pty. Ltd.
20-24 Apollo Drive
Hallam, Victoria 3803
Australia

NOTE: If you ever lose or misplace your copy of the manual or informational supplements, contact us and we will be glad to send you another copy immediately.



TWO-YEAR LIMITED WARRANTY

For two (2) years from the date of purchase, Pride will repair or replace at our option to the original purchaser any of the following parts found upon examination by an authorised representative of Pride to be defective in material and/or workmanship:

- Steel Frame
- Motor
- Hand Control
- Scissor Mechanisms
- Lift Mechanism

NOTE: Pride reserves the right to replace only the part of the steel lift mechanism that may be defective.

NOT COVERED UNDER WARRANTY

The following parts are classed as wear items, which may under normal wear and tear require replacing. These items are **not** therefore covered under warranty: all fabric. Warranty will also be refused if damage is deemed to have been caused through misuse or accident for which Pride Mobility Products Australia Pty. Ltd. cannot be deemed responsible.

NOTE: Pride Mobility Products Australia Pty. Ltd provides parts only under warranty. Your Pride Provider is responsible for labour and service. Please contact your Pride Provider for information about these services and for any applicable charges.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any foreseeable loss or damage. You are also entitled to have the goods repaired if the goods fail to be of acceptable quality and the failure does not amount to a major failure.



PRODUCT INSERT

LIFT CHAIR WEIGHT CAPACITIES

Your lift chair was developed for your convenience. See the below specification table for information on weight capacity. Refer to the owner's manual for any additional information concerning your lift chair.

MODEL	WEIGHT CAPACITY
C1, C11	113 kg (250 lbs.)
C6, 805, 806	136 kg (300 lbs.)
C15, 101	147 kg (325 lbs.)
C5, D30, SR525M, 550, 560, 660, T3	170 kg (375 lbs.)
358XL	227 kg (500 lbs.)



MANDATORY! Stay within the specified weight capacity of your lift chair. Exceeding the weight capacity voids your lift chair warranty.



PRODUCT INSERT

WALL MOUNT DIAGRAM

EXTERNAL TRANSFORMER OPTIONAL WALL MOUNT DIAGRAM

Cut out the diagram on the left and use it to mark an area of placement for the external transformer. If your external transformer has an oblong shape, simply hold the transformer on the wall where you intend to mount it, and either mark the holes or drill (as shown) to secure the transformer to the wall.

NOTE: Neither the manufacturer nor provider is responsible for damage resulting from improper installation.

