

### SAFETY

WELCOME to Pride Mobility Products (Pride). The product you have purchased combines state-of-the-art components with **safety**, comfort and styling in mind. We are confident the design features will provide you with the conveniences you expect during your daily activities. Understanding how to **safely** operate and care for this product should bring you years of trouble-free operation and service.

**Read** and **follow** all instructions, warnings and notes in this manual and all other accompanying literature before attempting to operate this product for the first time. In addition, your **safety** depends upon you, as well as your dealer, carer or healthcare professional in using good judgement.

If there is any information in this manual which you do not understand, or if you require additional assistance for setup or operation, please contact your authorised Pride Dealer. **Failure to follow the instructions, warnings and notes in this manual and those located on your Pride product can result in personal injury or product damage and will void Pride's product warranty.**

### PURCHASER'S AGREEMENT

By accepting delivery of this product, you promise that you will not change, alter or modify this product or remove or render inoperable or unsafe any guards, shields or other safety features of this product; fail, refuse or neglect to install any retrofit kits from time to time provided by Pride to enhance or preserve the safe use of this product.

Pride Mobility Products Ltd.  
32 Wedgwood Road  
Bicester, Oxfordshire OX26 4UL  
UK

***NOTE: If you ever lose or misplace your copy of this manual, contact us and we will be glad to send you a new one immediately.***



### TWO-YEAR LIMITED WARRANTY

For two (2) years from the date of purchase, Pride will repair or replace at our option to the original purchaser any of the following parts found upon examination by an authorised representative of Pride to be defective in material and/or workmanship:

- Steel Frame
- Motors
- Hand Control
- Scissor Mechanisms
- Steel Lift Mechanism

***NOTE: Pride reserves the right to replace only the part of the lift and recline mechanism or chair frame that may be defective.***

### NOT COVERED UNDER WARRANTY

The following parts are classed as wear items, which may under normal wear and tear require replacing. These items are **not** therefore covered under warranty: all fabric. Warranty will also be refused if damage is deemed to have been caused through misuse or accident for which Pride Mobility Products cannot be deemed responsible.

***NOTE: Pride Mobility Products provides parts only under warranty. Your Pride Dealer is responsible for labour and service. Please contact your Pride Dealer for information about these services and for any applicable changes.***



# PRODUCT INSERT

## WEIGHT CAPACITIES

Your lift chair was developed for your convenience. See the below specification table for information on weight capacity. Refer to the owner's manual for any additional information concerning your lift chair.

MODEL	WEIGHT CAPACITY
C1, C11	113 kg/17 stone 12 lbs. (250 lbs.)
C6	136 kg/21 stone 6 lbs. (300 lbs.)
LC101, LC107, C10KD, C15KD	147 kg/23 stone 3 lbs. (325 lbs.)
C5, C20KD, C30KD, D2KD, D30KD, 200KD, 300KD, 560KD, 570KD, 660KD, 670KD, L60, L70, 358M, 525M, T3	170 kg/26 stone 11 lbs. (375 lbs.)



**MANDATORY!** Stay within the specified weight capacity of your lift chair. Exceeding the weight capacity voids your lift chair warranty.



# PRODUCT INSERT

## EXTERNAL TRANSFORMER WALL MOUNT

### EXTERNAL TRANSFORMER OPTIONAL WALL MOUNT

Hold the transformer on the wall where you intend to mount it, and either mark the holes or drill to secure the transformer to the wall.

***NOTE: Neither the manufacturer nor provider is responsible for damage resulting from improper installation.***

