

# PRODUCT SPECIFICATIONS SHEET

## Model: Jazzy Select Series

### Document Information: INFINF B3151/Rev A/October 2018

Weight Capacity	300 lbs. (136 kg) 21 stone 6 lbs.
Maximum Speed <sup>1</sup>	Up to 4.3 mph (7 km/h)
Ground Clearance <sup>2</sup>	Jazzy Select: 0.75 in. (19 mm) at motor bracket Jazzy Select 6 and Jazzy Select 6 with Power Seat: 0.87 in. (22 mm) at motor bracket
Turning Radius <sup>2</sup>	Jazzy Select: 22.5 in. (571 mm) Jazzy Select 6 and Jazzy Select 6 with Power Seat: 21.5 in. (546 mm)
Overall Size <sup>2, 6</sup>	Jazzy Select: Length: 44 in. (1118 mm) Width: 23.74 in. (603 mm) Jazzy Select 6 and Jazzy Select 6 with Power Seat: Length: 43.25 in. (1099 mm) Width: 23.74 in. (603 mm)
Drive Wheels <sup>2</sup>	10 in. (254 mm), solid
Caster Wheels <sup>2</sup>	Jazzy Select: 5.75 in. (146 mm), solid, rear-mounted Jazzy Select 6 and Jazzy Select 6 with Power Seat: 5.75 in. (146 mm), solid, front- and rear-mounted
Anti-tip Wheels <sup>2</sup>	Jazzy Select only: 5 in. (129 mm), solid, front-mounted
Suspension	Jazzy Select: Active-Trac® Jazzy Select 6 and Jazzy Select 6 with Power Seat: Active-Trac® with Mid-Wheel 6 Technology
Range <sup>1, 4</sup>	Jazzy Select: Up to 19.4 miles (31.2 km) Jazzy Select 6 and Jazzy Select 6 with Power Seat: 13 miles (21 km)
Component Weights <sup>5, 2</sup>	Jazzy Select Base: 94 lbs. (42.58 kg) Jazzy Select 6 Base: 97.4 lbs. (44.18 kg) Jazzy Select 6 with Power Base: 114.4 lbs. (51.89 kg) Jazzy Select - Standard Seat (High Back) Weight: 51.4 lbs. (23.31 kg) Jazzy Select 6 and Jazzy Select 6 with Power Seat - Standard Seat Weight: 50.8 lbs. (23 kg) Batteries: 25.2 lbs. (11.42 kg) each
Drivetrain	Jazzy Select: Two motor, mid-wheel drive® Jazzy Select 6 and Jazzy Select 6 with Power Seat: Two motor, mid-wheel-6 drive®
Brakes	"Intelligent Braking" electronic regenerative, disc park brake
Batteries <sup>3</sup>	Two 12 volt, U-1 batteries
Battery Charger	3.5 amp, off-board
Freewheel Operational Force	Less than 60 N
Drive Control Operational Force	Less than 5 N
Maximum Rated Slope	6° (10.5%)
Maximum Obstacle Climbing Ability (Ascending and Descending)	1.7 in. (45 mm)
Measured Sound Level (A-weighted)	Less than 65 dB

<sup>1</sup> Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%.

<sup>2</sup> Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>3</sup> AGM or Gel-Cell type required. See "Batteries and Charging."


<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup> Battery weight may vary based on manufacturer.


<sup>6</sup> Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

**NOTE: This product conforms to all applicable ANSI/RESNA, ISO 7176 Series, AS/NZS 3695.2, and EN12184 standards. All specifications subject to change without notice.**





**Jazzy Select**



**Jazzy Select 6 and Jazzy Select 6 with Power Base**

