

NOTE: The following information replaces Appendix A in the Riser Recliner Series Owner's Manual that accompanied your Riser Recliner. It specifies the weight capacity of your NM415 Riser Recliner. You must read the Riser Recliner Series Owner's Manual and all supplemental information before use. Please retain this information for your records, and contact your authorised Dealer with any questions.

RISER RECLINER WEIGHT CAPACITY	
Model Number	Rated Weight Capacity
LC101, LC107, C10KD, C15KD	147 kg/23 stone (325 lbs.)
200KD, 300KD, 560KD, 570KD, 660KD, 670KD, C20KD, C30KD, CL500, CL600, CL800, D2KD, D30KD, L60, L70, NM415	170 kg/27 stone (375 lbs.)

RISER RECLINER MODEL NUMBERS
200KD, 300KD, 560KD, 570KD, 660KD, 670KD, C10KD, C15KD, C20KD, C30KD, CL500, CL600, CL800, D2KD, D30KD, L60, L70, LC101, LC107, NM415



MANDATORY! Stay within the specified weight capacity of your riser recliner. Exceeding the weight capacity voids your Pride Riser Recliner warranty.

