

## PRODUCT SPECIFICATIONS SHEET

**Model:** Commander 400

**Manufacturer:** Pride Mobility Products Corporation

**Internet Address:** www.pridemobility.com

**Document Information:** INFINF3146/Rev A/November 2012

<b>Overall Weight:</b>	109 lbs. (49.4 kg)
<b>Lifting Capacity:</b>	400 lbs. (181 kg)
<b>Fuse Rating:</b>	15 amp
<b>Circuit Breaker Rating:</b>	12 amp
<b>Minimum Height:</b>	31.75 in. (80.65 cm)
<b>Maximum Height:</b>	38.125 in. (96.84 cm)
<b>Arm Length:<sup>1</sup></b>	26 in. (66 cm) retracted 38 in. (96.5 cm) fully extended
<b>Angle Adjustment:</b>	0°-45°

<sup>1</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

**NOTE:** All specifications subject to change without notice.

## SPECIFICATIONS DE PRODUIT

**Modèle:** Commander 400

**Fabricant:** Pride Mobility Products Company

**Adresse Internet:** www.pridemobility.com

**Information sur le document:** INFINF3146/Rev A/November 2012

<b>Poids total:</b>	49.4 kg (109 lb)
<b>Capacité de levage:</b>	181 kg (400 lb)
<b>Calibre du fusible:</b>	15 amp
<b>Calibre du disjoncteur:</b>	12 amp
<b>Hauteur minimum:</b>	80.65 cm (31.75 po)
<b>Hauteur maximum:</b>	96.84 cm (38.125 po)
<b>Longueur du bras:<sup>1</sup></b>	66 cm (26 po) rétracté 96.5 cm (38 po) pleine extension
<b>Réglage de l'angle:</b>	0°-45°

<sup>1</sup> Considérant les tolérances de fabrication et les améliorations apportées au produit, cette spécifications est sujette à une variation de (+ ou -) 3%.

**NOTE:** Toutes les spécifications sont sujettes aux changements sans préavis.