

PRODUCT SPECIFICATIONS SHEET

Model: Commander 450

Manufacturer: Pride Mobility Products, Corp.

Internet Address: www.pridemobility.com

Document Information: INFINF3824/Rev A/July 2017

Overall Weight:	109 lbs. (49.4 kg)
Lifting Capacity:	450 lbs. (204 kg)
Fuse Rating:	15 amp
Circuit Breaker Rating:	20 amp
Minimum Height:	31.75 in. (80.65 cm)
Maximum Height at 0°:	38.125 in. (96.84 cm)
Arm Length: ¹	26 in. (66 cm) retracted 38 in. (96.5 cm) fully extended
Angle Adjustment:	0°-45°

¹ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

NOTE: All specifications subject to change without notice.

SPECIFICATIONS DE PRODUIT

Modèle: Commander 450

Fabricant: Pride Mobility Products Company

Adresse Internet: www.pridemobility.com

Information sur le document: INFINF3824/Rev A/July 2017

Poids total:	49.4 kg (109 lb)
Capacité de levage:	204 kg (450 lb)
Calibre du fusible:	15 amp
Calibre du disjoncteur:	20 amp
Hauteur minimum à 0°:	80.65 cm (31.75 po)
Hauteur maximum:	96.84 cm (38.125 po)
Longueur du bras: ¹	66 cm (26 po) rétracté 96.5 cm (38 po) pleine extension
Réglage de l'angle:	0°-45°

¹ Considérant les tolérances de fabrication et les améliorations apportées au produit, cette spécifications est sujette à une variation de (+ ou -) 3%.

NOTE: Toutes les spécifications sont sujettes aux changements sans préavis.

